

Convoid Trekker - WM0001374

Canadian Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	
<div style="text-align: center;"><u>Trekker Standard Features:</u></div> <div style="display: flex; justify-content: space-between;"> <ul style="list-style-type: none"> Trekker Folding Mobility Frame Accepts TR12 & TR14 Seating Systems Tilt-in-Space (-10° to +45°) 180° Reversible Seating Configuration 90° Fixed Front Rigging 7.5" x 2" Front and 11.5" Rear Solid Tires One Piece Height Adjustable Push Handle (Inside Mount) <ul style="list-style-type: none"> Silver Powder Coated Frame 90° Fixed Footplate Quick Release Rear Wheels Standard Recline 80° to 110° Depth Adjustable Seat Pan </div> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div> <u>Standard Models:</u> <input type="checkbox"/> TR Seating System - C904768 <input type="checkbox"/> TR Mobility Frame - C904661 </div> <div style="text-align: center;"> <u>Price</u> \$0 \$3,686 </div> <div style="text-align: right;"> All prices in CAD </div> </div>	

<u>Trekker Seating System</u>		<u>Price</u>
<input type="checkbox"/> Seat Width 12"		N/C
<input type="checkbox"/> Back Canes Only		N/C
<input type="checkbox"/> TR12 Height Adjustable Solid Back 9" - 17"		\$808
<input type="checkbox"/> Seat Width 14"		N/C
<input type="checkbox"/> Back Canes Only		N/C
<input type="checkbox"/> TR14 Height Adjustable Solid Back 17" -26"		\$808
Back Additional Recline Option		
<input type="checkbox"/> Recline System (80 to 130 degrees) ¹		\$740
Snuggle Back Cushions:**		
<input type="checkbox"/> TR12 Short 12" x 18"		N/C
<input type="checkbox"/> TR14 Tall 14" x 26"		N/C
Comfort Side Support Wings		
<input type="checkbox"/> Comfort Side Support Wing - Left		N/C
<input type="checkbox"/> Comfort Side Support Wing - Right		N/C
Extra Back Cushion Options		
<input type="checkbox"/> Lumbar Support		\$55
Seating Depth Position (Select One - See Chart on Bottom of Order Form):		
<input type="checkbox"/> Seat Depth Position 1 (TR12: 5.5"- 8"/TR14: 8"-10")		N/C
<input type="checkbox"/> Seat Depth Position 2 (TR12: 7.5"-10"/TR14: 10"-12")		N/C
<input type="checkbox"/> Seat Depth Position 3 (TR12: 9.5"-12"/TR14: 12"-14")		N/C
<input type="checkbox"/> Seat Depth Position 4 (TR12: NA/TR14: 13"-15")		N/C
Planar Seat Cushion**		
<input type="checkbox"/> Snuggle Planar Seat Short Cushion - Standard Fabric (TR12: 12x12) / (TR14: 14x14)		\$454
<input type="checkbox"/> Snuggle Planar Seat Short Cushion - Spacer Fabric (TR12: 12x12) / (TR14: 14x14)		\$454
Contoured Seat Cushion**		
<input type="checkbox"/> Snuggle Contoured Seat Long Cushion - Standard Fabric (TR12: 12x14) / (TR14: 14x16)		\$500
<input type="checkbox"/> Snuggle Contoured Seat Long Cushion - Spacer Fabric (TR12: 12x14) / (TR14: 14x16)		\$500
Cushion Accent Colors		
<input type="checkbox"/> Panther Black		N/C
<input type="checkbox"/> Royal Blue		N/C
<input type="checkbox"/> Forest Green		N/C
<input type="checkbox"/> Sassy Purple		N/C
<input type="checkbox"/> Princess Pink		N/C
<input type="checkbox"/> Cherry Red		N/C

¹**Manufacturer Installation Required.**

²**Fast tab is only added to Cushion if 3-Pt Belt is added. New Cushion is required if 3-Pt Belt is added after-market.**

***All seat cushions must have fast tab for aftermarket 3-Pt Belt and H-Harness.**

****Back cushion dimensions are seat pan to top of back pan.**

<u>Frame Options</u>		<u>Price</u>
<input type="checkbox"/> Height Adjustable Push Handle ¹ (Outside Mount)		N/C
<input type="checkbox"/> Caster Locks ¹		\$144
<input type="checkbox"/> Curb Tipper (Cannot use with Rear Anti-Tip Tubes)		\$151
<input type="checkbox"/> Rear Anti-Tip Tubes (Cannot use with Curb Tipper)		\$177

¹**Manufacturer Installation Required.**

<u>Transit Required Items</u>		
<input type="checkbox"/> Transit Anchor Brackets**		\$308
<input type="checkbox"/> 2-Point Pelvic Belt		\$69
<input type="checkbox"/> Using other manufacturer's Head Support		\$0

¹**Manufacturer Installation Required.**

***Requires Head Support and Pelvic Belt.**

TR12	TR14
C905515	
C508446	
C904662	
	C905516
	C508447
	C904663
C905053	C905053
C904696	
	C904699
C904617	C904617
C904618	C904618
C904664	C904664
C905354	C905354
C905355	C905355
C905356	C905356
	C905360
C905412	C904710
C904708	C905425
C905415	C904783
C905429	C905430
C903854	C903854
C903850	C903850
C903857	C903857
C903851	C903851
C903759	C903759
C903855	C903855

TR12	TR14
C904688	C904688
C904689	C904689
C905310	C905310
C904649	C904649

TR12	TR14
C904576	C904576
C905386	C905386
C995681	C995681

<u>Transit Option</u>	<u>Price</u>
<input type="checkbox"/> Q'Straint Transit Lap Belt	\$285

<u>Head Supports (Required for Transit)*</u>	<u>Price</u>
<input type="checkbox"/> Adjustable Wing Head Support (Requires height adjustable solid back system).	\$409
<input type="checkbox"/> Integrated Adj Head Support Assembly (Option only on TR12 Seating System)	\$272
<input type="checkbox"/> R82 Adjustable Width Headrest, Black Cover - Short (Hardware Required)	\$469
<input type="checkbox"/> Turtle Neck Bar System, Short (Required)	\$354
<input type="checkbox"/> Headrest Universal Mount (Required)	\$84
R82 Anatomic Headrest with Black Cover (Hardware Required)	
<input type="checkbox"/> Small	\$393
<input type="checkbox"/> Medium	\$393
<input type="checkbox"/> Turtle Neck Bar System, Short (Required)	\$354
<input type="checkbox"/> Headrest Universal Mount (Required)	\$84

**The Transit model requires head support selection if other manufacturer's head support is not selected under "Transit Required Items."*

<u>Trunk Positioning</u>	<u>Price</u>
<input type="checkbox"/> H-Harness w/ Padded Covers ^{^^}	\$151
Torso Vest by R82	
<input type="checkbox"/> LT Shoulder to RT Hip: 7.75" - 10.75" (Size XS)	\$308
<input type="checkbox"/> LT Shoulder to RT Hip: 10.75" - 12.5" (Size S)	\$308
<input type="checkbox"/> LT Shoulder to RT Hip: 12.5" - 16" (Size M)	\$308
<input type="checkbox"/> LT Shoulder to RT Hip: 16" - 18.5" (Size L)	\$308
<input type="checkbox"/> LT Shoulder to RT Hip: 18.5" - 20.5" (Size XL)	\$308
<input type="checkbox"/> Camlock Kit for Vest (Required)	N/C
Cross Vest by R82	
<input type="checkbox"/> LT Shoulder to RT Hip: 7.75" - 12.5" (Size XS)	\$292
<input type="checkbox"/> LT Shoulder to RT Hip: 10.25" - 13.5" (Size S)	\$292
<input type="checkbox"/> LT Shoulder to RT Hip: 13.5" - 17.75" (Size M)	\$292
<input type="checkbox"/> LT Shoulder to RT Hip: 17.75" - 20.75" (Size L)	\$292
<input type="checkbox"/> LT Shoulder to RT Hip: 20.75" - 24.5" (Size XL)	\$292
<input type="checkbox"/> Camlock Kit for Vest (Required)	N/C

^{^^}Requires the Height Adjustable Solid Back System and 3-Point Pelvic Belt.

<u>Laterals</u>	<u>Price</u>
Fixed Laterals - Width, Height, Adj. Lateral Thoracic Support	
Left Lateral Padding Selection:	
<input type="checkbox"/> Size 1: 4" (D) x 4.5" (H)	\$225
<input type="checkbox"/> Size 2: 5" (D) x 4.5" (H)	\$225
Right Lateral Padding Selection:	
<input type="checkbox"/> Size 1: 4" (D) x 4.5" (H)	\$225
<input type="checkbox"/> Size 2: 5" (D) x 4.5" (H)	\$225
Swing-Away Laterals - Height Adj. Lateral Thoracic Support	
Left Lateral Padding Selection:	
<input type="checkbox"/> Size 1: 4" (D) x 2.75" (H)	\$181
<input type="checkbox"/> Size 2: 5.75" (D) x 4.25" (H)	\$181
<input type="checkbox"/> Swing Away Adj. Hardware (Required for S/A Laterals)	\$366
Right Lateral Padding Selection:	
<input type="checkbox"/> Size 1: 4" (D) x 2.75" (H)	\$181
<input type="checkbox"/> Size 2: 5.75" (D) x 4.25" (H)	\$181
<input type="checkbox"/> Swing Away Adj. Hardware (Required for S/A Laterals)	\$366
<input type="checkbox"/> 1" Thick Padded Lateral Support Cover [^]	\$76

[^]For Swing-Away Laterals only

<u>Pelvic Positioning (Required for Transit)*</u>	<u>Price</u>
<input type="checkbox"/> 2-Point Pelvic Belt	\$69
<input type="checkbox"/> 3-Point Pelvic Belt	\$75
4-Point Pelvic Belt by R82 (Cannot be paired with 2-Point Pelvic Belt)*	
<input type="checkbox"/> Hip Width: 6.75" - 12.5" (Size XS)	\$207
<input type="checkbox"/> Hip Width: 8.25" - 13.75" (Size S)	\$207
<input type="checkbox"/> Hip Width: 9.75" - 16.25" (Size M)	\$207
<input type="checkbox"/> Hip Width: 12.5" - 17.25" (Size L)	\$207

**The Transit model requires belt selection if 2-Point Belt is not selected under "Transit Required Items."*

<u>Hip Guides</u>	<u>Price</u>
Hip Guides without Armrest - Height Adj. Lateral Pelvic Support	
Left Size Selection:	
<input type="checkbox"/> Small - 4.5" - 5.5" (H)*	\$184
<input type="checkbox"/> Large - 5.25" - 6.25" (H)*	\$184
Right Size Selection:	
<input type="checkbox"/> Small - 4.5" - 5.5" (H)*	\$184
<input type="checkbox"/> Large - 5.25" - 6.25" (H)*	\$184
Hip Guides with Armrest - Height Adj. Lateral Pelvic Support (Cannot be used with Reversible Tray)	
Left Size Selection:	
<input type="checkbox"/> Small - 5"-6"(H)*	\$225
<input type="checkbox"/> Large - 7"-8" (H)*	\$225
Right Size Selection:	
<input type="checkbox"/> Small - 5"-6"(H)	\$225
<input type="checkbox"/> Large - 7"-8" (H)*	\$225
Medial Knee Abductor (Pommel)	
<input type="checkbox"/> Small (1.5" front width, 3.5" depth, 2" back width, 2.25" height)	\$249
<input type="checkbox"/> Medium (2.6" front width, 4.1" depth, 2" back width, 2.7" height)	\$249

**Effective Height is reduced by 1" when used with Planar Seat Cushion or 2" if used with Contoured Seat Cushion.*

TR12	TR14
C903094	C903094

TR12	TR14
C904608	C904608
C904615	
C905287	C905287
C904948	C904948
C904946	C904946
C904950	C904950
C904951	C904951
C904948	C904948
C904946	C904946

TR12	TR14
C904695	C904695
C905476	C905476
C905152	C905152
C905153	C905153
C905154	C905154
C905480	C905480
C905394	C905394
C905176	C905176
C905177	C905177
C905178	C905178
C905179	C905179
C905180	C905180
C905394	C905394

C904760	C904760
C904876	C904876
C904761	C904761
C904886	C904886
C904864	C904864
C904866	C904866
C904884	C904884
C904865	C904865
C904867	C904867
C904885	C904885
C905411	C905411

TR12	TR14
C905386	C905386
C904017	C904017
C905148	C905148
C905149	C905149
C905150	C905150
C905151	C905151

TR12	TR14
C904935	C905042
C905040	C904937
C904936	C905043
C905041	C904938
C904611	C904613
	C904927
C904612	C904614
	C904928
C904557	C904557
C904591	C904591

<u>Foot & Leg Positioning</u>	<u>Price</u>
<input type="checkbox"/> Elevating Leg Rest (90° - 180°)*	\$384
<input type="checkbox"/> Padded Elevating Legrest Cover	\$144
<input type="checkbox"/> Calf Panel (E.I.) H: 3.5"	\$95
<input type="checkbox"/> Calf Panel (Standard) H: 6"	\$95
<input type="checkbox"/> Individual Angle Adjustable Footplates (Elevating Leg Rest Required)*	\$530
<input type="checkbox"/> One Piece Angle Adjustable Flip Up Footplate (-15° - +15°)	\$170
<input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required)	\$500
Foot Positioners	
<input type="checkbox"/> Small 7" - 10"	\$240
<input type="checkbox"/> Medium 8" - 11"	\$240
Therafin Shoe Holders w/Padded Straps	
<input type="checkbox"/> Small	\$289
<input type="checkbox"/> Medium	\$289

***Manufacturer Installation Required.**

<u>Tire Options</u>	<u>Price</u>
<input type="checkbox"/> 12.5" Rear Solid Knobby Tires	\$307
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tires	\$459

<u>Medical Support Accessories</u>	<u>Price</u>
<input type="checkbox"/> Medical Necessity Storage Basket	\$454
<input type="checkbox"/> IV Pole (Height Adjustable and Collapsible)*	\$294
<input type="checkbox"/> Adj. Swing-Away Hardware (Required for IV Pole)	\$212
<input type="checkbox"/> Vertical O2 Tank Holder**	\$269
<input type="checkbox"/> LTV Ventilator Holder*	\$693
<input type="checkbox"/> Ventilator/Suction/Hard Tray	\$693
<input type="checkbox"/> Extended Ventilator/Suction/Hard Tray	\$693

***Manufacturer Installation Required.**

***Requires Rear Anti-Tip Tubes.**

***Cannot be used during transit.**

<u>Accessories</u>	<u>Price</u>
<input type="checkbox"/> Headrest Cover (Canopy)	\$351
<input type="checkbox"/> Extended Headrest Cover (Canopy)	\$384
<input type="checkbox"/> Upper Extremity Support Tray (Clear Tray)*^	\$171
<input type="checkbox"/> Adj Removable Hardware (Required for Clear Tray)	\$350
<input type="checkbox"/> Padded Tray Cover* (Clear Tray Required)	\$101
<input type="checkbox"/> Under Seat Storage Basket	\$176
<input type="checkbox"/> Utility Bag	\$65
<input type="checkbox"/> Rain Cover (Canopy Required)	\$232
<input type="checkbox"/> Mosquito Net (Canopy Required)	\$232

***Cannot be used during transit.**

^Not available with Hip Guides with Armrest.

<u>R82 High-Low:xo Frame Options</u>	<u>Price</u>
Manual Lifts	
<input type="checkbox"/> High-Low:xo Frame, Manual Lift, Indoor, Size 1	\$3,002
<input type="checkbox"/> Rear Anti-Tip for High-Low:xo Frame	\$308

The Trekker Transit option includes four manufacturer installed black wheelchair transportation anchors. The Trekker Transit model has been crash tested and conforms to RESNA WC19 standards and ISO 7176-19:2022.

Do not use the foot positioners on the Convoid Trekker during transit in a motor vehicle.

If the Trekker with back canes only option is selected it is the providers responsibility to ensure the seating components are WC20 compatible, mounted on the Trekker Mobility Base that conforms to RESNA WC19 and ISO 7176-19:2022 standards.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard shipping rates apply • Expedited shipping available, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 1/5/2026

Convoid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF0701-ENG-REV.00

etac[®]
Creating Possibilities

TR12	TR14
C905046	C905047
C904619	C904619
C905370	C905370
C904580	C904580
C905049	C905050
C904568	C904596
C904955	C904956
C903415	C903415
C903416	C903416
C904679	C904679
C904680	C904680

TR12	TR14
C904117	C904117
C904118	C904118


TR12	TR14
C904606	C904606
C904566	C904566
C905612	C905612
C904565	C904565
C904671	C904671
C904595	C904595
C905416	C905416

TR12	TR14
C904572	C904599
C904690	C904691
C904556	C904590
C905361	C905361
C904586	C904610
C904605	C904605
C904531	C904531
C904790	C904791
C904788	C904789

TR12	TR14
C905632	C905632
C905633	C905633

Seat Depth Adjustable Positions:

Seat Depth Growth				
Postion:	1	2	3	4
TR12:	5.5" - 8"	7.5" - 10"	9.5" - 12"	N/A
TR14:	8" - 10"	10" - 12"	12" - 14"	13" - 15"



Trekker Top Vests Belt Slot (Shoulder) Height from Seat Pan

	Minimum	Maximum
TR12 Integrated Headrest	9"	13"
TR12 Separate Headrest	13"	17"
TR14 Separate Headrest	17"	26"

Seating Module Specifications:**Convoid Trekker Seating Module Specifications**

Seat Width with Accessories	TR 12 (Size 1)	TR 14 (Size 2)
Trunk Width between Fixed Laterals (With Comfort Wings/Without Comfort Wings)	7.5" - 8" / 7.5" - 12"	9.5" - 10" / 9.5" - 14"
Trunk Width between Swing-Away Laterals (With Comfort Wings/Without Comfort Wings)	6" - 8" / 6" - 10"	7" - 9.5" / 7" - 14"
Hip Width between Hip Guides Without Armrests	5.5" - 11.5"	6.5" - 14"
Hip Width between Hip Guides With Armrests	5.5" - 11.5"	6.5" - 14"

Convoid Trekker Seat to Footplate Dimensions

Seat to Footplate Dimensions	TR 12 (Size 1)	TR 14 (Size 2)
Fixed 90° with one piece footplate	5.5" - 11.5"	5.5" - 11.5"
Fixed 90° with one piece footplate (flipped)	3.25" - 9.25"	3.25" - 9.25"
ELR with one piece footplate	4.75" - 10.75"	6.25" - 12.25"
ELR with one piece footplate (flipped)	2" - 8"	4.25" - 10.25"
ELR with individual footplates	5" - 11"	7" - 13"
ELR with individual footplates (flipped)	2" - 5"	4" - 7"