

Convoid Carrot 3 Booster Seat - Child Restraint System

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Carrot 3 Booster Seat Standard Features:

- | | |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> Carrot 3 Backrest Recline Adjustment (90° - 180°) Head Support: 8.5" x 10.5" (H x L) 4 Point Positioning Harness H-Harness Attaching Bracket | <ul style="list-style-type: none"> Booster Seat Base Lower Anchors Connector System 2" Backrest Extension Shoulder Protector Wings Tether Strap (Black) |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

Standard Models:

<input type="checkbox"/> Carrot 3 Booster Seat CRT3BLACK - CCRT3BS-10	Price
<input type="checkbox"/> Carrot 3 Booster Seat CRT3BLUE - CCRT3BS-20	\$3,158
<input type="checkbox"/> Carrot 3 Booster Seat CRT3RED - CCRT3BS-30	\$3,158

Carrot Booster Seat Recommended Item

<input type="checkbox"/> Top Securement Strap*	Price
<input type="checkbox"/> Lower Securement Strap*	\$58
*For vehicles not equipped with anchor points.	
	\$73

Accessories

<input type="checkbox"/> Booster Support Tray (Crash Tested)	Price
<input type="checkbox"/> 3D Body Balance Pads L/R Set	\$323
<input type="checkbox"/> Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°)	\$168
<input type="checkbox"/> Backrest Extension - 4" (Extends back height to 32.8")	\$126
<input type="checkbox"/> Large Crotch Belt	\$101
	\$264

Black	Blue	Red
CCT3HSF-10	CCT3HSF-10	CCT3HSF-10
CCT3LFS-10	CCT3LFS-10	CCT3LFS-10

Black	Blue	Red
CCRT3BT-19	CCRT3BT-29	CCRT3BT-39
CCT3BBP-00	CCT3BBP-00	CCT3BBP-00
CCT3SA0-10	CCT3SA0-20	CCT3SA0-30
CCT3BE4-10	CCT3BE4-20	CCT3BE4-30
CCT3CBX-00	CCT3CBX-00	CCT3CBX-00

Standard domestic shipping rates apply • For expedited shipping rates, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.


Effective: 1/5/2026

Convoid
 2830 California Street
 Torrance, CA 90503
 1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF1100-ENG-REV.00

Convoid
 by Etac

CLIENT MEASUREMENTS:



A = _____ in.

B = _____ in.

C = _____ in.

D = _____ in.

E = _____ in.

F = _____ in.

A = Seat back height
B = Seat depth
C = Seat to floor

D = Shoulder width
E = Chest width
F = Hip width

Carrot 3 Base Model Specifications:

Seat Width:	16.5"
Seat Depth:	15" - 16.5" *
Backrest Height Adjustment:	28" - 33"
Head Support Height:	8.5"
Backrest Height:	29" - 34"
Backrest Extension Height (Standard):	2"
Backrest Extension Height (Option):	4"
Shoulder Wing Height (Standard):	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"
Product Weight:	21 lbs.
Recline Adjustment:	90° - 180°
Product Overall Width:	20" - 22"

User Measurement Guide:

Minimum User Weight:	79 lbs.
Maximum User Weight:	165 lbs.
Minimum User Height:	55"
Maximum User Height:	70"

*Factory Setting

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166