

Convoid Trekker - HCPCS E1234

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Trekker Standard Features:

- Trekker Folding Mobility Frame Accepts TR12 & TR14 Seating Systems
- Tilt-in-Space (-10° to +45°)
- 180° Reversible Seating Configuration
- 90° Fixed Front Rigging
- 7.5" x 2" Front and 11.5" Rear Solid Tires
- One Piece Height Adjustable Push Handle (Inside Mount)

Standard Models:

- ☐ TR Seating System - C904768
- ☐ TR Mobility Frame - C904661

Price

\$0
\$2,663

Transit Model Additional Features:

- Transit Anchor Brackets

Transit Included Models:

- ☐ TR Seating System - C904768
- ☐ TR Mobility Frame Transit - C905362

Price

\$0
\$2,886

<u>Trekker Seating System</u>			<u>TR12</u>	<u>TR14</u>
<input type="checkbox"/> Seat Width 12"	Price	HCPCS	C905515	
<input type="checkbox"/> Back Canes Only	N/C		C508446	
<input type="checkbox"/> TR12 Height Adjustable Solid Back 9" - 17"	\$584	E2293	C904662	
<input type="checkbox"/> Seat Width 14"	N/C			C905516
<input type="checkbox"/> Back Canes Only	N/C			C508447
<input type="checkbox"/> TR14 Height Adjustable Solid Back 17" -26"	\$584	E2293		C904663
Back Additional Recline Option				
<input type="checkbox"/> Recline Option (80 to 170 degrees) ¹	\$535	E1226	C905053	C905053
Snuggle Back Cushions:**				
<input type="checkbox"/> TR12 Short	N/C		C904696	
<input type="checkbox"/> TR14 Tall	N/C			C904699
Comfort Side Support Wings				
<input type="checkbox"/> Comfort Side Support Wing - Left	N/C		C904617	C904617
<input type="checkbox"/> Comfort Side Support Wing - Right	N/C		C904618	C904618
Extra Back Cushion Options				
<input type="checkbox"/> Lumbar Support	\$40	K0108	C904664	C904664
Seating Depth Position (Select One - See Chart on Bottom of Order Form):				
<input type="checkbox"/> Seat Depth Position 1 (TR12: 5.5"-8"/TR14: 8"-10")	N/C		C905354	C905354
<input type="checkbox"/> Seat Depth Position 2 (TR12: 7.5"-10"/TR14: 10"-12")	N/C		C905355	C905355
<input type="checkbox"/> Seat Depth Position 3 (TR12: 9.5"-12"/TR14: 12"-14")	N/C		C905356	C905356
<input type="checkbox"/> Seat Depth Position 4 (TR12: NA/TR14: 13"-15")	N/C			C905360
Planar Seat Cushion**				
<input type="checkbox"/> Snuggle Planar Seat Short Cushion - Standard Fabric (TR12: 12x12) / (TR14: 14x14)	\$328	E2292	C905412	C904710
<input type="checkbox"/> Snuggle Planar Seat Short Cushion - Spacer Fabric (TR12: 12x12) / (TR14: 14x14)	\$328	E2292	C904708	C905425
Contoured Seat Cushion**				
<input type="checkbox"/> Snuggle Contoured Seat Long Cushion - Standard Fabric (TR12: 12x14) / (TR14: 14x16)	\$361	E2294	C905415	C904783
<input type="checkbox"/> Snuggle Contoured Seat Long Cushion - Spacer Fabric (TR12: 12x14) / (TR14: 14x16)	\$361	E2294	C905429	C905430
Cushion Accent Colors				
<input type="checkbox"/> Panther Black	N/C		C903854	C903854
<input type="checkbox"/> Royal Blue	N/C		C903850	C903850
<input type="checkbox"/> Forest Green	N/C		C903857	C903857
<input type="checkbox"/> Sassy Purple	N/C		C903851	C903851
<input type="checkbox"/> Princess Pink	N/C		C903759	C903759
<input type="checkbox"/> Cherry Red	N/C		C903855	C903855

¹Manufacturer Installation Required.

²Fast tab is only added to Cushion if 3-Pt Belt is added. New Cushion is required if 3-Pt Belt is added after-market.

³All seat cushions must have fast tab for aftermarket 3-Pt Belt and H-Harness.

⁴Back cushion dimensions are seat pan to top of back pan.

<u>Frame Options</u>			<u>TR12</u>	<u>TR14</u>
<input type="checkbox"/> Height Adjustable Push Handle ¹ (Outside Mount)	Price	HCPCS	C904688	C904688
<input type="checkbox"/> Caster Locks ¹	N/C		C904689	C904689
<input type="checkbox"/> Curb Tipper (Cannot use with Rear Anti-Tip Tubes)	\$104	K0073	C905310	C905310
<input type="checkbox"/> Rear Anti-Tip Tubes (Cannot use with Curb Tipper)	\$109		C904649	C904649
<input type="checkbox"/> Rear Anti-Tip Tubes (Cannot use with Curb Tipper)	\$128	E0971		
¹Manufacturer Installation Required.				
<u>Transit Required Items</u>			<u>TR12</u>	<u>TR14</u>
<input type="checkbox"/> Transit Anchor Brackets ^{1A*}	Price	HCPCS	C904576	C904576
<input type="checkbox"/> 2-Point Pelvic Belt	\$223	E1022	C905386	C905386
<input type="checkbox"/> Using other manufacturer's Head Support	\$50	E0978	C995681	C995681
<input type="checkbox"/> Using other manufacturer's Head Support	\$0			
¹Manufacturer Installation Required.				
²The Transit Included Model excludes the price of the transit required items except the Transit Anchor Brackets.				
³Requires Head Support and Pelvic Belt.				
<u>Transit Option</u>			<u>TR12</u>	<u>TR14</u>
<input type="checkbox"/> Q'Straint Transit Lap Belt	Price	HCPCS	C903094	C903094
	\$206	A9270		

	Price	HCPCS
Head Supports (Required for Transit)*		
<input type="checkbox"/> Adjustable Wing Head Support (Requires height adjustable solid back system).	\$296	E0955
<input type="checkbox"/> Integrated Adj Head Support Assembly (Option only on TR12 Seating System)	\$197	E0955
<input type="checkbox"/> R82 Adjustable Width Headrest, Black Cover - Short (Hardware Required)	\$339	E0955
<input type="checkbox"/> Turtle Neck Bar System, Short (Required)	\$256	E1028
<input type="checkbox"/> Headrest Universal Mount (Required)	\$60	K0108
R82 Anatomic Headrest with Black Cover (Hardware Required)		
<input type="checkbox"/> Small	\$284	E0955
<input type="checkbox"/> Medium	\$284	E0955
<input type="checkbox"/> Turtle Neck Bar System, Short (Required)	\$256	E1028
<input type="checkbox"/> Headrest Universal Mount (Required)	\$60	K0108

***The Transit model requires head support selection if other manufacturer's head support is not selected under "Transit Required Items."**

TR12	TR14
C904608	C904608
C904615	
C905287	C905287
C904948	C904948
C904946	C904946
C904950	C904950
C904951	C904951
C904948	C904948
C904946	C904946

	Price	HCPCS
Trunk Positioning		
<input type="checkbox"/> H-Harness w/ Padded Covers ^{^^}	\$109	E0960
Torso Vest by R82		
<input type="checkbox"/> LT Shoulder to RT Hip: 7.75" - 10.75" (Size XS)	\$223	E0960
<input type="checkbox"/> LT Shoulder to RT Hip: 10.75" - 12.5" (Size S)	\$223	E0960
<input type="checkbox"/> LT Shoulder to RT Hip: 12.5" - 16" (Size M)	\$223	E0960
<input type="checkbox"/> LT Shoulder to RT Hip: 16" - 18.5" (Size L)	\$223	E0960
<input type="checkbox"/> LT Shoulder to RT Hip: 18.5" - 20.5" (Size XL)	\$223	E0960
<input type="checkbox"/> Camlock Kit for Vest (Required)	N/C	
Cross Vest by R82		
<input type="checkbox"/> LT Shoulder to RT Hip: 7.75" - 12.5" (Size XS)	\$211	E0960
<input type="checkbox"/> LT Shoulder to RT Hip: 10.25" - 13.5" (Size S)	\$211	E0960
<input type="checkbox"/> LT Shoulder to RT Hip: 13.5" - 17.75" (Size M)	\$211	E0960
<input type="checkbox"/> LT Shoulder to RT Hip: 17.75" - 20.75" (Size L)	\$211	E0960
<input type="checkbox"/> LT Shoulder to RT Hip: 20.75" - 24.5" (Size XL)	\$211	E0960
<input type="checkbox"/> Camlock Kit for Vest (Required)	N/C	

^^Requires the Height Adjustable Solid Back System and 3-Point Pelvic Belt.

TR12	TR14
C904695	C904695
C905476	C905476
C905152	C905152
C905153	C905153
C905154	C905154
C905480	C905480
C905394	C905394
C905176	C905176
C905177	C905177
C905178	C905178
C905179	C905179
C905180	C905180
C905394	C905394

	Price	HCPCS
Laterals		
Fixed Laterals - Width, Height, Adj. Lateral Thoracic Support		
Left Lateral Padding Selection:		
<input type="checkbox"/> Size 1: 4" (D) x 4.5" (H)	\$162	E0956
<input type="checkbox"/> Size 2: 5" (D) x 4.5" (H)	\$162	E0956
Right Lateral Padding Selection:		
<input type="checkbox"/> Size 1: 4" (D) x 4.5" (H)	\$162	E0956
<input type="checkbox"/> Size 2: 5" (D) x 4.5" (H)	\$162	E0956
Swing-Away Laterals - Height Adj. Lateral Thoracic Support		
Left Lateral Padding Selection:		
<input type="checkbox"/> Size 1: 4" (D) x 2.75" (H)	\$131	E0956
<input type="checkbox"/> Size 2: 5.75" (D) x 4.25" (H)	\$131	E0956
<input type="checkbox"/> Swing Away Adj. Hardware (Required for S/A Laterals)	\$264	E1028
Right Lateral Padding Selection:		
<input type="checkbox"/> Size 1: 4" (D) x 2.75" (H)	\$131	E0956
<input type="checkbox"/> Size 2: 5.75" (D) x 4.25" (H)	\$131	E0956
<input type="checkbox"/> Swing Away Adj. Hardware (Required for S/A Laterals)	\$264	E1028
<input type="checkbox"/> 1" Thick Padded Lateral Support Cover [^]	\$55	E0960

^For Swing-Away Laterals only

C904760	C904760
C904876	C904876
C904761	C904761
C904886	C904886
C904864	C904864
C904866	C904866
C904884	C904884
C904865	C904865
C904867	C904867
C904885	C904885
C905411	C905411

	Price	HCPCS
Pelvic Positioning (Required for Transit)*		
<input type="checkbox"/> 2-Point Pelvic Belt	\$50	E0978
<input type="checkbox"/> 3-Point Pelvic Belt	\$54	E0978
4-Point Pelvic Belt by R82 (Cannot be paired with 2-Point Pelvic Belt)*		
<input type="checkbox"/> Hip Width: 6.75" - 12.5" (Size XS)	\$150	K0108
<input type="checkbox"/> Hip Width: 8.25" - 13.75" (Size S)	\$150	K0108
<input type="checkbox"/> Hip Width: 9.75" - 16.25" (Size M)	\$150	K0108
<input type="checkbox"/> Hip Width: 12.5" - 17.25" (Size L)	\$150	K0108

***The Transit model requires belt selection if 2-Point Belt is not selected under "Transit Required Items."**

TR12	TR14
C905386	C905386
C904017	C904017
C905148	C905148
C905149	C905149
C905150	C905150
C905151	C905151

	Price	HCPCS
Hip Guides		
Hip Guides without Armrest - Height Adj. Lateral Pelvic Support		
Left Size Selection:		
<input type="checkbox"/> Small - 4.5" - 5.5" (H)*	\$133	E0956
<input type="checkbox"/> Large - 5.25" - 6.25" (H)*	\$133	E0956
Right Size Selection:		
<input type="checkbox"/> Small - 4.5" - 5.5" (H)*	\$133	E0956
<input type="checkbox"/> Large - 5.25" - 6.25" (H)*	\$133	E0956
Hip Guides with Armrest - Height Adj. Lateral Pelvic Support (Cannot be used with Reversible Tray)		
Left Size Selection:		
<input type="checkbox"/> Small - 5"-6"(H)*	\$162	E0956
<input type="checkbox"/> Large - 7"-8" (H)*	\$162	E0956
Right Size Selection:		
<input type="checkbox"/> Small - 5"-6"(H)	\$162	E0956
<input type="checkbox"/> Large - 7"-8" (H)*	\$162	E0956
Medial Knee Abductor (Pommel)		
<input type="checkbox"/> Small (1.5" front width, 3.5" depth, 2" back width, 2.25" height)	\$180	E0957
<input type="checkbox"/> Medium (2.6" front width, 4.1" depth, 2" back width, 2.7" height)	\$180	E0957

***Effective Height is reduced by 1" when used with Planar Seat Cushion or 2" if used with Contoured Seat Cushion.**

TR12	TR14
C904935	C905042
C905040	C904937
C904936	C905043
C905041	C904938
C904611	C904613
	C904927
C904612	C904614
	C904928
C904557	C904557
C904591	C904591

Retail Price List / Order Form

Convaid Trekker - HCPCS E1234

- Foot & Leg Positioning**
- ☐ Elevating Leg Rest (90° - 180°)¹
 - ☐ Padded Elevating Legrest Cover
 - ☐ Calf Panel (E.I.) H: 3.5"
 - ☐ Calf Panel (Standard) H: 6"
 - ☐ Individual Angle Adjustable Footplates (Elevating Leg Rest Required)¹
 - ☐ One Piece Angle Adjustable Flip Up Footplate (-15° - +15°)
 - ☐ Full Padded Footbox (One-Piece Footplate Required)

Price	HCPCS
\$278	E0990
\$104	K0108
\$69	E0995
\$69	E0995
\$383	K0040
\$123	K0040
\$361	K0108

Foot Positioners

- ☐ Small 7" - 10"
- ☐ Medium 8" - 11"

\$174	K0108
\$174	K0108

Therafin Shoe Holders w/Padded Straps

- ☐ Small
- ☐ Medium

\$209	K0108
\$209	K0108

¹Manufacturer Installation Required.

TR12	TR14
C905046	C905047
C904619	C904619
C905370	C905370
C904580	C904580
C905049	C905050
C904568	C904596
C904955	C904956
C903415	C903415
C903416	C903416
C904679	C904679
C904680	C904680

- Tire Options**
- ☐ 12.5" Rear Solid Knobby Tires
 - ☐ 7.5" Front, 12.5" Rear Pneumatic Knobby Tires

Price	HCPCS
\$222	
\$332	

TR12	TR14
C904117	C904117
C904118	C904118

- Medical Support Accessories**
- ☐ Medical Necessity Storage Basket
 - ☐ IV Pole (Height Adjustable and Collapsible)^o
 - ☐ Adj. Swing-Away Hardware (Required for IV Pole)
 - ☐ Vertical O2 Tank Holder**
 - ☐ LTV Ventilator Holder¹
 - ☐ Ventilator/Suction/Hard Tray
 - ☐ Extended Ventilator/Suction/Hard Tray

Price	HCPCS
\$328	
\$212	K0105
\$153	E1028
\$195	E2208
\$500	E1029
\$500	E1029
\$500	E1029

TR12	TR14
C904606	C904606
C904566	C904566
C905612	C905612
C904565	C904565
C904671	C904671
C904595	C904595
C905416	C905416

¹Manufacturer Installation Required.^oRequires Rear Anti-Tip Tubes.

*Cannot be used during transit.

- Accessories**
- ☐ Headrest Cover (Canopy)
 - ☐ Extended Headrest Cover (Canopy)
 - ☐ Upper Extremity Support Tray (Clear Tray)**
 - ☐ Adj Removable Hardware (Required for Clear Tray)
 - ☐ Padded Tray Cover* (Clear Tray Required)
 - ☐ Under Seat Storage Basket
 - ☐ Utility Bag
 - ☐ Rain Cover (Canopy Required)
 - ☐ Mosquito Net (Canopy Required)

Price	HCPCS
\$254	K0108
\$278	K0108
\$124	E0950
\$253	
\$73	K0108
\$127	
\$47	
\$168	
\$168	

TR12	TR14
C904572	C904599
C904690	C904691
C904556	C904590
C905361	C905361
C904586	C904610
C904605	C904605
C904531	C904531
C904790	C904791
C904788	C904789

*Cannot be used during transit.

^Not available with Hip Guides with Armrest.

- R82 High-Low:xo Frame Options**
- Manual Lifts**
- ☐ High-Low:xo Frame, Manual Lift, Indoor, Size 1
 - ☐ Rear Anti-Tip for High-Low:xo Frame

Price
\$2,169
\$223

TR12	TR14
C905632	C905632
C905633	C905633

The Trekker Transit option includes four manufacturer installed black wheelchair transportation anchors. The Trekker Transit model has been crash tested and conforms to RESNA WC19 standards and ISO 7176-19:2022.

Do not use the foot positioners on the Convaid Trekker during transit in a motor vehicle.

If the Trekker with back canes only option is selected it is the providers responsibility to ensure the seating components are WC20 compatible, mounted on the Trekker Mobility Base that conforms to RESNA WC19 and ISO 7176-19:2022 standards.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard domestic shipping rates apply • For expedited shipping rates, please call 1-844-876-6245 or email convaidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.


Effective: 1/5/2026

OF0700-ENG-REV.00

Convaid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

Seat Depth Adjustable Positions:

Seat Depth Growth				
Postion:	1	2	3	4
TR12:	5.5" - 8"	7.5" - 10"	9.5" - 12"	N/A
TR14:	8" - 10"	10" - 12"	12" - 14"	13" - 15"



Trekker Top Vests Belt Slot (Shoulder) Height from Seat Pan

	Minimum	Maximum
TR12 Integrated Headrest	9"	13"
TR12 Separate Headrest	13"	17"
TR14 Separate Headrest	17"	26"

Seating Module Specifications:**Convoid Trekker Seating Module Specifications**

Seat Width with Accessories	TR 12 (Size 1)	TR 14 (Size 2)
Trunk Width between Fixed Laterals (With Comfort Wings/Without Comfort Wings)	7.5" - 8" / 7.5" - 12"	9.5" - 10" / 9.5" - 14"
Trunk Width between Swing-Away Laterals (With Comfort Wings/Without Comfort Wings)	6" - 8" / 6" - 10"	7" - 9.5" / 7" - 14"
Hip Width between Hip Guides Without Armrests	5.5" - 11.5"	6.5" - 14"
Hip Width between Hip Guides With Armrests	5.5" - 11.5"	6.5" - 14"

Convoid Trekker Seat to Footplate Dimensions

Seat to Footplate Dimensions	TR 12 (Size 1)	TR 14 (Size 2)
Fixed 90° with one piece footplate	5.5" - 11.5"	5.5" - 11.5"
Fixed 90° with one piece footplate (flipped)	3.25" - 9.25"	3.25" - 9.25"
ELR with one piece footplate	4.75" - 10.75"	6.25" - 12.25"
ELR with one piece footplate (flipped)	2" - 8"	4.25" - 10.25"
ELR with individual footplates	5" - 11"	7" - 13"
ELR with individual footplates (flipped)	2" - 5"	4" - 7"