

# Convaid Cruiser New Style - HCPCS E1236

## Retail Price List / Order Form

Account #: \_\_\_\_\_

Purchase Order #: \_\_\_\_\_

Date: \_\_\_\_\_

Contact: \_\_\_\_\_

Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

### Cruiser (New Style) Standard Features:

- Removable, Height Adjustable Swing-Away Footplates
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire (CX10/12/14/16)
- 8" Front, 11" Rear Solid Tire (CX18)
- Quick Release Rolko Rear Wheels (CX10/12/14/16)
- Foot Operated Wheel Lock (CX10/12/14/16)
- Toggle Wheel Locks (CX18)
- Zippered Pockets in Seat & Back Upholstery (Cordura Only: CX10/12/14/16)
- One Piece Height Adjustable Push Handle (CX10/12/14/16)
- Two Piece Push Handles (CX18)
- Tool-less Adjustable Seat Depth Growth of 5"
- Tool-less Upholstery Removal
- 2-Point Positioning Belt
- 30° Fixed Tilt
- Silver Powder Coated Frame
- Self-Tensio®
- Tension Adjustable Back
- Made in the USA

### Standard Models:

	Price
<input type="checkbox"/> Cruiser CX10 - 903314	\$1,932
<input type="checkbox"/> Cruiser CX12 - 902845	\$1,932
<input type="checkbox"/> Cruiser CX14 - 900490	\$1,932
<input type="checkbox"/> Cruiser CX16 - 900145	\$1,932
<input type="checkbox"/> Cruiser CX18 - 902594	\$1,932

### Upholstery Colors

	Price
<input type="checkbox"/> Cordura Upholstery	\$0
<input type="checkbox"/> Royal Blue	\$0
<input type="checkbox"/> Forest Green	\$0
<input type="checkbox"/> Cherry Red	\$0
<input type="checkbox"/> Sassy Purple	\$0
<input type="checkbox"/> Princess Pink	\$0
<input type="checkbox"/> Panther Black	\$0
<input type="checkbox"/> Chocolate Brown	\$0
<input type="checkbox"/> Aqua Blue	\$0
<input type="checkbox"/> Apple Green	\$0
<input type="checkbox"/> Textilene Upholstery	\$0
<input type="checkbox"/> Navy Blue	\$0
<input type="checkbox"/> Violet Purple	\$0
<input type="checkbox"/> Teal Green	\$0
<input type="checkbox"/> Wild Raspberry	\$0

CX10	CX12	CX14	CX16	CX18
904498	904499	904500	904501	904502
903850	903850	903850	903850	903850
903857	903857	903857	903857	903857
903855	903855	903855	903855	903855
903851	903851	903851	903851	903851
903759	903759	903759	903759	903759
903854	903854	903854	903854	903854
904189	904189	904189	904189	904189
904753	904753	904753	904753	904753
904754	904754	904754	904754	904754
904480	904481	904482	904483	904484
903464	903464	903464	903464	903464
903465	903465	903465	903465	903465
903466	903466	903466	903466	903466
903467	903467	903467	903467	903467

### Push Handle Options

<input type="checkbox"/> Two Piece Push Handle (Manufacturer Installation Required) (Req for CX18)	Price \$0
--	--------------

CX10	CX12	CX14	CX16	CX18
904124	904125	904126	904127	904461

### Transit Required Items

<input type="checkbox"/> Transit Option <sup>1</sup>	Price \$205
Head Support - Select Style Below:	
<input type="checkbox"/> Headrest Extension - Curved - Cordura ^	\$133
<input type="checkbox"/> Headrest Extension - Rectangular - Textilene **	\$133
<input type="checkbox"/> Contoured Firm Headrest (Planar Option Only/Hrdw Req)	\$250
<input type="checkbox"/> Adj. Removable Hardware (Req for Contoured Headrest)	\$253
Trunk Positioning - Select Style Below:	
<input type="checkbox"/> 5-Point Positioning/Restraint Harness (CX10/12 only) * <sup>1</sup>	\$169
<input type="checkbox"/> H-Harness w/ Padded Covers (3-Pt Belt Req)	\$147
<input type="checkbox"/> 3-Pt Positioning Belt	\$50
<input type="checkbox"/> Full Torso Support Vest - Cordura	\$270
<input type="checkbox"/> Full Torso Support Vest - Textilene	\$270
Foot Positioners - Select Size Below:	
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	\$160
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	\$160
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	\$160

### HCPCS

	CX10	CX12	CX14	CX16	CX18
	903316	902847	900491	900173	902596
<b>E0966</b>	903332	902863	900484	900138	902619
<b>E0966</b>	904514	904515	904516	904517	904518
<b>E0955</b>	904032	904032	904033	904034	904034
<b>E1028</b>	904035	904035	904035	904035	904035
<b>E0960</b>	903682	903683			
<b>E0960</b>	903988	903988	903988	903988	903988
<b>E0978</b>	904017	904017	904017	904017	904017
	903330	902861	900680	900339	902616
	905436	905437	905438	905439	902616
	903415	903415	903415	903415	903415
	903416	903416	903416	903416	903416
	903417	903417	903417	903417	903417

<sup>^</sup> Not available with Planar Option.

\* Crash Tested - WC19 and FMVSS213 compliant. Not available with Reinforced Upholstery, Anatomic Back, Cushion, Planar Option, or H-Harness/3-Pt Belt.

\*\*Required when ordering the Contoured firm headrest and flat firm backrest extension.

<sup>1</sup>Manufacturer Installation Required.

### Transit Options

<input type="checkbox"/> Q'Straint Transit Lap Belt	Price \$92	HCPCS <b>E0978</b>
<input type="checkbox"/> Sure-Lok Transit Lap Belt	\$92	<b>E0978</b>

	CX10	CX12	CX14	CX16	CX18
	903094	903094	903094	903094	903094
	903232	903232	903232	903232	903232

### Scout Option

<input type="checkbox"/> Scout Option (Manufacturer Installation Required)*	Price \$614	HCPCS
• Anatomic Back	\$0	<b>K0015</b>
• Flip-Up Armrests	\$0	
• 8" x 2" Front Ribbed, 12 1/2" Rear Pneumatic Knobby Tire	\$0	
• Attendant Hand Brakes	\$0	
Upholstery Style - Choose One		
<input type="checkbox"/> Heavy-Duty Reinforced Upholstery - Cordura	\$0	
<input type="checkbox"/> Heavy-Duty Reinforced Upholstery - Textilene	\$0	

	CX10	CX12	CX14	CX16	CX18
		904532	904533	904534	904535
		902932B	902934B	902936B	902938B
		902933B	902933B	902933B	902933B
		903083B	903083B	903083B	903083B
		904536	904537	904538	904539
		902853B	900443B	900139B	902601B
		904486B	904487B	904488B	904489B

\* All parts listed are required when ordering the Scout Option. Not available with Planar Option.

**Retail Price List / Order Form**

**Convaid Cruiser New Style - HCPCS E1236**

<u>Scout Accessories</u>	<u>Price</u>
<input type="checkbox"/> Extended Handles - 4" *^1	<b>\$185</b>
<input type="checkbox"/> 8" x 2" Front Ribbed, 12 1/2" Rear Solid Knobby Tire *	<b>\$119</b>
<input type="checkbox"/> Scout Curb Tipper <sup>1</sup>	<b>\$100</b>

\* Requires Scout Option.  
 ^ Requires 2-Piece Push Handle Package.  
<sup>1</sup>Manufacturer Installation Required.

HCPCS	CX10	CX12	CX14	CX16	CX18
		902515	902541		
		903224	903224	903224	903224
		903916	903916	903916	903916

<u>Planar Option *</u>	<u>Price</u>
<input type="checkbox"/> Planar Seat Modification	<b>\$353</b>
<input type="checkbox"/> Planar Back with Attachment (Incl. 2 Backrest Ext. w/ Hrdw.)	<b>\$327</b>
Cushion - Select One Cushion Below:	
<input type="checkbox"/> Support - General Use Cushion (Solid Seat Insert Included) <sup>1</sup>	<b>\$247</b>
<input type="checkbox"/> Position - Medial Thigh Support/Anti Thrust (Solid Seat Insert Included) <sup>1</sup>	<b>\$478</b>
<input type="checkbox"/> Align - Lateral Pelvis/Lateral Thigh Support (Solid Seat Insert Included) <sup>1</sup>	<b>\$478</b>
Planar Head Support Options:	
<input type="checkbox"/> Contoured Firm Headrest (Hardware Required)	<b>\$250</b>
<input type="checkbox"/> Adj. Removable Hardware (Req for Contoured Firm Headrest)	<b>\$253</b>
<input type="checkbox"/> Flat Firm Backrest Extension (Hardware Req)	<b>\$234</b>
<input type="checkbox"/> Adj. Removable Hardware (Req for Flat Firm Backrest Ext)	<b>\$229</b>

\* Planar Option must be ordered with Planar Seat, Planar Back, Cushion, Planar Head Support, and Textilene Upholstery.  
<sup>1</sup>Manufacturer Installation Required.

HCPCS	CX10	CX12	CX14	CX16	CX18
	904019	904020	904021	904022	904023
<b>E2291</b>	904024	904025	904026	904027	904028
<b>E2601</b>	903639	903642	903645	903648	903651
<b>E2607</b>	903641	903644	903647	903650	903653
<b>E2607</b>	903640	903643	903646	903649	903652
<b>E0955</b>	904032	904032	904033	904034	904034
<b>E1028</b>	904035	904035	904035	904035	904035
<b>E0955</b>	904089	904090	904091	904092	904093
<b>E1028</b>	904057	904057	904057	904057	904057

<u>Accessories for Planar Option</u>	<u>Price</u>
Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)	
<input type="checkbox"/> Small = 3"x 5"	<b>\$140</b>
<input type="checkbox"/> Medium = 3.75"x 5.5"	<b>\$140</b>
<input type="checkbox"/> Large = 4.75"x 6"	<b>\$140</b>
Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)	
<input type="checkbox"/> Small = 3"x 5"	<b>\$140</b>
<input type="checkbox"/> Medium = 3.75"x 5.5"	<b>\$140</b>
<input type="checkbox"/> Large = 4.75"x 6"	<b>\$140</b>
<input type="checkbox"/> Right Adj. Removable Hardware (Req for Lateral Support)	<b>\$182</b>
<input type="checkbox"/> Left Adj. Removable Hardware (Req for Lateral Support)	<b>\$182</b>

HCPCS	CX10	CX12	CX14	CX16	CX18
<b>E0956</b>	902965	902965	902965	902965	902965
<b>E0956</b>	902969	902969	902969	902969	902969
<b>E0956</b>	902973	902973	902973	902973	902973
<b>E0956</b>	904462	904462	904462	904462	904462
<b>E0956</b>	904463	904463	904463	904463	904463
<b>E0956</b>	904464	904464	904464	904464	904464
<b>E1028</b>	904030	904030	904030	904031	904031
<b>E1028</b>	904030	904030	904030	904031	904031

<u>Pelvic Positioning</u>	<u>Price</u>
<input type="checkbox"/> 3-Pt Positioning Belt	<b>\$50</b>
<input type="checkbox"/> Heavy Duty Reinforced Upholstery - Cordura ^	<b>\$275</b>
<input type="checkbox"/> Heavy Duty Reinforced Upholstery - Textilene ^	<b>\$275</b>
<input type="checkbox"/> Custom Seat Depth Extension 3" ^1	<b>\$96</b>
<input type="checkbox"/> Support - General Use Cushion (Solid Seat Insert Included) *1	<b>\$247</b>
<input type="checkbox"/> Position - Medial Thigh Support/Anti Thrust (Solid Seat Insert Included) *1	<b>\$478</b>
<input type="checkbox"/> Align - Lateral Pelvis/Lateral Thigh Support (Solid Seat Insert Included) *1	<b>\$478</b>
<input type="checkbox"/> Incontinence Cushion Liner (Cushion Req) <sup>1</sup>	<b>\$74</b>

^ Not available with Planar Option.  
 \*Not available with Seat Depth Extension or Cordura Upholstery. Abd/Add Different Design.  
<sup>1</sup>Manufacturer Installation Required.

HCPCS	CX10	CX12	CX14	CX16	CX18
<b>E0978</b>	904017	904017	904017	904017	904017
	903322	902853	900443	900139	902601
	904485	904486	904487	904488	904489
<b>E1297</b>	903995	903997	903999	904001	
<b>E2601</b>	903639	903642	903645	903648	903651
<b>E2607</b>	903641	903644	903647	903650	903653
<b>E2607</b>	903640	903643	903646	903649	903652
<b>E2619</b>	902324	902324	902324	902324	902324

<u>Foot &amp; Leg Positioning</u>	<u>Price</u>
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	<b>\$160</b>
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	<b>\$160</b>
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	<b>\$160</b>
<input type="checkbox"/> Angle Adjustable Footplates (pair)	<b>\$183</b>
<input type="checkbox"/> Footplate Securement Strap	<b>\$31</b>
<input type="checkbox"/> Medial Thigh Support-Abductor (N/A w/ Position/Align Cushion) <sup>1</sup>	<b>\$166</b>
<input type="checkbox"/> Lateral Thigh Support-Adductor (N/A w/ Position/Align Cushion) <sup>1</sup>	<b>\$166</b>
<input type="checkbox"/> Calf Panel *	<b>\$63</b>

\* Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.  
<sup>1</sup>Manufacturer Installation Required.

HCPCS	CX10	CX12	CX14	CX16	CX18
	903415	903415	903415	903415	903415
	903416	903416	903416	903416	903416
	903417	903417	903417	903417	903417
<b>K0040</b>	903445	903268	903268	903269	903269
	903584	903584	903584	903584	903584
<b>E0957</b>	903326	902857	900480	900137	902607
<b>E0956</b>	903327	902858	900482	900188	902608
<b>E0995</b>	904580	904580	904580	904580	904580

<u>Trunk Positioning</u>	<u>Price</u>
<input type="checkbox"/> H-Harness with Padded Covers (3-Pt Positioning Belt Req)	<b>\$147</b>
<input type="checkbox"/> Soft Adjustable Lateral Support Single Flap w/ Scoli Strap - Cordura ^	<b>\$140</b>
<input type="checkbox"/> Soft Adjustable Lateral Support Single Flap w/ Scoli Strap - Textilene ^	<b>\$140</b>
<input type="checkbox"/> Soft Adjustable Lateral Support Double Flap w/ Scoli Strap - Cordura ^	<b>\$140</b>
<input type="checkbox"/> Soft Adjustable Lateral Support Double Flap w/ Scoli Strap - Textilene ^	<b>\$140</b>
<input type="checkbox"/> Full Torso Support Vest - Cordura*	<b>\$270</b>
<input type="checkbox"/> Full Torso Support Vest - Textilene*	<b>\$270</b>
<input type="checkbox"/> Reducer Seat Insert (Reduces Seat Width 2" & Depth 1") ^	<b>\$362</b>

^ Not available with Planar Option.  
 \* Not available with Single or Double Flap Lateral Support.

HCPCS	CX10	CX12	CX14	CX16	CX18
<b>E0960</b>	903988	903988	903988	903988	903988
<b>E0956</b>	905446	905447	905448	905449	
<b>E0956</b>	903328	902859	900678	900140	902614
<b>E0956</b>	905450	905451	905452	905453	
<b>E0956</b>	903329	902860	900679	900204	902615
	903330	902861	900680	900339	902616
	905436	905437	905438	905439	902616
	904211	904211	904110	904110	904110

<u>Head Supports</u>	<u>Price</u>
<input type="checkbox"/> Padded Headwings ^	<b>\$250</b>
<input type="checkbox"/> Occi Headwings ^	<b>\$234</b>
<input type="checkbox"/> Headrest Extension - Curved - Cordura ^	<b>\$133</b>
<input type="checkbox"/> Headrest Extension - Rectangular - Textilene ^	<b>\$133</b>

^ Not available with Planar Option.

HCPCS	CX10	CX12	CX14	CX16	CX18
<b>E0955</b>	903331	902862	900681	900142	902618
<b>E0955</b>	905270	905270	905271	905271	905272
<b>E0966</b>	903332	902863	900484	900138	902619
<b>E0966</b>	904514	904515	904516	904517	904518

<u>Frame Option</u>	<u>Price</u>
<input type="checkbox"/> Anatomic Back (Heavy Duty Reinforced Upholstery Required) *1	<b>\$302</b>
<input type="checkbox"/> Two Piece Push Handle <sup>1</sup>	<b>\$0</b>

\*Not available with Planar Option.  
<sup>1</sup>Manufacturer Installation Required.

HCPCS	CX10	CX12	CX14	CX16	CX18
	904129	904129	904130	904131	903096
	904124	904125	904126	904127	904461

<u>Tire Options</u>	<u>Price</u>
<input type="checkbox"/> 12.5" Rear Solid Knobby Tire *	<b>\$204</b>
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire *	<b>\$306</b>
<input type="checkbox"/> 8" x 2" Front Ribbed, 12.5" Rear Solid Knobby Tire	<b>\$306</b>
<input type="checkbox"/> 8" x 2" Front Ribbed, 12.5" Rear Pneumatic Knobby Tire	<b>\$306</b>

\* Not available with Scout Option.

HCPCS	CX10	CX12	CX14	CX16	CX18
	904117	904117	904117	904117	
	904118	904118	904118	904118	
					903084
					903083

<u>Accessories</u>	<u>Price</u>
<input type="checkbox"/> Basic Headrest Cover (Canopy)	\$177
<input type="checkbox"/> Headrest Cover (Canopy)	\$234
<input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy)	\$256
<input type="checkbox"/> Extended Headrest Cover - with Windows (Canopy)	\$302
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray) ^	\$114
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$233
<input type="checkbox"/> Under Seat Storage Basket - Fixed Tilt	\$70
<input type="checkbox"/> Wheelchair Travel Bag	\$211
<input type="checkbox"/> Utility Bag	\$39
<input type="checkbox"/> Saddle Bags (pair)	\$122
<input type="checkbox"/> Curb Tipper (Elevates front of chair to go up curbs)****	\$100
<input type="checkbox"/> Caster Lock for 2" Wide Front Wheels <sup>1</sup>	\$96
<input type="checkbox"/> Attendant Hand Brake ** <sup>1</sup>	\$242
<input type="checkbox"/> Rear Anti-Tip Tubes **^	\$118
<input type="checkbox"/> Rain Cover Size 1 (Canopy Req)	\$154
<input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (CX14 Extended Canopy Req)	\$154
<input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req)	\$154
<input type="checkbox"/> Mosquito Net Size 1 (Canopy req)	\$154
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (CX14 Extended Canopy Req)	\$154
<input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req)	\$154

<u>HCPCS</u>	<u>CX10</u>	<u>CX12</u>	<u>CX14</u>	<u>CX16</u>	<u>CX18</u>
	905016	905017	905018	905019	905020
	903336	902867	900162	900184	902623
	904720	904721	904722	904723	904724
	904736	904737	904738	904739	904740
<b>E0950</b>	903078	903385	903386	903387	903388
<b>E1028</b>	903352	903353	903354	903355	903356
	904046	904046	904047	904047	904048
	902022	902022	902022	902022	902022
	904531	904531	904531	904531	904531
	904094	904094	904094	904094	904094
	904128	904128	904128	904128	903916
<b>K0073</b>		902920	902920	902920	902920
	904134	904134	904134	904134	
<b>E0971</b>	904520	904520	904540	904540	
	904939	904939	904939		
			904940	904940	904940
				905283	905283
	904655	904655	904655		
			904656	904656	904656
				905319	905319

\* Highly recommended for users weighing 50lbs or less.  
 ^ Not available with Planar Option or any Cushion.  
 \*\* Not available with 2-Piece Push Handle Package.  
 ^^Not available with Scout Option.  
 \*\*\*Cannot be combined with Rear Anti-Tip Tubes but can still be ordered.  
<sup>1</sup>Manufacturer Installation Required.

The Cruiser transit option includes four manufacturer installed wheelchair transportation anchors. The Cruiser transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are CX10 and CX12: 66 lbs./30kg; CX14: 100 lbs./45.5kg; CX16 and CX18: 170 lbs./77kg.

The Convoid Cruiser was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convoid is unknown, and no warranty or claim about their performance is being, or can be, made.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

**Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245**

Prices in USD & subject to change, ex works Torrance, California.

Effective: 5/1/2021

OF1900-ENG-REV.03

Convoid  
 2830 California Street  
 Torrance, CA 90503  
 1 844-US-MOBILITY (844-876-6245)  
 www.etac.us.com

