

Convaid Rodeo Tilt - HCPCS E1232/E1161

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

Rodeo Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- One Piece Height Adjustable Push Handle
- One Piece Angle Adjustable Footplate
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire
- Quick Release Rolko Rear Tires
- Single Foot Activated Wheel Lock Bar
- Zippered Pockets in Seat Bottom Upholstery
- 2-Point Positioning Belt
- Seat Depth Growth of 5"
- Elevating Legrest
- HCPCS code E1161 is only valid on size 16 (RD16)
- Tool-less Upholstery Removal
- Made in the USA

Standard Models:

	<u>Price</u>
<input type="checkbox"/> Rodeo RD10 - 904340	\$2,530
<input type="checkbox"/> Rodeo RD12 - 902689	\$2,530
<input type="checkbox"/> Rodeo RD14 - 902691	\$2,530
<input type="checkbox"/> Rodeo RD16 - 903134	\$2,530

Upholstery Colors

	<u>Price</u>
<input type="checkbox"/> Royal Blue	\$0
<input type="checkbox"/> Forest Green	\$0
<input type="checkbox"/> Cherry Red	\$0
<input type="checkbox"/> Sassy Purple	\$0
<input type="checkbox"/> Princess Pink	\$0
<input type="checkbox"/> Panther Black	\$0
<input type="checkbox"/> Chocolate Brown	\$0
<input type="checkbox"/> Aqua Blue	\$0
<input type="checkbox"/> Apple Green	\$0

RD10	RD12	RD14	RD16
903850	903850	903850	903850
903857	903857	903857	903857
903855	903855	903855	903855
903851	903851	903851	903851
903759	903759	903759	903759
903854	903854	903854	903854
904189	904189	904189	904189
904753	904753	904753	904753
904754	904754	904754	904754

Transit Required Items

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Transit Option (Manufacturer Installation Required)	\$205	
<input type="checkbox"/> Headrest Extension	\$133	E0966
<input type="checkbox"/> Headrest Extension Hardware (Required for Headrest Extension)	N/C	
<input type="checkbox"/> H-Harness with Padded Covers	\$147	E0960
<input type="checkbox"/> 3-Pt Positioning Belt	\$50	E0978
Foot Positioners - Select Size Below:		
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	\$160	K0108
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	\$160	K0108
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	\$160	K0108

RD10	RD12	RD14	RD16
904341	903167	903168	903169
903332	902863	900484	900138
905462	905462	905462	905462
903988	903988	903988	903988
904017	904017	904017	904017
903415	903415	903415	903415
903416	903416	903416	903416
903417	903417	903417	903417

Transit Options

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Q'Straint Transit Lap Belt	\$92	E0978
<input type="checkbox"/> Sure-Lok Transit Lap Belt	\$92	E0978
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Small ^	\$155	E0950
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Medium ^	\$155	E0950

^ Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.

RD10	RD12	RD14	RD16
903094	903094	903094	903094
903232	903232	903232	903232
903669	903669		
		903670	903670

Pelvic Positioning

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 3-Pt Positioning Belt	\$50	E0978
<input type="checkbox"/> Support - General Use Cushion *1	\$247	E2601
<input type="checkbox"/> Position Cushion - Medial Thigh Support/Anti Thrust *1	\$478	E2607
<input type="checkbox"/> Align Cushion - Lateral Pelvis/Lateral Thigh Support *1	\$478	E2607
<input type="checkbox"/> Incontinence Cushion Liner (Cushion req) 1	\$74	E2619

RD10	RD12	RD14	RD16
904017	904017	904017	904017
904347	904350	904353	904356
904348	904351	904354	904357
904349	904352	904355	904358
902324	902324	902324	902324

Foot & Leg Positioning

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Individual Angle Adjustable Footplates 1	\$183	K0040
<input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required)	\$333	K0108
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	\$160	K0108
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	\$160	K0108
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	\$160	K0108
<input type="checkbox"/> Legrest Extension - 2" *1	\$56	
<input type="checkbox"/> Medial Thigh Support (Abductor)	\$166	E0957
<input type="checkbox"/> Lateral Thigh Support (Adductor)	\$166	E0956
<input type="checkbox"/> Calf Panel	\$63	E0995
<input type="checkbox"/> Padded Elevating Legrest Cover	\$15	

RD10	RD12	RD14	RD16
	904957	904958	904959
904953	904954	904955	904956
903415	903415	903415	903415
903416	903416	903416	903416
903417	903417	903417	903417
904095	904095	904096	904097
904359	902667	902647	903140
904360	902668	902648	903141
904361	904361	904361	904361
904619	904619	904619	904619

* 2" Legrest Extension interferes with caster swivel when seat is in upright position.
 Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.
 1 Manufacturer installation required if model was made prior to June 2019.

<u>Trunk Positioning</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> H-Harness with Padded Covers (3-Pt Positioning Belt Req)	\$147	E0960
Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)		
<input type="checkbox"/> X Small= 2.75"x 3.75" *	\$121	E0956
<input type="checkbox"/> Small = 3"x 5"	\$121	E0956
<input type="checkbox"/> Medium = 3.75"x 5.5"	\$121	E0956
<input type="checkbox"/> Large = 4.75"x 6"	\$121	E0956
Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)		
<input type="checkbox"/> X Small= 2.75"x 3.75" *	\$121	E0956
<input type="checkbox"/> Small = 3"x 5"	\$121	E0956
<input type="checkbox"/> Medium = 3.75"x 5.5"	\$121	E0956
<input type="checkbox"/> Large = 4.75"x 6"	\$121	E0956
<input type="checkbox"/> Right S/A Adj. Hardware for Lateral Trunk Support	\$244	E1028
<input type="checkbox"/> Left S/A Adj. Hardware for Lateral Trunk Support	\$244	E1028
<input type="checkbox"/> Full Torso Support Vest	\$244	
<input type="checkbox"/> Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")	\$362	

* RD10 requires X Small laterals.

RD10	RD12	RD14	RD16
903988	903988	903988	903988
904385	904385	904385	904385
	904071	904071	904071
	904072	904072	904072
	904073	904073	904073
904465	904465	904465	904465
	904466	904466	904466
	904467	904467	904467
	904468	904468	904468
904362	904362	904362	904362
904363	904363	904363	904363
904364	903261	903260	903259
904211	904211	904110	904110

<u>Head Supports</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Padded Headwings	\$250	E0955
<input type="checkbox"/> Occi Headwings	\$234	E0955
<input type="checkbox"/> Headrest Extension	\$133	E0966
<input type="checkbox"/> Headrest Extension Hardware (Required for Headrest Extension)	N/C	

RD10	RD12	RD14	RD16
904365	902673	902653	903147
905270	905270	905271	905271
903332	902863	900484	900138
905462	905462	905462	905462

<u>Frame Option</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Recline Option (Requires Manufacturer Installation)	\$473	E1225
<input type="checkbox"/> Height Adjustable Flip-Up Armrests (pair)	\$234	E0973
<input type="checkbox"/> Solid Side Panels	\$216	
<input type="checkbox"/> Mesh Side Panels	\$152	

RD10	RD12	RD14	RD16
904342	904343	904344	904345
904346	904346	904346	904346
904388	902672	902652	903145
904389	902680	902681	903146

<u>Tire Options</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 12.5" Rear Solid Knobby Tire	\$204	
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire	\$306	

RD10	RD12	RD14	RD16
904117	904117	904117	904117
904118	904118	904118	904118

<u>Accessories</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Headrest Cover (Canopy)	\$234	
<input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy)	\$256	
<input type="checkbox"/> Extended Headrest Cover - with Windows (Canopy)	\$302	
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray) ^	\$114	E0950
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$233	E1028
<input type="checkbox"/> Medical Necessity Storage Basket	\$302	
<input type="checkbox"/> Under Seat Storage Basket - Tilt-In-Space	\$77	
<input type="checkbox"/> IV Pole (height adjustable and collapsible) ²	\$154	K0105
<input type="checkbox"/> Oxygen Tank Bag	\$183	E2208
<input type="checkbox"/> Utility Bag	\$39	
<input type="checkbox"/> Attendant Hand Brakes (pair) ¹	\$242	
<input type="checkbox"/> Curb Tipper	\$100	
<input type="checkbox"/> Caster Lock for 2" Wide Front Tires ¹	\$96	K0073
<input type="checkbox"/> Rear Anti-Tip Tubes (Pair)*	\$118	E0971
<input type="checkbox"/> Rain Cover Size 1 (Canopy Req)	\$154	
<input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (RD14 Extended Canopy Req)	\$154	
<input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req)	\$154	
<input type="checkbox"/> Mosquito Net Size 1 (Canopy req)	\$154	
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (RD14 Extended Canopy Req)	\$154	
<input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req)	\$154	

RD10	RD12	RD14	RD16
904366	902679	902659	903152
904728	904729	904730	904731
904744	904745	904746	904747
904367	903405	903401	903402
903369	903369	903369	903369
904476	904476	904477	904477
904529	904529	904530	904530
904566	904566	904566	904566
902870	902870	902870	902870
904531	904531	904531	904531
		904370	904371
905310	905310	905310	905310
904650	904650	904650	904650
904649	904649	904649	904649
904939	904939	904939	
		904940	904940
			905283
904655	904655	904655	
		904656	904656
			905319

^ Flip-Up Armrests and Hardware required.

*Highly recommended for users weighing 50lbs or less.

¹ Manufacturer Installation Required.

² Anti-tip Tubes Required

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo mobility base was tested with Rodeo seating module equipped with head support, H-Harness with padded Covers, 3-point restraint system, foot positioners, and 4 transit anchor brackets.

It is Convaid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specially designed for the wheelchair to be secured in a forward facing configuration and not in any other direction.

For more information on RESNA WC-19 and WC-20 standards please visit <http://www.ercwts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166

Prices in USD & subject to change, ex works Torrance, California.

Convaid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

Effective: 5/1/2021

OF0100-ENG-REV.05

