

R82 Cricket

CAD Retail Price List/ Order Form

Account #: _____ Purchase Order No.: _____ Date: _____


Contact: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

US Configuration


Effective: 05.01.21 - Rev 5

- Cricket
- Cushion
- Arm rests
- Head support
- Tray
- Grip bar
- Folding top
- Transport fittings



- A (Back Height) _____
- B (Arm Height) _____
- C (Lower Leg Length) _____
- D (Seat Depth) _____
- E (Chest Width) _____
- F (Seat Width) _____
- G (Height from Seat to Top of Head) _____
- H (Shoulder Width) _____
- Age _____
- Weight _____

Measurements



Size 1 Inches

Effective seat width (A)	13½"
Effective seat depth (B)	8½"-12¼"
Backrest height (C1)	19¾"
Seat surface height at front edge (E)	19¾"-21"
Footrest to seat distance (D)	7-14¼"
Leg to seat surface angle	360°
Backrest width, outside (F1)	13½"
Backrest width, inside (F)	10¾"
Overall width (G)	20¾"
Length (H)	29½"
Overall height (I)	39¼"
Folded width	14"
Folded height	19"
Folding length	30"
Backrest angle	110°-125°
Minimum turning radius	27½"
Total mass	27.3 lbs.
Max load	77 lbs.
Max load at transportation	59½ lbs.
Required width of right angled corridor	32"
Static stability (sideways, downhill, uphill)	10°

Size 2 Inches

Effective seat width (A)	15½"
Effective seat depth (B)	12¼"-17"
Backrest height (C1)	22½"
Seat surface height at front edge (E)	21¼"-22¼"
Footrest to seat distance (D)	8¾"-14½"
Leg to seat surface angle	360°
Backrest width, outside (F1)	15½"
Backrest width, inside (F)	13
Overall width (G)	22¾"
Length (H)	29½"
Overall height (I)	41¼"
Folded width	14"
Folded height	19"
Folding length	30"
Backrest angle	110°-125°
Minimum turning radius	30"
Total mass	29½ lbs.
Max load	110 lbs.
Max load at transportation	94½ lbs.
Required width of right angled corridor	32"
Static stability (sideways, downhill, uphill)	10°


Notes:

Cricket Standard Features



- Back Recline
- Adjustable Seat Depth
- Height and Angle Adjustable Foot Support
- Armrest, Fixed (Pr.)

Description	Part No.	CAD Price
<input type="checkbox"/> Cricket sz 1	87043100	\$2,527
<input type="checkbox"/> Cricket sz 2	87043200	\$2,443



Description	Size	Part No.	CAD Price
Required, must select a color			
<input type="checkbox"/> Blue = 1		<input type="checkbox"/> Red = 2	
<input type="checkbox"/> Pink = 5		<input type="checkbox"/> Grey = 4	

<input type="checkbox"/> Cushion	1	8704311x	\$0
<input type="checkbox"/> Cushion	2	8704321x	\$0
◆ Colors Green, Orange & Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item			


<input type="checkbox"/> Seat Side Parts ◆	One Size	87043130	\$67
◆ Required with Seat Width Reduction Kit			



Description	Measurements	Size	Part No.	CAD Price
<input type="checkbox"/> Transit Option		One Size	87043188	\$256




Description	Measurements	Size	Part No.	CAD Price
<input type="checkbox"/> Armrest, Fixed (Pr.) ◆		1	87043120-1	\$0
<input type="checkbox"/> Armrest, Fixed (Pr.) ◆		2	87043120-2	\$0
◆ Required with Grip Bar, In Armrests				




Description	Measurements	Size	Part No.	CAD Price
<input type="checkbox"/> Seat Width Reduction Kit, Wedges (Pr.) ◆			87043130-1	\$92
<input type="checkbox"/> Seat Width Reduction Kit, Wedges (Pr.) ◆			87043130-2	\$92
◆ Requires Seat Side Parts (87043130)				
















<input type="checkbox"/> Wedges for Shaping Back (Pr.)	4½" x 6" x 11"	Type 1	87043140-1	\$92
<input type="checkbox"/> Wedges for Shaping Back (Pr.)	5¾" x 6" x 11"	Type 2	87043140-2	\$92








<input type="checkbox"/> Chest Belt, Upholstery (Ea.)	A:11¼ x B:7½ x C:9¼-12"	1	85428-1	\$345
<input type="checkbox"/> Chest Belt, Upholstery (Ea.)	A:14 x B:10 x C:10-12¾"	2	85428-2	\$345



<input type="checkbox"/> Torso vest (Ea.)	S-P: 7¾-10¾"	XS	85705-XS	\$256
<input type="checkbox"/> Torso vest (Ea.)	S-P: 10¾-12½"	S	85705-S	\$256
<input type="checkbox"/> Torso vest (Ea.)	S-P: 12½-16"	M	85705-M	\$256
<input type="checkbox"/> Torso vest (Ea.)	S-P: 16-18½"	L	85705-L	\$262

Accessories	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> 2-Point Pelvic Belt, Black	L: 13-14½"	1	85424-1	\$124
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 6¾-12½"	XS	85725-XS	\$173
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 8¼-13¾"	S	85725-S	\$173
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 9¾-16¼"	M	85725-M	\$172
	<input type="checkbox"/> Fix-Locks, Black, Bag of 6 Pcs ♦ ♦Required for use of Vests		One Size	862109-35	\$54
<hr/>					
	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Pelvic Harness	A:6½ x B:7½"	0	89108-0	\$314
	<input type="checkbox"/> Pelvic Harness	A:8½ x B:8¾"	1	89108-1	\$314
	<input type="checkbox"/> Pelvic Harness	A:10½ x B:9¾"	2	89108-2	\$314
	<input type="checkbox"/> Heel Straps (Pr.)		One Size	87043185	\$69
	<input type="checkbox"/> Upholstered Ankle Straps	A: 5½-7 x B:8½-10"	1	85426-1	\$236
	<input type="checkbox"/> Upholstered Ankle Straps	A: 6¼-7¾ x B:10-11¾"	2	85426-2	\$236
Other types of vests and belts available, please see Vest and Belts Order Form					
<hr/>					
	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Headrest Cushion, Black (Ea.)		One Size	87043150	\$256
	<input type="checkbox"/> Sculpted Headrest Cushion, Black (Ea.)		One Size	87043151	\$262
	<input type="checkbox"/> Backrest Height Extension (Ea.)		1	87043160-1	\$70
	<input type="checkbox"/> Backrest Height Extension (Ea.)		2	87043160-2	\$70
<hr/>					
	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Calf Support Strap (Ea.)		1	87043190-1	\$78
	<input type="checkbox"/> Calf Support Strap (Ea.)		2	87043190-2	\$78
	<input type="checkbox"/> Additional Footplate Locking Mechanism (Ea.)		One Size	87043186	\$38
	<input type="checkbox"/> Grip Bar, In Armrests ♦		1	87043181-1	\$220
	<input type="checkbox"/> Grip Bar, In Armrests ♦		2	87043181-2	\$220
♦Requires Armrest					
	<input type="checkbox"/> Grip Bar, In Frame		1	87043182-1	\$220
	<input type="checkbox"/> Grip Bar, In Frame		2	87043182-2	\$220
	<input type="checkbox"/> Caster Locks		One Size	87043187	\$120

Accessories	Description	Measurements	Size	Part No.	CAD Price
	Required, must select a color				
	x = Color				
	<input type="checkbox"/> Blue = 1	<input type="checkbox"/> Red = 2			
	<input type="checkbox"/> Pink = 5	<input type="checkbox"/> Grey = 4			
	<input type="checkbox"/> Folding Top (Canopy) (Ea.) ♦		1	8704317x	\$292
	<input type="checkbox"/> Folding Top (Canopy) (Ea.) ♦		2	8704327x	\$292
	♦ Required w/ Rain Cover				
	<input type="checkbox"/> Underseat Storage Basket ♦		1	87043184-1	\$89
	<input type="checkbox"/> Underseat Storage Basket ♦		2	87043184-2	\$89
	♦ 22 lbs. weight limit				
	<input type="checkbox"/> Tray ♦		1	87043180-1	\$432
	<input type="checkbox"/> Tray ♦		2	87043180-2	\$432
	♦ Must select same size as armrest				
	<input type="checkbox"/> Storage Strap		One Size	85100	\$38
	<input type="checkbox"/> Rain Cape with Hood, Clear	A:11 x B:29½"	4	87465-410	\$179
	<input type="checkbox"/> Rain Cape with Hood, Clear	A:11½ x B:34½"	5	87465-510	\$179
	<input type="checkbox"/> Rain Cover, Clear (Ea.) ♦		2	87467-210	\$147
	<input type="checkbox"/> Rain Cover, Clear (Ea.) ♦		3	87467-310	\$147
	♦ Requires Folding Top (Canopy)				

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

The Cricket complies with ISO7176-19, for safe transportation in motor vehicles. The R82 Cricket was tested with the Cricket back rest height extension and Y-hip belt and R82 transit anchors. The performance and safety of accessories and/or other seating options not tested by R82 is unknown, and no warranty or claim about their performance is or can be made.

When being transported in a vehicle, it is R82's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair; therefore R82 cannot recommend a transit system.

Prices in CAD are subject to change without notice
For Quotes and Orders please send completed Order Form to: sales.us@etac.com or fax: 1-800-760-0537

Etac
13137 Bleinheim Lane
Matthews, North Carolina 28105
+1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF5102-ENG-REV.05

Creating Possibilities