



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

# Heron

Company:	Ship To:
Account Number:	Date:
Quote: PO Number:	Address:
ATP:	City: State:
Marked For:	Zip: Phone:
Associated Quote/Order#:	Email/Fax:
	Effective: 3/13/2026

## Heron

### Heron Base Frame



#### Heron Base Frame

Heron Base Frame (880505-11)	\$5728
Heron Power Frame (880505-11EL)	\$6995
<i>Charger and Cable for Charger required for only Heron Power</i>	

#### Heron Power Charger and Cable

*Only compatible with Heron Power. Required selection for Heron Power.*

Heron Charger (880510-01)	N/C
Cable for Charger (880510-UL)	N/C

## Armrests and Trays

### Armrests Options



#### Flip Back Arms

*Required Selection*

Armrests, Flip Back H: 9.75" x L: 14.5" x W: 2" (Pair) (X880508-03)	\$625
---	-------



#### Armrest Mounted Tray

Tray, 17¼ x 7½" (880508-04)	\$425
-----------------------------	-------

## Seat and Back Options

### Seat and Back Padding and Seating Options



#### Upholstery

Pads, Complete (880507-01) (Set)	\$328
Pad, Back Only (101204) (Each)	\$163
Pad, Seat Only (117005) (Each)	\$181



#### Seat Options

*Reduces depth by 2"*

Seat Depth Reduction (880508-06)	\$57
----------------------------------	------

## Headrest and Hardware

### Headrests

#### Headrest Pad

Headrest, Grey, Size 2 A: 4" x B: 6.24" x C: 4.25" (6172000)	\$228
Headrest, Grey, Size 3, A: 4" x B: 5" x C: 5.75 x D: 3" (6173000)	\$228
Headrest, Grey, Size 4, A: 4" x B: 5.75" x C: 2.75 (6174000)	\$228

#### Headrest Hardware

Stainless Turtle Bar System, Long (98619-28)	\$256
--	-------

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0467.1 Pricing and specifications subject to change without notice.



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

# Heron

## Trunk and Hip Supports

### Shoulder and Hip Supports



Upholstery for Sides/Wings (Pair) (880507-04) \_\_\_\_\_ \$110



Sides/Wings (Pair) (880508-01) \_\_\_\_\_ \$93



Hip Supports, 8 $\frac{1}{4}$ "W x 4"H (Pair) (880508-02) \_\_\_\_\_ \$351



Hip Support Spacer Blocks (880508-02LB) \_\_\_\_\_ \$85

## Lateral Support Cushion and Fittings

### Swing Away Lateral Support

*Swing Away Bracket, Bars, and Curved Cushion required to make complete Swing Away Lateral Supports*

#### Swing Away Lateral Bracket - Left

Swing-Away Lateral Support Bracket, Left (Each) (E1028) (83000-L) \_\_\_\_\_ \$267

#### Swing Away Lateral Bracket - Right

Swing-Away Lateral Support Bracket, Right (Each) (E1028) (83000-R) \_\_\_\_\_ \$267

### Pad Mounting Bar

*Required selection for Swing Away Laterals. Quantity 2 required to mount left and right sides.*

#### Straight Bars

Bar, Straight (Each) (82648) \_\_\_\_\_ N/C

#### Offset Bars

Offset Bar 1" Stainless (150238) (Each) \_\_\_\_\_ N/C

### Fixed Side Supports

Fixed Side Support Pad 6" x 6" (Pair) (880508-09) \_\_\_\_\_ \$250

### Lateral Cushions

#### Curved Lateral Cushion

Curved Lateral Pad, 3.5" x 7.5" (Each) (E0956) (6175407) \_\_\_\_\_ \$125

Curved Lateral Pad, 2.5" x 6.75" (Each) (6175408) (E0956) \_\_\_\_\_ \$125

Curved Lateral Pad, 2.5" x 6" (Each) (6175409) (E0956) \_\_\_\_\_ \$125

#### Protective Cushion Cover

*Covers the Swing Away Hardware Bar*

Protective Cushion for Laterals, Black (89151) (Set) \_\_\_\_\_ \$112

## Foot Supports and Accessories

### Foot Supports



Foot Support, Right-Angled Knee, 12"-15" Short (880506-03) (Pair) \_\_\_\_\_ \$625

Foot Support, Right-Angled Knee, 15.4" - 21.7" Long (880506-01) (Pair) \_\_\_\_\_ \$625

Foot Support, Open Knee Angle, 14" - 21" (880506-02) (Pair) \_\_\_\_\_ \$625



Padding for Foot Supports (880506-12) (Pair) \_\_\_\_\_ \$85

### Foot Support Accessories



Calf Support, Strap, 4" (880509-01) (Each) \_\_\_\_\_ \$90



Heel Stops, Straps (880509-02) (Pair) \_\_\_\_\_ \$131



Phone (877) 227-0097 | Fax (704) 882-0751  
 sales.us@etac.com


# Heron

## Vest and Belts

### Vests and Hip-Belts



Vest, Size 3, 2.75" x 4.5" x 11" x 5.75" x 18"  
 (89111-340) \_\_\_\_\_ \$219  
 Vest, Size 4, 2.75" x 6.5" x 13.5" x 7.5" x 19"  
 (89111-440) \_\_\_\_\_ \$219



Fix-Locks, White, Bag of 4 Pieces (890505-36) \_\_\_\_\_ \$46



Cross Vest, Water-Proof Mesh, Size 3 S-P:  
 17¼-20¼" (89625-340) \_\_\_\_\_ \$220  
 Cross Vest, Water-Proof Mesh, Size 4 S-P: 20¼ -  
 24 ½" (89625-440) \_\_\_\_\_ \$220



Neoprene Hip Belt, 1", Size 2 12 x 3 x 18"  
 (85423-2) \_\_\_\_\_ \$107

## Accessories

### Bathroom Accessories



Commode Pan (880130) \_\_\_\_\_ \$57



Lid for Commode Pan (880130-L) \_\_\_\_\_ \$73



Splash Guard (880508-05) (One Size) \_\_\_\_\_ \$79  
 Splash Guard, Low (880508-07) \_\_\_\_\_ \$79



Potty Ring for Heron Power (128623) \_\_\_\_\_ \$51  
 Potty Ring, Standard (109072) \_\_\_\_\_ \$51

### Accessories

Shock Absorber, White (116924-KIT) \_\_\_\_\_ \$43



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

# Heron

## Measurements

mm (inch) kg (lb)	Hydraulic	Power
 Overall width	670 mm (26½")	670 mm (26½")
 Length	710 mm (28")	710 mm (28")
 Overall height	1080–1350 mm (44½–53¼")	1130–1350 mm (44½–53¼")
 Effective seat width	550 mm (21¾")	550 mm (21¾")
Seat width between hip supports (accessories)	350–490 mm (13¾ - 19¼")	350–490 mm (13¾–19¼")
 Effective seat depth	450 mm (17¾")	450 mm (17¾")
 Back width, top of backrest	290 mm (11½")	290 mm (11½")
 Backrest height	555–625 mm (21¾–24½")	555–625 mm (21¾–24½")
 Seat surface height at front edge	510–780 mm (20–30¾")	560–780 mm (22–30¾")
 Seat plane angle	(-5°)–(25°)	(-5°)–(25°)
 Backrest angle	(-10°)–(30°)	(-10°)–(30°)
 Lower leg length	400–560 mm (15¾–22")	400–560 mm (15¾–22")
 Minimum turning diameter	900 mm (35½")	900 mm (35½")
 Mass device	28 kg (61.8 lb)	31 kg (68.3 lb)
 Maximum user mass	100 kg (220 lb)	130 kg (286.6 lb)
 Maximum rated load	130 kg (286.6 lb)	163 kg (359.3 lb)
 Static stability sideways	5°	5°
 Static stability uphill and downhill	10°	10°