

# R82 Stingray - E1233

## USD Retail Price List/ Order Form

Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

### US Configuration


**Effective: 05.27.21 - Rev 8**

- ☐ Frame (Stingray, Trekker, High-Low:v)
- ☐ Stingray Seat
- ☐ Cushions
- ☐ Head Support
- ☐ Side Support
- ☐ Hip Support
- ☐ Pommel
- ☐ Vests and Belts
- ☐ Tray
- ☐ Transport Fittings



A (Back Height)	
B (Arm Height)	
C (Lower Leg Length)	
D (Seat Depth)	
E (Chest Width)	
F (Seat Width)	
G (Height: Seat to Top of Head)	
H (Shoulder Width)	
Age	
Weight	

### Measurements



	Size 1 Inches	Size 2 Inches
Effective Seat Width (A)	7-11 $\frac{3}{4}$ "	9-13 $\frac{3}{4}$ "
Seat Width Between Knee Supports	11"	12 $\frac{3}{4}$ "
Seat Width Between Hip Supports (Accessories)	7", 8 $\frac{1}{2}$ ", 10 $\frac{1}{4}$ "	9", 10 $\frac{1}{2}$ ", 12"
Effective Seat Depth (B)	7 - 11 $\frac{3}{4}$ "	9 $\frac{1}{4}$ - 15 $\frac{1}{2}$ "
Back Height at Shoulders (C)	10 $\frac{1}{2}$ - 16"	14 - 19 $\frac{1}{4}$ "
Backrest Height (C1)	16 - 21 $\frac{1}{2}$ "	19 $\frac{1}{2}$ - 25 $\frac{1}{4}$ "
Footrest to Seat Distance (D)	5 $\frac{1}{2}$ - 17 $\frac{1}{2}$ "	5 $\frac{1}{2}$ - 17 $\frac{1}{2}$ "
Seat Surface Height at Front Edge (E)	20 $\frac{1}{4}$ "	20 $\frac{1}{4}$ "
Back Width, Inside (F)	11 $\frac{3}{4}$ "	14"
Back Width, Outside (F1)	17 $\frac{1}{2}$ "	20 $\frac{1}{4}$ "
Overall Width (G)	23"	23"
Length with Turnable Front Wheels (H1)	28 $\frac{3}{4}$ "	28 $\frac{3}{4}$ "
Length with Large Front Wheels	32 $\frac{3}{4}$ "	32 $\frac{3}{4}$ "
Overall Height, Push Brace 15° (I)	39", 40", 41 $\frac{1}{4}$ "	39", 40", 41 $\frac{1}{4}$ "
Foot Plate, Width x Depth	12" x 7 $\frac{3}{4}$ "	12" x 7 $\frac{3}{4}$ "
Folded Length	35"	35"
Folded Height	19"	19"
Total Mass, Frame	29 $\frac{1}{2}$ lbs.	29 $\frac{1}{2}$ lbs.
Total Mass, Seat	12 lbs.	14 $\frac{1}{2}$ lbs.
Max Load	88 lbs.	88 lbs.
Back Recline / Seat Plane Angle	Back: 170° / Seat: 45°	Back: 170° / Seat: 45°
Approx. Minimum User Height	29 $\frac{1}{2}$ "	35 $\frac{1}{2}$ "
Approx. Maximum User Height	36"	47 $\frac{1}{2}$ "

NOTES:

Stingray Standard Features



- Tilt-in-Space - 45°
- Back Recline - 170°
- Turnable Seat - 180°
- One-Size Frame
- Accepts Two Seat Sizes for Maximum Growth Potential
- Adjustable Seat Depth
- Adjustable Back Height
- Suspension

Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> Solid Tires, Individual Wheel Locks	8724212	\$2,506	E1233
<input type="checkbox"/> Solid Tires, Individual Wheel Locks & Hand Brakes ••	8724211	\$2,891	E1233
<input type="checkbox"/> Solid Tires, Central Wheel Locks	8724222	\$2,609	E1233
<input type="checkbox"/> Solid Tires, Central Wheel Locks & Hand Brakes	8724221	\$2,993	E1233
<input type="checkbox"/> Solid Tires, Individual Wheel Locks & Quick Release Rear Axels	8722212	\$2,762	E1233
••Custom order; Approx. lead time: 6-8 weeks. Non-returnable item			

Trekker Frame Standard Features



- Tilt-in-Space (-10° to +45°)
- 180° Reversible Seating Configuration
- Single Foot Activated Wheel Lock Bar
- 7½" x 2" Front and 11½" x 2½" Rear Solid Tires
- Quick Release Rear Wheels
- Made in the USA
- Transit Option - **Transit Ready Frame Only**


Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> Trekker Frame (TR)	C904661	\$2,516	E1234
<input type="checkbox"/> Trekker Frame - Transit Ready (TR)	C905362	\$2,726	E1234
<input type="checkbox"/> Stingray Seat Module - Interface Adaptor Kit for Trekker Base ♦	C905363	N/C	
♦Required with Trekker Frame			


Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> Transit Option (TR)	C904576	\$210	
<input type="checkbox"/> Using other manufacturer's Head Support	C995681	N/C	
♦Requires Head Support & Pelvic Positioning Belt			

Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> Height Adjustable Push Handle w/Outside Mount (TR)	C904688	N/C	
^Trekker Frame comes standard with push buttons mounted inside Height Adjustable Push Handle			

Description	Part No.	USD Price	HCPCS
<input type="checkbox"/> 12½ x 2" Rear Solid Knobby Tire (TR) ••	C904117	\$209	
<input type="checkbox"/> 7½ x 2" Front & 11½ x 2½" Rear Pneumatic Knobby Tires (TR) ••	C904118	\$314	
••Custom order; Approx. lead time: 6-8 weeks			

Trekker Frame Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> IV Pole (height adj. & collapsible) (TR) Anti-Tipper Required		One Size	C904566	\$210	K0105
	<input type="checkbox"/> O2 Tank Holder (TR) *°		One Size	C904565	\$184	E2208
	<input type="checkbox"/> LTV Ventilator Holder (TR) °		One Size	C904671	\$473	E1029
	<input type="checkbox"/> Ventilator/Suction/Hard Tray (TR)		One Size	C904595	\$473	E1029
	<input type="checkbox"/> Curb Tipper (TR)		One Size	C905310	\$103	
	<input type="checkbox"/> Rear Anti-Tipper (TR)		One Size	C904649	\$121	E0971
	<input type="checkbox"/> Under Seat Storage Basket (TR)		One Size	C904605	\$126	
	<input type="checkbox"/> Medical Necessity Storage Basket (TR)		One Size	C904606	\$310	
	* <b>Cannot be used during transit</b>					
	° <b>Manufacturer Installation Required</b>					

High-Low Frame (Trekker)	Description	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> High-Low:x sz 1, Manual Lift	1,2	8911603-02	\$1,620	
	<input type="checkbox"/> Push Brace - Fixed Angle	One Size	8910545	\$164	
	<input type="checkbox"/> Push Brace - Angle Adjustable	1,2	8910555	\$395	











High-Low Frame	Description	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> High-Low:v Indoor Frame	1	8911871-20	\$1,814	
	<input type="checkbox"/> High-Low:v Indoor Frame	1, 2	8911872-20	\$2,248	
	<input type="checkbox"/> Push Brace - Fixed Angle	One Size	8910545	\$164	
	<input type="checkbox"/> Push Brace, Angle Adjustable	One Size	8910555	\$395	


Stingray Seat	Description	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Seat and Back size 1	1	8711120		
	<input type="checkbox"/> Seat size 1	1	x8711000	\$420	E2292
	<input type="checkbox"/> Back size 1	1	x8712000	\$518	E2291
	<input type="checkbox"/> Seat and Back size 2	2	8712220		
	<input type="checkbox"/> Seat size 2	2	x8721000	\$420	E2292
	<input type="checkbox"/> Back size 2	2	x8722000	\$518	E2291
	<input checked="" type="checkbox"/> Fully Reclining Back (Required)	One Size	x8702000	\$415	E1014

*Adjustable seat depth and height adjustable back*


Description	Size	Part No.	USD Price	HCPCS
<b>Required, must select a color</b>				
<b>X = Color</b>				
<input type="checkbox"/> Blue = 1	<input type="checkbox"/> Red = 2			
<input type="checkbox"/> Pink = 5	<input type="checkbox"/> Grey = 4			
<input type="checkbox"/> Seat and Back Cushion (Pr.)	1	875051x	N/C	
<input type="checkbox"/> Seat and Back Cushion (Pr.)	2	875052x	N/C	
*Colors Green, Orange & Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item				

Description	Measurements	Size	Part No.	USD Price	HCPCS
<input type="checkbox"/> Footrest Rod, Elevating ♦	8¼"	Short	x8711100	\$113	K0108
<input type="checkbox"/> Footrest Rod, Elevating ♦	11"	Long	x8712200	\$113	K0108
♦ <b>Required with Angle Adjustable Footplate (x8710010)</b>					
<input type="checkbox"/> Angle Adjustable Footplate ♦		One Size	x8710010	\$97	K0040
♦ <b>Requires Footrest Rod Short or Long (x8711100 or x8712200)</b>					








Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Transit Option		One Size	8720902	\$205	
	<b>Description</b>	<b>Measurements</b>	<b>Size</b>	<b>Part No.</b>	<b>USD Price</b>	<b>HCPCS</b>
	<input type="checkbox"/> Hip Supports, Black (Pr.)	6½" x 7"	1	8740230	\$231	E0956
	<input type="checkbox"/> Hip Supports, Black (Pr.)	7" x 8½"	2	8740240	\$231	E0956
	<input type="checkbox"/> Hip Supports, Pad, Black (Pr.) ♦		1	8740250	\$231	E0956
	<input type="checkbox"/> Multi-Adjustable, Removable Hardware (Pr.) (Required) ♦		1	x8740250	\$231	E1028
	<input type="checkbox"/> Hip Supports, Pad, Black (Pr.) ♦		2	8740260	\$231	E0956
	<input type="checkbox"/> Multi-Adjustable, Removable Hardware (Pr.) (Required) ♦		2	x8740260	\$231	E1028
	♦ Hip Support and Multi-Adjustable Hardware must be ordered together					
	<input type="checkbox"/> Side Supports, Multi Adjustable (Pr.)	4" x 4"	1	8720301	\$226	E0956
	<input type="checkbox"/> Side Supports, Multi Adjustable (Pr.)	5" x 5¾"	2	8720302	\$226	E0956
	<b>Description</b>	<b>Measurements</b>	<b>Size</b>	<b>Part No.</b>	<b>USD Price</b>	<b>HCPCS</b>
	<input type="checkbox"/> Covers for Knobs on Side Supports (Pr.) ●●		One Size	8720360	\$31	
	●● Custom order; Approx. lead time: 6-8 weeks. Non-returnable item					
	<input type="checkbox"/> Medial Thigh Support (Ea.) ♦	3½" x 1½" x 2¾"	1	8720501	\$138	E0957
	<input type="checkbox"/> Adjustable, Removable Hardware (Required) ♦		1	x8720501	\$169	E1028
	<input type="checkbox"/> Medial Thigh Support (Ea.) **	3½" x 3½" x 3½"	2	8720502	\$138	E0957
	<input type="checkbox"/> Adjustable, Removable Hardware (Required) ♦		2	x8720502	\$169	E1028
	<input type="checkbox"/> Medial Thigh Support (Ea.) **	3½" x 3½" x 4¾"	3	8720503	\$138	E0957
	<input type="checkbox"/> Adjustable, Removable Hardware (Required) ♦		3	x8720503	\$169	E1028
	♦ Medial Thigh Support and Adjustable, Removable Hardware must be ordered together					
	<input type="checkbox"/> Anti-Thrust Wedge, Black (Pr.)	4" x 6"	1, 2, 3	8910405-3	\$47	K0108
	<b>Description</b>	<b>Measurements</b>	<b>Size</b>	<b>Part No.</b>	<b>USD Price</b>	<b>HCPCS</b>
	<input type="checkbox"/> Contoured Headrest, Black		1	8740810	\$205	E0955
	<input type="checkbox"/> Contoured Headrest, Black		2	8740820	\$205	E0955
	<input type="checkbox"/> Sculpted Head Support ①	9½" x 4½"	Small	8720831	\$210	E0955
	<input type="checkbox"/> Sculpted Head Support ①	11" x 5¾"	Large	8720832	\$210	E0955
	① Preferred option for Stingray seating					
	<input type="checkbox"/> Anatomic Head Support ♦	A:5½ x B:5½ x C:9½"	Small	5099001	\$261	E0955
	<input type="checkbox"/> Head Support Cover (Required)		Small	9923-83	N/C	
	<input type="checkbox"/> Anatomic Head Support ♦	A:6 x B:5½ x C:10"	Medium	5099002	\$261	E0955
	<input type="checkbox"/> Head Support Cover (Required)		Medium	9945-83	N/C	
	<input type="checkbox"/> Anatomic Head Support ♦	A:6½ x B:5½ x C:12½"	Large	5099003	\$261	E0955
	<input type="checkbox"/> Head Support Cover (Required)		Large	9967-83	N/C	
	♦ Turtle Bar System Short or Long Required; 98619-25 or 98619-26					
	<b>Swing-away Mounting Brackets - Height, Depth, Angle and Lateral Adjustable</b>					
	<input type="checkbox"/> Turtle Bar System ♦	A: 5", B: 4-8½"	Short	98619-25	\$236	E1028
	<input type="checkbox"/> Turtle Bar System ♦	A: 10", B: 4-8½"	Long	98619-26	\$236	E1028
	♦ Required with Anatomic Head Support; 5099001, 5099002 or 5099003					

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> 2-Part Headrest		One Size	8720870	\$277	E0955
	<input type="checkbox"/> Headrest Fittings ♦ ♦ <b>Required with Turtle Bar System or Universal Headrest Mount</b>		One Size	8720850	\$52	K0108
	<input type="checkbox"/> Universal Headrest Mount ♦ <b>Required when using other manufactures Head Support</b> <b>Allow most Whitmyer, Stealth, Otto Bock headrest, I2I to mount to our products (x:panda, Serval, etc.). Please contact Customer Service for confirmation</b>		One Size	87590	\$63	E1028







Other types of headrest available, please see Head Supports Order Form




Description	Measurements	Size	Part No.	USD Price	HCPCS
					




■ XX = 60 (Black)

	<input type="checkbox"/> Cross Vest (Ea.)	S-P: 7 <sup>3</sup> / <sub>4</sub> –12 <sup>1</sup> / <sub>2</sub> "	XS	85720-XS	\$195	E0960
	<input type="checkbox"/> Cross Vest (Ea.)	S-P: 10 <sup>1</sup> / <sub>4</sub> –13 <sup>1</sup> / <sub>2</sub> "	S	85720-S	\$195	E0960
	<input type="checkbox"/> Cross Vest (Ea.)	S-P: 13 <sup>1</sup> / <sub>2</sub> - 17 <sup>3</sup> / <sub>4</sub> "	M	85720-M	\$195	E0960
	<input type="checkbox"/> Cross Vest (Ea.)	S-P: 17 <sup>3</sup> / <sub>4</sub> –20 <sup>3</sup> / <sub>4</sub> "	L	85720-L	\$195	E0960
	<input type="checkbox"/> Cross Vest (Ea.)	S-P: 20 <sup>3</sup> / <sub>4</sub> - 24 <sup>1</sup> / <sub>2</sub> "	XL	85720-XL	\$195	E0960
	<input type="checkbox"/> Vest (Ea.)	1 <sup>1</sup> / <sub>4</sub> " x 2 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 3 <sup>1</sup> / <sub>2</sub> " x 11"	0	89111-0XX	\$205	E0960
	<input type="checkbox"/> Vest (Ea.)	1 <sup>1</sup> / <sub>4</sub> " x 2 <sup>1</sup> / <sub>2</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 4 <sup>1</sup> / <sub>2</sub> " x 12"	<sup>3</sup> / <sub>4</sub>	89111-075XX	\$205	E0960
	<input type="checkbox"/> Vest (Ea.)	2 <sup>1</sup> / <sub>4</sub> " x 2 <sup>1</sup> / <sub>2</sub> " x 7" x 5" x 13 <sup>1</sup> / <sub>4</sub> "	1	89111-1XX	\$205	E0960
	<input type="checkbox"/> Vest (Ea.)	2 <sup>1</sup> / <sub>2</sub> " x 3 <sup>1</sup> / <sub>2</sub> " x 9 <sup>1</sup> / <sub>4</sub> " x 5 <sup>1</sup> / <sub>2</sub> " x 15 <sup>1</sup> / <sub>4</sub> "	2	89111-2XX	\$205	E0960
	<input type="checkbox"/> Vest (Ea.)	2 <sup>3</sup> / <sub>4</sub> " x 4 <sup>3</sup> / <sub>4</sub> " x 11 x 5 <sup>3</sup> / <sub>4</sub> " x 18 <sup>1</sup> / <sub>4</sub> "	3	89111-3XX	\$205	E0960
	<input type="checkbox"/> Vest (Ea.)	2 <sup>3</sup> / <sub>4</sub> " x 6 <sup>3</sup> / <sub>4</sub> " x 13 <sup>3</sup> / <sub>4</sub> " x 7 <sup>1</sup> / <sub>2</sub> " x 19"	4	89111-4XX	\$205	E0960
	<input type="checkbox"/> Torso vest (Ea.)	S-P: 7 <sup>3</sup> / <sub>4</sub> –10 <sup>3</sup> / <sub>4</sub> "	XS	85705-XS	\$205	E0960
	<input type="checkbox"/> Torso vest (Ea.)	S-P: 10 <sup>3</sup> / <sub>4</sub> –12 <sup>1</sup> / <sub>2</sub> "	S	85705-S	\$205	E0960
	<input type="checkbox"/> Torso vest (Ea.)	S-P: 12 <sup>1</sup> / <sub>2</sub> –16"	M	85705-M	\$205	E0960
	<input type="checkbox"/> Torso vest (Ea.)	S-P: 16–18 <sup>1</sup> / <sub>2</sub> "	L	85705-L	\$205	E0960
	<input type="checkbox"/> Torso vest (Ea.)	S-P: 18 <sup>1</sup> / <sub>2</sub> –20 <sup>1</sup> / <sub>2</sub> "	XL	85705-XL	\$205	E0960
	<input type="checkbox"/> Torso vest (Ea.)	S-P: 20 <sup>1</sup> / <sub>2</sub> –23 <sup>1</sup> / <sub>2</sub> "	XXL	85705-XXL	\$205	E0960
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 6 <sup>3</sup> / <sub>4</sub> –12 <sup>1</sup> / <sub>2</sub> "	XS	85725-XS	\$138	E0978
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 8 <sup>1</sup> / <sub>4</sub> –13 <sup>3</sup> / <sub>4</sub> "	S	85725-S	\$138	E0978
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 9 <sup>3</sup> / <sub>4</sub> –16 <sup>1</sup> / <sub>4</sub> "	M	85725-M	\$138	E0978
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 12 <sup>1</sup> / <sub>2</sub> –17 <sup>1</sup> / <sub>4</sub> "	L	85725-L	\$138	E0978
	<input type="checkbox"/> 4-Point Padded Hip Belt, 1"	H-H: 16 <sup>1</sup> / <sub>4</sub> –20 <sup>1</sup> / <sub>2</sub> "	XL	85725-XL	\$138	E0978
	<input type="checkbox"/> 2-Point Pelvic Belt, Black	L: 13-14 <sup>1</sup> / <sub>2</sub> "	1	85424-1	\$99	E0978
	<input type="checkbox"/> Fix-Locks, Black, Bag of 6 Pcs ♦ ♦ <b>Required for use of Vests</b>		One Size	862109-35	\$43	
	<input type="checkbox"/> Foot Straps, Plastic Lock (Pr.)	Max: 22", A: 13 <sup>1</sup> / <sub>2</sub> -20 <sup>1</sup> / <sub>2</sub> "	One Size	85096	\$37	K0040
	<input type="checkbox"/> Foot Straps, Velcro (Pr.)		One Size	85102	\$52	K0040

Other types of vests and belts available, please see Vest and Belts Order Form

Accessories	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Folding Top (Canopy), Black		1	8740110	\$338	
	<input type="checkbox"/> Folding Top (Canopy), Black		2	8740120	\$338	
	<input type="checkbox"/> Foam Back Plate ●●	5" x 12"	1	8776325	\$24	
	<input type="checkbox"/> Foam Back Plate ●●	5½" x 13¾"	2	8776326	\$24	
	<b>Future Growth - Fills out the gap between seat and the back</b>					
	<b>●●Custom order; Approx. lead time: 6-8 weeks. Non-returnable item</b>					
	<input type="checkbox"/> Tray with Edge, Grey ♦	17.3" x 13.3"	One Size	8776350	\$103	E0950
	<input type="checkbox"/> Adjustable Removable Tray Hardware ♦		One Size	x8776350	\$169	E1028
	<b>♦Tray and Adjustable Removable Hardware must be ordered together</b>					
	<input type="checkbox"/> Rain Cover ♦		1	8720141	\$57	
	<input type="checkbox"/> Rain Cover ♦		2	8720142	\$57	
	<b>♦Requires Folding Top (Canopy)</b>					
	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Handle Grip	14"	1	8720401	\$205	
	<input type="checkbox"/> Handle Grip	16"	2	8720402	\$205	
	<input type="checkbox"/> Backpack, Black ●●		Small	87043183-1	\$57	
	<b>●●Custom order; Approx. lead time: 6-8 weeks. Non-returnable item</b>					
	<input type="checkbox"/> Medical Essential Carrier		One Size	8720600	\$99	
	<input type="checkbox"/> Vent Tray ♦*	13" x 17"	One Size	8776290	\$472	E1029
	<b>♦Requires Central Wheel Locks</b>					
	<b>*44 lbs. weight capacity</b>					
	<input type="checkbox"/> O2 Tank Holder, Horizontal Mounted *			8776295	\$359	E2208
	<b>*Cannot be used during transit and manufacturer installation is required</b>					
	<input type="checkbox"/> O2 Tank Holder, Vertical Mounted *			8776296	\$105	E2208
	<b>*Cannot be used during transit and manufacturer installation is required</b>					
	<input type="checkbox"/> IV Pole, Height Adjustable and Collapsible		One Size	C904566	\$205	K0105
	<input type="checkbox"/> Mounting bracket for Stingray backrest*		One Size	8776291	\$144	
	<b>*Required with IV Pole #C904566</b>					

Spare Parts	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<b>Required, must select a color</b>					
	<b>X = Color</b>					
	<input type="checkbox"/> Blue = 1	<input type="checkbox"/> Red = 2				
	<input type="checkbox"/> Pink = 5	<input type="checkbox"/> Grey = 4				
	<input type="checkbox"/> Seat and Back Cushion (Pr.) ●●		1	875051x	\$444	
	<input type="checkbox"/> Seat and Back Cushion (Pr.) ●●		2	875052x	\$444	
	<b>●●Colors Green, Orange &amp; Yellow available as a custom order; Approx. lead time: 6-8 weeks. Non-returnable item</b>					
	<input type="checkbox"/> Adductor Covers, Black (Pr.)		1	8750410	\$57	
	<input type="checkbox"/> Adductor Covers, Black (Pr.)		2	8750420	\$57	
	<input type="checkbox"/> Hardware Kit for Hip Supports, Fixed		One Size	872023-KIT	\$36	
	<input type="checkbox"/> Directional Arrow, Replacement		One Size	149889	\$6	

Spare Parts	Description	Measurements	Size	Part No.	USD Price	HCPCS
	<input type="checkbox"/> Release Strap for Back Angling		1	87102250	\$47	
	<input type="checkbox"/> Release Strap for Back Angling		2	87102255	\$47	
	<input type="checkbox"/> Footrest Rod, Elevating	8 1/4"	Short	8775317	\$113	
	<input type="checkbox"/> Footrest Rod, Elevating	11"	Long	8775318	\$113	
	<input type="checkbox"/> Angle Adjustable Footplate		One Size	8775010	\$267	

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

When being transported in a vehicle, it is R82's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair; therefore R82 cannot recommend a transit system.

The R82 Stingray mobility base was tested with a footrest and R82 transit anchors. The performance and safety of accessories and/or other seating options not tested by R82 is unknown, and no warranty or claim about their performance is or can be made.

Prices in USD & is subject to change without notice

**For Quotes and Orders please send completed Order Form to: [sales.us@etac.com](mailto:sales.us@etac.com) or fax: 1-800-760-0537**

Etac United States  
13137 Bleinheim Lane  
Matthews, North Carolina 28105  
+1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF6300-ENG-REV.07

  
Creating Possibilities