


Tetraknapp till 22", 24", 26" drivhjul

Quadraplegic wheel release for 22", 24", 26" Rear wheel

Tetraknöpfe für 22", 24", 26" Hinterräder

74011.2 2022-08-11

1



\* = Infästningspunkter Tetraknapp  
 \* = Housing points for the Quadraplegic wheel release  
 \* = Punkte für Festung der Tetraknöpfe

Pos.	Benämning	Specification	Spezifikation
1.	Tetra knapp	Quadraplegic wheel release	Tetraknöpfe
2.	Skruv	Screw	Schraube
3.	Eker	Spoke	Speichen
4.	Nav	Hub	Nabe
5.	Bricka	Washer	Dichtung
6.	Fjäder	Spring	Feder
7.	Mutter	Nut	Mutter