

Etac Cross 5 prescription form

1A15715.1 2020-03-18

1 Model						
Short seat depth			Long seat depth			Item. no.
35 cm		45 cm	35 cm		45 cm	
37.5 cm		47.5 cm	37.5 cm		47.5 cm	
40 cm		50 cm	40 cm		50 cm	
42.5 cm			42.5 cm			

2 Frame Active
No
Yes

3 Frame colour	
Black -01	
White -92	
Dark blue -93	

4 Back support system	Back support height		Seat to back angle	Lumbar angle	Cover		
	Std	High			Standard	Comfort	Tidy
Cross back 3A			o	o			
	cm	cm					
Cross back 2B (with fixed push handles)			o		Standard		

5 Push handles/Bow handle	
Push handles, height adjustable	
Bow handle, height adjustable	
Push handles, fixed	

7 Attendant brake	
Without attendant brake	
With attendant brake	

If visible, required field. Please, select.

8 Seat	
Seat upholstery	
Solid seat, detachable	

6 User brakes		
Model	R	L
Standard angled black		
Standard angled red		
With extension arm (foldable)		
Single-hand brake right or left mounted		

10 Leg support		
Model	R	L
Standard, 70° angle		
Standard, 70° angle with locking device		
Narrow, 90° angle		
Narrow, 90° angle with locking device		
Short narrow, 90° angle		
Short narrow, 90° angle with locking device		
Elevating, foldable calf support		
Elevating, padded calf support adjustable laterally		
Elevating, padded, flip-up calf support adjustable laterally		

11 Other leg supports		
Amputee support, plush cover		
Amputee support, hygiene cover		
Amputee weight		
Amputee support+weight, plush cover		
Amputee support+weight, hygiene cover		
Plaster leg support, 0° knee angle		

12 Foot support		
Model	R	L
Depth 13 cm, medium attachment tube		
Depth 13 cm, short attachment tube		
Depth 13 cm, long attachment tube		
One-piece foot support		

Seat height and balance position

1

Camber angle	
0°	2°



2

Rear seat height and rear wheel

Rear wheel	20"	22"	24"	24" Lightweight wheel (carbon)	26"
G	34 cm	36,5 cm	39 cm	39 cm	
F	36 cm	38,5 cm	41 cm	41 cm	43 cm*
E	38 cm	40,5 cm	43 cm	43 cm	45 cm*
D	40 cm	42,5 cm	45 cm	45 cm	47 cm
C	42 cm	44,5 cm	47 cm	47 cm	49 cm
B	44 cm	46,5 cm	49 cm	49 cm	51 cm
A	46 cm	48,5 cm	51 cm	51 cm	53 cm

* Not with user brake in this position

3 Balance position



4

Front seat height and castor

Castor	Front fork extra short	Front fork short	Front fork Medium	Front fork Long	Castor housing	
4"	cm	-	-	-	Low widened castor housing	
5"	-	cm	-	-		
6½"	-	cm	-	-		
5"	-	cm	cm	-	Standard castor housing	Widened castor housing
6"	-	cm	cm	-		
6½"	-	cm	cm	cm		
7"	-	-	cm	cm		
8"	-	-	cm	cm		

4 mm fork thickness

5

Hand rim

Aluminium	
Stainless	
Cellular rubber	
Vinyl coated	
Titanium	

7

Single hand drive

Without single hand drive	
Mounted on right side	
Mounted on left side	

6

Tyre rear wheel

Pneumatic high pressure 1"	
Pneumatic 1 3/8"	
Solid	

8

Tyre castor 8"

8" (200X35) semipneumatic	
8" (200X50) semipneumatic	
8" (200X30) pneumatic	

For correct setting of the front fork attachment: Max. backward inclination 6 cm higher front than rear (9°), max. forward inclination 2 cm lower front than rear (3°).
Measures in brackets are theoretical. Hover over the editable text fields to see possible front seat height.

Accessories

Arm support/mudguard		R	L	Item. no.
Solid moulded plate, height adjustable	Length 25			25172-1 25172-2
	Length 25, with finger screw			26379-1 26379-2
	Length 38			25173-1 25173-2
	Length 38, with finger screw			26380-1 26380-2
Detachable arm support plate, height adjustable	Length 25			27230-1 27230-2
	Length 25, with finger screw			27349-1 27349-2
	Length 35			27231-1 27231-2
	Length 35, with finger screw			27350-1 27350-2
Mudguard, carbon fibre				27397-1 27397-2
Without arm support/mudguard				
Arm support cushion dark grey plush, detachable	Foam	25 cm		62282-60
		38 cm		62283-60
	Gel	25 cm		26240-60
		38 cm		26241-60
Arm support cushion Hemi				
Arm support cushion, wide and soft				26873
Without arm support cushion				
Thigh support				
Soft pad for thigh support				82488-1 82488-2
Without thigh support				
Widening kit				27533
Without widening kit				

				Item. no.
Tray, to be mounted on long arm supports				
Tray Hemi	For short arm support			
	For long arm support			
Non-slip cover for tray Hemi				81633
Tray Hemi communication right/left				27352-1 27352-2
Without tray				
Seat extension Cross back 3A				
Without seat extension				
Seat cushion, plush cover				
Comfort wedge				
Seat cushion Tidy, hygiene cover				
Without seat cushion				
		Ri	Le	
Calf support cover, for elevating leg support	Plush			62762-60
	Hygiene			62762-63
Without calf support cover				
Calf strap				
Padding for calf strap				
Calf support wide				
Without calf strap				
Foot support extension 21 cm		Ri	Le	
Without foot support extension				
Foot support padding				
For foot support depth 13 cm				
For one-piece foot support				
Without foot support padding				

Heel strap			24458
Without heel strap			
Positioning belt		pcs	25668

Accessories

			Item. no.
Head support	Dark grey plush		28534-60
	Grey Hygiene		28534-65
Head support with side support	Grey Hygiene		28536-65
Without head support			
Head support attachment for Cross back 3A			
Head support attachment for bow handle			25795
Without head support attachment			
Back support brace, foldable			
Without back support brace			
Back support wedges	H:2, W:9, L:13 cm	pair	28030
	H:2, W:11, L:18.5 cm	pair	28139
	H:2.5, W:17.5, L:20 cm	pair	28140
Back support wedge plush	Short 15 cm	pcs	63012-60
	Long 30 cm	pcs	63013-60
Side cushion			
Side cushion			63171-60-1 63171-60-2
Without side cushion			
Side stop	Straight	pcs	28136
	Widening 15 mm	pcs	28137
	Widening 30 mm	pcs	28173

	Ri	Le	
Trunk support, swing-away	0°		28141-1 28141-2
	-15°		28180-1 28180-2
	-45°		28181-1 28181-2
Trunk support, detachable	+15°		28143
	0°		28142
	-15°		28144
	-45°		28145
Trunk support, soft			27395
Without trunk support			
Anti-tipper			28457
Anti-tipper for frame Active. Flip-up	Settable in height and length		28470
	Settable in height and length		28471

	Ri	Le	Item. no.
Cane holder universal			27327
Without cane holder			
Spoke guard	Grey		
	Yellow		
Without spoke guard			
Locking strap for folded wheelchair			27151
Without locking strap			
High mounting kit for solid seat			27330
Without high mounting kit			
Tool kit			25174
Without tool kit			

Reference

User

Phone

Date

Comments