

M100

Prescription form

Name: _____
 Reference: _____
 Standard version is marked with* Date: _____

M100 Seat width (cm) Seat depth 40 cm	37,5 cm 13400103 qty_____	40 cm 13400105 qty_____	42,5 cm 13400107 qty_____
	45 cm 13400109 qty_____	47,5 cm 13400111 qty_____	50 cm 13400113 qty_____
M100 Long Seat width (cm) Seat depth 46 cm	37,5 cm 13400104 qty_____	40 cm 13400106 qty_____	42,5 cm 13400108 qty_____
	45 cm 13400110 qty_____	47,5 cm 13400112 qty_____	50 cm 13400114 qty_____
Frame colour	Grey -42 <input type="checkbox"/>		
Back/Push handles	Push handles fixed <input checked="" type="checkbox"/>	Height adjustable <input type="checkbox"/>	Bow handle <input type="checkbox"/>
Backrest upholstery	Fixed <input checked="" type="checkbox"/> Adjustable <input type="checkbox"/>		
Brakes	Standard Right <input checked="" type="checkbox"/> Left <input checked="" type="checkbox"/>	With brake lever extension Right <input type="checkbox"/> Left <input type="checkbox"/>	
Legrests Lockable, detachable/swing away	Standard Right <input checked="" type="checkbox"/> Left <input checked="" type="checkbox"/>	Narrow (only with 5" and 6" castor) Right <input type="checkbox"/> Left <input type="checkbox"/>	Narrow, short (only with 5" and 6" castor) Right <input type="checkbox"/> Left <input type="checkbox"/>
	Elevating legrest Flip-up calf support Right <input type="checkbox"/> Left <input type="checkbox"/>	Amputee support adjustable in height, angle, depth and sideways Right <input type="checkbox"/> Left <input type="checkbox"/>	Grey plush <input type="checkbox"/> Black Dartex <input type="checkbox"/>
Footrest	Flip-up standard Adjustable in, angle and height Right <input checked="" type="checkbox"/> Left <input checked="" type="checkbox"/>	Long length tube (+10 cm) Right <input type="checkbox"/> Left <input type="checkbox"/>	
Rear wheel (alum handrim, Quick release)		<input checked="" type="checkbox"/> Rear wheel 24"	
Castor, seat height rear (seat angle 3°)		Seat height	Castor
Massive 1 3/8" <input checked="" type="checkbox"/>	Balance position front rear	<input type="checkbox"/> 40	<input type="checkbox"/> 6"
Drum brake <input type="checkbox"/>	<input checked="" type="checkbox"/> front <input type="checkbox"/> rear	<input type="checkbox"/> 42,5	<input type="checkbox"/> 6"
		<input checked="" type="checkbox"/> 45	<input type="checkbox"/> 6" <input checked="" type="checkbox"/> 8" (200x35) <input type="checkbox"/> 8" (200x50)
		<input type="checkbox"/> 47,5	<input type="checkbox"/> 8" (200x35) <input type="checkbox"/> 8" (200x50)
		<input type="checkbox"/> 50	<input type="checkbox"/> 8" (200x35) <input type="checkbox"/> 8" (200x50)
<input type="checkbox"/> Rear wheel 20"		<input type="checkbox"/> Rear wheel 22"	
Seat height	Castor	Seat height	Castor
<input type="checkbox"/> 36	<input type="checkbox"/> 5"	<input type="checkbox"/> 38	<input type="checkbox"/> 5" <input type="checkbox"/> 6"
<input type="checkbox"/> 38,5	<input type="checkbox"/> 5" <input type="checkbox"/> 6"	<input type="checkbox"/> 40,5	<input type="checkbox"/> 5" <input type="checkbox"/> 6"
<input type="checkbox"/> 41	<input type="checkbox"/> 5" <input type="checkbox"/> 6"	<input type="checkbox"/> 43	<input type="checkbox"/> 6"
<input type="checkbox"/> 43,5	<input type="checkbox"/> 6"	<input type="checkbox"/> 45,5	<input type="checkbox"/> 6" <input type="checkbox"/> 8" (200x35) <input type="checkbox"/> 8" (200x50)
<input type="checkbox"/> 46	<input type="checkbox"/> 6" <input type="checkbox"/> 8" (200x35) <input type="checkbox"/> 8" (200x50)	<input type="checkbox"/> 48	<input type="checkbox"/> 8" (200x35) <input type="checkbox"/> 8" (200x50)




74992E

18-08-07













Etac, Box 203, SE-334 24 Anderstorp

Tel. +46/371-58 73 00, Fax. +46/371-58 73 90

www.etac.com

Armrest with massive plate	 Length 25 cm <input checked="" type="checkbox"/>	 Length 35 cm <input type="checkbox"/>
Anti-tip	Tilter <input checked="" type="checkbox"/>	Anti-tip 27234 <input type="checkbox"/> 

Accessories M100

	Backrest brace only together with height adjustable push handles	seat width cm 37,5 25611 <input type="checkbox"/>	40 25612 <input type="checkbox"/>	42,5 25613 <input type="checkbox"/>	45 25614 <input type="checkbox"/>	47,5 25739 <input type="checkbox"/>	50 25615 <input type="checkbox"/>
	Tray assembly on long armrests	seat width cm 37,5 (47,5cm*) 25187 <input type="checkbox"/>	40 (50cm*) 24418 <input type="checkbox"/>	42,5 (52,5cm*) 25188 <input type="checkbox"/>	45 (55cm*) 24419 <input type="checkbox"/>	47,5 (57,5cm*) 25731 <input type="checkbox"/>	50 (60cm*) 24601 <input type="checkbox"/>
	*=table width						
	Seat cushion Universal dark grey plush, length 56 cm, anatomically shaped, cut to individual seat depth	seat width cm 37,5 62231-60 <input type="checkbox"/>	40 62232-60 <input type="checkbox"/>	42,5 62233-60 <input type="checkbox"/>	45 62234-60 <input type="checkbox"/>	47,5 62374-60 <input type="checkbox"/>	50 62266-60 <input type="checkbox"/>
	Calf strap	37,5 cm 25356 <input type="checkbox"/>	40-42,5 cm 25357 <input type="checkbox"/>	45-47,5 cm 25358 <input type="checkbox"/>	50 cm 25359 <input type="checkbox"/>		
	Heel strap footrest, adjustable in length		Right 24458 <input type="checkbox"/>	Left 24458 <input type="checkbox"/>			
	Positioning belt with snap lock, to be mounted on wheelchair frame					80209440 <input type="checkbox"/>	
	Headrest mod. 1 incl. tube, to be mounted on bow handle, adjustable in height, depth and angle		Grey plush 24462-60 <input type="checkbox"/>	Black Dartex 24462-63 <input type="checkbox"/>			
	Headrest bracket, for bow handle					25795 <input type="checkbox"/>	
	Spokeguard	grey with logo 20" 25502-28 <input type="checkbox"/>	22" 25503-28 <input type="checkbox"/>	24" 25504-28 <input type="checkbox"/>	yellow with logo 20" 25502-29 <input type="checkbox"/>	22" 25503-29 <input type="checkbox"/>	24" 25504-29 <input type="checkbox"/>
	Cane holder		Right 27235 <input type="checkbox"/>	Left 27235 <input type="checkbox"/>			
	IV stand					28324 <input type="checkbox"/>	
	Tool kit					82362 <input type="checkbox"/>	
	Gloves For alternatives look at www.etac.com Item no. _____						