

Etac Prio prescription form

75326F 2017-08-24

③ Back support system		Cover		
		Grey plush	Grey hygiene	Green hygiene
Height settable without dynamic pelvic support. Height 55-70 cm				
Height settable with dynamic pelvic support. Height 55-70 cm				
3D without dynamic pelvic support. Height settable 45-60 cm				
3D with dynamic pelvic support. Height settable 45-60 cm				

④ Head support		Cover		
		Grey plush	Grey hygiene	Green hygiene
Soft				
Multi adjustable				

⑤ Comfort tilt lock	
Without Comfort tilt lock	
With Comfort tilt lock, preset at 0°	
With Comfort tilt lock preset at 4°	
With Comfort tilt lock preset at 8°	
With Comfort tilt lock preset at 12°	

⑥ Arm support		
Model	R	L
Standard		
Hemi		

① Wheelchair		
Seat width	Item. no.	
40	13420205	
45	13420209	
50	13420213	
55	13420217	

② Seat depth	Seat cushion								Seat plate			
	Short	Long	Basic	Com-fort	Grey plush	Grey hygiene	Green hygiene	No seat cushion	Wood	Plas-tic	Wood w. Velcro	Plastic w. Velcro
Pos 1= 39-46 cm short cushion												
Pos 2= 41,5-48,5 cm short cushion												
Pos 3= 44-51 cm short or long cushion												



⑦ Leg support		
Model	R	L
Elevating		
Fixed angle		
Amputee plush		
Amputee Hygiene		
Amputee weight		
Amputee+weight plush		
Amputee+weight Hygiene		
Leg support Cross		
Solid footplate (not with 8" castor)		

⑧ Foot supports		
Model	R	L
Model 1 long attachment tube (38.5-55 cm from seat plate)		
Model 1 med. attachment tube (31-47.5 cm from seat plate)		
Model 2 long attachment tube (38.5-55 cm from seat plate)		
Model 2 med. attachment tube (31-47.5 cm from seat plate)		
Solid footplate (see Leg supports)		

Seat height and balance position

⑨

Seat height, rear wheels, castor wheels and position

			Front fork short						Front fork medium					
														
			S1		S2		S3		M1		M2		M3	
														
		Seat height	6"	8"	6"	8"	6"	8"	6"	8"	6"	8"	6"	8"
E	20"***	43,0	-	-	-	-	-	-	-	-	-	-	-	-
	22"***	45,5	-	-	-	-	-	-	-	-	-	-	-	-
	24"	47,5	-	-	-	-	-	-	*	-	-	-	-	-
D	20"***	45,5	-	-	-	-	-	-	-	-	-	-	-	-
	22"***	48,0	-	-	-	-	-	-	*	-	-	-	-	-
	24"	50,0	-	-	-	-	-	-	-	-	-	-	-	-
C	20"***	48,0	-	-	-	-	-	-	*	-	-	-	-	-
	22"***	50,5	-	-	-	-	-	-	-	-	-	-	-	-
	24"	52,5	-	-	-	-	-	-	-	-	-	-	-	-
B	16"	45,0	-	-	-	-	-	-	-	-	-	-	-	-
A	16"	47,5	-	-	-	-	-	-	*	-	-	-	-	-

⑪

Castor wheel 8"

8x13/8" (200x35 mm)	
8x2" (200x50 mm)	

* Only 8x13/8" (200x35 mm), not 8x2" (200x50 mm)

** In combination with user brake, only in balance position 1


*** In combination with user brake, only in balance position 1 and 2

⑫

Single hand drive

Without single hand drive	
Single hand drive mounted on right side	
Single hand drive mounted on left side	

⑩ Balance position



⑬

Brakes

Model	R	L
Standard black		
Standard red		
Extended handle		
Single hand brake, assembled on right or left		
Attendant brake (std with 16" rear wheel)		

⑭

Push handles/Bow handle

Push handles, height adjustable		
Bow handle, height adjustable		
Bow handle, height and angle adjustable		

Accessories

Accessories		
Trunk support		Item. no.
Detachable without tools	pcs	28020
Swing-away, right		28019-1
Swing-away, left		28019-2
Widening kit +2.5 cm		27892
Padding side support		28021
Tray		
40		27971
45		27972
50		27973
55		27974
Pommel		28022
Calf strap		
40		25357
45		25358
50		25359
55		26062
Padding for calf strap		
20 cm		26845
30 cm		26846
40 cm		26847
Positioning belt	pcs	27904

Accessories		
		Item. no.
Chest strap bracket		28186
Comfort wedges		
H:2, W:11, L:18.5 cm	pair	28139
H:2.5, W:17.5, L:20 cm	pair	28140
H:2, W:9, L:13 cm	pair	28030
Gas piston clips (10 pcs/set)		28163
Lock seat and back angle controls	pcs	28047
Oxygen holder/Transport basket		28162
IV stand		28160
Velcro kit seat cushion for separate mounting		28183
Spoke guards grey		
20"		25502-28
22"		25503-28
24"		25504-28
Spoke guards yellow		
20"		25502-29
22"		25503-29
24"		25504-29

Accessories		
		Item. no.
Cane holder	pcs	27327
Extended foot support (21.5 cm)		
40-42.5		50798
45-47.5		50799
50		50800

Reference

User

Phone

Date

Comments