



# Immedia E-Board Oval

Transfer board - for horizontal transfer, repositioning and application of a sling

## Sliding boards for multiple use

E-Board Oval are thin sliding boards used for bed to bed transfers, application of slings and repositioning.

The flexible boards have a smooth, low-friction surface that reduces the friction and makes it very easy to perform different transfers with very little effort for the carer.

The thin and flexible design of the boards makes it very easy to apply and remove.

The E-Board Oval comes in two different sizes - S and L.

The E-Board Oval S are used for:

Repositioning a user further back in a wheelchair

Application of a sling

Other small repositioning purposes

The E-Board Oval L are used for:

Horizontal transfer from a bed to another lying surface

Application of a sling in a wheelchair or in a bed

Other repositioning purposes

The handles make the boards easy to transport and hang on a wall. Both E-Board Oval S and L are always delivered in a set of 2.



## Features



### Thin and flexible

Immedia E-board oval are very thin and flexible boards which reduces the friction.



### Easy handling

The thin and smooth surface makes the board very easy to apply and remove. The low weight and the handles makes it easy to carry.






### Two sizes

Immedia E-Board Oval is available in two sizes - S and L.

Product	Dimensions	Item no.
Immedia E-Board Oval L	W850 x L390 x H3 mm	IM411
Immedia E-Board Oval S	W450 x L190 x H3 mm	IM412

**Material**  
Polyethylene PEHD

**CE**

This product can be wiped off with soapy water or cleaned with a solvent-free cleaning agent with a pH level between 5-9, or a 70% disinfection solution.