



Clean, Clean 24, Clean Height adj.

Shower- and toiletchair

Sparepart list



78299H 14-10-20







Clean

Clean 24"

Clean Height adj.

Accessories



Spareparts			Product			
Pos.	Description	Item no.	Clean	Clean 24	Clean Height adj.	Accessories
1	Armrest, pair, grey	8400 1270	●	●	●	
	Armrest, pair, green	8400 5035	●	●		
2	Armrest lock, pair	8400 5036	●	●	●	
3	Castors Clean 125 mm lockable, 2pcs, green	8400 5033	●	●		
	Castors Clean 125 mm lockable, 2pcs, grey	8400 1280	●	●		
	Castor Clean 125 mm, pair, green	8400 5034	●	●		
	Castor Clean 125 mm, pair, grey	8400 1281	●	●		
	Castors Clean height adj lockable compl	8400 1276			●	
4	Screw/ nut kit Clean, 4+4	8400 5042	●	●		
5	Screw/nut kit, 4+4	8400 1271			●	
6	Pin/nut kit Clean height adj, 2+2	8400 1272			●	
7	Tube plug Clean height adj, for castors 4 pcs	8400 1273			●	
8	End cap holder footrest stop, Clean height adj 2 pcs	8400 1274			●	
9	Footrest Clean, grey, grå	8400 1275	●	●	●	
	Footrest Clean, green, grön	8400 5037	●	●		
10	Seat Clean, grey	8400 1277	●	●	●	
	Seat Clean, green	8400 1278	●	●		
11	Backrest Clean, grey	8400 5085	●	●	●	
	Backrest Clean, green	8400 5086	●	●		
12	Rear wheel Clean pair, 2st	8400 5039		●		
13	Brake Clean right	8400 5040		●		
	Brake Clean left	8400 5041		●		
14	Handle for pan lid, grey	8400 1282				●
	Handle for pan lid, green	8400 5038				●
15	Snap fasteners Protect cushion	8400 5056				●
16	Attachment legrest right, Clean height adjustable	8400 1285				●
	Attachment legrest left, Clean height adjustable	8400 1286				●
17	Touch up paint white	8400 5052				●
18	Torxkey T 30, 2 pcs, Clean	8400 5043				●
19	Castor housing incl. screw – green	8400 5106	●	●		
	Castor housing incl. screw – grey	8400 5107	●	●		