



 **Convaid**

Coaster

Enhances family adventures, comfort, and safety while on the go!

Made for Families

The Coaster makes mobility easy for families with its lightweight, compact folding design with an intuitive one-step unfolding feature for quick setup. Families can travel confidently, knowing the crash-tested frame meets WC19 and ISO 7176-19 standards for safety.

Made for Adventures

Transform everyday outings into exciting family adventures with the Coaster, designed to keep your child comfortable by providing tailored positioning for maximum comfort, encouraging engagement with peers and family. Whether enjoying a day at the beach or exploring a theme park, the Coaster is an ideal companion for accessible family fun, creating cherished memories as you explore together.

Made to Last

Built for durability and safety, the Coaster is engineered to withstand the everyday rigors of busy family life, ensuring it remains a reliable mobility solution. Families can relax, knowing the Coaster delivers long-term performance and dependability for countless adventures together.



Coaster

The Coaster is thoughtfully designed with families and children in mind, providing easy mobility and enhancing comfort through adjustable fixed tilt and recline settings. This feature is particularly beneficial for children with limited body strength, empowering them to engage during family outings and actively create lasting memories together. The Coaster is not just a mobility solution; it prioritizes user comfort and accessibility, making it a must-have for enhancing the quality of life for families.

Key Features



Compact-folding, durable, adjustable fixed tilt and recline



Adjustable Fixed Tilt and Recline

With the simple replacement of two bolts, you can adjust the fixed tilt setting to 10°, 20°, or 30° to accommodate the user's needs. For extra positioning, the Coaster offers four fixed recline settings, dependent on the fixed tilt selected.



Zippered Pockets in Seat and Back

The upholstery with built in seat and back zippers allow additional padding and support for the user.



Height and Tension Adjustable Back Positioner

Provides increased lumbar support to encourage proper body positioning and alignment.



Compact Clamshell Folding

Like many commercial buggies, the Coaster can compactly fold into an easily transportable package. The rear wheels are also quick release to offer an even smaller option when loading this buggy.



Continued Convaid Durability

Convaid's reputation-building robustness is once again found in the Coaster, exceeding all industry standards for durability testing.

Accessories



Mosquito net
Canopy required



Rain cover
Canopy required



Headrest cover (canopy)



Extended Headrest cover (canopy)



Padded headwings



Full torso support vest



Height and tension adjustable back positioner



Calf panel



Foot positioners



Height adjustable flip-up armrests



Curb Tipper





















Rear anti-tip tubes

For full product offering visit our website: www.etac.com



Measurements

| | CTR 12 | CTR 14 | CTR 16 | CTR 18 |
|---|-------------------------------|--|-------------------------------|-------------------------------|
|  Seat width | 305mm | 356mm | 406mm | 457mm |
|  Seat depth (Standard) | 254-305mm | 279-330mm | 305-381mm | 368-445mm |
|  +76mm Seat depth (Option) | 330-381mm | 356-406mm | 381-457mm | 445-521mm |
|  Seat back height (10° / 20° / 30°) | 559/584/610mm | 635/660/686mm | 660/686/711mm | 711/711/737mm |
|  Seat to footplate | 152-406mm | 152-406mm | 152-508mm | 152-508mm |
|  Seat to floor (10° / 20° / 30°) | 483/508/533mm | 483/508/533mm | 533/559/584mm | 533/559/584mm |
|  Overall height | 1016mm | 1041mm | 1067mm | 1118mm |
|  Overall length | 762mm | 787mm | 813mm | 838mm |
|  Overall width | 559mm | 610mm | 660mm | 711mm |
|  Folded height | 432mm | 432mm | 432mm | 432mm |
|  Folded length | 813mm | 813mm | 889mm | 889mm |
|  Folded width | 559mm | 610mm | 660mm | 711mm |
|  Headrest extension height | 203mm | 203mm | 203mm | 203mm |
|  Shoulder strap height (10° / 20° / 30°) | 267-406 / 292-432 / 318-457mm | 330-483 / 356-508 / 381-533mm | 356-508 / 381-533 / 406-559mm | 419-610 / 445-635 / 470-660mm |
|  Push handle adjustments | 813-1168mm | 864-1194mm | 864-1270mm | 889-1295mm |
|  Weight capacity regular/transit | 34/30kg | 45/45kg | 77/77kg | 77/77kg |
|  Weight of buggy without front rigging | 13kg | 13.5kg | 14kg | 14.5kg |
|  Fixed tilt and recline settings | 10° fixed tilt | 100°, 105°, 110°, 115° recline options | | |
| | 20° fixed tilt | 90°, 95°, 100°, 105° recline options | | |
| | 30° fixed tilt | 80°, 85°, 90°, 95° recline options | | |



Coaster_Brochure_01_INT | Reserving our rights for printer's error discontinued products. Please find an updated edition on our website © 2025 Etac All rights reserved

Etac is a world-leading developer of innovative assistive devices and patient handling equipment. Since 1973, we have been committed to improving quality of life for the individual, family members and caregivers.

For the intended purpose of the products and latest information, please visit www.etac.com

Etac

+1 844-US-MOBILITY (844-876-6245)

international.convaid.us@etac.com

www.etac.com

etac[®]
Creating Possibilities