



R82[®]

R82 Crocodile

Supporting the motivation to move in all terrains

Product highlights

The Crocodile is the perfect companion for everyday adventures. A wide range of accessories can be added to meet your child's need for support while exploring the world.

Growth and support

The Crocodile adapts effortlessly, thanks to its easy, adjustment of the height and its two angle adjustable handles. Positioning the handles close to the body provides the child with better support and a better walking position. Additional accessories can be added guaranteeing that your child enjoys a stable walking position, fostering interaction and participation with both peers and family.

Early intervention

Choosing the Crocodile in a size 0 or size 1 with the Starter kit accessory, provides supported standing and walking, giving the opportunity for early intervention and facilitation of the child's social interaction with parents, carers, and the environment.

Safety

Our Crocodile is designed with comfort and safety in mind and comes equipped with an anti reverse wheel stop, preventing the Crocodile from rolling backwards. Additionally, hand brakes can be mounted for extra security ensuring that your child is safe to go on family adventures in all terrains.

Easy transportation

The Crocodile can be easily folded with the seat and accessories mounted, making it the perfect companion for you to bring on the go.

Outdoor adaptable

The Crocodile is equipped with front wheels with directional stabilizers to provide your child with a supported walk in all terrains. Choosing between the Crocodile with tension adjustable rear wheels or the Crocodile All-Terrain size variant with larger wheels and a deeper tread, gives you the opportunity to find the walker that best suit your child's needs.

Being active or taking a break?

Adding the Flip down seat to your child's Crocodile provides them with the opportunity to take a rest when needed.



The Crocodile family

The Crocodile family consists of seven sizes with two different variants; the Crocodile and the Crocodile All-Terrain. We have listed the most important features of the two size variants, so that you can make an informed choice, when choosing which walking aid is the right one for your child.



Crocodile:

The Crocodile is available in size 0, 1, 2 and 3

- Height adjustable frame
- Angle adjustable handles
- Anti-reverse wheel stop
- Slim wheels with an adjustable tension on the rear wheels
- Foldable frame
- Multiple accessories available for additional support
- Accommodating children with a user height ranging from 41.3" - 70.9" (1050 mm to 1800 mm)
- Possibility of mounting the Starter kit on size 1

Crocodile All-Terrain:

The Crocodile AT variant comes in size 1, 2 and 3

- Height adjustable frame
- Angle adjustable handles
- Anti-reverse wheel stop
- Bigger wheels with a deeper tread and better compound
- Possibility of dismounting the rear wheels for transportation
- Foldable frame
- Multiple accessories available for additional support
- Accommodating children with a user height ranging from 46.9" - 72.4" (1190 mm to 1840 mm)

Find out more about the Crocodile



Measurements

		Size 0	Size 1	Size 1 AT	Size 2	Size 2 AT	Size 3	Size 3 AT
	Width between handles	3.9 - 9.8" (100-250 mm)	7.5 - 13.8" (190-350 mm)	7.5 - 13.8" (190-350 mm)	9.4 - 16.5" (240-420 mm)	9.4 - 16.5" (240-420 mm)	13 - 18.9" (330-480 mm)	13 - 18.9" (330-480 mm)
	Height from floor to handles	11.8 - 24" (300-610 mm)	16.5 - 28" 15 - 27.6** (420-710 mm 380-700 mm*)	18.1 - 29.5" (460-750 mm)	22 - 33.5" (560-850 mm)	23.6 - 35" (600-890 mm)	29.1 - 40.4" (740-1025 mm)	30.7 - 41.9" (780-1065 mm)
	Max. user height	41.3" (1050 mm)	45.3" (1150 mm)	46.9" (1190 mm)	59" (1500 mm)	60.6" (1540 mm)	70.9" (1800 mm)	72.4" (1840 mm)
	Max. user mass	66 lbs (30 kg)	66 lbs (30 kg)	66 lbs (30 kg)	99 lbs (45 kg)	99 lbs (45 kg)	176 lbs (80 kg)	176 lbs (80 kg)
	Length	25.2" (640 mm)	28.3" (720 mm)	30.7" (780 mm)	33.1" (840 mm)	35.4" (900 mm)	34.3" (870 mm)	36.6" (930 mm)
	Overall width	22.8" (580 mm)	24.8" (630 mm)	26" (660 mm)	27.2" (690 mm)	27.6" (700 mm)	27.2" (690 mm)	28.3" (720 mm)
	Folded length	25.6" (650 mm)	28.3" (720 mm)	29.5" (750 mm)	33.9" (860 mm)	35" (890 mm)	39.4" (1000 mm)	40.6" (1030 mm)
	Folded width	22.8" (580 mm)	24.4" (620 mm)	26" (660 mm)	27.2" (690 mm)	27.6" (700 mm)	27.2" (690 mm)	28.3" (720 mm)
	Folded height	11.8" (300 mm)	11.8" (300 mm)	13.6" (345 mm)	10.6" (270 mm)	12.2" (310 mm)	11.4" (290 mm)	13" (330 mm)

*Mounted with the Starter kit accessory. Further measurements are available on our website.

Further measurements can be found online

Colors

