

# Rogue2

Solid og aktiv rullestol for optimal ytelse!

Beskrivelse	Lev. artnr	HMS art.nr.	Setebredde (cm)	Setedybde (cm)	Setehøyde uten bytte av deler (cm)	Setehøyde med bytte av deler (cm)	Rygghøyde (cm)	Totalbredde (cm)	Totalvekt (kg)	Transportvekt (kg)
Rogue2 sb31 sd36	KI13751214	317551	31	36	46-48	34-52	28-51 (standard 30-36)	55,5	10,4	6,6
Rogue2 sb31 sd38	KI13751215	317552	31	38	46-48	34-52	28-51 (standard 30-36)	55,5	10,5	6,7
Rogue2 sb31 sd41	KI13751216	317553	31	41	46-48	34-52	28-51 (standard 30-36)	55,5	10,6	6,8
Rogue2 sb31 sd43	KI13751217	317554	31	43	46-48	34-52	28-51 (standard 30-36)	55,5	10,6	6,8
Rogue2 sb31 sd46	KI13751218	317555	31	46	46-48	34-52	28-51 (standard 30-36)	55,5	10,7	6,9
Rogue2 sb33 sd36	KI13751314	317556	33	36	46-48	34-52	28-51 (standard 30-36)	57,5	10,6	6,8
Rogue2 sb33 sd38	KI13751315	317557	33	38	46-48	34-52	28-51 (standard 30-36)	57,5	10,7	6,9
Rogue2 sb33 sd41	KI13751316	317558	33	41	46-48	34-52	28-51 (standard 30-36)	57,5	10,7	6,9
Rogue2 sb33 sd43	KI13751317	317559	33	43	46-48	34-52	28-51 (standard 30-36)	57,5	10,8	7
Rogue2 sb33 sd46	KI13751318	317560	33	46	46-48	34-52	28-51 (standard 30-36)	57,5	10,8	7
Rogue2 sb33 sd48	KI13751319	317561	33	48	46-48	34-52	28-51 (standard 30-36)	57,5	10,9	7,1
Rogue2 sb36 sd36	KI13751414	317562	36	36	46-48	34-52	28-51 (standard 30-36)	60,5	10,9	7,1
Rogue2 sb36 sd38	KI13751415	317563	36	38	46-48	34-52	28-51 (standard 30-36)	60,5	10,9	7,1
Rogue2 sb36 sd41	KI13751416	301422	36	41	46-48	34-52	28-51 (standard 30-36)	60,5	11	7,2
Rogue2 sb36 sd43	KI13751417	317564	36	43	46-48	34-52	28-51 (standard 30-36)	60,5	11	7,2
Rogue2 sb36 sd46	KI13751418	317565	36	46	46-48	34-52	28-51 (standard 30-36)	60,5	11	7,2
Rogue2 sb36 sd48	KI13751419	317566	36	48	46-48	34-52	28-51 (standard 30-36)	60,5	11,1	7,3
Rogue2 sb36 sd51	KI13751420	317567	36	51	46-48	34-52	28-51 (standard 30-36)	60,5	11,1	7,3
Rogue2 sb38 sd36	KI13751514	317568	38	36	46-48	34-52	28-51 (standard 30-36)	62,5	10,9	7,1
Rogue2 sb38 sd38	KI13751515	317569	38	38	46-48	34-52	28-51 (standard 30-36)	62,5	11	7,2
Rogue2 sb38 sd41	KI13751516	317570	38	41	46-48	34-52	28-51 (standard 30-36)	62,5	11	7,2
Rogue2 sb38 sd43	KI13751517	317571	38	43	46-48	34-52	28-51 (standard 30-36)	62,5	11,1	7,3
Rogue2 sb38 sd46	KI13751518	317572	38	46	46-48	34-52	28-51 (standard 30-36)	62,5	11,2	7,4
Rogue2 sb38 sd48	KI13751519	317573	38	48	46-48	34-52	28-51 (standard 30-36)	62,5	11,2	7,4
Rogue2 sb38 sd51	KI13751520	317574	38	51	46-48	34-52	28-51 (standard 30-36)	62,5	11,2	7,4
Rogue2 sb41 sd36	KI13751614	317575	41	36	46-48	34-52	28-51 (standard 30-36)	65,5	11	7,2
Rogue2 sb41 sd38	KI13751615	317576	41	38	46-48	34-52	28-51 (standard 30-36)	65,5	11	7,2
Rogue2 sb41 sd41	KI13751616	317577	41	41	46-48	34-52	28-51 (standard 30-36)	65,5	11,1	7,3
Rogue2 sb41 sd43	KI13751617	317578	41	43	46-48	34-52	28-51 (standard 30-36)	65,5	11,1	7,3
Rogue2 sb41 sd46	KI13751618	317579	41	46	46-48	34-52	28-51 (standard 30-36)	65,5	11,2	7,4
Rogue2 sb41 sd48	KI13751619	317580	41	48	46-48	34-52	28-51 (standard 30-36)	65,5	11,2	7,4
Rogue2 sb41 sd51	KI13751620	317581	41	51	46-48	34-52	28-51 (standard 30-36)	65,5	11,3	7,5
Rogue2 sb43 sd36	KI13751714	317582	43	36	46-48	34-52	28-51 (standard 30-36)	67,5	11,1	7,3
Rogue2 sb43 sd38	KI13751715	317583	43	38	46-48	34-52	28-51 (standard 30-36)	67,5	11,2	7,4
Rogue2 sb43 sd41	KI13751716	317584	43	41	46-48	34-52	28-51 (standard 30-36)	67,5	11,2	7,4
Rogue2 sb43 sd43	KI13751717	301295	43	43	46-48	34-52	28-51 (standard 30-36)	67,5	11,2	7,4
Rogue2 sb43 sd46	KI13751718	317585	43	46	46-48	34-52	28-51 (standard 30-36)	67,5	11,3	7,5
Rogue2 sb43 sd48	KI13751719	317586	43	48	46-48	34-52	28-51 (standard 30-36)	67,5	11,3	7,5

# Rogue2

Solid og aktiv rullestol for optimal ytelse!

Beskrivelse	Lev. artnr	HMS art.nr.	Setebredde (cm)	Setedybde (cm)	Setehøyde uten bytte av deler (cm)	Setehøyde med bytte av deler (cm)	Rygghøyde (cm)	Totalbredde (cm)	Totalvekt (kg)	Transportvekt (kg)
Rogue2 sb43 sd51		317587	43	51	46-48	34-52	28-51 (standard 30-36)	67,5	11,4	7,6
Rogue2 sb46 sd38		317588	46	38	46-48	34-52	28-51 (standard 30-36)	70,5	11,2	7,4
Rogue2 sb46 sd41		317589	46	41	46-48	34-52	28-51 (standard 30-36)	70,5	11,2	7,4
Rogue2 sb46 sd43		317590	46	43	46-48	34-52	28-51 (standard 30-36)	70,5	11,3	7,5
Rogue2 sb46 sd46		317591	46	46	46-48	34-52	28-51 (standard 30-36)	70,5	11,3	7,5
Rogue2 sb46 sd48		317592	46	48	46-48	34-52	28-51 (standard 30-36)	70,5	11,4	7,6
Rogue2 sb46 sd51		317593	46	51	46-48	34-52	28-51 (standard 30-36)	70,5	11,4	7,6
Rogue2 sb48 sd41		317594	48	41	46-48	34-52	28-51 (standard 30-36)	72,5	11,2	7,4
Rogue2 sb48 sd43		317595	48	43	46-48	34-52	28-51 (standard 30-36)	72,5	11,3	7,5
Rogue2 sb48 sd46		317596	48	46	46-48	34-52	28-51 (standard 30-36)	72,5	11,3	7,5
Rogue2 sb48 sd48		317597	48	48	46-48	34-52	28-51 (standard 30-36)	72,5	11,3	7,5
Rogue2 sb48 sd51		317598	48	51	46-48	34-52	28-51 (standard 30-36)	72,5	11,4	7,6
Rogue2 sb51 sd41		317599	51	41	46-48	34-52	28-51 (standard 30-36)	75,5	11,3	7,5
Rogue2 sb51 sd43		317600	51	43	46-48	34-52	28-51 (standard 30-36)	75,5	11,3	7,5
Rogue2 sb51 sd46		317601	51	46	46-48	34-52	28-51 (standard 30-36)	75,5	11,4	7,6
Rogue2 sb51 sd48		317602	51	48	46-48	34-52	28-51 (standard 30-36)	75,5	11,4	7,6
Rogue2 sb51 sd51		317603	51	51	46-48	34-52	28-51 (standard 30-36)	75,5	11,5	7,7

