

# Combination agreement overview

Wheelchairs and seating



1A15211.16 26-04-08

This list is an overview of combinations between Ki devices and devices from other suppliers. The combinations are only applicable as long as the devices are compliant with MDR 2017/745 (unless other is specified). The IFU's of each device must be followed.

Product	Ki devices																					
	Axiom PB PL DL	Catalyst 4	Catalyst 4C	Catalyst 5	Catalyst 5 Ti	Catalyst 5 VX	Ethos	Focus CR	Liberty	Little Wave	Little Wave Cliik	Little Wave Cliik XP	Little Wave Flip	Little Wave XP	Rogue	Rogue ALX	Rogue TTL	Rogue XP	Rogue2	Spark	Tsunami	Tsunami ALX
<b>Alber</b>																						
e-fix E35/E36		✓		✓		✓	✓	✓	✓		✓		✓	✓	✓			✓	✓		✓	
e-motion M25		✓		✓		✓	✓	✓	✓		✓		✓	✓	✓			✓	✓		✓	
e-pilot P15							✓											✓	✓			
Scalamobil S35/S38		✓		✓																		
SMOOV 010		✓		✓			✓							✓				✓	✓		✓	
Twion T24		✓		✓																		
Via GO V24				✓																		
Viamobil V25		✓		✓																		
<b>AMI Autonomie</b>																						
SIROCCO Baroudeur				✓		✓	✓									✓		✓	✓			
SIROCCO Compact				✓		✓	✓									✓		✓	✓			
SIROCCO Extreme				✓		✓	✓									✓		✓	✓			
SIROCCO Voyageur				✓		✓	✓									✓		✓	✓			
<b>Batec</b>																						
Batec Mini 2			✓			✓	✓								✓	✓						✓
Batec Electric 2			✓			✓	✓								✓	✓						✓

# Combination agreement overview

Wheelchairs and seating



1A15211.16 26-04-08

Product	Ki devices																						
	Axiom PB PL DL	Catalyst 4	Catalyst 4C	Catalyst 5	Catalyst 5 TI	Catalyst 5 VX	Ethos	Focus CR	Liberty	Little Wave	Little Wave Cliik	Little Wave Cliik XP	Little Wave Flip	Little Wave XP	Rogue	Rogue ALX	Rogue TTL	Rogue XP	Rogue2	Spark	Tsunami	Tsunami ALX	
Batec Scrambler 2			✓			✓	✓								✓	✓						✓	
Batec Rapid 2			✓			✓	✓								✓	✓						✓	
Batec Hybrid 2						✓																	
Batec Hybrid 3			✓																				
Batec Hybrid 4							✓																
Batec Hybrid 5																✓						✓	
Batec Hybrid 6															✓								
<b>Decon</b>																							
Easywheel				✓	✓	✓	✓		✓		✓				✓			✓					
E-Drive				✓	✓	✓	✓		✓		✓				✓			✓					
E-Move				✓	✓	✓	✓		✓			✓			✓			✓					
<b>Etac</b>																							
Crissy SA	✓																						
Cross 5	✓																						
Cross 6	✓																						
<b>Indes</b>																							
Indes Wheeldrive											✓				✓								
<b>Klaxon</b>																							
Twist and Twist R				✓		✓										✓		✓	✓				

# Combination agreement overview

Wheelchairs and seating



1A15211.16 26-04-08

Product	Ki devices																						
	Axiom PB PL DL	Catalyst 4	Catalyst 4C	Catalyst 5	Catalyst 5 Ti	Catalyst 5 VX	Ethos	Focus CR	Liberty	Little Wave	Little Wave Cliik	Little Wave Cliik XP	Little Wave Flip	Little Wave XP	Rogue	Rogue ALX	Rogue TTL	Rogue XP	Rogue2	Spark	Tsunami	Tsunami ALX	
<b>Pandhora</b>																							
Ep3 Mini, Ep3Cross, Ep3 Smart, Ep3 Standard, Ep3 Plus, Ep3 Enduro, Ep30ffRoad		✓	✓	✓	✓	✓	✓				✓			✓		✓		✓	✓	✓			✓
<b>Rehasense</b>																							
PAWS		✓		✓	✓	✓	✓			✓	✓				✓			✓		✓	✓		
<b>Sunrise Medical</b>																							
Empulse F55							✓								✓							✓	
Jay Jay3 backs															✓		✓						
Nitrum, Helium, Argon, Life T	✓																						
Quickie Wheeldrive															✓								
Quickie/Sopur Attitude															✓		✓						
<b>TriRide</b>																							
Base (12" e 14") Kids, Pieghevole, Special Light ( 12" e 14 ") Special, Compact, L14, HP16, Mad Max, T-Rocks				✓		✓	✓				✓				✓							✓	
<b>Yamaha</b>																							
Navione				✓	✓	✓	✓		✓		✓			✓	✓								