

## Certificate Regarding Air Transportation by Passengers with Walking Disabilities

We hereby certify that the product(s):

**Klaxon Klick Power**  
48 V / 5,8 Ah / 288Wh  
[11,0 kg Total Gross weight / 2.00 kg Net weight]

**Klaxon Klick Race**  
48 V / 5,8 Ah / 288Wh  
[15,0 kg Total Gross weight / 2.00 kg Net weight]

**Klaxon Klick Monster**  
48 V / 5,8 Ah / 288Wh  
[16,0 kg Total Gross weight / 2.00 kg Net weight]

**Klaxon Klick Electric**  
48 V / 5,8 Ah / 288Wh  
[11,0 kg Total Gross weight / 2.00 kg Net weight]

**Klaxon Klick Light**  
48 V / 5,8 Ah / 288Wh  
[9,0 kg Total Gross weight / 2.00 kg Net weight]

**Klaxon Klick Power Plus**  
48 V / 5,8 Ah / 288Wh  
[14,0 kg Total Gross weight / 2.00 kg Net weight]

**Klaxon Klick Scrambler**  
48 V / 5,8 Ah / 288Wh  
[13,0 kg Total Gross weight / 2.00 kg Net weight]

**Klaxon Twist**  
36 V / 4,0 Ah / 144Wh  
[8,0 kg Total Gross weight / 1,00 kg Net weight]

**Klaxon Klick Power**  
48 V / 2,9 Ah / 144Wh  
[11,0 kg Total Gross weight / 1.00 kg Net weight]

**Klaxon Klick Race**  
48 V / 2,9 Ah / 144Wh  
[15,0 kg Total Gross weight / 1.00 kg Net weight]

**Klaxon Klick Monster**  
48 V / 2,9 Ah / 144Wh  
[16,0 kg Total Gross weight / 1.00 kg Net weight]

**Klaxon Klick Electric**  
48 V / 2,9 Ah / 144Wh  
[11,0 kg Total Gross weight / 1.00 kg Net weight]

**Klaxon Klick Light**  
48 V / 2,9 Ah / 144Wh  
[9,0 kg Total Gross weight / 1.00 kg Net weight]

**Klaxon Klick Power Plus**  
48 V / 2,9 Ah / 144Wh  
[14,0 kg Total Gross weight / 1.00 kg Net weight]

**Klaxon Klick Scrambler**  
48 V / 2,9 Ah / 144Wh  
[13,0 kg Total Gross weight / 1.00 kg Net weight]

**Klaxon Twist R**  
36 V / 4,0 Ah / 144Wh  
[8,0 kg Total Gross weight / 1,00 kg Net weight]

including Lithium-ion batteries of type **UN3481, P.I. 967 Section II** (which meets the requirements of UN *Manual of Tests and Criteria*, Part III subsection 38.3), may be carried on passenger aircraft without restriction according to **IATA 2.3.2.4**.

The conditions necessary to meet the above listed regulation are:

- The drive unit must be turned off with the power switch and removed from the wheelchair.
- The battery unit *and/or* drive unit must be carried in the passenger cabin and stored in a manner to prevent powering on and damage by the movement of baggage, mail or other cargo.

- The pilot-in-command must be informed of the location of the lithium battery unit in the cabin.
- A passenger may carry a maximum of one spare lithium ion battery not exceeding 300 Wh or two spare batteries each not exceeding 160 Wh.

**It is strongly recommended that you make advanced arrangements with each airline operator and receive approval prior to departure date. This certificate does not apply to damaged or defective batteries.**

KLAXON MOBILITY GmbH  
Industriestraße 1  
9601 Arnoldstein, Austria  
FN 442221a, UID-Nr.: ATU69961437

Arnoldstein, 01/01/2025

