

GMDN codes

GMDN code	GMDN text	Medical device
13088	<p>Portable bathing bath</p> <p>A bathing (washing) receptacle designed to be easily moved (e.g., it may be readily dismantled and packed for transport). It is typically filled with warm water to accommodate the occupant's full body for bathing in a seated or reclined position. The device is typically used in settings where there is no access to conventional bathing facilities (e.g., field military hospitals, isolated dwellings).</p>	Orca Penguin
31092	<p>Bathing stretcher</p> <p>A manually-operated device consisting of a platform mounted on a frame designed to support a recumbent patient for bathing. It is designed to be submerged in an appropriate bath. The device is used primarily for a patient incapable of supporting them self in an upright position for bathing.</p>	Manatee
31119	<p>Basic walker, non-foldable</p> <p>An all- or partially-wheeled mobility aid in the form of a waist-high, non-foldable framework with handgrips and four or three legs designed to provide a stable support to a person with a disability, a geriatric or an infirm person when standing or ambulating; it may also be used during rehabilitation. It is constructed of lightweight metal, plastic, and rubber materials in the shape of a mobile support structure which the user holds and pushes in front of them as they walk. It is non-foldable and includes a braking mechanism; it may be width-adjustable and some types may include a seat on which the user can rest.</p>	Pony
31162	<p>Assistive pushchair</p> <p>A wheeled mobility device manually propelled by an attendant that consists of an open framework with one or more handle(s) designed to carry one or more person(s) with a disability in a seated or semi-lying (semi-recumbent) position.</p>	Cricket Stingray
37201	<p>Stand-up bed</p> <p>A mains electricity (AC-powered), motorized bed used for a person with a disability or a bedridden patient to enable them to achieve a standing position whilst being held in the bed. It can be raised from a horizontal (flat) to a vertical (standing) position and helps prevent undesirable effects, e.g., pressure sores; contracture in the hip, knee and ankle joints; or deterioration of cardio-circulatory and respiratory performance.</p>	Caribou
37951	<p>Basic walker, foldable</p> <p>An all- or partially-wheeled mobility aid in the form of a waist-high, foldable framework with handgrips and four or three legs designed to provide a stable support to a person with a disability, a geriatric or an infirm person when standing or ambulating; it may also be used during rehabilitation. It is constructed of lightweight metal, plastic, and rubber materials and in the form of a mobile support structure which the user holds and pushes in front of them as they walk. It is foldable for storage and transportation, and includes a braking mechanism; it may be width-adjustable and some types may include a seat on which the user can rest.</p>	Crocodile Mustang

GMDN code	GMDN text	Medical device
40539	<p>Commode</p> <p>A special kind of chair with a built-in receptacle for toileting away from a bathroom or fixed toilet and is designed for persons with a disability. Included in this group of devices are sanichairs and commode chairs used for showering. It may or may not be supplied with wheels (castors).</p>	<p>Flamingo Flamingo High-low Heron Swan</p>
41622	<p>Wheelchair, attendant/occupant-driven, rear-wheels-operated, collapsible</p> <p>A wheeled personal mobility device that incorporates a seat-support system for a person with a disability or a person without the full capacity to walk (not bariatric) designed to be manually propelled by the user while seated in the device or by an attendant. The occupant moves the device by rotating both rear wheels. An attendant moves the device by pushing or pulling handles on the device. The device may be disassembled or folded for transport.</p>	<p>Cougar Kudu Combi Frame:x Multi Frame Multi Frame:x</p>
58827	<p>Assistive ergonomic chair</p> <p>A manually-operated device intended to support and position a patient (disabled and/or discomforted) in a sitting posture while performing activities of daily living, including work-related activities. It is designed with ergonomic adjustments to provide improved personal body posture and/or to help prevent/reduce physical problems associated with poor mobility or sedentary lifestyle. It typically includes a variety of mechanical operator functions (e.g., height and position adjustments), often pneumatically-powered, and has wheels for ease of mobility. It is not a wheelchair.</p>	<p>High-low High-low:x Strong Base Wombat Living Wombat Solo Scallop Panda Futura Panda Futura 5 x:panda</p>
61700	<p>Manual stander</p> <p>A manually-powered, mobile device designed to support and lift a disabled (paraplegic) patient to an upright standing position to provide the health benefits associated with standing (e.g., improved breathing, circulation, and flexibility); it is not intended to facilitate gait. It typically consists of a seat, sling and/or platform with appropriate supports (e.g., footplates, kneepads, chest pad), wheels, and a manual mechanical lifting mechanism whereby the patient is raised from a supine, prone, and/or seated position by themselves or with the assistance of a healthcare provider. It includes wheels for reorientation and may include a manual mechanism to allow a degree of mobility.</p>	<p>Gazelle Meerkat Rabbit Up Toucan</p>