

Compliance statement

Company name:	Etac A/S
Postal address:	Parallevej 3
Postcode:	DK-8751
City:	Gedved
Telephone number:	+45 796 858 33

Statement Regarding Air Transportation

Etac A/S declare that the compliance statement is issued under our sole responsibility and belongs to the following product:

Device(s)	<p>Molift Smart 150 Molift Mover 180 Molift Mover 205 Molift Mover 300 Molift Partner 255</p>
Statement	<p>The above mentioned models may contain Lithium-Ion batteries.</p> <p>Batteries is tested in accordance to the mentioned standards with positive results and may be carried-on passenger aircraft without restriction according to IATA 2.3.2.4</p> <p>Type:</p> <p style="padding-left: 40px;">Lithium-Ion Battery Pack Mass: 0,800 kg Watt-hour rating: 73,25 Wh</p> <p>Standards:</p> <p style="padding-left: 40px;">UN Manual of Tests and Criteria Part III, Subsection 38.3 Lithium Metal and Lithium Ion Batteries, UN ST/SG/AC.10/11/Rev.8</p> <p style="padding-left: 40px;">EN IEC 62281:2019+A1:2021; IEC 2281:2019+A1:2021</p> <p>The conditions necessary to meet the above listed regulations are:</p> <ul style="list-style-type: none"> - The battery unit and/or drive unit must be carried in the passenger cabin and stored in a manner to prevent powering on and damage by the movement of baggage, mail or other cargo. - The pilot-in-command must be informed of the location of the lithium battery unit in the cabin. - A passenger may carry a maximum of one spare lithium ion battery not exceeding 300 Wh or two spare batteries each not exceeding 160 Wh. <p>It is strongly recommended that the customer make advanced arrangements with each airline operator and receive approval prior to departure date. This statement does not apply to damaged or defective batteries.</p>

Approval:

Quality Specialist: Sampavi Rajkumar

Date: 16/2/2026

Signature: 