



# Upright EZ Planar


## Standard Specifications

### MEASUREMENTS:

INCHES

HCPCS:E1236

\* - Plus or minus .5 inches due to fabric  
 Δ - Weight of chair is without front rigging



	EZ 12	EZ 14	EZ 16	EZ 18
Seat Width*	12	14	16	18
Seat Depth	8 – 13	9 – 14	10 – 15	13 – 18
Seat Back Height	16.5	19	20	28
Seat to Footplate*	6 - 20	5 – 22	6 – 24	4.5 – 24.5
Seat to Floor*	21.5	23	23	24.5
Back Angle	100°	100°	100°	100°
Seat Angle	10°	10°	10°	10°
Overall Height	38	38	36	48
Overall Length	43	46	49	47.5
Overall Width	21.5	23	26	26
Folded Length	45.5	45.5	45.5	45.5
Folded Width	14	13.5	14	14.5
Folded Height	14.5	15.5	16	15.5
Headrest Extension	4.5	5	6	6
Shoulder Strap Height	12 & 21	13 & 24	14.5 & 25	17 & 28
Lateral Support Width Adjustment	5 – 10.75	7 – 12.75	9 – 14.75	11 – 16.75
Weight Capacity Regular/Transit	75/66 lbs	100/100 lbs	170/170 lbs	250/170 lbs
Weight of Chair Δ	32 lbs	35 lbs	37 lbs	49 lbs