


Y-hip belt

When considering using vests and belts for a child, it is important to choose the correct size in order to ensure the correct positioning and fit. When measuring the child, the measurements should be taken with a flexible tape measure in order to follow the curves of the body.



The y-hip belt maintains a stable pelvic position by providing support at two angles, down and back. It aims to maintain a neutral pelvic position and prevent the child from tilting the pelvis out of position and sliding forward in the seat.

For the y-hip belt the measurement should be taken from the hip bone (H) to hip bone (H).

	XS	S	M	L	XL
Hip (H) to hip (H)	6¾" – 12½"	8¼" – 13¾"	9¾" – 16¼"	12½" – 17¼"	16¼" – 20½"