

# Convoid EZ Rider Transit- HCPCS E1236

## Retail Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ ATP/RTS: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

### EZ Rider Standard Features:

- 10° Fixed Tilt
- Silver Powder Coated Frame
- Removable, Height Adjustable Swing-Away Footplates
- Tool-less Adjustable Seat Depth Growth of 5"
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire (EZ12/14/16)
- 8" Front, 11" Rear Solid Tire (EZ18)
- Quick Release Rolko Rear Wheels (EZ12/14/16)
- One Piece Height Adjustable Push Handle (EZ12/14/16)
- Foot Operated Wheel Lock (EZ12/14/16)
- Zippered Pockets in Seat & Back Upholstery (Cordura Only EZ12/14/16)

- Tool-less Upholstery Removal
- Toggle Wheel Locks (EZ18)
- H-Harness with Padded Covers
- 3-Pt Positioning Belt
- Foot Positioners
- Tension Adjustable Back
- Self-Tensio®
- Two Piece Push Handles (EZ18)
- Made in the USA

### EZ Rider Transit Features:

- Transit Anchor Brackets
- Headrest Extension
- H-Harness w/ Padded Covers & 3-PT Positioning Belt
- Foot Positioners

### EZ Rider Convertible Features:

- Mag Solid Tires
- Quick Release Axles
- Axle Attachment Brackets

### Standard Models:

- EZ Rider Transit EZ12T - 904856
- EZ Rider Transit EZ14T - 904857
- EZ Rider Transit EZ16T - 904858
- EZ Rider Transit EZ18T - 904859

### Price

- \$2,731**
- \$2,731**
- \$2,731**
- \$2,731**

### Convertible Self-Propel Transit Model

- EZ Rider Transit EZ12T - C905506
- EZ Rider Transit EZ14T - C905507
- EZ Rider Transit EZ16T - C905508
- EZ Rider Transit EZ18T - C905509

### Price

- \$3,356**
- \$3,356**
- \$3,356**
- \$3,356**

### Upholstery Colors

- Cordura Upholstery
  - Panther Black
  - Royal Blue
  - Sassy Purple
  - Princess Pink
  - Forest Green
  - Cherry Red
  - Chocolate Brown
  - Aqua Blue
  - Apple Green
- Textilene Upholstery
  - Navy Blue
  - Violet Purple
  - Wild Raspberry
  - Teal Green

### Price

- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$0
- \$280
- \$280

	EZ12	EZ14	EZ16	EZ18
C903854	C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851	C903851
C903759	C903759	C903759	C903759	C903759
C903857	C903857	C903857	C903857	C903857
C903855	C903855	C903855	C903855	C903855
C904189	C904189	C904189	C904189	C904189
C904753	C904753	C904753	C904753	C904753
C904754	C904754	C904754	C904754	C904754
C903464	C903464	C903464	C903464	C903464
C903465	C903465	C903465	C903465	C903465
C903467	C903467	C903467	C903467	C903467
C903466	C903466	C903466	C903466	C903466
C900835	C900321	C901031	C900392	C900392
C904494	C904495	C904496	C904497	C904497

### Frame Option

- Height Adjustable Flip-Up Armrests (pair)\*
- Two Piece Push Handle<sup>1</sup>
- Caster Lock For 2" Wide Front Wheels<sup>1</sup>
- Rear Anti-Tip Tubes (pair)

### Price

- \$212
- \$0
- \$98
- \$109

### HCPCS

- E0973
- 
- K0073
- E0971

	EZ12	EZ14	EZ16	EZ18
C904523	C904523	C904524	C904524	C904524
C904510	C904511	C904512	C904513	C904513
C902920	C902920	C902920	C902920	C902920
C900377	C900377	C900378	C900378	C900378

<sup>1</sup>Manufacturer Installation Required.

\*See EZ Rider Convertible Options for item compatible with the Convertible Self-Propel Package.

### Transit Required Items

- Head Support - Select Style Below:
- Headrest Extension - Rectangular - Textilene\*\*
  - Headrest Extension - Curved - Cordura (Not available in Planar Option)
  - Contoured Firm Headrest (Planar Option Only/Hrdw Req)
  - Adj. Removable Hardware (Req for Contoured Headrest)
- Trunk Positioning - Select Style Below:
- 5-Point Harness with Padded Covers (EZ12 only)\*<sup>1</sup>
  - H-Harness with Padded Covers (3-Pt Belt Required)
  - 3-Pt Positioning Belt
  - Full Torso Support Vest - Cordura (Not available in Planar Option)
  - Full Torso Support Vest - Textilene
- Foot Positioners - Select Size Below:
- Foot Positioners - Small 7" - 10" (pair)
  - Foot Positioners - Medium 8" - 11" (pair)
  - Foot Positioners - Large 9" - 12" (pair)

### Price

- \$0
- \$0
- \$254
- \$257
- \$150
- \$0
- \$0
- \$275
- \$275
- \$0
- \$0
- \$0

### HCPCS

- E0966
- E0966
- E0955
- E1028
- E0960
- E0960
- E0978
- K0108
- K0108
- 
- 
- 

	EZ12	EZ14	EZ16	EZ18
C900850B	C900685B	C901029B	C900390B	C900390B
C902863B	C900484B	C900138B	C902619B	C902619B
C904032	C904033	C904034	C904034	C904034
C904035	C904035	C904035	C904035	C904035
C903702				
C903988B	C903988B	C903988B	C903988B	C903988B
C904017B	C904017B	C904017B	C904017B	C904017B
C900848	C900338	C901020	C905501	C905501
C905437	C905438	C905439	C900387	C900387
C903415B	C903415B	C903415B	C903415B	C903415B
C903416B	C903416B	C903416B	C903416B	C903416B
C903417B	C903417B	C903417B	C903417B	C903417B

<sup>1</sup>Manufacturer Installation Required.

\*Crash Tested - WC19 Compliant. Not available with Reinforced Upholstery, Cushion, Planar Option or H-Harness and 3 Pt. Belt.

\*\*Required when ordering the Contoured firm headrest and flat firm backrest extension.

<u>Planar Option*</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Planar Seat Modification	\$359	
<input type="checkbox"/> Planar Back with Attachment (Incl. 1 Backrest Ext. w/ Hrdw.)	\$333	E2291
Cushion - select one cushion below:		
<input type="checkbox"/> Support - General Use Cushion (Solid Seat Insert Included) <sup>1</sup>	\$251	E2601
<input type="checkbox"/> Position - Medial Thigh Support/Anti Thrust (Solid Seat Insert Included) <sup>1</sup>	\$486	E2607
<input type="checkbox"/> Align - Lateral Pelvis/Lateral Thigh Support (Solid Seat Insert Included) <sup>1</sup>	\$486	E2607
Planar Head Support Options - Must Select One		
<input type="checkbox"/> Contoured Firm Headrest (Hardware Required) <sup>1</sup>	\$254	E0955
<input type="checkbox"/> Adj. Removable Hardware (Req for Contoured Headrest) <sup>1</sup>	\$257	E1028
<input type="checkbox"/> Flat Firm Backrest Extension (Hardware Required) <sup>1</sup>	\$238	E0955
<input type="checkbox"/> Adj. Removable Hardware (Req for Flat Firm Backrest Ext.) <sup>1</sup>	\$233	E1028

<sup>1</sup>Manufacturer Installation Required.

\*Planar Option must be ordered with Planar Seat, Planar Back, Cushion, Planar Head Support, and Textilene Upholstery.

EZ12	EZ14	EZ16	EZ18
C904049	C904050	C904051	C904052
C904053	C904054	C904055	C904056
C903684	C903687	C903690	C903693
C903686	C903688	C903692	C903695
C903685	C903689	C903691	C903694
C904032	C904033	C904034	C904034
C904035	C904035	C904035	C904035
C903928	C903582	C902151	C903674
C904057	C904057	C904057	C904057

<u>Tire Options</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 12 1/2" Solid Knobby Rear Tires*	\$208	
<input type="checkbox"/> 7 1/2" x 2" Front, 12 1/2" Rear Pneumatic Knobby Tires*	\$311	
<input type="checkbox"/> 8" x 2" Front Ribbed, 12 1/2" Rear Solid Knobby Tires*	\$311	
<input type="checkbox"/> 8" x 2" Front Ribbed, 12 1/2" Rear Pneumatic Knobby Tires*	\$311	
<input type="checkbox"/> 8" x 2" Solid Front Tires	\$162	
Convertible Self-Propel Package		
<input type="checkbox"/> 20" Mag Solid Tires	\$0	
<input type="checkbox"/> 22" Mag Solid Tires	\$0	
<input type="checkbox"/> 24" Mag Solid Tires	\$0	
<input type="checkbox"/> Quick Release Axles	\$0	
<input type="checkbox"/> Axle Attachment Brackets	\$0	

\*Not available with the Convertible Self-Propel Package.

^With this option, standard tires are 6" Front, 8" Rear Solid Tires (EZ12/14) or 8" Front, 8" Rear Solid Tires (EZ16/18).

EZ12	EZ14	EZ16	EZ18
C904117	C904117	C904117	
C904118	C904118	C904118	
			C903084
			C903083
			C903080
C902940B			
	C902946B		
		C902946B	C902946B
C902942B	C902942B	C902942B	C902942B
C900377	C900377	C900378	C900378

<u>Head Supports</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Padded Headwings^	\$254	E0955
<input type="checkbox"/> Occi Headwings^	\$238	E0955

^Not available with Planar Option.

EZ12	EZ14	EZ16	EZ18
C900849	C900682	C900683	C900684
C905271	C905271	C905271	C905272

<u>Trunk Positioning</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Soft Adjustable Lateral Support Single Flap w/ Scoli Strap - Cordura^	\$142	E0956
<input type="checkbox"/> Soft Adjustable Lateral Support Double Flap w/ Scoli Strap - Cordura^	\$142	E0956
<input type="checkbox"/> Soft Adjustable Lateral Support Single Flap w/ Scoli Strap - Textilene^	\$142	E0956
<input type="checkbox"/> Soft Adjustable Lateral Support Double Flap w/ Scoli Strap - Textilene^	\$142	E0956
<input type="checkbox"/> Full Torso Support Vest - Cordura*	\$275	K0108
<input type="checkbox"/> Full Torso Support Vest - Textilene*	\$275	K0108
<input type="checkbox"/> Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")^	\$368	
Planar Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req) <sup>1</sup>		
<input type="checkbox"/> Small = 3"x 5"	\$142	E0956
<input type="checkbox"/> Medium = 3.75"x 5.5"	\$142	E0956
<input type="checkbox"/> Large = 4.75"x 6"	\$142	E0956
Planar Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req) <sup>1</sup>		
<input type="checkbox"/> Small = 3"x 5"	\$142	E0956
<input type="checkbox"/> Medium = 3.75"x 5.5"	\$142	E0956
<input type="checkbox"/> Large = 4.75"x 6"	\$142	E0956
<input type="checkbox"/> Planar Right Adj. Removable Hardware (Req for Lateral Support) <sup>1</sup>	\$185	E1028
<input type="checkbox"/> Planar Left Adj. Removable Hardware (Req for Lateral Support) <sup>1</sup>	\$185	E1028

<sup>1</sup>Manufacturer Installation Required.

^Not available with Planar Option.

\*Not available with Single or Double Flap Lateral Support.

EZ12	EZ14	EZ16	EZ18
C900846	C900327	C901018	C900386
C900847	C900650	C900651	C900652
C902859	C900678	C900140	C902614
C902860	C900679	C900204	C902615
C900848	C900338	C901020	C905501
C905437	C905438	C905439	C900387
C904211	C904110	C904110	C904110
C902965	C902965	C902965	C902965
C902969	C902969	C902969	C902969
C902973	C902973	C902973	C902973
C904462	C904462	C904462	C904462
C904463	C904463	C904463	C904463
C904464	C904464	C904464	C904464
C904030	C904030	C904031	C904031
C904030	C904030	C904031	C904031

<u>Pelvic Positioning</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Support - General Use Cushion (Solid Seat Insert Included) <sup>1</sup> *	\$251	E2601
<input type="checkbox"/> Position - Medial Thigh Support/Anti Thrust (Solid Seat Insert Included) <sup>1</sup> *	\$486	E2607
<input type="checkbox"/> Align - Lateral Pelvis/Lateral Thigh Support (Solid Seat Insert Included) <sup>1</sup> *	\$486	E2607
<input type="checkbox"/> Incontinence Cushion Liner (Cushion req) <sup>1</sup>	\$75	E2619
<input type="checkbox"/> Custom Seat Depth Extension 3"^(1)	\$98	E1297

<sup>1</sup>Manufacturer Installation Required.

^Not available with Planar Option.

\*Not available with Seat Depth Extension or Cordura Upholstery.

EZ12	EZ14	EZ16	EZ18
C903684	C903687	C903690	C903693
C903686	C903688	C903692	C903695
C903685	C903689	C903691	C903694
C902324	C902324	C902324	C902324
C904003	C904005	C904007	C904009

<u>Foot &amp; Leg Positioning</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Medial Thigh Support-Abductor (N/A Position/Align Cushion)	\$169	E0957
<input type="checkbox"/> Lateral Thigh Support-Adductor (N/A Position/Align Cushion)	\$169	E0956
<input type="checkbox"/> Calf Panel*	\$64	E0995
<input type="checkbox"/> Angle Adjustable Footplates (pair)	\$186	K0040
<input type="checkbox"/> Footplate Securement Strap	\$32	

\*Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.

EZ12	EZ14	EZ16	EZ18
C900843	C900329	C901024	C900388
C900844	C900330	C901026	C900389
C904580	C904580	C904580	C904580
C903268	C903268	C903269	C903269
C903584	C903584	C903584	C903584

<u>EZ Rider Convertible Options</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 8" x 2" Solid Front Tires	\$162	
<input type="checkbox"/> Height Adjustable Flip-Up Armrests (pair)	\$212	E0973
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray)^	\$116	E0950
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$237	E1028

\*Requires Convertible Self-Propel Package.

^Flip-Up Armrests and Hardware required.

EZ12	EZ14	EZ16	EZ18
		C903080	C903080
C904044	C904044	C904045	C904045
C903393	C903394	C903395	C903396
C903361	C903361	C903361	C903361

- Transit Options**
- Q'Straint Transit Lap Belt **Price \$94 HCPCS E0978**
  - Sure-Lok Transit Lap Belt **Price \$94 HCPCS E0978**
- ^Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.**

EZ12	EZ14	EZ16	EZ18
C903094	C903094	C903094	C903094
C903232	C903232	C903232	C903232

- Accessories**
- Basic Headrest Cover (Canopy) **Price \$180 HCPCS**
  - Headrest Cover (Canopy) **Price \$238**
  - Extended Headrest Cover - No Windows (Canopy) **Price \$260**
  - Extended Headrest Cover - with Windows (Canopy) **Price \$307**
  - Upper Extremity Support Surface (Clear Tray)\*^ **Price \$116 HCPCS E0950**
  - Adj. Removable Hardware (Required for Clear Tray)\* **Price \$237 HCPCS E1028**
  - Upper Extremity Support Surface (Foam Tray) Small^ **Price \$158 HCPCS E0950**
  - Upper Extremity Support Surface (Foam Tray) Medium^ **Price \$158 HCPCS E0950**
  - Upper Extremity Support Surface (Foam Tray) Large^ **Price \$158 HCPCS E0950**
  - Wheelchair Travel Bag **Price \$215**
  - Travel Bag for Large Wheels\* **Price \$135**
  - Utility Bag **Price \$40**
  - Rain Cover Size 1 (Canopy Req) **Price \$157**
  - Rain Cover Size 2 (Canopy Req) (EZ14 Extended Canopy Req) **Price \$157**
  - Rain Cover Size 3 (Extended Canopy Req) **Price \$157**
  - Mosquito Net Size 1 (Canopy req) **Price \$157**
  - Mosquito Net Size 2 (Canopy req) (EZ14 Extended Canopy Req) **Price \$157**
  - Mosquito Net Size 3 (Extended Canopy req) **Price \$157**

EZ12	EZ14	EZ16	EZ18
C905021	C905022	C905023	C905024
C900853	C901210	C901205	C901211
C904721	C904725	C904726	C904727
C904737	C904741	C904742	C904743
C904525	C904526	C904527	C904528
C903369	C903369	C903369	C903369
C903669			
	C903670	C903670	
			C903671
C902022	C902022	C902022	C902022
C902484	C902485	C902486	C902487
C904531	C904531	C904531	C904531
C904939	C904939		
	C904940	C904940	C904940
			C905283
C904655	C904655	C904655	
	C904656	C904656	C904656
			C905319

**<sup>1</sup>Manufacturer Installation Required.**  
**\*See EZ Rider Convertible Options for item compatible with the Convertible Self-Propel Package.**  
**^Flip-Up Armrests and Hardware required.**

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

The EZ Rider transit option includes four manufacturer installed wheelchair transportation anchors. The EZ Rider transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are EZ12: 66 lbs./30kg.; EZ14: 100 lbs./45.5kg; EZ16 and EZ18: 170 lbs./77kg. The Convaid EZ Rider was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convaid is unknown, and no warranty or claim about their performance is being, or can be, made.

**Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email [convaidsales.us@etac.com](mailto:convaidsales.us@etac.com)**

Prices in USD & subject to change, ex works Torrance, California.

Effective: 6/1/2022

Convaid  
 2830 California Street  
 Torrance, CA 90503  
 1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

