

Convoid Safari Tilt - HCPCS E1232

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Safari Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire (SFT16)
- 6" Front, 11" Rear Solid Tires (SFT18)
- Quick Release Rear Wheels (SFT16)
- Wheel Locks
- 2-Point Positioning Belt
- Heavy-Duty Back Upholstery
- One Piece, Non-Folding/Non-Elevating Footplate
- Support Cushion with Solid Seat Insert
- Made in the USA

Standard Models:

	Price
<input type="checkbox"/> Safari SFT16 - C900526	\$3,050
<input type="checkbox"/> Safari SFT18 - C902432	\$3,050

Upholstery Colors

	Price
<input type="checkbox"/> Panther Black	\$0
<input type="checkbox"/> Royal Blue	\$0
<input type="checkbox"/> Sassy Purple	\$0

Frame Option

	Price	HCPCS
<input type="checkbox"/> Anatomic Back	\$441	
<input type="checkbox"/> Recline Option	\$608	E1225
<input type="checkbox"/> Height Adjustable Flip-Up Armrests (pair)	\$301	E0973
<input type="checkbox"/> Attendant Hand Brakes (pair)*	\$311	
<input type="checkbox"/> Solid Side Panels (N/A w/ Recline Option)	\$278	
<input type="checkbox"/> Curb Tipper (elevates front of chair to go up curbs)	\$105	

**Available with Knobby Tires only. Not available with Recline Option.*

Transit Required Items

	Price	HCPCS
<input type="checkbox"/> Transit Option (Manufacturer Installation Required)	\$264	
<input type="checkbox"/> Headrest Extension	\$138	E0966
<input type="checkbox"/> H-Harness with Padded Covers	\$105	E0960
<input type="checkbox"/> 3-Pt Positioning Belt	\$52	E0978
<input type="checkbox"/> Elevating Legrests with Swing-Away Footplates (pair)	\$356	E0990
Foot Positioners - Select Size Below:		
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	\$167	
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	\$167	
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	\$167	

Transit Options

	Price	HCPCS
<input type="checkbox"/> Q'Straint Transit Lap Belt	\$198	E0978

Head Supports

	Price	HCPCS
<input type="checkbox"/> Padded Headwings	\$260	E0955
<input type="checkbox"/> Occi Headwings	\$244	E0955
<input type="checkbox"/> Headrest Extension	\$138	E0966

Trunk Positioning

	Price	HCPCS
<input type="checkbox"/> H-Harness with Padded Covers (3-Pt Positioning Belt Req)	\$105	E0960
<input type="checkbox"/> Full Torso Support Vest*	\$348	
<input type="checkbox"/> Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")	\$405	

Pelvic Positioning

	Price	HCPCS
<input type="checkbox"/> 3-Pt Positioning Belt	\$52	E0978
<input type="checkbox"/> Position Cushion - Medial Thigh Support/Anti Thrust*	\$458	E2607
<input type="checkbox"/> Align Cushion - Lateral Pelvis/Lateral Thigh Support*	\$458	E2607
<input type="checkbox"/> Incontinence Cushion Liner (Cushion and Solid Seat req)	\$83	E2619
<input type="checkbox"/> Solid Seat Insert (Required with Cushion)	\$157	E0992

**Solid Seat Insert required with any Cushion.*

SFT16	SFT18
C903854	C903854
C903850	C903850
C903851	C903851

SFT16	SFT18
C901319	C902476
C901321	
C901322	C902477
C904064	C904064
C904392	C904393
C904063	C904063

SFT16	SFT18
C900527	C902435
C901310	C902471
C903988	C903988
C904017	C904017
C901299	C902460
C903415	C903415
C903416	C903416
C903417	C903417

SFT16	SFT18
C903095	C903095

SFT16	SFT18
C901309	C902470
C905271	C905273
C901310	C902471

SFT16	SFT18
C903988	C903988
C901308	C902469
C904110	C904110

SFT16	SFT18
C904017	C904017
C903590	C903591
C903594	C903595
C902324	C902324
C904029	C904029

Foot & Leg Positioning

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Medial Thigh Support (Abductor)	\$213	E0957
<input type="checkbox"/> Lateral Thigh Support (Adductor)	\$213	E0956
<input type="checkbox"/> Elevating Legrests w/ Swing-Away Footplates (pair)	\$356	E0990
<input type="checkbox"/> Legrest Extension - 3"*	\$73	
<input type="checkbox"/> Calf Panel	\$66	E0995
<input type="checkbox"/> Angle Adj. Footplates - Adult Size (pair)^	\$235	K0040
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)^	\$167	
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)^	\$167	
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)^	\$167	
<input type="checkbox"/> Footplate Securement Strap	\$33	

***3" Legrest Extension interferes with caster swivel when seat is in upright position.**

Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.

^Elevating Legrests with Swing-Away Footplates required.

SFT16	SFT18
C901297	C902458
C901298	C902459
C901299	C902460
C901324	C902478
C904580	C904580
C903271	C903271
C903415	C903415
C903416	C903416
C903417	C903417
C903584	C903584

Tire Options

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 12.5" Rear Solid Knobby Tire	\$229	
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire	\$342	
<input type="checkbox"/> 12.5" x 2" Rear Solid Knobby Tire	\$230	

SFT16	SFT18
C904117	
C904118	
	C903420

Accessories

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Headrest Cover (Canopy)	\$262	
<input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy)	\$286	
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray)^	\$128	E0950
<input type="checkbox"/> Adj. Removeable Hardware (Required for Clear Tray)	\$261	E1028
<input type="checkbox"/> Wheelchair Travel Bag	\$237	
<input type="checkbox"/> Utility Bag	\$45	
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req)	\$173	

SFT16	SFT18
C901313	C902475
C904734	C904735
C903403	C903404
C904519	C904519
C902022	C902022
C904531	C904531
C904656	C904656

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convaidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.

Effective: 1/10/2024

Convoid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF0700-ENG-REV.09

Creating Possibilities