

Standard configuration

- Rear seat height: 40-49 cm (preset 45 cm).
- Back support: Cross back 3A, height 32-45 cm (preset 40 cm).
- Back support angle: Standard 2°. Settable -5° - +20° (std +2°).
- Back support upholstery: Plush.
- Rear wheels: 24" high pressure tyre.
- Handrims: Aluminium.
- Camber angle: 2°.
- Castors: 165x30 mm, solid (6,5").
- Front fork: Medium fork with steel axle.
- Front fork attachment: Settable in height and angle.
- Leg support: Black, wide angle.
- Foot support: Flip-up.
- Arm support: Length 25 cm, settable in height.
- User brakes: High mounted, with red brake handles.
- Anti tip: Yes.
- Push handle: Height adjustable.



Define size

Seat width (cm)	35	37.5	40	42.5	45	47.5	50	52.5
-----------------	----	------	----	------	----	------	----	------

Frame length	Standard	Long	Item no
--------------	----------	------	---------

Widening kit (15 mm)	Mounted	Order as extra	Item no
		pcs	27533

Frame colour	Grey -78	Black -01	White -92
--------------	----------	-----------	-----------

You have now decided your wheelchair. Press "Ready" to save predefined configuration...



...or **continue** to "Most frequent adaptations", such as colour, leg support angle, camber angle, rear wheels and rear seat height.

## Most frequent adaptations

<b>Leg supports</b>	<b>Standard 70°</b>		Narrow 90°		Narrow 90° short with short att. tube		Lockable		Without leg support		<b>Castor housing</b>	<b>Standard</b>	Widened
	R	L	R	L	R	L	No	Yes	R	L			
<b>Rear/front seat height 20" wheel</b>	32/32	32/34	34/34	34/36	36/36	36/38	38/38	38/40	40/40	40/42			
<b>Rear/front seat height 22" wheel</b>	34.5/34.5	34.5/36.5	36.5/36.5	36.5/38.5	38.5/38.5	38.5/40.5	40.5/40.5	40.5/42.5	42.5/42.5	42.5/44.5			
<b>Rear/front seat height 24" wheel</b>	39/39	39/41	41/41	<b>41/43</b>	43/43	43/45	45/45	45/47	47/49	49/51	<b>24" wheel type</b>	<b>Std</b>	Carbon
<b>Front fork thickness</b>	<b>Standard (3 mm)</b>		4 mm										
<b>Balance position</b>	Active		<b>Medium</b>		Passive								
	1	2	3	4	5	6							
<b>Handrims</b>	<b>Aluminium</b>		Stainless steel		Cellular rubber		Vinyl coated		Titanium				
<b>Rear wheel tyre</b>	<b>High pressure</b>		Solid		Schwalbe Rightrun		Schwalbe Marathon plus						

**!** Do you need more advanced seat height/wheel options? Please, download the document "Advanced seat height options" from [www.etac.com](http://www.etac.com)

You have now completed the second step in configuring the wheelchair. You can choose to stop here if you are satisfied or...



More options and accessories.

# Prescription form Etac Cross 6

## Seat adaptions/Accessories

<b>Seat</b>	<b>Upholstery seat</b>		Solid seat	<b>Seat extender</b>	Yes	<b>No</b>	Item no	<b>Seat cushion</b>	Plush	Tidy	<b>No</b>	Item no
<b>Comfort wedge</b>	Yes	<b>No</b>	Item no	<b>Positioning belt</b>	L 100 mm	L 128 mm	<b>No</b>	Item no				

## Back support/head support adaptions/Accessories

<b>Back support</b>	<b>3A standard</b>		3A high	Back support angle (if other than than preset, 2°)				<b>Back support cover</b>	<b>Standard</b>	Comfort	Tidy	Airmesh			
				°											
<b>Back wedges</b>	13x9.5x2	Item no	19x12x2	Item no	20x17x2.5	Item no	Plush 15 cm	Item no	Plush 30 cm	Item no					
	pcs	28030	pcs	28139	pcs	28140	pcs	63012-60	pcs	63013-60					
<b>Side stop</b>	Straight	Item no	Angled 30 mm	Item no											
	pcs	28136	pcs	28173											
<b>Lateral support swing away</b>		R	L	Item no	<b>Lateral support detachable</b>		R	L	Item no	<b>Lateral support others</b>		R	L	Item no	
	0 mm			28181-1 28181-2		0 mm			28142		Soft			27395	
	-15 mm			28180-1 28180-2		-15 mm			28144		Large			27917	
	-45 mm			28141-1 28141-2		-45 mm			28145		Swing-away Medium			27915-1 27915-2	
	+15 mm			28143	Grey plush			28647	Swing-away Large				27916-1 27916-2		
					<b>Cover</b>	<b>No cover</b>									
											<b>No lateral support</b>				

<b>Head support</b>	Traditional				With side support		<b>No</b>	Item no
	One handle		Three handles		One handle	Three handles		
	Plush	Hygiene	Plush	Hygiene	Hygiene	Hygiene		
<b>Head support attachment</b>	Mounted on back support frame		Mounted on bow handle		<b>No attachment</b>		Item no	

**Leg support adaptations/Accessories**

**Elevating leg supports\***

Model 1		Model 3	
R	L	R	L

**Calf support cover, elevating calf support**

Plush	Hygiene	No	Item no

**Other leg supports\*\***

Amputee leg support		Amputee leg support with weight	
R	L	R	L

**Plush Hygiene**

Amputee weight		Plaster leg support	
R	L	R	L

\*(selection replaces std or narrow leg support. To deselect, go to [Leg support options](#))

\*\* (selection replaces std, narrow or elevating leg support. To deselect, select Elevating leg supports or go to [Leg support options](#))

**Foot support adaptations/Accessories**

**Foot support**

Std attachment tube	Long attachment tube	One-piece foot support

**Padding for foot plate**

Right	Left	Both	One piece	No	Item no

**Calf strap**

Yes	No	Item no

**Padding for calf strap**

Yes	No	Item no

**Calf strap, wide**

Yes	No	Item no

**Foot plate extension**

Right	Left	No	Item no

**Foot plate extension Add-ons**

Right	Left	Both	No	Item no	With side protection
					28772

**Footbox**

Yes	No	Item no

**Heel strap**

Yes	No	Item no
		24458

**Arm support adaptations/Accessories**

Arm supports  
/  
Side supports

	R	L	Item no
<b>Arm support L 25 cm</b>			25172-1 25172-2
Arm support L 25 cm with finger screw			26379-1 26379-2
Arm support L 38 cm			25173-1 25173-2
Arm support L 38 cm with finger screw			26380-1 26380-2
Arm support L 25 cm, detachable pad			28780-1 28780-2
Arm support L 25 cm, detachable pad, finger screw			28782-1 28782-2
Arm support L 35 cm, detachable pad			28781-1 28781-2
Arm support L 35 cm, detachable pad, finger screw			28783-1 28783-2
Mud guard			27397-1 27397-2

**Arm support cushion**

	R	L	Item no
Foam, dark grey plush			L=25 62282-60
			L=38 62283-60
Foam, hygiene			L=25 62282-63
			L=38 62283-63
Gel, dark grey plush			L=25 26240-60
			L=38 26241-60
Gel, hygiene			L=25 26240-63
			L=38 26241-63
Hemi			25770-60-1 25770-60-2
Wide and soft			26873
<b>No cushion</b>			
<b>Thigh support</b>			
Basic			
With soft pad			82488-1 82488-2
Knee support			28777-1 28777-2
<b>No thigh support</b>			

**Other adaptations/Accessories**

**Tray**

	R	L	Item no
Tray, for long arm supports			
Tray hemi for short arm support			
Tray hemi for long arm support			
Tray Hemi, foldable			27351-1 27351-2
<b>No tray</b>			

**Non slip cover for Hemi tray**

Item no  
81633

**Push handle/  
Bow handle**

<b>Push handle</b>	Bow handle	Fixed push handle

**User brake**

<b>With red handle</b>	With black handle	With extended handle	Single hand brake	
			Right	Left

**Cross-brace**

Yes	<b>No</b>	Item no

**Drip stand**

Yes	<b>No</b>	Item no
		28323

Other:

**Cane holder**

Yes	<b>No</b>	Item no
		27327

**Single hand drive**

<b>No</b>	Right	Left

**Attendant brake**

Yes	<b>No</b>

**Tool kit**

Yes	<b>No</b>	Item no
		25174

**Spoke guards**

Grey	Yellow	<b>No</b>	Item no

**Tetra quick release**

Yes	<b>No</b>	Item no
		24051-01

You have now completed the last step in configuring the wheelchair.  
Save your choices or...



Return to "Most frequent adaptations".