

# Etac Prio 3A prescription form

1a17652.5 2023-09-29

Applies to Etac Prio with serial no beginning with 38D

③ Seat tilt and back recline function	Manual		Electric
	Std	User operated seat tilt	
Manual or electric		Right	
		Left	

④ Back support system	With ext. 52-68	Cushion/cover			Material	
		Cush-ion	Cush-ion 3D	Cover	Grey plush	Grey hygiene
Height settable without dynamic pelvic support. Height 46-62 cm						
Height settable with dynamic pelvic support. Height 46-62 cm						

⑤ Back support locking	Standard	QR-pin
Standard or quick-release		

⑥ Head support	Knobs		Cover	
	One knob	Three knobs	Grey plush	Grey hygiene
Multi adjustable				
Multi adjustable with side support				
Anatomic medium				
Anatomic large				
Without head support				

⑦ Comfort tilt lock	
Without Comfort tilt lock	
With Comfort tilt lock preset at 0°	
With Comfort tilt lock preset at 4°	
With Comfort tilt lock preset at 8°	
With Comfort tilt lock preset at 12°	

⑧ Arm support		
Model	R	L
D-shaped 49x8.5 depth adj		
D-shaped 49x8.5 depth adj with finger screw		
Solid plate length 38 cm		
S-shaped 38x8 depth adj		

## ① Wheelchair

Seat width	Item. no.	
40	13420405	
45	13420409	
50	13420413	
55	13420417	

## ② Seat depth

	Seat cushion						Seat plate				
	Short	Long	Basic	Com-fort	Grey plush	Grey hygiene	No seat cushion	Wood	Plastic	Wood w. Velcro	Plastic w. Velcro
Pos 1= 39-46 cm short cushion											
Pos 2= 41,5-48,5 cm short cushion											
Pos 3= 44-51 cm short or long cushion											

## ⑨ Leg support


Model	R	L
Elevating		
Fixed angle		
Amputee plush		
Amputee Hygiene		
Amputee weight		
Amputee+weight plush		
Amputee+weight Hygiene		

## ⑩ Foot supports

Model	R	L
Depth 13 cm, long attachment tube (38.5-55 cm from seat plate)		
Depth 13 cm, medium attachment tube (31-47.5 cm from seat plate)		
Solid footplate angle adjustable		
Solid footplate angle adjustable Add-on		


# Seat height and balance position

11 Seat height, rear wheels, castor wheels and position



		Front fork short								Front fork medium								
		S1		S2		S3		M1		M2		M3						
		6"	6"	6.5"	7"	8"	6"	6.5"	7"	8"	6"	6.5"	7"	8"	6"	6.5"	7"	8"
		6"	6"	6.5"	7"	8"	6"	6.5"	7"	8"	6"	6.5"	7"	8"	6"	6.5"	7"	8"
		Seat height																
E	20"	38,0	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	22"	40,5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	24"	42,5	-	-	-	-	-	-	-	*	-	-	-	-	-	-	-	-
D	20"	40,5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	22"	43,0	-	-	-	-	-	-	-	*	-	-	-	-	-	-	-	-
	24"	45,0	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
C	20"	43,0	-	-	-	-	-	-	-	*	-	-	-	-	-	-	-	-
	22"	45,5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	24"	47,5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
B	20"	45,5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	22"	48,0	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	24"	50,0	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	16"	40,0	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
A	20"	48,0	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	22"	50,5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	24"	52,5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	16"	42,5	-	-	-	-	-	-	-	*	-	-	-	-	-	-	-	-
a	16"	45,0	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

12 Balance position



13 Tyre

Solid	
High pressure	

14 Handrim

Aluminium	
Plastic coated	
Stainless steel	

15 Front fork attachment

Standard	
Widened	

16 Castor wheel 8"

8x13/8" (200x35 mm)	
8x2" (200x50 mm)	

17 Single hand drive

Without single hand drive	
With single hand drive right side	
With single hand drive left side	

18 Brakes

Model	R	L
Standard black		
Standard red		
Extended handle		
Single hand brakes right or left mounted		
Attendant brake (std when selecting 16" rear wheel)		
Without attendant brake		

19 Push handle/Bow handle

Push handle, height adjustable	
Bow handle, height adjustable	
Bow handle, height and angle adjustable	

\*Only 8x13/8" (200x35 mm), not 8x2" (200x50 mm).

# Accessories

	R	L	Item. no.
Trunk support detachable without tools 0 mm			28142
Trunk support detachable without tools -15 mm			28144
Trunk support detachable without tools -45 mm			28145
Trunk support Swing-away 0 mm			28181-1 28181-2
Trunk support Swing-away -15 mm			28180-1 28180-2
Trunk support Swing-away -45 mm			28141-1 28141-2
Trunk support Swing-away +15 mm			28143
Trunk support soft			27395
Swing-away medium			27915-1 27915-2
Swing-away large			27916-1 27916-2
Without trunk support			
Cover for trunk support			28647
Without cover			
Side support cushion			63171-60-1 63171-60-2
Without side support cushion			
Arm support padding std			63193-60-1 63193-60-2
Arm support padding solid	Foam		62283-60
	Gel		26241-60
Hemi cushion			25770-60-1 25770-60-2
Cushion wide and soft			28640
Without arm support padding			

Side stop straight 40 mm	pcs	28136
Side stop angled 30 mm	pcs	28173
Back support widening		28786
Without back support widening		
Pressure distributor	pcs	28637
Widening kit +2.5 cm		27892
Without widening kit		
Side support padding		63271-1 63271-2
Without side support padding		
Tray		
Without tray		

# Accessories

		Item. no.
Pommel		28022
Without Pommel		
<b>Calf support</b>		
Without calf support		
<b>Calf support padding</b>		
20 cm		26845
30 cm		26846
40 cm		26847
Without calf support padding		
<b>Calf support wide</b>		
70 cm		63186-60
90 cm		63187-60
Without calf support wide		
<b>Positioning belt</b>		
100 cm	pcs	27904
128 cm	pcs	28544
<b>Attachment positioning belt front</b>		
Without attachment positioning belt front		28641
<b>Chest strap attachment</b>		
Without chest strap attachment		28186
<b>Comfort wedges</b>		
H:2, W:11, L:18.5 cm	pair	28139
H:2.5, W:17.5, L:20 cm	pair	28140
H:2, W:9, L:13 cm	pair	28030
Short plush 150 mm	st	63012-60
Long plush 300 mm	st	63013-60

		Item. no.
Gas piston clips (10 pcs/set)		28163
Without gas piston clips		
<b>Lock seat and back angle controls</b>		
	pcs	28047
<b>Back angle control</b>		
Without back angle control		28573
<b>Basket for oxygen tank</b>		
Without basket		28316
<b>IV stand</b>		
Without IV stand		28341

# Accessories

		Item. no.
Velcro kit seat cushion for separate mounting		28183
Without Velcro kit		
Spoke guards grey		
Spoke guards yellow		
Without spoke guards		
Cane holder	st	27327
	R	L
Extended foot support (21.5 cm)		
Without extended foot support		
Extended footplate Add-on	S	28774
	M	28775-1 28775-2
	L	28776-1 28776-2
Without extended foot support Add-on		
Side protection extended foot support Add-on		28772
Without side protection Add-on		
Foot support padding		
For foot supports model 2		
For solid foot support, angle adjustable		
Without padding		

Heel strap for foot support depth 13 cm		24458
Without heel strap		
Knee support		28777-1 28777-2
Without knee support		
Footbox		
Without footbox		
Memory function head support		28564
Without memory function head support		
Back cushion extender		28787-65
Without back cushion extender		
Extension back support upholstery		28657
Without Extension back support upholstery		
Safety edge for 16" wheel		28415
Without safety edge for 16" wheel		

Reference

User

Phone

Date

Comments