






Next Comfort

Prescription form



Name: _____
 Reference: _____
 Standard version is marked with * Date: _____

Short frame Seat width (cm) Seat depth 35-41 cm (Measured from the back support tube)	35 13202001 ☐ 37,5 13202003 ☐ 40 13202005 ☐ 42,5 13202007 ☐	45 13202009 ☐ 47,5 13202011 ☐ 50 13202013 ☐									
Long frame Seat width (cm) Seat depth 41-47 cm (Measured from the back support tube)	35 13202002 ☐ 37,5 13202004 ☐ 40 13202006 ☐ 42,5 13202008 ☐	45 13202010 ☐ 47,5 13202012 ☐ 50 13202014 ☐									
Frame colour	Blueberry blue -48 ☒ Grey blast -53 ☐										
Back support angle	Adjustable -5°-+15° Standard +2°=92° between back/seat Indicate angle ____°										
Back support	<table border="1"> <tr> <td> Back support cushion ☒ 59 cm adjustable from 59-69 cm  </td> <td> Back upholstery adjustable ☐ (incl. brace) Stepless adjustable 30-45 cm Standard 35 cm Indicate height: ____cm </td> <td> Back support high ☐ (incl. brace) Stepless adjustable 38-55 cm Standard 45 cm Indicate height: ____cm </td> </tr> <tr> <td> Back support cushion ☐ 48 cm adjustable from 48 - 62 cm  </td> <td> Back support cover ☐ Standard dark grey plush  </td> <td></td> </tr> <tr> <td> Black Dartex ☒ Dark grey plush ☐ </td> <td> Comfort ☐  </td> <td></td> </tr> </table>		Back support cushion ☒ 59 cm adjustable from 59-69 cm 	Back upholstery adjustable ☐ (incl. brace) Stepless adjustable 30-45 cm Standard 35 cm Indicate height: ____cm	Back support high ☐ (incl. brace) Stepless adjustable 38-55 cm Standard 45 cm Indicate height: ____cm	Back support cushion ☐ 48 cm adjustable from 48 - 62 cm 	Back support cover ☐ Standard dark grey plush 		Black Dartex ☒ Dark grey plush ☐	Comfort ☐ 	
Back support cushion ☒ 59 cm adjustable from 59-69 cm 	Back upholstery adjustable ☐ (incl. brace) Stepless adjustable 30-45 cm Standard 35 cm Indicate height: ____cm	Back support high ☐ (incl. brace) Stepless adjustable 38-55 cm Standard 45 cm Indicate height: ____cm									
Back support cushion ☐ 48 cm adjustable from 48 - 62 cm 	Back support cover ☐ Standard dark grey plush 										
Black Dartex ☒ Dark grey plush ☐	Comfort ☐ 										
Seat	Solid seat ☒ 										
Seat cushion	<table border="1"> <tr> <td>  </td> <td> Flat PE ☐ Flat memory foam ☐ </td> </tr> <tr> <td>  </td> <td> Formed PE ☒ Formed memory foam ☐ </td> </tr> </table>			Flat PE ☐ Flat memory foam ☐		Formed PE ☒ Formed memory foam ☐					
	Flat PE ☐ Flat memory foam ☐										
	Formed PE ☒ Formed memory foam ☐										
Seat cover	Black Dartex ☒ Dark grey plush ☐										
	Seat upholstery adjustable ☐ 										











74844H









08-04-2014





Etac, Box 203, SE-334 24 Anderstorp









Tel. +46/371-58 73 00, Fax. +46/371-58 73 90



www.etac.com

Push handles height adjustable	Push handles ☒  Bow handle ☐  Bow handle angle adj. ☐ 																				
Brakes	Standard Ri ☒ Le ☒ Single hand brake Ri ☐ Le ☐																				
Legrests alt 1 Detachable/swing away With locking device Calf support Foot support	Standard Ri ☒ Le ☒  Narrow Ri ☐ Le ☐  Short /Narrow Ri ☐ Le ☐  Elevating Ri ☐ Le ☐  Elevating with foldable calf support Ri ☐ Le ☐  Support console excl calf support Ri ☐ Le ☐  Calf support Ri ☐ Le ☐  Foldable. Adjustable in depth, angle and height.  Tube-standard length ☒ Tube-short length (-8cm) Only for short/narrow legrest ☐ Tube-long length (+10 cm) ☐ Foldable. Adjustable in angle.  Foot support one piece ☐																				
Legrests alt 2 Detachable/swing away self-extending incl attachment pipe Calf support Foot support	Standard 41-50 cm Ri ☐ Le ☐  Short 34-37 cm Ri ☐ Le ☐ Long 50-70 cm Ri ☐ Le ☐ Standard formad Ri ☐ Le ☐ With edge Ri ☐ Le ☐  Shark fin shaped Ri ☐ Le ☐  Wave shaped Ri ☐ Le ☐  Black plastic assembeled with heel strap as support Ri ☐ Le ☐ 																				
Legrests Others	Amputee support adjustable in height Ri ☐ Le ☐  Amputee weight compensating compl. Ri ☐ Le ☐  Foot support one-piece children compl. 35/37.5 adjustable in height, angle and depth, detachable ☐  Plaster legrest Ri ☐ Le ☐ 																				
Rear wheels with quick release Handrims	Pneumatic high pressure 1" Pneumatic 1 3/8" Solid Solid inner tube 20" --- ☐ ☐ ☐ ☐ 22" --- ☐ ☐ ☐ ☐ 24" ☒ --- ☐ ☐ ☐  Aluminium (not with solid inner tube) ☐ Stainless steel ☐ Plastic coated --- Cellular rubber ☐ Assistant brake ☐ 																				
Castors	5" (125x30 mm) solid ☐ 6" (150x30 mm) solid ☒ 7" (175x45 mm) semipneumatic (only with 24" rear wheels) ☐																				
Camber angle	0° 2° ☐ ☒																				
Front seat height (cm) with seat tilted backwards +3°	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>Rear wheel</th><th>Pos</th><th>Size 20"</th><th>22"</th><th>24"</th></tr> </thead> <tbody> <tr> <td>D</td><td></td><td>☐ 40 cm</td><td>☐ 42 cm</td><td>☐ 44 cm</td></tr> <tr> <td>C</td><td></td><td>☐ 42,5 cm</td><td>☐ 44,5 cm</td><td>☒ 46,5 cm</td></tr> <tr> <td>B</td><td></td><td>☐ 45 cm</td><td>☐ 47 cm</td><td>☐ 49 cm</td></tr> </tbody> </table>	Rear wheel	Pos	Size 20"	22"	24"	D		☐ 40 cm	☐ 42 cm	☐ 44 cm	C		☐ 42,5 cm	☐ 44,5 cm	☒ 46,5 cm	B		☐ 45 cm	☐ 47 cm	☐ 49 cm
Rear wheel	Pos	Size 20"	22"	24"																	
D		☐ 40 cm	☐ 42 cm	☐ 44 cm																	
C		☐ 42,5 cm	☐ 44,5 cm	☒ 46,5 cm																	
B		☐ 45 cm	☐ 47 cm	☐ 49 cm																	

	ARM SUPPORT LONG	lockable, height-adjustable 38 cm long, solid top (included in standard features)	right left	81203 81204	☒ ☒
	ARM SUPPORT SHORT	lockable, height-adjustable 25 cm long, solid top	right left	81201 81202	☐ ☐
	ARM SUPPORT LOW	lockable, height-adjustable 25 cm long, solid top lockable, height-adjustable 38 cm long, solid top	right left right left	26289-1 26289-2 26290-1 26290-2	☐ ☐ ☐ ☐
	ARM SUPPORT COVER detachable, 25 or 38 cm long dark grey plush	padded short long gel short long	quantity____ quantity____	62282-60 62283-60	☐ ☐
			quantity____ quantity____	26240-60 26241-60	☐ ☐
	detachable, 25 or 38 cm long black Dartex	padded short long gel short long	quantity____ quantity____	62282-63 62283-63	☐ ☐
			quantity____ quantity____	26240-63 26241-63	☐ ☐
	COVER ARM SUPPORT wide and soft, to be mounted on existing arm support		size 8 cm x 40 cm	26873	☐
	ARM SUPPORT HEMI detachable cover Only for long arm support!		right left	25770-60-1 25770-60-2	☐ ☐
	TRAY to be mounted on long arm support seat width cm	35 24417 37,5 25187 40 24418	☐ ☐ ☐ ☐	42,5 25188 45 24419 47,5 25731 50 24601	☐ ☐ ☐ ☐
	TABLE HEMI incl. cover with velcro lock seat width cm	arm support short 35 26556 37,5 26557 40 26558 42,5 26559 45 26560 47,5 26561 50 26562	☐ ☐ ☐ ☐ ☐ ☐ ☐	arm support long 35 26563 37,5 26564 40 26565 42,5 26566 45 26567 47,5 26568 50 26569	☐ ☐ ☐ ☐ ☐ ☐ ☐
	ANTI GLIDE TABLE HEMI	35-40 cm 42,5-50 cm	quantity____ quantity____	81633 81634	☐ ☐
	CUSHION grey plush, 56 cm long anatomically contoured, adapted for sling seat to be cut to individual seat depth	seat width cm 35 62230-60 37,5 62231-60 40 62232-60 42,5 62233-60 45 62234-60 47,5 62374-60 50 62266-60	☐ ☐ ☐ ☐ ☐ ☐ ☐		
	COMFORT WEDGE increases the anatomically shape (To be placed in the upholstery pocket)	seat width cm 35-37,5 80714 40-45 80715 47,5-50 80779	☐ ☐ ☐		
	ATTACHMENT KIT SEAT FLAT ECS assembly kit that will elevate the seat to the same height as the seat pipes			27330	☐

	CALF STRAP detachable, adjustable in length	seat width cm	35-37,5	25356	☐
			40-42,5	25357	☐
	PADDING for calf strap		45-47,5	25358	☐
			50	25359	☐
			20 cm	26845	☐
	EXTENDED FOOT SUPPORT to be mounted on foot supports	seat width cm	30 cm	26846	☐
			40 cm	26847	☐
	HEEL STRAP (Included in standard features)	adjustable in length	40-42,5	quantity _____	50798 ☐
			45-47,5	quantity _____	50799 ☐
			50	quantity _____	50800 ☐
	HEEL STRAP (Included in standard features)	adjustable in length	right	24458	☒
			left	24458	☒
	COVER black Dartex Only fits calf support for elevating legrest with foldable calf support		right	400220	☐
			left	400221	☐
	CALF SUPPORT BRACKET side swing away for legrest alt 2		long	quantity _____	400072 ☐
			short	quantity _____	400097 ☐
	LEGREST COVER black Dartex for legrest alt 2		right	400216	☐
			left	400216	☐
	KNEE SUPPORT for legrest alt 2		right	400216	☐
			left	400216	☐
	BACK SUPPORT WEDGE darkgrey plush	short 15 cm long 30 cm	quantity _____	62346-60	☐
			quantity _____	62347-60	☐
	POSITIONING BELT two-part with snap-lock, attachment points on wheelchair frame			26273	☐
	HEAD SUPPORT MODEL 1 adjustable in height, depth and angle (Included in standard features)		black Dartex	24462-63	☒
			dark grey plush	24462-60	☐
	ATTACHMENT BRACKET for head support (Included in standard features) In combination with adjustable back upholstery, choose bow handle and attachment bracket			406256	☒
				25795	☐
	HEAD SUPPORT adjustable tube 45°	Head support model 2 With side support 320(240inv)x170 mm	Height adjustable tube	45°	☐
				406255	☐
	HEAD SUPPORT adjustable tube 45°	Head support model 3 Curved 250x140 mm	Height adjustable tube	90°	☐
				400260	☐
	HEAD SUPPORT adjustable tube 45°	Head support model 2 with sideways position	Height adjustable tube	curved	☐
				400262	☐

	BACKREST BRACE seat width cm											
			35	25610	☐	47,5	25739	☐				
	with snap-lock		37,5	25611	☐	50	25615	☐				
			40	25612	☐	52,5	26030	☐				
			42,5	25613	☐	55	26031	☐				
			45	25614	☐							
	TRUNK SUPPORT adjustable in depth and angle, black Dartex, can be fold up, incl. attachment											
	Small curved	ri	405802	☐	Ear-shaped	ri	405800	☐	Large curved	ri	405804	☐
	140mm x 100mm	le	405803	☐	100mm x 200mm	le	405801	☐	150mm x 120mm	le	405805	☐
												
	Upholstery trunk support, black Dartex											
							quantity ____	410319	☐			
	ANTI-TIPS flip-up (Included in standard features) for low seat height 20" wheels											
							quantity 2	25969	☒			
	ANTI-TIPS telescoping For 20" to 24" rear wheels											
							quantity ____	26895	☐			
	CANE HOLDER universal											
							right	27327	☐			
							left	27327	☐			
	SPOKEGUARD with grey logo		20"	25502-28	☐	with yellow logo	20"	25502-29	☐			
			22"	25503-28	☐		22"	25503-29	☐			
			24"	25504-28	☐		24"	25504-29	☐			
	Other spokeguards on see www.etac.se											
	TETRA QUICK RELEASE ADAPTER											
							quantity ____st	24051-01	☐			
	TRANSPORT ATTACHMENT complete set											
								26236	☐			
	LOCKING of folded chassis											
								27151	☐			
	ANGLE LIMITATION +20° - 0°											
								27791	☐			
	TOOLKIT											
								26204	☐			

	BAGS Backpack Mobile holder	92009	☐		RAINCAP Helly Hansen Pro For alternatives see www.etac.com Art.no: _____
	COCOON Helly Hansen Pro For alternatives see www.etac.com Art.no: _____				LEG WARMER Helly Hansen Pro For alternatives see www.etac.com Art.no: _____
	BAGS Case Logic For alternatives see www.etac.com				Art.no: _____