

Cross 5 XL




















Prescription form

75188L 2023-05-08

Reference:	Name:
Tel:	Date:

Standard version is marked with
 The measurements given can vary $\pm 2\%$

Short frame Seat width (cm) <small>Seat depth 36-42 (-49*)cm Measured from back post (*Triangular back support)</small>	47,5	13123111	<input type="checkbox"/>	55	13123117	<input type="checkbox"/>		
	50	13123113	<input type="checkbox"/>	57,5	13123119	<input type="checkbox"/>		
	52,5	13123115	<input type="checkbox"/>	60	13123121	<input type="checkbox"/>		
Long frame Seat width (cm) <small>Seat depth 42-48 (-55*)cm Measured from back post (*Triangular back support)</small>	47,5	13123112	<input type="checkbox"/>	55	13123118	<input type="checkbox"/>		
	50	13123114	<input type="checkbox"/>	57,5	13123120	<input type="checkbox"/>		
	52,5	13123116	<input type="checkbox"/>	60	13123122	<input type="checkbox"/>		
Frame colour	Black	-01	<input type="checkbox"/>					
Back support Cross back 3A	Type	Lumbar angle standard $+7^\circ$						
	Triangular	<input type="checkbox"/>	-8 - +18°					°
	Back Support Height Standard	32-45 cm (std 40 cm)	Indicate height					cm
	Back Support Height High	38-51 cm (std 46 cm)	Indicate height					cm
Back Support Angle	-5° - +20° (std +2°)	Indicate angle			°			
Cover	Standard <input type="checkbox"/>	Comfort <input type="checkbox"/>	Tidy <input type="checkbox"/>					
	Back support Cross back 2	Continuously adjustable 30-45 cm. Standard 35 cm	<input type="checkbox"/>	Indicate height: _____cm	Cover standard	<input type="checkbox"/>		
	Back support Cross back 2, widened	Continuously adjustable 30-45 cm. Standard 35 cm Extends the width 5 cm	<input type="checkbox"/>	Indicate height: _____cm	Cover comfort	<input type="checkbox"/>		
Back support tubes for hard shell back support, ECS		<input type="checkbox"/>	Continuously adjustable 30-45 cm Standard 35 cm	(for seat and back support system ECS, see separate prescription form at www.etac.com)				
Back support tubes for hard shell widened back support, ECS (max chair width 55 cm)		<input type="checkbox"/>	Continuously adjustable 37-42 cm Standard 37 cm Extends the back support width by 5 cm					
Back support angle	-5° - +15°	Standard $+2^\circ = 92^\circ$ between back and seat	Indicate angle _____°					
Push handles	Push handles <input checked="" type="checkbox"/> Detachable, height adjustable	Bow handle <input type="checkbox"/> Detachable, height adjustable						
	Seat upholstery <input checked="" type="checkbox"/> Adjustable seat depth 6 cm	Solid seat, detachable <input type="checkbox"/>						


Brakes		With angled brake handle Ri <input checked="" type="checkbox"/> Le <input checked="" type="checkbox"/> Red <input checked="" type="checkbox"/> Black <input type="checkbox"/>		With extension arm, foldable Ri <input type="checkbox"/> Le <input type="checkbox"/>		Single-hand brake Ri <input type="checkbox"/> Le <input type="checkbox"/> With angled brake handle <input type="checkbox"/>
		Standard (not for bow handle) <input type="checkbox"/>		With larger grip <input type="checkbox"/>		
Leg supports Black Detachable/swing-away		Standard Ri <input checked="" type="checkbox"/> Le <input checked="" type="checkbox"/>		Narrow Ri <input type="checkbox"/> Le <input type="checkbox"/>		
	with locking device	Yes <input type="checkbox"/> No <input checked="" type="checkbox"/>	Yes <input type="checkbox"/> No <input type="checkbox"/>			
Elevating leg supports		Mod 1 without calf support Ri <input type="checkbox"/> Le <input type="checkbox"/>		Mod 3 without calf support Ri <input type="checkbox"/> Le <input type="checkbox"/>		Mod 2 with calf support foldable, adjustable in angle, depth and sideways Ri <input type="checkbox"/> Le <input type="checkbox"/>
		Mod 3 with calf support foldable sideways in height two positions. From 40 cm seat width (with two leg supports) Ri <input type="checkbox"/> Le <input type="checkbox"/>				
Options for Mod 1 and Mod 3 without calf support		Attachment for calf support Ri <input type="checkbox"/> Le <input type="checkbox"/>		Calf support Ri <input type="checkbox"/> Le <input type="checkbox"/>		
		Attachment for calf-support, foldable Ri <input type="checkbox"/> Le <input type="checkbox"/>				
Foot supports		Flip-up, adjustable	Tube, standard length <input checked="" type="checkbox"/> Tube, short length(-8 cm) only for short/narrow leg support <input type="checkbox"/> Tube, long length (+10 cm) <input type="checkbox"/>	Foot support one piece <input type="checkbox"/> Flip-up, adjustable in angle <input type="checkbox"/>		
		Amputee support adjustable in height, angle, depth and sideways Ri <input type="checkbox"/> Le <input type="checkbox"/>		Amputation weight, can be combined with Amputee support Ri <input type="checkbox"/> Le <input type="checkbox"/>		Plaster leg support with fixed 0° knee angle. Support plates are adjustable Ri <input type="checkbox"/> Le <input type="checkbox"/>
Material	Plush <input type="checkbox"/> Grey hygiene <input type="checkbox"/>					
Rear wheels with quick release	20" Pneumatic high pressure 1" --- <input type="checkbox"/> Pneumatic 1 3/8" (not with drum brakes) <input type="checkbox"/> Solid <input type="checkbox"/>	22" Pneumatic high pressure 1" <input type="checkbox"/> Pneumatic 1 3/8" (not with drum brakes) --- <input type="checkbox"/> Solid <input type="checkbox"/>	24" Pneumatic high pressure 1" <input checked="" type="checkbox"/> Pneumatic 1 3/8" (not with drum brakes) --- <input type="checkbox"/> Solid <input type="checkbox"/>	26" Pneumatic high pressure 1" <input type="checkbox"/> Pneumatic 1 3/8" (not with drum brakes) --- <input type="checkbox"/> Solid --- <input type="checkbox"/>		
	Aluminium <input type="checkbox"/> Stainless steel <input type="checkbox"/> Cellular rubber <input type="checkbox"/> Vinyl coated (not with solid tyre) <input type="checkbox"/>	Aluminium <input type="checkbox"/> Stainless steel <input type="checkbox"/> Cellular rubber <input type="checkbox"/> Vinyl coated <input type="checkbox"/>	Aluminium <input checked="" type="checkbox"/> Stainless steel <input type="checkbox"/> Cellular rubber <input type="checkbox"/> Vinyl coated <input type="checkbox"/>	Aluminium --- <input type="checkbox"/> Stainless steel <input type="checkbox"/> Cellular rubber --- <input type="checkbox"/> Vinyl coated --- <input type="checkbox"/>		
		20" <input type="checkbox"/> 22" <input type="checkbox"/> 24" <input type="checkbox"/>	Handrims on Right <input type="checkbox"/> Left <input type="checkbox"/>			

Castors	4" (100x30 mm) solid <input type="checkbox"/>	7" (175x45 mm) semipneumatic <input type="checkbox"/>
	5" (125x30 mm) solid (max 135 kg user weight) <input type="checkbox"/>	8" (200x35 mm) semipneumatic <input type="checkbox"/>
	6" (150x30 mm) solid <input type="checkbox"/>	8" (200x50 mm) semipneumatic <input type="checkbox"/>
	6,5" (165x30 mm) solid <input checked="" type="checkbox"/>	8" (200x30 mm) pneumatic <input type="checkbox"/>

Front fork attachment	Standard <input checked="" type="checkbox"/>	Widened (moves castor 3 cm to the side for foot propulsion) <input type="checkbox"/>
-----------------------	--	--



Camber angle	0° <input type="checkbox"/>	2° <input checked="" type="checkbox"/>
--------------	-----------------------------	--

Rear seat height (cm)	Rear 20"		22"		24"		26"	
	Pos A	46	<input type="checkbox"/>	48,5	<input type="checkbox"/>	51	<input type="checkbox"/>	53
B	44	<input type="checkbox"/>	46,5	<input type="checkbox"/>	49	<input type="checkbox"/>	51	<input type="checkbox"/>
C	42	<input type="checkbox"/>	44,5	<input type="checkbox"/>	47	<input checked="" type="checkbox"/>	49	<input type="checkbox"/>
D	40	<input type="checkbox"/>	42,5	<input type="checkbox"/>	45	<input type="checkbox"/>	47	<input type="checkbox"/>
E	38	<input type="checkbox"/>	40,5	<input type="checkbox"/>	43	<input type="checkbox"/>	45	<input type="checkbox"/>
F	36	<input type="checkbox"/>	38,5	<input type="checkbox"/>	41	<input type="checkbox"/>	43	<input type="checkbox"/>
G	34	<input type="checkbox"/>	36,5	<input type="checkbox"/>	39	<input type="checkbox"/>	-	---

Front seat height (cm)		Castor Extra short front fork Short front fork Medium front fork Long front fork	4" * <input type="checkbox"/>	34-39,5 cm	Indicate height	cm	* Max 100 kg user weight ** Max 135 kg user weight
			5" ** <input type="checkbox"/>	36-44,5 cm		cm	
			6" <input type="checkbox"/>	40-45,5 cm		cm	
			6,5" <input type="checkbox"/>	41-46,5 cm		cm	
			5" ** <input type="checkbox"/>	43-49 cm		cm	
			6" <input type="checkbox"/>	45-50 cm		cm	
			6,5" <input checked="" type="checkbox"/>	46-51 cm		cm	
			8" <input type="checkbox"/>	48-52 cm		cm	
			8"x1" (pneumatic) <input type="checkbox"/>	46,5-50,5 cm		cm	
			6,5" <input type="checkbox"/>	50,5-55,5 cm		cm	
	8" <input type="checkbox"/>	52,5-56,5 cm		cm			
Fork thickness		3 mm <input type="checkbox"/>	4 mm <input checked="" type="checkbox"/>				

Balance position


For alternative seat height/seat angles, see Manual Cross































0° Cambervinkel















Indicate balance position 1-6 (standard pos 4)

Accessories

Arm support		With detachable arm support plate, length 25 cm, height adjustable	right	27230-1	<input type="checkbox"/>
			left	27230-2	<input type="checkbox"/>
		Adjustable with fingerscrew	right	27349-1	<input type="checkbox"/>
			left	27349-2	<input type="checkbox"/>

Arm support		Massive moulded plate, length 25 cm, height adjustable	right	25172-1	<input type="checkbox"/>			
			left	25172-2	<input type="checkbox"/>			
		Adjustable with fingerscrew	left	26379-1	<input type="checkbox"/>			
			left	26379-2	<input type="checkbox"/>			
Arm support		With detachable arm support plate, length 35 cm, height adjustable	right	27231-1	<input type="checkbox"/>			
			left	27231-2	<input type="checkbox"/>			
		Adjustable with fingerscrew	right	27350-1	<input type="checkbox"/>			
			left	27350-2	<input type="checkbox"/>			
Arm support		Massive moulded plate, length 38 cm, height adjustable	right	25173-1	<input type="checkbox"/>			
			left	25173-2	<input type="checkbox"/>			
		Adjustable with fingerscrew	right	26380-1	<input type="checkbox"/>			
			left	26380-2	<input type="checkbox"/>			
Cover arm support <small>Fits arm support with massive moulded plate</small>		Dark grey plush, detachable, length 25 or 38 cm	Padded	25 cm	Quantity	pcs	62282-60	<input type="checkbox"/>
			Padded	38 cm	Quantity	pcs	62283-60	<input type="checkbox"/>
			Gel	25 cm	Quantity	pcs	26240-60	<input type="checkbox"/>
			Gel	38 cm	Quantity	pcs	26241-60	<input type="checkbox"/>
Arm support Hemi		Detachable cover, dark grey plush, sideways angle adjustable Fits with long arm supports 27231-1/2 and 27350-1/2	right	27313-60-1	<input type="checkbox"/>			
			left	27313-60-2	<input type="checkbox"/>			
		Detachable cover, dark grey plush Fits with long arm supports 25173-1/2 and 26380-1/2	right	25770-60-1	<input type="checkbox"/>			
			left	25770-60-2	<input type="checkbox"/>			
Arm support pad wide and soft <small>Fits arm support with massive moulded plate</small>		8x40 cm, to be mounted on arm supports		pcs	26873	<input type="checkbox"/>		
Side support		Black, adjustable in height	right	25189-1	<input type="checkbox"/>			
			left	25189-2	<input type="checkbox"/>			
Thigh support		To be mounted on short or long arm supports, front attachment is included. Limitations: 24" wheel, not in B1-2, C1-2, D1-2, E1-2, F1-2, G1-2 22" wheel, not in C1-2, D1-2, E1-2, F1-2, G1 20" wheel, not in E1-2, F1-2, G1	Short	right	27391-1	<input type="checkbox"/>		
			Short	left	27391-2	<input type="checkbox"/>		
			Long	right	27392-1	<input type="checkbox"/>		
			Long	left	27392-2	<input type="checkbox"/>		
Soft pad		For thigh support, velcro attachment	right	82488-1	<input type="checkbox"/>			
			left	82488-2	<input type="checkbox"/>			
Mudguard		Carbon fibre, adjustable in height and depth to fit rear wheel position	right	27397-1	<input type="checkbox"/>			
			left	27397-2	<input type="checkbox"/>			
Width expansion kit		Extends rear wheels, brakes and arm supports 5 mm on each side (not together with drum brakes)			27533	<input type="checkbox"/>		
Tray		To be mounted on long arm supports, transparent	Seat width cm	45	24419	<input type="checkbox"/>		
			47,5	25731	<input type="checkbox"/>			
			50	24601	<input type="checkbox"/>			
			52,5	26020	<input type="checkbox"/>			
			55	26021	<input type="checkbox"/>			
			57,5	26022	<input type="checkbox"/>			
Tray Hemi/Communication		Tip up right/left	right	27352-1	<input type="checkbox"/>			
			left	27352-2	<input type="checkbox"/>			
Seat extension Cross back 3A		Extends the seat depth up to 7 cm at rear, in combination with triangular back support	Seat width 47.5-60 cm		62977	<input type="checkbox"/>		

Seat cushion		Dark grey plush, 56 cm long, anatomically shaped, cut to individual seat depth	Seat width cm		45	62234-60	<input type="checkbox"/>		
			47,5	62374-60	<input type="checkbox"/>				
			50	62266-60	<input type="checkbox"/>				
			52,5	62420-60	<input type="checkbox"/>				
			55	62421-60	<input type="checkbox"/>				
			57,5-60	62422-60	<input type="checkbox"/>				
Comfort wedge		Increases the contour. To be placed in the pocket of the inner front upholstery	Seat width cm		45	80715	<input type="checkbox"/>		
			47,5-60	80779	<input type="checkbox"/>				
Seat cushion Tidy		Wipeable	Seat width cm		45	63102-63	<input type="checkbox"/>		
			47,5	63103-63	<input type="checkbox"/>				
			50	63104-63	<input type="checkbox"/>				
			52,5	63105-63	<input type="checkbox"/>				
			55	63106-63	<input type="checkbox"/>				
57,5-60	63107-63	<input type="checkbox"/>							
ECS seat- and back-support		The seat- and back support system offers products and covers in different kinds of upholstery and fabrics (see separate prescription form at www.etac.com)							
Calf support cover		Dark grey plush For elevating leg support model 3	Quantity	pcs	62762-60	<input type="checkbox"/>			
		Black Dartex For elevating leg support model 3	Quantity	pcs	62762-63	<input type="checkbox"/>			
Calf strap		Detachable, adjustable in length	Seat width cm		45-47,5	25358	<input type="checkbox"/>		
			50-52,5	25359	<input type="checkbox"/>				
			55-60	26062	<input type="checkbox"/>				
Padding for calf strap					20 cm	26845	<input type="checkbox"/>		
					30 cm	26846	<input type="checkbox"/>		
					40 cm	26847	<input type="checkbox"/>		
Foot support extension		Mounted on foot support	Seat width cm		45-47,5	Quantity	pcs	50799	<input type="checkbox"/>
					50-57,5	Quantity	pcs	50800	<input type="checkbox"/>
Heel strap		Adjustable in length			right	24458	<input type="checkbox"/>		
					left	24458	<input type="checkbox"/>		
Positioning belt		Two-part with snap-lock. 2 alt. attachment points on wheelchair frame				25668	<input type="checkbox"/>		
Head support		Curved, padded. Adjustable in height, depth, side-ways and angle.			grey hygiene	28223-65	<input type="checkbox"/>		
					dark grey plush	28223-60	<input type="checkbox"/>		
Head support attachment		For Cross back 3A, mod 4	45	28129	<input type="checkbox"/>	55	28133	<input type="checkbox"/>	
			47,5	28130	<input type="checkbox"/>	57,5	28134	<input type="checkbox"/>	
			50	28131	<input type="checkbox"/>	60	28135	<input type="checkbox"/>	
			52,5	28132	<input type="checkbox"/>				
Head support bracket		For bow handle				25795	<input type="checkbox"/>		
Back support wedge		Dark grey plush	Short 15 cm	Quantity	pcs	63012-60	<input type="checkbox"/>		
			Long 30 cm	Quantity	pcs	63013-60	<input type="checkbox"/>		

Side cushion		Provides gentle lateral support. Used when the trunk support is too stabilising. (Combined with arm support and, if necessary, side stop)	Right	63171-60-1	<input type="checkbox"/>			
			Left	63171-60-2	<input type="checkbox"/>			
Side stop		Straight or angled. Combined with wedge or side support cushion. Creates a more distinct stop and holds the wedge or side support cushion in place.	Straight	pcs	28136	<input type="checkbox"/>		
			Widening 15 mm	pcs	28137	<input type="checkbox"/>		
			Widening 30 mm	pcs	28173	<input type="checkbox"/>		
Trunk support		Adjustable, swing-away and detachable	0°	Right	28141-1	<input type="checkbox"/>		
			0°	Left	28141-2	<input type="checkbox"/>		
			-15°	Right	28180-1	<input type="checkbox"/>		
			-15°	Left	28180-2	<input type="checkbox"/>		
			-45°	Right	28181-1	<input type="checkbox"/>		
			-45°	Left	28181-2	<input type="checkbox"/>		
Trunk support		Settable and detachable	0°	pcs	28142	<input type="checkbox"/>		
			+15°	pcs	28143	<input type="checkbox"/>		
			-15°	pcs	28144	<input type="checkbox"/>		
			-45°	pcs	28145	<input type="checkbox"/>		
Trunk support		Dartex, fixed, adjustable in height, depth and width	Quantity	pcs	27395	<input type="checkbox"/>		
		PUR, fixed, adjustable in height, depth and width	Quantity	pcs	27917	<input type="checkbox"/>		
Trunk support		PUR, swing away, adjustable in height and depth	Medium	right	27915-1	<input type="checkbox"/>		
			Medium	left	27915-2	<input type="checkbox"/>		
			Large	right	27916-1	<input type="checkbox"/>		
			Large	left	27916-2	<input type="checkbox"/>		
Anti-tip		Telescoping, flip-up	Quantity	pcs	27526	<input type="checkbox"/>		
Cane holder				right	27327	<input type="checkbox"/>		
				left	27327	<input type="checkbox"/>		
Spoke guard		Other spoke guards see www.etac.com	Grey with logo		Yellow with logo			
			20"	25502-28	<input type="checkbox"/>	20"	25502-29	<input type="checkbox"/>
			22"	25503-28	<input type="checkbox"/>	22"	25503-29	<input type="checkbox"/>
			24"	25504-28	<input type="checkbox"/>	24"	25504-29	<input type="checkbox"/>
Tetra quick release adapter			Quantity	st	24051-01	<input type="checkbox"/>		
Locking strap		Of folded chassis			27151	<input type="checkbox"/>		
Toolkit					25174	<input type="checkbox"/>		
Gloves		For sizes see www.etac.com	Bags		Cellular phone holder	92009	<input type="checkbox"/>	
		Item no				92010	<input type="checkbox"/>	

Etac
Box 203, 334 24 Anderstorp
Sweden
Tel +46 371 58 73 00
Fax +46 371 58 73 90
www.etac.com

