

Etac Prio prescription form

75326P 2021-06-23

Applies to Etac Prio with serial no beginning with 38C

③ Seat tilt and back recline function	Manual	Electric
Manual or electric		

④ Back support system		Cushion/cover		Cover material	
		Cush-ion	Cover	Grey plush	Grey hygiene
Height settable without dynamic pelvic support. Height 55-70 cm					
Height settable with dynamic pelvic support. Height 55-70 cm					
3D without dynamic pelvic support. Height settable 45-60 cm					
3D with dynamic pelvic support. Height settable 45-60 cm					

⑤ Back support locking	Standard	QR-pin
Standard or quick-release		

⑥ Head support		Cover	
		Grey plush	Grey hygiene
Multi adjustable			
Multi adjustable with side support			
Anatomic medium			
Anatomic large			
Without head support			

⑦ Comfort tilt lock	
Without Comfort tilt lock	
With Comfort tilt lock preset at 0°	
With Comfort tilt lock preset at 4°	
With Comfort tilt lock preset at 8°	
With Comfort tilt lock preset at 12°	

⑧ Arm support			
Model	R	L	
Standard			
Solid plate length 38 cm			

① Wheelchair

Seat width	Item. no.	
40	13420305	
45	13420309	
50	13420313	
55	13420317	

② Seat depth

	Seat cushion							Seat plate			
	Short	Long	Basic	Com-fort	Grey plush	Grey hygiene	No seat cushion	Wood	Plastic	Wood w. Velcro	Plastic w. Velcro
Pos 1= 39-46 cm short cushion											
Pos 2= 41,5-48,5 cm short cushion											
Pos 3= 44-51 cm short or long cushion											

⑨ Leg support

Model	R	L
Elevating		
Fixed angle		
Amputee plush		
Amputee Hygiene		
Amputee weight		
Amputee+weight plush		
Amputee+weight Hygiene		
Leg support Cross		


⑩ Foot supports

Model	R	L
Depth 15 cm, long attachment tube (38.5-55 cm from seat plate)		
Depth 15 cm, medium attachment tube (31-47.5 cm from seat plate)		
Depth 13 cm, long attachment tube (38.5-55 cm from seat plate)		
Depth 13 cm, medium attachment tube (31-47.5 cm from seat plate)		
Solid footplate angle adjustable		
Solid footplate angle adjustable Add-on		


Seat height and balance position

⑪

Seat height, rear wheels, castor wheels and position



Front fork short



Front fork medium

S1

S2

S3

M1

M2

M3

6"

6.5"

8"

6"

6.5"

8"

6"

6.5"

8"

6"

6.5"

8"

6"

6.5"

8"

6"

6.5"

8"

6"

6.5"

8"

6"

6.5"

8"



Seat height

20" **

38,0

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

22"

40,5

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

24"

42,5

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

20" **

40,5

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

22"

43,0

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

24"

45,0

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

-

20" **

43,0

-

-

Accessories

Accessories			
	R	L	Item. no.
Trunk support detachable without tools			28020
Trunk support Swing-away			28019-1 28019-2
Without trunk support			
Cover for trunk support			28647
Without cover			
Adapter kit for trunk support cushion R82			28240
Without adapter kit			
Side support cushion			63171-60-1 63171-60-2
Without side support cushion			
Arm support padding std			63193-60-1 63193-60-2
Arm support padding solid	Foam		62283-60
	Gel		26241-60
Hemi cushion			25770-60-1 25770-60-2
Cushion wide and soft			28640
Without arm support padding			
Widening kit +2.5 cm			27892
Without widening kit			
Side support padding			63271-1 63271-2
Without side support padding			
Tray			
Without tray			

Accessories			
			Item. no.
Pommel			28022
Without Pommel			
Calf support			
Without calf support			
Calf support padding			
20 cm			26845
30 cm			26846
40 cm			26847
Without calf support padding			
Calf support wide	70 cm		63186-60
	90 cm		63187-60
Without calf support wide			
Positioning belt	100 cm	pcs	27904
	128 cm	pcs	28544
Attachment positioning belt			28641
Without attachment positioning belt			
Chest strap bracket			28186
Without chest strap bracket			
Comfort wedges			
H:2, W:11, L:18.5 cm	pair		28139
H:2.5, W:17.5, L:20 cm	pair		28140
H:2, W:9, L:13 cm	pair		28030

Accessories		
		Item. no.
Gas piston clips (10 pcs/set)		28163
Without gas piston clips		
Lock seat and back angle controls	pcs	28047
Back angle control		28573
Without back angle control		
Oxygen holder/Transport basket		28316
Without oxygen holder		
Drip stand		28341
Without drip stand		

Accessories

Accessories			
			Item. no.
Velcro kit seat cushion for separate mounting			28183
Without Velcro kit			
Spoke guards grey			
Spoke guards yellow			
Without spoke guards			
Cane holder	st		27327
	R	L	
Extended foot support (21.5 cm)			
Without extended foot support			
Extended footplate Add-on	S		28774
	M		28775-1 28775-2
	L		28776-1 28776-2
Without extended foot support Add-on			
Side protection extended foot support Add-on			28772
Without side protection Add-on			
Foot support padding			
For foot supports model 2			
For solid foot support, angle adjustable			
Without padding			

Heel strap for foot support depth 13 cm			24458
Without heel strap			
Footbox			
Without footbox			
Memory function head support			28564
Without memory function head support			
Extension back support upholstery			28657
Without Extension back support upholstery			
Pressure distributor	pcs		28637
Safety edge for 16" wheel			28415
Without safety edge for 16" wheel			

Reference

User

Phone

Date

Comments