

Cross 5XL

Therapist		User
Phone no		Order date
Purchase Order no		Customer Ref
Base price		€XXX
		Standard option included in base price
Colour	01 Black	




Item no	Seat width (sw)	Frame length	Seat depth
13123111	47.5 cm	Standard	36-49 cm
13123112	47.5 cm	Long	42-55 cm
13123113	50 cm	Standard	36-49 cm
13123114	50 cm	Long	42-55 cm
13123115	52.5 cm	Standard	36-49 cm
13123116	52.5 cm	Long	42-55 cm
13123117	55 cm	Standard	36-49 cm
13123118	55 cm	Long	42-55 cm
13123119	57.5 cm	Standard	36-49 cm
13123120	57.5 cm	Long	42-55 cm
13123121	60 cm	Standard	36-49 cm
13123122	60 cm	Long	42-55 cm

OPTIONS

Back supports & Covers				Config no	Item no
3A back support 	Back height 32-45 cm		inkl cross bar & adjustable strap back item:		
	recline -5° - +20°	sw 47.50-50 cm	62654-60	43438	28845
	lumbar angle -16° - +16°	sw 52.50-60 cm	62655-60	43439	28846
	Pre-set back height				40 cm
	Pre-set recline angle		22041		2°
	Back height 38-51 cm		inkl back 32-45 and 3A back extension item:		
	recline -5° - +20°	sw 47.50-50 cm	28845 + 28659	43441	
	lumbar angle -16° - +16°	sw 52.50-60 cm	28846 + 28660	43442	
	Pre-set back height				46 cm
	Pre-set recline angle		22041		2°

2B back support 	Back height Low, 35-40 cm		with fixed push handles		
	sw 47.5-50 cm		43255		28382
	sw 52.5-60 cm		43256		28383
	Back height High, 42.5-47.5 cm		with fixed push handles		
	sw 47.5-50 cm		43258		28385
	sw 52.5-60 cm		43259		28386

Back supports & Covers					Config no	Item no
Covers	Plush					
	Comfort					
	Airmesh					
	Tidy/hygiene					
						
	Back height 32-45 cm					
		Plush	Comfort	Airmesh	Tidy/hygiene	
	sw 47.5 cm	62946-60	62961-60	63209-68	63092-63	
	sw 50 cm	62947-60	62962-60	63210-68	63093-63	
	sw 52.5 cm	62948-60	62963-60	63211-68	63094-63	
sw 55 cm	62949-60	62964-60	63212-68	63095-63		
sw 57.5 cm	62950-60	62965-60	63213-68	63096-63		
sw 60 cm	62951-60	62966-60	63214-68	63097-63		
Back height 38-51 cm						
sw 47.5 cm	63025-60	63038-60	-	-		
sw 50 cm	63026-60	63039-60	-	-		
sw 52.5 cm	63027-60	63040-60	-	-		
sw 55 cm	63028-60	63041-60	-	-		
sw 57.5 cm	63029-60	63042-60	-	-		
sw 60 cm	63030-60	63043-60	-	-		

Brakes					Config no	Item no	
Knee lever brake	 	Red angled handle	Right	Left	45678	27528-1/-2	
		Black angled handle	Right	Left		27532-1/-2	
		Extended foldable handle	Right	Left		27530-1/-2	
One hand brake		Extended foldable handle	Right	Left		27531-1/-2	
Drum brake							
Select rear wheels with drumbrake and height adjustable push handles with attendant brake below							
Leg supports (hangers)					Config no	Item no	
Fixed angle	 	Swing away, detachable	lower leg length 31-49 cm				
		Wide angle 70°	Right	Left		25574-01	
		Narrow angle 80°	Right	Left		25570-01	
		Lockable, swing away, detachable	lower leg length 31-49 cm				
		Wide angle 70°	Right	Left		25575-01	
Elevated		Prio incl calf support, hygiene	lower leg length 31-49 cm	Right	Left		
		sw 47.5-50 cm			43392/43393	27858-65-1/-2	
		sw 52.5-55 cm			43392/43393	27859-65-1/-2	
		sw 57.5-60 cm			43392/43393	27860-65-1/-2	
		Mod 1	lower leg length 31-49 cm	Right	Left		25660-01
		add on:					
		Attachment for calf support				24544-01	
		Attachment for calf support flip-up				27322-01	
		Calf support, plush item:					
		sw 47.5 cm			62484-60		
		sw 50-60 cm			62451-60		
		Mod 3	lower leg length 31-49 cm	Right	Left		
		excl calf support				27149	
		incl calf support				26100	
	Plaster		3 cushions	Hygiene	Right	Left	27417-1/-2

Amputation



Amputation support	Right	Left
Plush	26210-60-1/-2	
Hygiene	26210-63-1/-2	
Amputation weight	26876	

Foot supports

Item no

Divided



Flip up	attachment tube:	
standard	Right	Left
+10 cm	Right	Left
-8 cm	Right	Left
Heel strap	for divided foot supports	
	24458	

One piece

Flip up				
				
Divided standard	Divided +10 cm	Divided -8 cm	One-piece	
47.5 cm	25736-1/2	25737-1/2	25738-1/2	25937
50 cm	25646-1/2	25652-1/2	25658-1/2	25937
52.5-60 cm	25646-1/2 + 51170	25652-1/2 + 51170	25658-1/2 + 51170	26853

Rear wheels

Config no

Item no

C6




20"	Aluminium	Stainless steel	Plastic coated
Solid	83885	83910	83911
High pressure, Speedlite	83908	83909	83884
22"			
Solid	83882	83907	83883
High pressure, Speedlite	83880	83906	83881
24"			
Solid	83876	83877	83878
High pressure, Speedlite	83873	83874	83875

Maxx


25"			
High pressure, Trailblazer	28964-1/-2 (43701)	-	28965-1/-2 (43702)

C5 Drum brake



20" DB	Aluminium	Stainless steel	Plastic coated
Solid	81337 (23608)	-	-
22" DB			
Solid	81339 (23602)	81341 (23601)	81340 (23603)
High pressure	81342 (23605)	81344 (23604)	81343 (23606)
24" DB			
Solid	81346 (23595)	81348 (23594)	81347 (23596)
High pressure	28625 (23598)	28624 (23597)	81350 (23599)
26" DB			
Solid innertube	-	28112 (23992)	-

Rear seat height



Heavy to propel and hard to tip. Efficient to propel and easy to tip.



	20"	22"	24"	25"	26"
G	34 cm	36.5 cm	39 cm	40.5 cm	41.5 cm
F	36 cm	38.5 cm	41 cm	42.5 cm	43.5 cm
E	38 cm	40.5 cm	43 cm	44.5 cm	45.5 cm
D	40 cm	42.5 cm	45 cm	46.5 cm	47.5 cm
C	42 cm	44.5 cm	47 cm	48.5 cm	49.5 cm
B	44 cm	46.5 cm	49 cm	50.5 cm	51.5 cm
A	46 cm	48.5 cm	51 cm	52.5 cm	53.5 cm

Seating

Config no


Item no

Upholstery

Depth adjustable

Hard seat

Height adjustable



	Upholstery basic std	Upholstery basic long	Hard seat std	Hard seat long
sw 47.5 cm	62978-60	62984-60	26722	26727
sw 50 cm	62979-60	62985-60	26723	26728
sw 52.5 cm	62980-60	62986-60	26723	26728
sw 55 cm	62981-60	62987-60	26724	26729
sw 57.5 cm	62982-60	62988-60	26724	26729
sw 60 cm	62983-60	62989-60	26660	26662

Seat depth extension



sw 47.5-50 cm	63292
sw 52.5-60 cm	63293

Push handles

Config no

Item no

Push handles



Height adjustable

100 mm

27589-01

incl Cross bar

sw 47.5 cm	43178	28695
sw 50 cm	43179	28696
sw 52.5 cm	43180	28697
sw 55 cm	43181	28698
sw 57.5 cm	43182	28699
sw 60 cm	43183	28700

Push handles with drum brake



Height adjustable	100 mm handle and 100 mm grip		27626-01
incl Cross bar	sw 47.5 cm	43170	28695
	sw 50 cm	43171	28696
	sw 52.5 cm	43172	28697
	sw 55 cm	43173	28698
	sw 57.5 cm	43174	28699
	sw 60 cm	43175	28700

Bow handle



Height adjustable			
	sw 47.5 cm		27434-01
	sw 50 cm		26063-01
	sw 52.5 cm		26064-01
	sw 55 cm		26068-01
	sw 57.5 cm		26069-01
	sw 60 cm		26070-01

Arm supports and mudguards

Config no

Item no

Arm support



Solid pad				
	25 cm	Right	Left	25172-1/-2
	38 cm	Right	Left	25173-1/-2
w fingerscrew for height adjustment	25 cm	Right	Left	23910/23911 26379-1/-2
w fingerscrew for height adjustment	38 cm	Right	Left	23912/23913 26380-1/-2



Removable pad				
	25 cm	Right	Left	28780-1/-2
	35 cm	Right	Left	28781-1/-2
w fingerscrew for height adjustment	25 cm	Right	Left	43557/43558 28782-1/-2
w fingerscrew for height adjustment	35 cm	Right	Left	43559/43560 28783-1/-2

Mudguard



		Right	Left	27397-1/-2
--	--	-------	------	-------------------

Anti tips

Config no

Item no

Anti tips multi cc59



	applicable for all seat heights	43427	28457
--	---------------------------------	-------	--------------

Seat cushions

Universal

Plush cover

Tidy cover (hygiene)




	Plush	Tidy
sw 47.5 cm	62374-60	63135-63
sw 50 cm	62266-60	63136-63
sw 52.5 cm	62420-60	63137-63
sw 55 cm	62421-60	63138-63
sw 57.5 cm	62422-60	63139-63
sw 60 cm	62422-60	63139-63

Widening kit

Config no

Item no

Widening kit



sw +1 cm	Not with drum brake	43080	27533
----------	---------------------	-------	--------------

ACCESSORIES

Arm supports				Pieces	Item no	
Hemi arm support		for item 25173-1/-2	Right	Left	25770-60-1/-2	
Wide and soft arm pad		8 x 40 cm	for item 25172-1/-2 and 25173-1/-2		226873	
Cover		25 cm	for item 25172-1/-2			
		padded plush		2	62282-60	
		padded hygiene		2	62282-63	
		gel plush		2	26240-60	
		gel hygiene		2	26240-63	
		38 cm	for item 25173-1/-2			
		padded plush		2	62283-60	
		padded hygiene		2	62283-63	
		gel plush		2	26241-60	
		gel hygiene		2	26241-63	
Side guard pad		Cross, standard shape	reduce the seat width 30mm/pair	Right	Left	63373-1/-2
Thigh support	 	Short 50 cm	for item 25172-1/-2 and 25173-1/-2	Right	Left	27391-1/-2
		Long 57 cm	for item 25172-1/-2 and 25173-1/-3	Right	Left	27392-1/-2
		Soft pad for thigh support		Right	Left	82488-1/-2
Attachments for medical devices				Pieces	Item no	
Crutch and cane holder		Universal		1	27327	
IV stand		For two IV bags complete with attachment		1	28323	
Oxygen holder		For oxygen tubes, maximum load 5 kg		1	28937	
Fixations and brackets					Item no	
Locking strap		to lock the wheelchair when folded			27151	
Head supports					Item no	
Attachment	 	Assembly on 3A back support				
		sw 47.5 cm			28130	
		sw 50 cm			28131	
		sw 52.5 cm			28132	
		sw 55 cm			28133	
		sw 57.5 cm			28134	
		sw 60 cm			28135	
		Assembly on bow handle				
			25795			

Head support		Traditional	26x14 cm incl turtle bar 350 mm			
			with plush cover		28534-60	
			with hygiene cover		28534-65	
Memory function		With side support	incl turtle bar 250 mm			
			with plush cover		28536-60	
			with hygiene cover		28536-65	
		Memory function	for head support			28564

Leg and foot supports						Pieces	Item no
Cover		to cover attachment when not using leg support				2	51213
Foot support extension		Extension, plastic	for right & left foldable foot support				
			sw 47.5 cm		2	50799	
			sw 50-60 cm		2	50800	
Add-on extension		Extension, metal	for right & left foldable foot support				
		Medium	sw 47.5 cm	Right	Left	28775-1/-2	
		Large	sw 50-60 cm	Right	Left	28776-1/-2	
		Guide to link right & left extension into one piece					28773
		Side protection	for extension, metal		2	28772	
Foot box		Foot box	Soft padded				
		Medium	sw 47.5 cm			28642	
		large	sw 50-60 cm			28643	
Thigh support		Distal thigh support	for fixed leg supports,	Right	Left	28777-1/-2	
Cover foot support		For flip up foot plates	padded plush	Right	Left		
			sw 47.5-50 cm			63189-60-1/-2	
			sw 52.5-60 cm			63192-60-1/-2	
		For one piece foot plate	padded plush				
		sw 47.5-50 cm			63191-60		
Calf support		Cover, hygiene	for elevated leg support item 26100			62762-63	
		Calf strap	adjustable in length				
		sw 47.5 cm			25358		
		sw 50-52.5 cm			25359		
			sw 55-60 cm			26062	
		Padding for calf strap	30 cm			26846	
		40 cm			26847		
		Panel	Plush				
		75 cm (sw 47.5 cm)			63186-60		
		90 cm (sw 50-60 cm)			63187-60		

Postural supports			Pieces	Item no
Positioning belt, 2 point		plastic push release		
		100 cm		25668
		128 cm		26163
Wedge		push button buckle		
		170 cm		62190
		Soft		
		Plush		
		Short, 15 cm		63012-60
		Long, 30 cm		63013-60
Lateral support		Cell foam		
		Polyethylene		
		90 x 20 x 130 mm		28030
		110 x 20 x 185 mm		28139
		175 x 25 x 200 mm		28140
		Side cushion	Right Left	63171-60-1/-2
Back width extension		Side stop, straight		
		gives extra support to side cushion		28136
		Width extension, 30 mm	combine with wedges	28173
Trunk support 3A		Pressure distributor		
		padded cover with flexible insert		28637
		Swing away & detachable		
		Adjustable in height, depth and width		
		0 mm (straight profile)	Right Left	28181-1/-2
		-15 mm	Right Left	28180-1/-2
Cover for trunk support		Swing away R82 pad		
		Adjustable in height and depth		
		Medium	Right Left	27915-1/-2
		Large	Right Left	27916-1/-2
		Soft		
		Adjustable in height		27395
Attachment for harness		Soft cover		
		Fits item 28142, 28143, 28144, 28145, 28181-1/-2, 28180-1/-2, 28141-1/-2		28647
		Upper attachment		
		to be mounted on head support stem 15x15 mm		28186

Seating				
Seat wedge		Comfort wedge		
		Provides extra seat contour for seat cushion Universal.		
		sw 47.5-60 cm		80779

Trays

Tray



Transparent	for 38 cm arm support 25173-1/-2		
	sw 47.5 cm		25731
	sw 50 cm		24601
	sw 52.5 cm		26020
	sw 55 cm		26021
	sw 57.5 cm		26022
	sw 60 cm		26907

Hemi tray



Transparent	for 25 cm arm support 25173-1/-2		
	sw 47.5 cm		26561
	sw 50 cm		26562
Transparent	for 38 cm arm support 25173-1/-2		
	sw 47.5 cm		26568
	sw 50 cm		26569
Anti glide	for hemi tray		81633

Half tray



Foldable	for arm support 25172-1/-2 & 25173-1/-2	Right	Left	27352-1/-2
----------	---	-------	------	------------

Attachment for tray

Round tubes	fits tubes up to ø 18mm	Right	Left	27956-1/-2
-------------	-------------------------	-------	------	------------

Wheels

Pieces

Item no

Spokeguard



Transparent			
	20"	2	28934
	22"	2	28935
	24"	2	28936



Transparent	for CS DB		
	20"	2	28816-28
	22"	2	28817-28
	24"	2	28818-28

Additional information to order:

Etac Supply Center AB
 Phone: +46 371 58 73 30
 Mail: info@etac.com
 Website: www.etac.com


 Creating Possibilities