



# Cross 5XL

Therapist	User	
Phone no	Order date	
Purchase Order no	Customer Ref	
Base price	€XXX	Standard option included in base price
Colour	01 Black	



Item no	Seat width (sw)	Frame length	Seat depth
13123111	47.5 cm	Standard	36-49 cm
13123112	47.5 cm	Long	42-55 cm
13123113	50 cm	Standard	36-49 cm
13123114	50 cm	Long	42-55 cm
13123115	52.5 cm	Standard	36-49 cm
13123116	52.5 cm	Long	42-55 cm
13123117	55 cm	Standard	36-49 cm
13123118	55 cm	Long	42-55 cm
13123119	57.5 cm	Standard	36-49 cm
13123120	57.5 cm	Long	42-55 cm
13123121	60 cm	Standard	36-49 cm
13123122	60 cm	Long	42-55 cm

## OPTIONS

### Back supports & Covers Config no    Item no

#### 3A back support



Back height 32-45 cm		inkl cross bar & adjustable strap back item:	
recline -5° - +20°	sw 47.50-50 cm	62654-60	43438 <b>28845</b>
lumbar angle -16° - +16°	sw52.50-60 cm	62655-60	43439 <b>28846</b>

Pre-set back height	<b>40 cm</b>
Pre-set recline angle	22041 <b>2°</b>

Back height 38-51 cm		inkl back 32-45 and 3A back extension item:	
recline -5° - +20°	sw 47.50-50 cm	28845 + 28659	43441
lumbar angle -16° - +16°	sw52.50-60 cm	28846 + 28660	43442

Pre-set back height	<b>46 cm</b>
Pre-set recline angle	22041 <b>2°</b>

**Back supports & Covers** Config no **Item no**

**Covers**

Plush

---

Comfort

---

Airmesh

---

Tidy/hygiene


---



	Plush	Comfort	Airmesh	Tidy/hygiene
<b>Back height 32-45 cm</b>				
sw 47.5 cm	62946-60	62961-60	63209-68	63092-63
sw 50 cm	62947-60	62962-60	63210-68	63093-63
sw 52.5 cm	62948-60	62963-60	63211-68	63094-63
sw 55 cm	62949-60	62964-60	63212-68	63095-63
sw 57.5 cm	62950-60	62965-60	63213-68	63096-63
sw 60 cm	62951-60	62966-60	63214-68	63097-63
<b>Back height 38-51 cm</b>				
sw 47.5 cm	63025-60	63038-60	-	-
sw 50 cm	63026-60	63039-60	-	-
sw 52.5 cm	63027-60	63040-60	-	-
sw 55 cm	63028-60	63041-60	-	-
sw 57.5 cm	63029-60	63042-60	-	-
sw 60 cm	63030-60	63043-60	-	-


**Brakes** Config no **Item no**

**Knee lever brake**



Red angled handle	Right	Left	45678	<b>27528-1/-2</b>
Black angled handle	Right	Left		<b>27532-1/-2</b>
Extended foldable handle	Right	Left		<b>27530-1/-2</b>

**One hand brake**




Single hand brake	Right	Left		<b>27531-1/-2</b>
-------------------	-------	------	--	-------------------

**Drum brake** Select rear wheels with drumbrake and height adjustable push handles with attendant brake below


**Leg supports (hangers)** Config no **Item no**

**Fixed angle**




Swing away, detachable	lower leg length 31-49 cm			
Wide angle 70°		Right	Left	<b>25574-01</b>
Narrow angle 80°		Right	Left	<b>25570-01</b>
Lockable, swing away, detachable	lower leg length 31-49 cm			
Wide angle 70°		Right	Left	<b>25575-01</b>

**Elevated**



Prio incl calf support, hygiene	lower leg length 31-49 cm	Right	Left	
sw 47.5-50 cm				43392/43393 <b>27858-65-1/-2</b>
sw 52.5-55 cm				43392/43393 <b>27859-65-1/-2</b>
sw 57.5-60 cm				43392/43393 <b>27860-65-1/-2</b>
Mod 1	lower leg length 31-49 cm	Right	Left	<b>25680-01</b>
add on:	Attachment for calf support			<b>24544-01</b>
	Attachment for calf support flip-up			<b>27322-01</b>
	Calf support, plush item:			
	sw 47.5 cm			62484-60
	sw 50-60 cm			62451-60
Mod 3	lower leg length 31-49 cm	Right	Left	
	excl calf support			<b>27149</b>
	incl calf support			<b>26100</b>

**Plaster**



3 cushions	Hygiene	Right	Left	<b>27417-1/-2</b>
------------	---------	-------	------	-------------------

**Amputation**



Amputation support	Right	Left
Plush	<b>26210-60-1/-2</b>	
Hygiene	<b>26210-63-1/-2</b>	
Amputation weight	<b>26876</b>	

**Foot supports**

**Item no**

**Divided**



Flip up	attachment tube:	
standard	Right	Left
+10 cm	Right	Left
-8 cm	Right	Left
Heel strap	for divided foot supports	
	<b>24458</b>	

**One piece**

Flip up				
				
	Divided standard	Divided +10 cm	Divided -8 cm	One-piece
47.5 cm	25736-1/2	25737-1/2	25738-1/2	25937
50 cm	25646-1/2	25652-1/2	25658-1/2	25937
52.5-60 cm	25646-1/2 + 51170	25652-1/2 + 51170	25658-1/2 + 51170	26853

**Rear wheels**

*Config no*

**Item no**

**C6**



20"	Aluminium	Stainless steel	Plastic coated
Solid	83885	83910	83911
High pressure, Speedlite	83908	83909	83884
22"			
Solid	83882	83907	83883
High pressure, Speedlite	83880	83906	83881
24"			
Solid	83876	83877	83878
High pressure, Speedlite	83873	83874	83875
25"			
High pressure, Trailblazer	28964-1/-2 (43701)	-	28965-1/-2 (43702)

**Maxx**

**C5 Drum brake**



20" DB	Aluminium	Stainless steel	Plastic coated
Solid	81337B (23608)	-	-
22" DB			
Solid	81339B (23602)	81341B (23601)	81340B (23603)
High pressure	81342B (23605)	81344B (23604)	81343B (23606)
24" DB			
Solid	81346B (23595)	81348B (23594)	81347B (23596)
High pressure	28625B (23598)	28624B (23597)	81350B (23599)
26" DB			
Solid innertube	-	28112B (23992)	-

**QR-axle**



Axle QR D12,7x110	for C6 and Maxx	43727	<b>83872</b>
-------------------	-----------------	-------	--------------



Axle QR Drum brake	for C5 Drum brake		<b>81325</b>
--------------------	-------------------	--	--------------

**One arm drive**



Drive unit	Right	Left	
20"			<b>26534</b>
22"			<b>26533</b>
24"			<b>26522</b>
<b>Telescopic shaft</b>			
sw 47,5-50			27319
sw 52,5-57,5			27320
sw 60			27321

**Castors** (see seat height chart below)

Config no

4"		100 x 34	seat height 34-39.5	max user weight 100kg	43690	<b>28919</b>
5"		125 x 34	seat height 34.5-48.5	max user weight 135kg	43691	<b>28921</b>
6"		150 x 34	seat height 38.5-50		43692	<b>28922</b>
8"		200 x 38	seat height 43.5-56.5		43693	<b>28526</b>

**Front fork housing** (see seat height chart below)

Config no

<b>Standard</b>		H100	seat height 36-52.5		22019	<b>28147</b>
<b>Widened</b>		H100	seat height 36-52.5		22017	<b>25763</b>
<b>Low Widened</b>		H100/H75	seat height 34-39.5		23478	<b>28633</b>
<b>High</b>		H142	seat height 50-56.5		43694	<b>28968</b>

**Front forks** (see seat height chart below)

Config no

<b>XS</b>		Aluminium	seat height 34-39.5			<b>28923</b>
<b>S</b>		Steel	seat height 34.5-47.5			<b>24038</b>
<b>M</b>		Steel	seat height 41-56.5			<b>28969</b>

**Camber**


Config no

Item no

2°		outside/inside (2x2 pcs)	23476	<b>82567/82568</b>
0°		outside/inside (2x2 pcs)	23477	<b>82569/82570</b>

**Seat height**



**Front seat height**



Front fork position

1	2	3	4
---	---	---	---

Rear seat height



	20*	22*	24*	25*	26*
G	34 cm	36.5 cm	39 cm	40.5 cm	41.5 cm
F	36 cm	38.5 cm	41 cm	42.5 cm	43.5 cm
E	38 cm	40.5 cm	43 cm	44.5 cm	45.5 cm
D	40 cm	42.5 cm	45 cm	46.5 cm	47.5 cm
C	42 cm	44.5 cm	47 cm	48.5 cm	49.5 cm
B	44 cm	46.5 cm	49 cm	50.5 cm	51.5 cm
A	46 cm	48.5 cm	51 cm	52.5 cm	53.5 cm


Seating

Config no

Item no

Upholstery  
Hard seat

Depth adjustable  
Height adjustable



	Upholstery basic std	Upholstery basic long	Hard seat std	Hard seat long
sw 47.5 cm	62978-60	62984-60	26722	26727
sw 50 cm	62979-60	62985-60	26723	26728
sw 52.5 cm	62980-60	62986-60	26723	26728
sw 55 cm	62981-60	62987-60	26724	26729
sw 57.5 cm	62982-60	62988-60	26724	26729
sw 60 cm	62983-60	62989-60	26660	26662

Seat depth extension



sw 47.5-50 cm	63292
sw 52.5-60 cm	63293

Push handles

Config no

Item no

Push handles



	Height adjustable	100 mm	27589-01
incl Cross bar	sw 47.5 cm	43178	28695
	sw 50 cm	43179	28696
	sw 52.5 cm	43180	28697
	sw 55 cm	43181	28698
	sw 57.5 cm	43182	28699
	sw 60 cm	43183	28700

**Push handles with drum brake**



Height adjustable	100 mm handle and 100 mm grip			<b>27626-01</b>
incl Cross bar	sw 47.5 cm		43170	<b>28695</b>
	sw 50 cm		43171	<b>28696</b>
	sw 52.5 cm		43172	<b>28697</b>
	sw 55 cm		43173	<b>28698</b>
	sw 57.5 cm		43174	<b>28699</b>
	sw 60 cm		43175	<b>28700</b>

**Bow handle**



Height adjustable				
	sw 47.5 cm			<b>27434-01</b>
	sw 50 cm			<b>26063-01</b>
	sw 52.5 cm			<b>26064-01</b>
	sw 55 cm			<b>26068-01</b>
	sw 57.5 cm			<b>26069-01</b>
	sw 60 cm			<b>26070-01</b>


**Arm supports and mudguards**

Config no Item no

**Arm support**



Solid pad				
	25 cm	Right	Left	<b>25172-1/-2</b>
	38 cm	Right	Left	<b>25173-1/-2</b>
w fingerscrew for height adjustment	25 cm	Right	Left	23910/23911 <b>26379-1/-2</b>
w fingerscrew for height adjustment	38 cm	Right	Left	23912/23913 <b>26380-1/-2</b>



Removable pad				
	25 cm	Right	Left	<b>28780-1/-2</b>
	35 cm	Right	Left	<b>28781-1/-2</b>
w fingerscrew for height adjustment	25 cm	Right	Left	43557/43558 <b>28782-1/-2</b>
w fingerscrew for height adjustment	35 cm	Right	Left	43559/43560 <b>28783-1/-2</b>

**Mudguard**



Mudguard profile 220		Right	Left	
	S			<b>24094-1/-2</b>
	M			<b>24095-1/-2</b>
	L			<b>24096-1/-2</b>

	20°	22°	24°	25°	26°
H	L				
G	M/L	L			
F	S/M	M/L	L	L	
E	S	M	M/L	L	L
D	S*	S	M	M	M/L
C	S**	S*	S	S/M	M
B	S***	S**	S*	S	S
A	S***	S***	S**	S*	S*

Note: conflict may occur, check manual.

**Anti tips**

Config no Item no

**Anti tips multi cc59**



applicable for all seat heights				43427	<b>28457</b>
---------------------------------	--	--	--	-------	--------------

**Seat cushions**

**Universal**

Plush cover  
Tidy cover (hygiene)




	<b>Plush</b>	<b>Tidy</b>
sw 47.5 cm	62374-60	63135-63
sw 50 cm	62266-60	63136-63
sw 52.5 cm	62420-60	63137-63
sw 55 cm	62421-60	63138-63
sw 57.5 cm	62422-60	63139-63
sw 60 cm	62422-60	63139-63

**Widening kit**

Config no Item no

**Widening kit**



sw +1 cm	Not with drum brake		43080	<b>27533</b>
----------	---------------------	--	-------	--------------

## ACCESSORIES

### Arm supports

Pieces

Item no

#### Hemi arm support



for item 25173-1/-2

Right

Left

25770-60-1/-2

#### Wide and soft arm pad



8 x 40 cm

for item 25172-1/-2 and 25173-1/-2

2

26873

#### Cover



25 cm

for item 25172-1/-2

padded plush

2

62282-60

padded hygiene

2

62282-63

gel plush

2

26240-60

gel hygiene

2

26240-63

38 cm

for item 25173-1/-2

padded plush

2

62283-60

padded hygiene

2

62283-63

gel plush

2

26241-60

gel hygiene

2

26241-63

#### Side guard pad



Cross, standard shape

reduce the seat width 30mm/pair

Right

Left

63373-1/-2

#### Thigh support

Not with 22" in pos C1

Not with 24" in pos C1 or B1



Short 50 cm

for item 25172-1/-2 and 25173-1/-2

Right

Left

27391-1/-2

Long 57 cm

for item 25172-1/-2 and 25173-1/-3

Right

Left

27392-1/-2

Soft pad for thigh support

Right

Left

82488-1/-2

### Attachments for medical devices

Pieces

Item no

#### Crutch and cane holder

Universal

1

27327

#### IV stand



For two IV bags complete with attachment

1

28323

#### Oxygen holder



For oxygen tubes, maximum load 5 kg

1

28937

### Head supports

Item no

#### Attachment



Assembly on 3A back support

sw 47.5 cm

28130

sw 50 cm

28131

sw 52.5 cm

28132

sw 55 cm


28133

sw 57.5 cm

28134

sw 60 cm

28135



Assembly on bow handle

25795

**Head support**



Traditional	26x14 cm incl turtle bar 350 mm		
	with plush cover		28534-60
	with hygiene cover		28534-65



With side support	incl turtle bar 250 mm		
	with plush cover		28536-60
	with hygiene cover		28536-65

**Extension**



Extension head support	for head support		29012
------------------------	------------------	--	-------

Only extension, not complete head support.

**Memory function**



Memory function	for head support		28564
-----------------	------------------	--	-------

**Leg and foot supports**

Pieces

Item no

**Cover**



to cover attachment when not using leg support		2	51213
--	--	---	-------

**Foot support extension**



Extension, plastic	for right & left foldable foot support		
	sw 47.5 cm	2	50799
	sw 50-60 cm	2	50800

**Add-on extension**



Extension, metal	for right & left foldable foot support		
Medium	sw 47.5 cm	Right	Left
			28775-1/-2
Large	sw 50-60 cm	Right	Left
			28776-1/-2

Guide to link right & left extension into one piece			28773
---	--	--	-------

**Foot box**



Side protection	for extension, metal	2	28772
-----------------	----------------------	---	-------



Foot box	Soft padded		
Medium	sw 47.5 cm		28642
large	sw 50-60 cm		28643

**Thigh support**

Not with lockable or elevated leg supports



Distal thigh support	for fixed leg supports,	Right	Left
			28777-1/-2

**Cover foot support**



For flip up foot plates	padded plush	Right	Left
	sw 47.5-50 cm		
	sw 52.5-60 cm		
			63189-60-1/-2



For one piece foot plate	padded plush		
	sw 47.5-50 cm		

**Calf support**



Cover, hygiene	for elevated leg support item 26100		62762-63
----------------	-------------------------------------	--	----------



Calf strap	adjustable in length		
	sw 47.5 cm		25358
	sw 50-52.5 cm		25359
	sw 55-60 cm		26062



Padding for calf strap			
	30 cm		26846
	40 cm		26847



Panel	Plush		
	75 cm (sw 47.5 cm)		63186-60
	90 cm (sw 50-60 cm)		63187-60

**Postural supports**

Pieces

Item no

**Positioning belt, 2 point**



plastic push release

100 cm	25668
128 cm	26163

push button buckle

170 cm	62190
--------	-------

**Wedge**



Soft

Plush		pcs
Short, 15 cm		63012-60
Long, 30 cm		63013-60

Cell foam

Polyethylene		pair
90 x 20 x 130 mm		28030
110 x 20 x 185 mm		28139
175 x 25 x 200 mm		28140

**Lateral support**



Side cushion

Right	Left	63171-60-1/-2
-------	------	---------------

Side stop, straight

gives extra support to side cushion	28136
-------------------------------------	-------

**Back width extension**



Width extension, 30 mm

combine with wedges	28173
---------------------	-------

Widening pad

padded cover with flexible insert	28637
-----------------------------------	-------

**Trunk support 3A**



Swing away & detachable

Adjustable in height, depth and width			
0 mm (straight profile)	Right	Left	28181-1/-2
-15 mm	Right	Left	28180-1/-2
-45 mm	Right1	Left	28141-1/-2



Detachable

Adjustable in height, depth and width			
0 mm (straight profile)			2 28142
+15 mm			2 28143
-15 mm			2 28144
-45 mm			2 28145



Soft

Adjustable in height	27395
----------------------	-------



Swing away R82 pad

Adjustable in height and depth			
Medium	Right	Left	27915-1/-2
Large	Right	Left	27916-1/-2

**Cover for trunk support**



Soft cover

Fits item 28142, 28143, 28144, 28145, 28181-1/-2, 28180-1/-2, 28141-1/-2	28647
--	-------

**Attachment for harness**



Upper attachment

to be mounted on head support stem 15x15 mm	28186
---	-------

**Seating**

**Seat wedge**



Comfort wedge

Provides extra seat contour for seat cushion Universal.	
sw 47.5-60 cm	80779

**Trays**

**Tray**



Transparent	for 38 cm arm support 25173-1/-2		
	sw 47.5 cm		<b>25731</b>
	sw 50 cm		<b>24601</b>
	sw 52.5 cm		<b>26020</b>
	sw 55 cm		<b>26021</b>
	sw 57.5 cm		<b>26022</b>
	sw 60 cm		<b>26907</b>

**Hemi tray**



Transparent	for 25 cm arm support 25172-1/-2		
	sw 47.5 cm		<b>26561</b>
	sw 50 cm		<b>26562</b>

Transparent	for 38 cm arm support 25173-1/-2		
	sw 47.5 cm		<b>26568</b>
	sw 50 cm		<b>26569</b>

Anti glide	for hemi tray		
			<b>81633</b>

**Half tray**



Foldable	for arm support 25172-1/-2 & 25173-1/-2	Right	Left	
				<b>27352-1/-2</b>

**Attachment for tray**

Round tubes	fits tubes up to ø 18mm	Right	Left	
				<b>27956-1/-2</b>

**Wheels**

Pieces

Item no

**Spokeguard**



Transparent			
	20"	<b>2</b>	<b>28934</b>
	22"	<b>2</b>	<b>28935</b>
	24"	<b>2</b>	<b>28936</b>



Transparent	for C5 DB		
	20"	<b>2</b>	<b>28816-28</b>
	22"	<b>2</b>	<b>28817-28</b>
	24"	<b>2</b>	<b>28818-28</b>

**Additional information to order:**

---



---



---



---



---



---

**Etac Supply Center AB**  
 Phone: +46 371 58 73 30  
 Mail: [Info@etac.com](mailto:Info@etac.com)  
 Website: [www.etac.com](http://www.etac.com)



**etac**<sup>®</sup>  
 Creating Possibilities