

# Cross 6



<b>Therapist</b>	<b>User</b>	
<b>Phone no</b>	<b>Order date</b>	
<b>Purchase Order no</b>	<b>Customer Ref</b>	
<b>Base price</b>	€XXX	Standard option included in base price
<b>Colour</b>	01 Black	
	78 Grey	
	92 White	

Item no	Seat width (sw)	Frame length	Seat depth
13240101	35 cm	Standard	36-49 cm
13240102	35 cm	Long	42-55 cm
13240103	37.5 cm	Standard	36-49 cm
13240104	37.5 cm	Long	42-55 cm
13240105	40 cm	Standard	36-49 cm
13240106	40 cm	Long	42-55 cm
13240107	42.5 cm	Standard	36-49 cm
13240108	42.5 cm	Long	42-55 cm
13240109	45 cm	Standard	36-49 cm
13240110	45 cm	Long	42-55 cm
13240111	47.5 cm	Standard	36-49 cm
13240112	47.5 cm	Long	42-55 cm
13240113	50 cm	Standard	36-49 cm
13240114	50 cm	Long	42-55 cm
13240115	52.5 cm	Standard	36-49 cm
13240116	52.5 cm	Long	42-55 cm

## OPTIONS

### Back supports & Covers Config no    Item no

#### 3A back support



Back height 32-45 cm	inkl adjustable strap back item:		
recline -5° - +20°	sw 35-40 cm	62953-60	43449 <b>28586</b>
lumbar angle -16° - +16°	sw 42.5-50 cm	62654-60	43450 <b>28587</b>
	sw 52.5 cm	62955-60	43451 <b>28588</b>

Pre-set back height	
Pre-set recline angle	22041

Back height 38-51 cm	inkl back 32-45 and 3A back extension item:		
recline -5° - +20°	sw 35-40 cm	28586 + 28658	43452
lumbar angle -16° - +16°	sw 42.5-50 cm	28587 + 28659	43453
	sw 52.5 cm	28588 + 28660	43454

Pre-set back height	
Pre-set recline angle	22041

#### 2B back support



Back height Low, 35-40 cm	with fixed push handles		
	sw 35-40 cm		43504 <b>28381</b>
	sw 42.5-47.5 cm		43505 <b>28382</b>
	sw 50-55 cm		43506 <b>28383</b>

with attachment for height adjustable push handles			
	sw 35-40 cm		43534 <b>28597</b>
	sw 42.5-47.5 cm		43535 <b>28598</b>
	sw 50-55 cm		43536 <b>28599</b>



Back height High, 42.5-47.5 cm	with fixed push handles		
	sw 35-40 cm		43507 <b>28384</b>
	sw 42.5-47.5 cm		43508 <b>28385</b>
	sw 50-55 cm		43509 <b>28386</b>

with attachment for height adjustable push handles			
	sw 35-40 cm		43537 <b>28594</b>
	sw 42.5-47.5 cm		43538 <b>28595</b>
	sw 50-55 cm		43539 <b>28596</b>

**Back supports & Covers**
*Config no*
**Item no**
**Covers**

Plush

Comfort

Airmesh

Tidy/hygiene



Back height 32-45 cm	Plush	Comfort	Airmesh	Tidy/hygiene
sw 35 cm	62941-60	62956-60	63204-68	63087-63
sw 37.5 cm	62942-60	62957-60	63205-68	63088-63
sw 40 cm	62943-60	62958-60	63206-68	63089-63
sw 42.5 cm	62944-60	62959-60	63207-68	63090-63
sw 45 cm	62945-60	62960-60	63208-68	63091-63
sw 47.5 cm	62946-60	62961-60	63209-68	63092-63
sw 50 cm	62947-60	62962-60	63210-68	63093-63
sw 52.5 cm	62948-60	62963-60	63211-68	63094-63
<b>Back height 38-51 cm</b>				
sw 35 cm	63020-60	63033-60	-	-
sw 37.5 cm	63021-60	63034-60	-	-
sw 40 cm	63022-60	63035-60	-	-
sw 42.5 cm	63023-60	63036-60	-	-
sw 45 cm	63024-60	63037-60	-	-
sw 47.5 cm	63025-60	63038-60	-	-
sw 50 cm	63026-60	63039-60	-	-
sw 52.5 cm	63027-60	63040-60	-	-

**Brakes**
*Config no*
**Item no**
**Knee lever brake**


Red angled handle	Right	Left	45678	<b>27528-1/-2</b>
Black angled handle	Right	Left		<b>27532-1/-2</b>
Extended foldable handle	Right	Left		<b>27530-1/-2</b>

**One hand brake**


Single hand brake	Right	Left		<b>27531-1/-2</b>
-------------------	-------	------	--	-------------------

**Drum brake**

Select rear wheels with drumbrake and height adjustable push handles with attendant brake below

**Leg supports (hangers)**
*Config no*
**Item no**
**Fixed angle**


Swing away, detachable	lower leg length 31-49 cm			
	Wide angle 70°	Right	Left	<b>25574-01</b>
	Narrow angle 80°	Right	Left	<b>25570-01</b>
Narrow angle 80° short (23-33 cm)		Right	Left	<b>25571-01</b>
	Lockable, swing away, detachable	lower leg length 31-49 cm		
		Wide angle 70°	Right	Left
Narrow angle 80°		Right	Left	<b>25573-01</b>
Narrow angle 80° short (23-33 cm)		Right	Left	<b>25572-01</b>

**Elevated**



Prio incl calf support, hygiene	lower leg length 31-49 cm	Right	Left	
	sw 40-45 cm			43392/43393
	sw 47.5-50 cm			43392/43393
	sw 52.5 cm			43392/43393




Mod 1	lower leg length 31-49 cm	Right	Left	<b>25680-01</b>
	add on:			
	Attachment for calf support			<b>24544-01</b>
	Attachment for calf support flip-up			<b>27322-01</b>
	Calf support, plush item:			
	sw 35-37.5 cm			62105-60
	sw 40-42.5 cm			62104-60
	sw 45-47.5 cm			62484-60
	sw 50-52.5 cm			62451-60








Mod 3	lower leg length 31-49 cm	Right	Left	
	excl calf support			<b>27149</b>
	incl calf support			<b>28100</b>

<b>Plaster</b>		3 cushions	Hygiene	Right	Left	27417-1/-2
----------------	---	------------	---------	-------	------	------------


<b>Amputation</b>		Amputation support				
			Plush			25516-60
			Hygiene			25516-63
		Amputation weight				26876


**Foot supports** Item no


<b>Divided</b>		Flip up		attachment tube:			
			standard	Right	Left		
			+10 cm	Right	Left		
			-8 cm	Right	Left		
		Heel strap				for divided foot supports	24458


<b>One piece</b>	Flip up				
					
	<b>Divided standard</b>	<b>Divided +10 cm</b>	<b>Divided -8 cm</b>	<b>One-piece</b>	
	sw 35 cm	25641-1/2	25647-1/2	25653-1/2	25936
	sw 37.5 cm	25642-1/2	25648-1/2	25654-1/2	25936
	sw 40 cm	25643-1/2	25649-1/2	25655-1/2	25937
	sw 42.5 cm	25644-1/2	25650-1/2	25656-1/2	25937
	sw 45 cm	25645-1/2	25651-1/2	25657-1/2	25937
	sw 47.5 cm	25736-1/2	25737-1/2	25738-1/2	25937
	sw 50 cm	25646-1/2	25652-1/2	25658-1/2	25937
sw 52.5 cm	25646-1/3	25652-1/3	25658-1/3	26853	


**Rear wheels** Config no Item no


<b>C6</b>		20"	Aluminium	Stainless steel	Plastic coated
		Solid	83885	83910	83911
		High pressure, Speedlite	83908	83909	83884
		22"			
		Solid	83882	83907	83883
		High pressure, Speedlite	83880	83906	83881
24"					
Solid	83876	83877	83878		
High pressure, Speedlite	83873	83874	83875		


<b>Maxx</b>		22"	Hard anodised aluminium	Stainless steel	Plastic coated
		High pressure, Speedlite	28956	-	28957
		High pressure, Trailblazer	24052-1/-2 (43703)	-	24053-1/-2 (43704)
		24"			
		High pressure, Speedlite	28950	-	28951
		High pressure, Trailblazer	28952-1/-2 (43699)	-	28953-1/-2 (43700)
25"					
High pressure, Trailblazer	28964-1/-2 (43701)	-	28965-1/-2 (43702)		


<b>C5 Drum brake</b>		20"	Aluminium	Stainless steel	Plastic coated
		Solid	81337B (23608)	-	-
		22"			
		Solid	81339B (23602)	81341B (23601)	81340B (23603)
		High pressure	81342B (23605)	81344B (23604)	81343B (23606)
		24"			
		Solid	81346B (23595)	81348B (23594)	81347B (23596)
		High pressure	28625B (23598)	28624B (23597)	28626B (43461)
		26"			
		Solid innertube	-	28112B (23992)	-


<b>QR-axle</b>		Axle QR D12,7x110	for C6 and Maxx	43727	<b>83872</b>
		Axle QR Drum brake	for C5 Drum brake		<b>81325</b>

<b>One arm drive</b>		Drive unit	Right	Left	
			20"		<b>26534</b>
			22"		<b>26533</b>
			24"		<b>26522</b>
		Telescopic shaft			
			sw 35-37,5 cm		27317
	sw 40-42,5 cm		27318		
	sw 45-50 cm		27319		
	sw 52,5 cm		27320		

<b>Castors</b> (see seat height chart below)				Config no	Item no	
<b>4"</b>		100 x 34	seat height 32-38.5	max user weight 100kg	43690	<b>28919</b>
		125 x 34	seat height 34.5-50		43691	<b>28921</b>
		150 x 34	seat height 39.5-51.5		43692	<b>28922</b>
		200 x 38	seat height 44.5-57.5		43693	<b>28526</b>


<b>Front fork housings</b> (see seat height chart below)				Config no	Item no	
<b>Standard Widened Low Widened High</b>		H100	seat height 37.5-53.5		22019	<b>28147</b>
		H100	seat height 37.5-53.5		22017	<b>25763</b>
		H100/75	seat height 32-38.5		23478	<b>28633</b>
		H142	seat height 46-57.5		43694	<b>28968</b>

<b>Front forks</b> (see seat height chart below)				Config no	Item no	
<b>XS S M</b>		Aluminium	seat height 32-36			<b>28923</b>
		Steel	seat height 34.5-48.5			<b>24038</b>
		Steel	seat height 42.5-57.5			

<b>Camber</b>				Config no	Item no
<b>2° 0°</b>		outside/inside (2x2 pcs)		23476	<b>82567/82568</b>
		outside/inside (2x2 pcs)		23477	<b>82569/82570</b>


**Seat height**

**Front seat height**



Front fork position	1	2	3	4
---------------------	---	---	---	---

Rear seat height



	20*	22*	24*	25*	26*
H	32 cm	34.5 cm	37 cm	38.5 cm	39.5 cm
G	34 cm	36.5 cm	39 cm	40.5 cm	41.5 cm
F	36 cm	38.5 cm	41 cm	42.5 cm	43.5 cm
E	38 cm	40.5 cm	43 cm	44.5 cm	45.5 cm
D	40 cm	42.5 cm	45 cm	46.5 cm	47.5 cm
C	42 cm	44.5 cm	47 cm	48.5 cm	49.5 cm
B	44 cm	46.5 cm	49 cm	50.5 cm	51.5 cm
A	46 cm	48.5 cm	51 cm	52.5 cm	53.5 cm

Seating


Config no

Item no

Upholstery  
Hard seat

Depth adjustable

Height adjustable



	Upholstery basic std	Upholstery basic long	Hard seat std	Hard seat long
sw 35 cm	62718-60	62728-60	26720	26725
sw 37.5 cm	62719-60	62729-60	26720	26725
sw 40 cm	62720-60	62730-60	26721	26726
sw 42.5 cm	62721-60	62731-60	26721	26726
sw 45 cm	62722-60	62732-60	26722	26727
sw 47.5 cm	62978-60	62984-60	26722	26727
sw 50 cm	62979-60	62985-60	26723	26728
sw 52.5 cm	62980-60	62986-60	26723	26728

Seat depth extension



sw 35-40 cm	63291
sw 42.5-50 cm	63292
sw 52.5 cm	63293

Push handles

Config no

Item no

Push handles

Height adjustable



100 mm	23479	<b>27588-01</b>
75 mm	22018	<b>27338-01</b>
Active		<b>28873</b>
Swing away	43733	<b>24591</b>



Fixed

**28875**

Push handles with drum brake



Height adjustable

100 mm handle and 100 mm grip

23540

**27625-01**

**Bow handle**

Height adjustable

Height & angle adjustable



Height adjustable



Height & angle adjustable

sw 35 cm	27433-01	-
sw 37.5 cm	25578-01	27161
sw 40 cm	25579-01	27162
sw 42.5 cm	25580-01	27163
sw 45 cm	25581-01	27164
sw 47.5 cm	25582-01	27165
sw 50 cm	25730-01	27166
sw 52.5 cm	25583-01	-

**Arm supports and mudguards**

Config no

Item no

**Arm support**



Solid pad

	25 cm	Right	Left		<b>25172-1/-2</b>
	38 cm	Right	Left		<b>25173-1/-2</b>
w fingerscrew for height adjustment	25 cm	Right	Left	23910/23911	<b>25172-1/-2</b>
w fingerscrew for height adjustment	38 cm	Right	Left	23912/23913	<b>25173-1/-2</b>



Removable pad

	25 cm	Right	Left		<b>28780-1/-2</b>
	35 cm	Right	Left		<b>28781-1/-2</b>
w fingerscrew for height adjustment	25 cm	Right	Left	43557/43558	<b>28782-1/-2</b>
w fingerscrew for height adjustment	35 cm	Right	Left	43559/43560	<b>28783-1/-2</b>



Flip-up

	25 cm	Right	Left		<b>28842-1/-2</b>
	38 cm	Right	Left		<b>28843-1/-2</b>

**Mudguard**



	Right		Left		
S					<b>24094-1/-2</b>
M					<b>24095-1/-2</b>
L					<b>24096-1/-2</b>

	20°	22°	24°	25°	26°
H	L				
G	M/L	L			
F	S/M	M/L	L	L	
E	S	M	M/L	L	L
D	S*	S	M	M	M/L
C	S**	S*	S	S/M	M
B	S***	S**	S*	S	S
A	S***	S**	S*	S**	S*

S					<b>24094-1/-2</b>
M					<b>24095-1/-2</b>
L					<b>24096-1/-2</b>

Note: conflict may occur, check manual.

**Anti tips**

Config no

Item no

**Anti tips multi cc59**



applicable for all seat heights

43427

**28457**

**Seat cushions**

**Universal**

Plush cover

Tidy cover (hygiene)



Plush



Tidy

sw 35 cm	62230-60	63130-63
sw 37.5 cm	62231-60	63131-63
sw 40 cm	62232-60	63132-63
sw 42.5 cm	62233-60	63133-63
sw 45 cm	62234-60	63134-63
sw 47.5 cm	62374-60	63135-63
sw 50 cm	62266-60	63136-63
sw 52.5 cm	62420-60	63137-63

**Widening kit**

Config no

Item no

**Widening kit**












sw +1 cm

Not with drum brake

43080

**27533**

## ACCESSORIES

Arm supports		Pieces	Item no
<b>Hemi arm support</b>		for item 25173-1/-2	Right Left 25770-60-1/-2
		8 x 40 cm	for item 25172-1/-2 and 25173-1/-2 2 26873
<b>Wide and soft arm pad</b>		25 cm	for item 25172-1/-2
			padded plush 2 62282-60 padded hygiene 2 62282-63 gel plush 2 26240-60 gel hygiene 2 26240-63
<b>Cover</b>		38 cm	for item 25173-1/-2
			padded plush 2 62283-60 padded hygiene 2 62283-63 gel plush 2 26241-60 gel hygiene 2 26241-63
<b>Side guard</b>		Active (to be mounted), pair	for item 25172-1/-2 and 25173-1/-2 Right Left 28857
		Active	Right Left 24781-1/-2
<b>Side guard pad</b>		Cross, standard shape	reduce the seat width 30mm/pair Right Left 63373-1/-2
		Active	reduce the seat width 30mm/pair Right Left 63372-1/-2
<b>Thigh support</b> Not with 22" in pos C1 Not with 24" in pos C1 or B1		Short 50 cm	for item 25172-1/-2 and 25173-1/-2 Right Left 27391-1/-2
		Long 57 cm	for item 25172-1/-2 and 25173-1/-3 Right Left 27392-1/-2
		Soft pad for thigh support	Right Left 82488-1/-2
Attachments for medical devices		Pieces	Item no
<b>Crutch and cane holder</b>		Universal	1 27327
<b>IV stand</b>		For two IV bags complete with attachment	1 28323
<b>Oxygen holder</b>		For oxygen tubes, maximum load 5 kg	1 28937
Back supports			Item no
<b>Axiom back</b>		Assembly kit DEX 3A Axiom	24046
		Axiom PL	width x height
		36 x 26	KI124045
		36 x 31	KI124046
		36 x 36	KI124047
		36 x 41	KI124048
		36 x 46	KI124049
		41 x 26	KI124050
		41 x 31	KI124051
		41 x 36	KI124052
		41 x 41	KI124053
		41 x 46	KI124054
		46 x 26	KI124057
		46 x 31	KI124058
		46 x 36	KI124059
		46 x 41	KI124060
		46 x 46	KI124061
		51 x 26	KI124064
		51 x 31	KI124065
		51 x 36	KI124066
		51 x 41	KI124067
		51 x 46	KI124068

**Fixations and brackets**

Item no

**Cross brace**


with snaplock for 3A back support

sw 35 cm					28114
sw 37.5 cm					28115
sw 40 cm					28116
sw 42.5 cm					28117
sw 45 cm					28118
sw 47.5 cm					28695
sw 50 cm					28696
sw 52.5 cm					28697

**Head supports**

Item no

**Attachment**


Assembly on 3A back support

sw 35 cm					28125
sw 37.5 cm					28126
sw 40 cm					28127
sw 42.5 cm					28128
sw 45 cm					28129
sw 47.5 cm					28130
sw 50 cm					28131
sw 52.5 cm					28132

**Head support**


Assembly on bow handle

25795



Traditional

26x14 cm incl turtle bar 350 mm

with plush cover					28534-60
with hygiene cover					28534-65



With side support

incl turtle bar 250 mm

with plush cover					28536-60
with hygiene cover					28536-65

**Extension**


Extension head support

for head support

29012

Only extension, not complete head support.

**Memory function**


Memory function

for head support

28564

**Leg and foot supports**

Pieces

Item no

**Cover**


to cover attachment when not using leg support

2

51213

**Foot support extension**


Extension, plastic

for right &amp; left foldable foot support

sw 40-42.5 cm					2	50798
sw 45-47.5 cm					2	50799
sw 50-52.5 cm					2	50800

**Add-on extension**


Extension, metal

for right &amp; left foldable foot support

Small, one piece	sw 35-50 cm					28774
Medium	sw 40-47.5 cm	Right	Left			28775-1/-2
Large	sw 50-52.5 cm	Right	Left			28776-1/-2

Guide to link right &amp; left extension into one piece

28773

**Foot box**


Foot box

Soft padded

Medium	sw 40-47.5 cm					28642
large	sw 50-52.5 cm					28643

**Thigh support**

Not with lockable or elevated leg supports

























Distal thigh support

for fixed leg supports,


Right

Left


28771-1/-2


Cover foot support		For flip up foot plates	padded plush	Right	Left			
			sw 35-40 cm			63188-60-1/-2		
			sw 42.5-50 cm			63189-60-1/-2		
			sw 52.5-60 cm			63192-60-1/-2		
Calf support		For one piece foot plate	padded plush					
			sw 35-37.5 cm			63234-60		
			sw 40-45 cm			63190-60		
			sw 45-50 cm			63191-60		
Calf support	   	Cover, hygiene	for elevated leg support item 26100			2	62762-63	
		Calf strap	adjustable in length					
			sw 35-37.5 cm				25356	
			sw 40-42.5 cm				25357	
			sw 45-47.5 cm				25358	
			sw 50-52.5 cm				25359	
		Padding for calf strap						
			20 cm					26845
			30 cm					26846
			40 cm					26847
Panel		Panel	Plush					
			75 cm ( sw 35-47.5)				63186-60	
			90 cm (sw 50-52.5)					63187-60
<b>Postural supports</b>						Pieces	Item no	
Positioning belt		2 point, plastic push release						
			100 cm				27904	
			128 cm					28544
Wedge	 	2 point, push button buckle						
			170 cm				62190	
		Soft	Plush					pcs
			Short, 15 cm			2		63012-60
	Long, 30 cm			2		63013-60		
Cell foam		Cell foam	Polyethylene				pair	
			90 x 20 x 130 mm				28030	
			110 x 20 x 185 mm					28139
	175 x 25 x 200 mm					28140		
Lateral support	 	Side cushion		Right	Left		63171-60-1/-2	
		Side stop, straight	gives extra support to side cushion				2	28136
Back width extension	 	Width extension, 30 mm	combine with wedges				2	28173
		Widening pad	padded cover with flexible insert					2
Trunk support 3A	  	Swing away & detachable	Adjustable in height, depth and width					
			0 mm (straight profile)	Right	Left		28181-1/-2	
			-15 mm	Right	Left			28180-1/-2
			-45 mm	Right	Left			28141-1/-2
		Detachable	Adjustable in height, depth and width					
			0 mm (straight profile)					28142
	+15 mm					28143		
	-15 mm					28144		
	-45 mm					28145		
Soft	 	Soft	Adjustable in height				2	27395
		Swing away R82 pad	Adjustable in height and depth					
			Medium	Right	Left			27915-1/-2
	Large	Right	Left			27916-1/-2		
Cover for trunk support		Soft cover	Fits item 28142, 28143, 28144, 28145, 28181-1/-2, 28180-1/-2, 28141-1/-2					28647
Attachment for harness		Upper attachment	to be mounted on head support stem 15x15 mm					28186

**Seating** Pieces **Item no**

<b>Seat wedge</b>		<b>Comfort wedge</b>	Provides extra seat contour for seat cushion Universal.	
			sw 35-37.5 cm	<b>80714</b>
			sw 40-45 cm	<b>80715</b>
			sw 47.5-52.5 cm	<b>80779</b>

**Trays** Pieces **Item no**

<b>Tray</b>		<b>Transparent</b>	for 38 cm arm support 25173-1/-2	
			sw 35 cm	<b>24417</b>
			sw 37.5 cm	<b>25187</b>
			sw 40 cm	<b>24418</b>
			sw 42.5 cm	<b>25188</b>
			sw 45 cm	<b>24419</b>
			sw 47.5 cm	<b>25731</b>
			sw 50 cm	<b>24601</b>
sw 52.5 cm	<b>26020</b>			

<b>Hemi tray</b>		<b>Transparent</b>	for 25 cm arm support 25172-1/-2	
			sw 35 cm	<b>26556</b>
			sw 37.5 cm	<b>26557</b>
			sw 40 cm	<b>26558</b>
			sw 42.5 cm	<b>26559</b>
			sw 45 cm	<b>26560</b>
			sw 47.5 cm	<b>26561</b>


<b>Transparent</b>	for 38 cm arm support 25173-1/-2	
	sw 35 cm	<b>26563</b>
	sw 37.5 cm	<b>26564</b>
	sw 40 cm	<b>26565</b>
	sw 42.5 cm	<b>26566</b>
	sw 45 cm	<b>26567</b>
	sw 47.5 cm	<b>26568</b>
	sw 50 cm	<b>26569</b>

<b>Anti glide</b>	for hemi tray	<b>81633</b>
-------------------	---------------	--------------

<b>Half tray</b>		<b>Foldable</b>	for arm support 25172-1/-2 & 25173-1/-2		<b>27352-1/-2</b>
			Right	Left	

<b>Attachment for tray</b>		<b>Round tubes</b>	fits tubes up to ø 18mm		<b>27956-1/-2</b>
			Right	Left	

**Wheels** Pieces **Item no**

<b>Spokeguard</b>		<b>Transparent</b>	for C6 and Maxx	
			20"	<b>2</b> <b>28934</b>
			22"	<b>2</b> <b>28935</b>
			24"	<b>2</b> <b>28936</b>

<b>Transparent</b>	for C5 DB	
	20"	<b>2</b> <b>28816-28</b>
	22"	<b>2</b> <b>28817-28</b>
	24"	<b>2</b> <b>28818-28</b>

**Additional information to order:**

---



---



---



---



---

**Etac Supply Center AB**  
 Phone: +46 371 58 73 30  
 Mail: info@etac.com  
 Website: www.etac.com



Creating Possibilities