

# Xact



Therapist	User	
Phone no	Order date	
Purchase Order no	Customer Ref	
Base price	€XXX	Standard option included in base price
Colour	01 Black	

Item no	Seat width (sw)	Frame length	Seat depth
13250101	35 cm	Standard	38-46 cm
13250102	35 cm	Short	35-43 cm
13250103	37.5 cm	Standard	38-46 cm
13250104	37.5 cm	Short	35-43 cm
13250105	40 cm	Standard	38-46 cm
13250106	40 cm	Short	35-43 cm
13250107	42.5 cm	Standard	38-46 cm
13250108	42.5 cm	Short	35-43 cm
13250109	45 cm	Standard	38-46 cm
13250110	45 cm	Short	35-43 cm
13250111	47.5 cm	Standard	38-46 cm
13250112	47.5 cm	Short	35-43 cm
13250113	50 cm	Standard	38-46 cm
13250114	50 cm	Short	35-43 cm

## OPTIONS

### Seat height Config no    Item no

#### Seat height tabel

Seat angle is always 4°

- Decide seat height and check the box to the right of the row
- Choose castor size
- Choose rear wheel size

cm	cm	4°	5°	6°	20°	22°	24°	25°	Model	
37.5	40.5	1XS			E				1	43630
38.8	41.8	2XS			D				1	43631
40.0	43.0	3XS	1S		C	E			2	43632
41.3	44.3		2S		B	D			2	43633
42.5	45.5		3S	2S	A	C	E		3	43634
43.8	46.8		4S	3S		B	D	E	3	43635
45.0	48.0		3S+2S	4S		A	C	D	3	43636
46.3	49.3		4S+2S	3S+2S		B	C		4	43637
47.5	50.5			4S+2S		A	B		4	43638

### Rear wheels Config no    Item no

#### Maxx




	Hard anodised aluminium	Titanium	Plastic coated
22"			
High pressure, Speedlite	28956	24087	28957
High pressure, Trailblazer	24052-1/-2 (43703)	24115-1/-2 (43776)	24053-1/-2 (43704)
Solid	24123	24124	24125
24"			
High pressure, Speedlite	28950	24088	28951
High pressure, Trailblazer	28952-1/-2 (43699)	24116-1/-2 (43775)	28953-1/-2 (43700)
Solid	24126	24127	24128
25"			
High pressure, Trailblazer	28964-1/-2 (43701)	24117-1/-2 (43754)	28965-1/-2 (43702)
Solid	24129	24130	24131

#### C6



	Aluminium	Stainless steel	Plastic coated
20°			
Solid	83885	-	-
High pressure, Speedlite	83908	-	83884
22°			
Solid	83882	-	83883
24°			
Solid	83876	-	83878

#### C5 Drum brake



	Aluminium	Stainless steel	Plastic coated
20° DB			
Solid	81337B (43682)	-	-
22° DB			
Solid	81339B (43707)	81341B (43732)	-
High pressure	81342B (43681)	-	-

**C5 Drum brake**




24" DB	Aluminium	Stainless steel	Plastic coated
Solid	81346B (43706)	81348B (43731)	81347B (43755)
High pressure	28625B (43680)	-	-

**QR axle**

Axle QR D12,7x99	for C6 and Maxx	<b>83492</b>
Axle QR Drum brake	for C5 Drum brake	<b>81325</b>

**One arm drive**



**Drive unit QR**

20"	84101
22"	84102
24"	84103


**Rear axle attachment**

	Config no	
Standard D-E	43678, 43762 to 43767	<b>24004</b>
High A-C	43679, 43768 to 43773	<b>24019</b>

**Rear axle attachment - One arm drive QR**


	Config no	
Standard D-E	43734 to 43740	<b>24004</b>
High C (A+B not applicable)	43741 to 43747	<b>24019</b>

\*including assembly kit, crossbar and telescopic shaft



	Assembly kit	Crossbar one arm drive QR	Telescopic shaft
sw 35 cm	24071	24072	24171
sw 37.5 cm	24071	24073	24171
sw 40 cm	24071	24074	24172
sw 42.5 cm	24071	24075	24172
sw 45 cm	24071	24076	24173
sw 47.5 cm	24071	24077	24174
sw 50 cm	24071	24078	24174

**Castors**



**Castors**

		Config no	Item no
4" 100 x 34	max user weight 100 kg	83344	<b>28919</b>
5" 125 x 28	lightweight	83867	<b>27299</b>
5" 125 x 34		83345	<b>28921</b>
6" 150 x 34		83346	<b>28922</b>

**Front forks inkl axle**

	Seat height and castor size determines the required front fork.	
XS		<b>28923</b>
S		<b>28924</b>
S+25		<b>28925</b>

**Anti-tips**



**Without anti-tips**

Without anti-tips	22166	
With anti-tips		
Model 1	43641	<b>24163-1/-2</b>
Model 2	43642	<b>24164-1/-2</b>
Model 3	43643	<b>24165-1/-2</b>
Model 4	43644	<b>24166-1/-2</b>

**Brakes**

**Scissor brake**



	Config no	Item no
not with solid tyres	Right Left 289150-1/-2	<b>28915-1/-2</b>

**Knee lever brake**









Black angled handle	Right Left 289170-1/-2	<b>28917-1/-2</b>
Red angled handle	Right Left 289160-1/-3	<b>28916-1/-2</b>


**One hand brake**





Single hand brake	Right Left	<b>28918-1/-2</b>
-------------------	------------	-------------------

Leg supports (hangers)		Config no	Item no
		Without leg supports	Right Left 22148/22149
Fixed angle		Active, angle 80°, lightweight alu lower leg length 31-49 cm	Right Left <b>28888</b>
		Narrow short, angle 80° steel lower leg length 23-33 cm	Right Left <b>25571-01</b>
Elevated		Mod 3 lower leg length 31-49 cm	Right Left <b>26100</b>
Amputee		Amputee support	
		Plush	<b>25516-60</b>
		Hygiene	<b>25516-63</b>

Foot supports		Config no	Item no		
		Without foot supports	Right Left 22150/22151		
Divided		Flip up attachment tube:			
		short	Right Left		
		std +8 cm	Right Left		
One piece		Flip up			
					
		<b>Divided short</b>	<b>Divided std +8</b>		
		sw 35 cm	25653-1/2	25641-1/2	25936
		sw 37.5 cm	25654-1/2	25642-1/2	25936
		sw 40 cm	25655-1/2	25643-1/2	25937
		sw 42.5 cm	25656-1/2	25644-1/2	25937
		sw 45 cm	25657-1/2	25645-1/2	25937
		sw 47.5 cm	25738-1/2	25736-1/2	25937
		sw 50 cm	25658-1/2	25646-1/2	25937

Push handles		Config no	Item no
Push handles		Without push handles	22618
		Height adjustable	
		Active	43688 <b>28873</b>
		Swing away	43733 <b>24591</b>
		Fixed	43689 <b>28875</b>
Push handles with drumbrake		Height adjustable 100 mm handle and 100 mm grip	23540 <b>27625-01</b>
Bow handle		Height adjustable	
		sw 35 cm	<b>27439-01</b>
		Height & angle adjustable	
		sw 37.5 cm	<b>27161</b>
		sw 40 cm	<b>27162</b>
		sw 42.5 cm	<b>27163</b>
sw 45 cm	<b>27164</b>		
sw 47.5 cm	<b>27165</b>		
sw 50 cm	<b>27166</b>		

Arm supports and mudguards		Config no	Item no
		Without arm supports	Right Left 22178/22179
Arm support		Active	
		25 cm	Right Left 288550-1/-2 <b>28855-1/-2</b>
		38 cm	Right Left 288560-1/-2 <b>28856-1/-2</b>
Side guard		Active profile 200	Right Left 245920-1/-2 <b>24592-1/-2</b>

**Back supports & Covers**

Config no

Item no

**3A back support**



inkl tension adjustable back

back height 32 - 45 cm (extension +6 cm available as accessory)	Pre-set back height	<b>36 cm</b>
recline -14° - +11°	Pre-set recline angle	<b>- 4°</b>
lumbar angle -16° - +16°	Pre-set lumbar angle	<b>0°</b>
	sw 35-40 cm	<b>62953-60</b>
	sw 42.50-50 cm	<b>62954-60</b>

**Covers**

Airmesh

Plush

Plush high (for back height over 46-51 cm)

Without cover



	<b>Airmesh</b>	<b>Plush</b>	<b>Plush high</b>
sw 35 cm	63204-68	62941-60	63020-60
sw 37.5 cm	63205-68	62942-60	63021-60
sw 40 cm	63206-68	62943-60	63022-60
sw 42.5 cm	63207-68	62944-60	63023-60
sw 45 cm	63208-68	62945-60	63024-60
sw 47.5 cm	63209-68	62946-60	63025-60
sw 50 cm	63210-68	62947-60	63026-60

**Seating**

Config no

Item no

Seat upholstery

sw 35 cm	<b>63362-60</b>
sw 37.5 cm	<b>63363-60</b>
sw 40 cm	<b>63364-60</b>
sw 42.5 cm	<b>63365-60</b>
sw 45 cm	<b>63366-60</b>
sw 47.5 cm	<b>63367-60</b>
sw 50 cm	<b>63368-60</b>

**ACCESSORIES**

**Arm supports**

Pieces

Item no

**Mudguard**



Mudguard profile 155

	20°	22°	24°	25°
	E	S	M	L
	D	S	S	M
	C	S*	S	M
	B	S**	S*	S
	A	S**	S**	S

Right Left

S	<b>24091-1/-2</b>
M	<b>24092-1/-2</b>
L	<b>24093-1/-2</b>

Note: conflict may occur, check manual.

**Sideguard pad**



Active	reduce the seat width 30mm/pair	Right	Left	<b>63372-1/-2</b>
--------	---------------------------------	-------	------	-------------------

**Hemi arm support**



Hemi, horizontal rotation	for item 28856-1/-2	Right	Left	<b>25770-60-1/-2</b>
---------------------------	---------------------	-------	------	----------------------

**Wide and soft arm pad**



8 x 40 cm	for item 28855-1/-2 and 28856-1/-2		<b>2</b>	<b>26873</b>
-----------	------------------------------------	--	----------	--------------

**Cover**



25 cm	for item 28855-1/-2			
	padded plush		<b>2</b>	<b>62282-60</b>
	padded hygiene		<b>2</b>	<b>62282-63</b>
	gel plush		<b>2</b>	<b>26240-60</b>
	gel hygiene		<b>2</b>	<b>26240-63</b>
38 cm	for item 28856-1/-2			
	padded plush		<b>2</b>	<b>62283-60</b>
	padded hygiene		<b>2</b>	<b>62283-63</b>
	gel plush		<b>2</b>	<b>26241-60</b>
	gel hygiene		<b>2</b>	<b>26241-63</b>

**Attachments for medical devices**

Pieces

Item no

**Crutch and cane holder**



				<b>28932</b>
--	--	--	--	--------------

**IV stand**



For two IV bags complete with attachment			<b>1</b>	<b>28323</b>
--	--	--	----------	--------------

**Back supports** Item no

**Back extension 3A**



Adds 6 cm to the back, for 38-51 cm adjustment

sw 35-40 cm	28658
sw 42.5-50 cm	28659

**Axiom back**



Assembly kit DEX 3A Axiom

24046

Axiom PL width x height

36 x 26	KI124045
36 x 31	KI124046
36 x 36	KI124047
36 x 41	KI124048
36 x 46	KI124049
41 x 26	KI124050
41 x 31	KI124051
41 x 36	KI124052
41 x 41	KI124053
41 x 46	KI124054
46 x 26	KI124057
46 x 31	KI124058
46 x 36	KI124059
46 x 41	KI124060
46 x 46	KI124061
51 x 26	KI124064
51 x 31	KI124065
51 x 36	KI124066
51 x 41	KI124067
51 x 46	KI124068

**Covers**

Enter item number from below:

	Comfort	Comfort high	Tidy/hygiene
sw 35 cm	62956-60	63033-60	63020-60
sw 37.5 cm	62957-60	63034-60	63021-60
sw 40 cm	62958-60	63035-60	63022-60
sw 42.5 cm	62959-60	63036-60	63022-61
sw 45 cm	62960-60	63037-60	63022-62
sw 47.5 cm	62961-60	63038-60	63025-60
sw 50 cm	62962-60	63039-60	63026-60

**Brakes** Item no

**Extension**



Brake extension for knee lever brake

Right Left

27530-1/-2

**Fixations and brackets** Item no

**Cross brace**



with snaplock for 3A back support

sw 35 cm	28114
sw 37.5 cm	28115
sw 40 cm	28116
sw 42.5 cm	28117
sw 45 cm	28118
sw 47.5 cm	28695
sw 50 cm	28696

**Head supports** Item no

**Attachment**



Assembly on 3A back support

sw 35 cm	28125
sw 37.5 cm	28126
sw 40 cm	28127
sw 42.5 cm	28128
sw 45 cm	28129
sw 47.5 cm	28130
sw 50 cm	28131

**Head support**



Traditional

26x14 cm incl turtle bar 350 mm	
with plush cover	28534-60
with hygiene cover	28534-65

With side support

incl turtle bar 250 mm	
with plush cover	28536-60
with hygiene cover	28536-65

**Memory function**



Memory function

for head support	28564
------------------	-------

**Leg and foot supports**

Pieces

Item no

**Fixed angle**



Swing away, detachable	lower leg length 31-49 cm			
	Narrow angle 80°	Right	Left	25570-01
	Wide angle 70°	Right	Left	25574-01



Lockable, swing away, detachable	lower leg length 31-49 cm			
	Narrow angle 80°	Right	Left	25573-01
	Narrow angle 80° short (23-33 cm)	Right	Left	25572-01
	Wide angle 70°	Right	Left	25575-01

Mod 1	lower leg length 31-49 cm	Right	Left	25660-01
add on:	Attachment for calf support			24544-01
	Attachment for calf support flip-up			27322-01
	Calf support, plush item:			
	sw 35-37.5 cm			62105-60
	sw 40-42.5 cm			62104-60
	sw 45-47.5 cm			62484-60
	sw 50 cm			62451-60

**Foot support extension**



Extension, plastic	for right & left foldable foot support			
	sw 40-42.5 cm			2 50798
	sw 45-47.5 cm			2 50799
	sw 50-52.5 cm			2 50800

**Add-on extension**



Extension, metal	for right & left foldable foot support			
Small, one piece	sw 35-50 cm			28774
Medium	sw 40-47.5 cm	Right	Left	28775-1/-2
Large	sw 50 cm	Right	Left	28776-1/-2

Guide to link right & left extension into one piece				28773
---	--	--	--	-------



Side protection	for extension, metal			2 28772
-----------------	----------------------	--	--	---------

**Foot box**



Foot box	Soft padded			
Medium	sw 40-47.5 cm			28642
large	sw 50 cm			28643

**Thigh support**

Not with lockable or elevated leg supports



Distal thigh support	for fixed leg supports (not active)	Right	Left	28777-1/-2
----------------------	-------------------------------------	-------	------	------------

**Heel strap**



Heel strap	for divided foot supports	Right	Left	24458
------------	---------------------------	-------	------	-------

**Cover foot support**



For flip up foot plates	padded plush	Right	Left	
	sw 35-40 cm			63188-60-1/-2
	sw 42.5-50 cm			63189-60-1/-2



For one piece foot plate	padded plush			
	sw 35-37.5 cm			63234-60
	sw 40-45 cm			63190-60
	sw 45-50 cm			63191-60

**Calf support**



Cover, hygiene	for elevated leg support item 26100			2 62762-63
----------------	-------------------------------------	--	--	------------

















Calf strap	adjustable in length			
	sw 35-37.5 cm			25356
	sw 40-42.5 cm			25357
	sw 45-47.5 cm			25358
	sw 50 cm			25359







Padding for calf strap				
	20 cm			26845
	30 cm			26846
	40 cm			26847



Panel	Plush			
	75 cm ( sw 35-47.5)			63186-60
	90 cm (sw 50)			63187-60

Postural supports		Pieces	Item no			
<b>Attachment</b>		Attachment for assembly of 2- and 4-point positioning belts, set of 2	28933			
<b>Positioning belt</b>		2 point, plastic push release				
		100 cm	27904			
		128 cm	28544			
<b>Wedge</b>	 	Soft	Plush			
		Short, 15 cm	2	63012-60		
		Long, 30 cm	2	63013-60		
		Cell foam	Polyethylene	pair		
		90 x 20 x 130 mm		28030		
		110 x 20 x 185 mm		28139		
		175 x 25 x 200 mm		28140		
<b>Lateral support</b>	 	Side cushion	Right	Left	63171-60-1/2	
		Side stop, straight	gives extra support to side cushion		2	28136
<b>Back width extension</b>	 	Width extension, 30 mm	combine with wedges		2	28173
		Widening pad	padded cover with flexible insert		2	28637
<b>Trunk support 3A</b>	   	Swing away & detachable	Adjustable in height, depth and width			
		0 mm (straight profile)	Right	Left		28181-1/-2
		-15 mm	Right	Left		28180-1/-2
		-45 mm	Right	Left		28141-1/-2
		Detachable	Adjustable in height, depth and width			
		0 mm (straight profile)				28142
		+15 mm				28143
		-15 mm				28144
		-45 mm				28145
		Soft	Adjustable in height		2	27395
		Swing away R82 pad	Adjustable in height and depth			
		Medium	Right	Left		27915-1/-2
		Large	Right	Left		27916-1/-2
<b>Cover for trunk support</b>		Soft cover	Fits item 28142, 28143, 28144, 28145, 28181-1/-2, 28180-1/-2, 28141-1/-2			28647
<b>Attachment for harness</b>		Upper attachment	to be mounted on head support stem 15x15 mm			28186

Seating		Pieces	Item no			
<b>Universal seat cushion</b>		Plush cover				
		<b>Plush</b>				
		<b>Tidy</b>				
		sw 35 cm	62230-60	63130-63		
		sw 37.5 cm	62231-60	63131-63		
		sw 40 cm	62232-60	63132-63		
		sw 42.5 cm	62233-60	63133-63		
		sw 45 cm	62234-60	63134-63		
		sw 47.5 cm	62374-60	63135-63		
		sw 50 cm	62266-60	63136-63		
<b>Seat wedge</b>		Comfort wedge	Provides extra seat contour for seat cushion Universal.			
		sw 35-37.5 cm				80714
		sw 40-45 cm				80715
		sw 47.5-50 cm				80779
<b>Widening kit</b>		sw + 1.2 cm	Not with drum brake			24079

**Trays**

Pieces

Item no

**Tray**



Transparent	for 38 cm arm support 28856-1/-2		
	sw 35 cm		<b>24417</b>
	sw 37.5 cm		<b>25187</b>
	sw 40 cm		<b>24418</b>
	sw 42.5 cm		<b>25188</b>
	sw 45 cm		<b>24419</b>
	sw 47.5 cm		<b>25731</b>
	sw 50 cm		<b>24601</b>

**Hemi tray**



Transparent	for 25 cm arm support 28856-1/-2		
	sw 35 cm		<b>26556</b>
	sw 37.5 cm		<b>26557</b>
	sw 40 cm		<b>26558</b>
	sw 42.5 cm		<b>26559</b>
	sw 45 cm		<b>26560</b>
	sw 47.5 cm		<b>26561</b>
	sw 50 cm		<b>26562</b>

Transparent	for 38 cm arm support 28856-1/-2		
	sw 35 cm		<b>26563</b>
	sw 37.5 cm		<b>26564</b>
	sw 40 cm		<b>26565</b>
	sw 42.5 cm		<b>26566</b>
	sw 45 cm		<b>26567</b>
	sw 47.5 cm		<b>26568</b>
	sw 50 cm		<b>26569</b>

Anti glide	for hemi tray		
			<b>81633</b>

**Half tray**



Foldable	for arm support 28855-1/-2 & 28856-1/-2	Right	Left	
				<b>27352-1/-2</b>

**Attachment for tray**

Round tubes	fits tubes up to ø 18mm	Right	Left	
				<b>27956-1/-2</b>

**Wheels**

Pieces

Item no

**Spokeguard**



Transparent	for C6 & Maxx		
	20"	2	<b>28934</b>
	22"	2	<b>28935</b>
	24"	2	<b>28936</b>



Transparent	for C5 DB		
	20"	2	<b>28816-28</b>
	22"	2	<b>28817-28</b>
	24"	2	<b>28818-28</b>

**Additional information to order:**

---



---



---



---



---

Etac Supply Center AB  
 Phone: +46 371 58 73 30  
 Mail: orderexport@etac.com  
 Website: www.etac.com

