

Xact

Therapist		User	
Phone no		Order date	
Purchase Order no		Customer Ref	

☒ Standard option

 Colour 01 Black ☒


Item no	Seat width (sw)	Frame length	Seat depth	
13250101	35 cm	Standard	38-46 cm	☐
13250102	35 cm	Short	35-43 cm	☐
13250103	37.5 cm	Standard	38-46 cm	☐
13250104	37.5 cm	Short	35-43 cm	☐
13250105	40 cm	Standard	38-46 cm	☐
13250106	40 cm	Short	35-43 cm	☐
13250107	42.5 cm	Standard	38-46 cm	☐
13250108	42.5 cm	Short	35-43 cm	☐
13250109	45 cm	Standard	38-46 cm	☐
13250110	45 cm	Short	35-43 cm	☐
13250111	47.5 cm	Standard	38-46 cm	☐
13250112	47.5 cm	Short	35-43 cm	☐
13250113	50 cm	Standard	38-46 cm	☐
13250114	50 cm	Short	35-43 cm	☐

OPTIONS

Seat height	Config no	Item no
-------------	-----------	---------

Seat height tabel

Seat angle is always 4°

- Decide seat height and check the box to the right of the row
- Choose castor size
- Choose rear wheel size

cm	cm	4"	5"	6"	20"	22"	24"	25"	Model		
37.5	40.5	☐ 1XS			☐ E				1	43630	☐
38.8	41.8	☐ 2XS			☐ D				1	43631	☐
40.0	43.0	☐ 3XS	☐ 1S		☐ C	☐ E			2	43632	☐
41.3	44.3		☐ 2S		☐ B	☐ D			2	43633	☐
42.5	45.5		☒ 3S	☐ 2S	☐ A	☐ C	☒ E		3	43634	☒
43.8	46.8		☐ 4S	☐ 3S		☐ B	☐ D	☐ E	3	43635	☐
45.0	48.0		☐ 3S+2S	☐ 4S		☐ A	☐ C	☐ D	3	43636	☐
46.3	49.3		☐ 4S+2S	☐ 3S+2S			☐ B	☐ C	4	43637	☐
47.5	50.5			☐ 4S+2S			☐ A	☐ B	4	43638	☐

Rear wheels	Config no	Item no
-------------	-----------	---------

Maxx



22"	Hard anodised aluminium	Stainless steel	Plastic coated
High pressure, Speedlite	☐ 28956	-	☐ 28957
High pressure, Trailblazer	☐ 24052-1/-2 (43703)	-	☐ 24053-1/-2 (43704)

24"			
High pressure, Speedlite	☒ 28950	-	☐ 28951
High pressure, Trailblazer	☐ 28952-1/-2 (43699)	-	☐ 28953-1/-2 (43700)

25"			
High pressure, Trailblazer	☐ 28964-1/-2 (43701)	-	☐ 28965-1/-2 (43702)

C6





20"	Aluminium	Stainless steel	Plastic coated
Solid	☐ 83885	-	☐ 83911
High pressure, Speedlite	☐ 83908	-	☐ 83884

22"			
Solid	☐ 83882	-	☐ 83883

24"			
Solid	☐ 83876	-	☐ 83878

C5 Drum brake



20" DB	Aluminium	Stainless steel	Plastic coated
Solid	☐ 81337 (43682)	-	-


22" DB			
Solid	☐ 81339 (43707)	☐ 81341 (43732)	
High pressure	☐ 81342 (43681)		

24" DB			
Solid	☐ 81346 (43706)	☐ 81348 (43731)	☐ 81347 (43755)
High pressure	☐ 28625 (43680)		

QR axle			
Axle QR D12.7x99	for C6 and Maxx		83492 ☒
Axle QR Drum brake	for C5 Drum brake		81325

Rear axle attachment			
Standard D-E		43678	24004
High A-C		43679	24019

Castors



Castors			
4" 100 x 34	max user weight 100 kg	83344	28919 ☐
5" 125 x 28	lightweight	83867	27299 ☒
5" 125 x 34		83345	28921 ☐
6" 150 x 34		83346	28922 ☐

Front forks inkl axle	Seat height and castor size determines the required front fork.		
XS			28923
S			28924
S+25			28925

Anti-tips



Anti-tips			
Without anti-tips		22166	☐
With anti-tips			☒
Model 1		43641	24163-1/-2
Model 2		43642	24164-1/-2
Model 3		43643	24165-1/-2
Model 4		43644	24166-1/-2

Brakes

Scissor brake



not with solid tyres	☒ Right	☒ Left	289150-1/-2	28915-1/-2 ☒
----------------------	---	--	-------------	---

Knee lever brake



Black angled handle	☐ Right	☐ Left	289170-1/-2	28917-1/-2 ☐
Red angled handle	☐ Right	☐ Left	289160-1/-3	28916-1/-2 ☐

One hand brake



Extended foldable handle	☐ Right	☐ Left		28918-1/-2 ☐
--------------------------	--------------------------------	-------------------------------	--	--


Leg supports (hangers)

Fixed angle



Without leg supports	☐ Right	☐ Left	22148/22149	☐
Active, angle 80°, lightweight alu	lower leg length 31-49 cm	☒ Right	☒ Left	28888 ☒
Narrow short, angle 80° steel	lower leg length 23-33 cm	☐ Right	☐ Left	25571-01 ☐

Elevated



Mod 3	lower leg length 31-49 cm	☐ Right	☐ Left	26100 ☐
-------	---------------------------	--------------------------------	-------------------------------	---------------------------------------

Amputee



Amputee support

Plush	25516-80	☐
Hygiene	25516-63	☐

Foot supports

Item no

Without foot supports	☐ Right	☐ Left	22150/22151	☐
-----------------------	--------------------------------	-------------------------------	-------------	--------------------------

Divided

Flip up	attachment tube:			
short	☒ Right	☒ Left		☒
std +8 cm	☐ Right	☐ Left		☐

One piece

Flip up				☐
---------	--	--	--	--------------------------



Divided short



Divided std +8



One-piece

sw 35 cm	25653-1/2	25641-1/2	25936
sw 37.5 cm	25654-1/2	25642-1/2	25936
sw 40 cm	25655-1/2	25643-1/2	25937
sw 42.5 cm	25656-1/2	25644-1/2	25937
sw 45 cm	25657-1/2	25645-1/2	25937
sw 47.5 cm	25738-1/2	25736-1/2	25937
sw 50 cm	25658-1/2	25646-1/2	25937

Push handles

Config no

Item no

Push handles



Height adjustable			
Active	43688	28873	☐
Swing away	43733	24591	☐
Fixed	43689	28875	☐

Push handles with drumbrake



Height adjustable	100 mm handle and 100 mm grip	23540	27625-01	☐
-------------------	-------------------------------	-------	----------	--------------------------

Bow handle



Height adjustable	sw 35 cm	27433-01	☐
Height & angle adjustable	sw 37.5 cm	27161	☐
	sw 40 cm	27162	☐
	sw 42.5 cm	27163	☐
	sw 45 cm	27164	☐
	sw 47.5 cm	27165	☐
	sw 50 cm	27166	☐

Arm supports and mudguards

Config no

Item no


Arm support



Without arm supports	☐ Right	☐ Left	22178/22179	☐
----------------------	--------------------------------	-------------------------------	-------------	--------------------------

Active					
25 cm	☒ Right	☒ Left	288550-1/2	28855-1/2	☒
38 cm	☐ Right	☐ Left	288560-1/2	28856-1/2	☐

Side guard



Active profile 200	☐ Right	☐ Left	245920-1/2	24592-1/2	☐
--------------------	--------------------------------	-------------------------------	------------	-----------	--------------------------

Back supports & Covers

Config no

Item no

3A back support



inkl tension adjustable back				☒
back height 32 - 45 cm (extension +6 cm available as accessory)	Pre-set back height	36 cm		
recline -14° - +11°	Pre-set recline angle	- 4°		
lumbar angle -16° - +16°	Pre-set lumbar angle	0°		
sw 35-40 cm		62953-80		
sw 42.50-50 cm		62954-80		

Covers

Airmesh	☒		
Plush	☐		
Plush high (for back height over 46-51 cm)	☐		
Without cover	☐		
			
	Airmesh	Plush	Plush high
sw 35 cm	63204-68	62941-60	63020-60
sw 37.5 cm	63205-68	62942-60	63021-60
sw 40 cm	63206-68	62943-60	63022-60
sw 42.5 cm	63207-68	62944-60	63023-60
sw 45 cm	63208-68	62945-60	63024-60
sw 47.5 cm	63209-68	62946-60	63025-60
sw 50 cm	63210-68	62947-60	63026-60

Seating

Seat upholstery		☒
sw 35 cm	63362-60	
sw 37.5 cm	63363-60	
sw 40 cm	63364-60	
sw 42.5 cm	63365-60	
sw 45 cm	63366-60	
sw 47.5 cm	63367-60	
sw 50 cm	63368-60	

ACCESSORIES

Arm supports

Mudguard



Mudguard profile 155		☐ Right	☐ Left	
S				24091-1/-2 ☐
M				24092-1/-2 ☐
L				24093-1/-2 ☐



	20°	22°	24°	25°
E	S	M	M	L
D	S	S	M	L
C	S*	S	M	M
B	S**	S*	S	M
A	S**	S**	S	-

Sideguard pad



Active	reduce the seat width 30mm/pair	☐ Right	☐ Left	63372-1/-2 ☐
--------	---------------------------------	--------------------------------	-------------------------------	-------------------------------------

Hemi arm support



Hemi, horizontal rotation	for item 28855-1/-2	☐ Right	☐ Left	25770-60-1/-2 ☐
---------------------------	---------------------	--------------------------------	-------------------------------	--

Wide and soft arm pad



8 x 40 cm	for item 28855-1/-2 and 28856-1/-2	2	26873 ☐
-----------	------------------------------------	---	--------------------------------

Cover



25 cm	for item 28855-1/-2		
padded plush		2	62282-60 ☐
padded hygiene		2	62282-63 ☐
gel plush		2	26240-60 ☐
gel hygiene		2	26240-63 ☐
38 cm	for item 28856-1/-2		
padded plush		2	62283-60 ☐
padded hygiene		2	62283-63 ☐
gel plush		2	26241-60 ☐
gel hygiene		2	26241-63 ☐

Fixed angle



Swing away, detachable	lower leg length 31-49 cm			
	Narrow angle 80°	☒ Right	☒ Left	25570-01 ☐

Lockable, swing away, detachable	lower leg length 31-49 cm			
	Narrow angle 80°	☐ Right	☐ Left	25573-01 ☐
	Narrow angle 80° short (23-33 cm)	☐ Right	☐ Left	25572-01 ☐



Mod 1	lower leg length 31-49 cm	☐ Right	☐ Left	25660-01 ☐
add on:	Attachment for calf support			24544-01 ☐
	Attachment for calf support flip-up			27322-01 ☐
	Calf support, plush item:			
	sw 35-37.5 cm			62105-60 ☐
	sw 40-42.5 cm			62104-60 ☐
	sw 45-47.5 cm			62484-60 ☐
	sw 50 cm			62451-60 ☐

Foot support extension



Extension, plastic	for right & left foldable foot support			
	sw 40-42.5 cm		2	50798 ☐
	sw 45-47.5 cm		2	50799 ☐
	sw 50-52.5 cm		2	50800 ☐

Add-on extension



Extension, metal	for right & left foldable foot support			
Small, one piece	sw 35-50 cm			28774 ☐
Medium	sw 40-47.5 cm	☐ Right	☐ Left	28775-1/-2 ☐
Large	sw 50 cm	☐ Right	☐ Left	28776-1/-2 ☐

Guide to link right & left extension into one piece				28773 ☐
---	--	--	--	--------------------------------



Side protection	for extension, metal		2	28772 ☐
-----------------	----------------------	--	---	--------------------------------

Foot box



Foot box	Soft padded			
Medium	sw 40-47.5 cm			28842 ☐
large	sw 50 cm			28843 ☐

Thigh support

Not with lockable or elevated leg supports



Distal thigh support	for fixed leg supports (not active)	☐ Right	☐ Left	28777-1/-2 ☐
----------------------	-------------------------------------	--------------------------------	-------------------------------	-------------------------------------

Heel strap



Heel strap	for divided foot supports	☐ Right	☐ Left	24458 ☐
------------	---------------------------	--------------------------------	-------------------------------	--------------------------------

Cover foot support



For flip up foot plates	padded plush	☐ Right	☐ Left	
	sw 35-40 cm			63188-60-1/-2 ☐
	sw 42.5-50 cm			63189-60-1/-2 ☐



For one piece foot plate	padded plush			
	sw 35-37.5 cm			63234-60 ☐
	sw 40-45 cm			63190-60 ☐
	sw 45-50 cm			63191-60 ☐

Calf support



Cover, hygiene	for elevated leg support item 26100		2	62762-63 ☐
----------------	-------------------------------------	--	---	-----------------------------------

Calf strap	adjustable in length			
	sw 35-37.5 cm			25356 ☐
	sw 40-42.5 cm			25357 ☐
	sw 45-47.5 cm			25358 ☐
	sw 50 cm			25359 ☐

Padding for calf strap				
	20 cm			26845 ☐
	30 cm			26846 ☐
	40 cm			26847 ☐




Panel	Plush
	75 cm (sw 35-47.5) 63186-60 ☐
	90 cm (sw 50) 63187-60 ☐

Postural supports

Pieces

Item no

Attachment



Attachment for assembly of 2- and 4-point positioning belts, set of 2	28933 ☐
---	--------------------------------

Positioning belt




2 point, plastic push release		
	100 cm	27904 ☐
	128 cm	28544 ☐

Wedge



Soft	Plush		
	Short, 15 cm	2	63012-60 ☐
	Long, 30 cm	2	63013-60 ☐



Cell foam	Polyethylene		
	90 x 20 x 130 mm	28030	☐
	110 x 20 x 185 mm	28139	☐
	175 x 25 x 200 mm	28140	☐

Lateral support



Side cushion	☐ Right ☐ Left	63171-60-1/-2 ☐
--------------	--	--

Side stop, straight	gives extra support to side cushion	2	28136 ☐
---------------------	-------------------------------------	---	--------------------------------

Back width extension



Width extension, 30 mm	combine with wedges	2	28173 ☐
------------------------	---------------------	---	--------------------------------

Pressure distributor	padded cover with flexible insert	2	28637 ☐
----------------------	-----------------------------------	---	--------------------------------

Trunk support 3A



Swing away & detachable	Adjustable in height, depth and width					
0 mm (straight profile)		Right		Left	28181-1/-2	☐
-15 mm		Right		Left	28180-1/-2	☐
-45 mm		Right		Left	28141-1/-2	☐



Detachable	Adjustable in height, depth and width	
	0 mm (straight profile)	28142 ☐
	+15 mm	28143 ☐
	-15 mm	28144 ☐
	-45 mm	28145 ☐



Soft	Adjustable in height	2	27395 ☐
------	----------------------	---	--------------------------------

Swing away R82 pad	Adjustable in height and depth					
Medium	☐	Right	☐	Left	27915-1/-2	☐
Large	☐	Right	☐	Left	27916-1/-2	☐

Cover for trunk support



Soft cover	Fits item 28142, 28143, 28144, 28145, 28181-1/-2, 28180-1/-2, 28141-1/-2	28647 ☐
------------	--	--------------------------------

Attachment for harness



Upper attachment	to be mounted on head support stem 15x15 mm	28186 ☐
------------------	---	--------------------------------

Universal seat cushion

Plush cover ☐

Tidy cover (hygiene) ☐



Plush



Tidy

sw 35 cm	62230-60	63130-63
sw 37.5 cm	62231-60	63131-63
sw 40 cm	62232-60	63132-63
sw 42.5 cm	62233-60	63133-63
sw 45 cm	62234-60	63134-63
sw 47.5 cm	62374-60	63135-63
sw 50 cm	62266-60	63136-63

Seat wedge



Comfort wedge	Provides extra seat contour for seat cushion Universal.	
sw 35-37.5 cm		80714 ☐
sw 40-45 cm		80715 ☐
sw 47.5-50 cm		80779 ☐

Widening kit

sw + 1.2 cm Not with drum brake 24079 ☐

Trays	Pieces	Item no
-------	--------	---------

Tray



Transparent	for 38 cm arm support 28856-1/-2	
sw 35 cm		24417 ☐
sw 37.5 cm		25187 ☐
sw 40 cm		24418 ☐
sw 42.5 cm		25188 ☐
sw 45 cm		24419 ☐
sw 47.5 cm		25731 ☐
sw 50 cm		24601 ☐

Hemi tray



Transparent	for 25 cm arm support 28856-1/-2	
sw 35 cm		26556 ☐
sw 37.5 cm		26557 ☐
sw 40 cm		26558 ☐
sw 42.5 cm		26559 ☐
sw 45 cm		26560 ☐
sw 47.5 cm		26561 ☐
sw 50 cm		26562 ☐

Transparent	for 38 cm arm support 28856-1/-2	
sw 35 cm		26563 ☐
sw 37.5 cm		26564 ☐
sw 40 cm		26565 ☐
sw 42.5 cm		26566 ☐
sw 45 cm		26567 ☐
sw 47.5 cm		26568 ☐
sw 50 cm		26569 ☐

Anti glide for hemi tray 81633 ☐

Half tray



Foldable for arm support 28855-1/-2 & 28856-1/-2 ☐ Right ☐ Left 27352-1/-2 ☐

Attachment for tray

Round tubes fits tubes up to 18mm ☐ Right ☐ Left 27956-1/-2 ☐

One arm drive



Drive unit QR			
20"		84101	☐
22"		84102	☐
24"		84103	☐

Assembly kit one arm drive QR	24071	☐
-------------------------------	-------	--------------------------

Cross bar one arm drive QR			
sw 35 cm		24072	☐
sw 37.5 cm		24073	☐
sw 40 cm		24074	☐
sw 42.5 cm		24075	☐
sw 45 cm		24076	☐
sw 47.5 cm		24077	☐
sw 50 cm		24078	☐

Telescopic shaft			
sw 35-37.5 cm		24171	☐
sw 40-42.5 cm		24172	☐
sw 45 cm		24173	☐
sw 47.5-50 cm		24174	☐

Spokeguard



Transparent			
for C6 & Maxx			
20"	2	28934	☐
22"	2	28935	☐
24"	2	28936	☐

Transparent			
for C5 DB			
20"	2	28816-28	☐
22"	2	28817-28	☐
24"	2	28818-28	☐

Additional information to order: