



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

# Crocodile

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 7/21/2025

## Measurements and Fitting

		Size 0	Size 1	Size 1 AT	Size 2	Size 2 AT	Size 3	Size 3 AT
	Width between handles	100-250 mm (3.9" - 9.8")	190-350 mm (7.5" - 13.8")	190-350 mm (7.5" - 13.8")	240-420 mm (9.4" - 16.5")	240-420 mm (9.4" - 16.5")	330-480 mm (13" - 18.9")	330-480 mm (13" - 18.9")
	Height from floor to handles	300-610 mm (11.8" - 24")	420-710 mm (16.5" - 28") 380-700 mm* (15" - 27.6")	460-750 mm (18.1" - 29.5")	560-850 mm (22" - 33.5")	600-890 mm (23.6" - 35")	740-1025 mm (29.1" - 40.4")	780-1065 mm (30.7" - 41.9")
	Max. user height	1050 mm (41.3")	1150 mm (45.3")	1190 mm (46.9")	1500 mm (59")	1540 mm (60.6")	1800 mm (70.9")	1840 mm (72.4")
	Max. user mass	30 kg (66 lbs)	30 kg (66 lbs)	30 kg (66 lbs)	45 kg (99 lbs)	45 kg (99 lbs)	80 kg (176 lbs)	80 kg (176 lbs)
	Length	640 mm (25.2")	720 mm (28.3")	780 mm (30.7")	840 mm (33.1")	900 mm (35.4")	870 mm (34.3")	930 mm (36.6")
	Overall width	580 mm (22.8")	630 mm (24.8")	660 mm (26")	690 mm (27.2")	700 mm (27.6")	690 mm (27.2")	720 mm (28.3")
	Folded length	650 mm (25.6")	720 mm (28.3")	750 mm (29.5")	860 mm (33.9")	890 mm (35")	1000 mm (39.4")	1030 mm (40.6")
	Folded width	580 mm (22.8")	620 mm (24.4")	660 mm (26")	690 mm (27.2")	700 mm (27.6")	690 mm (27.2")	720 mm (28.3")
	Folded height	300 mm (11.8")	300 mm (11.8")	345 mm (13.6")	270 mm (10.6")	310 mm (12.2")	290 mm (11.4")	330 mm (13")

\*Mounted with the Starter kit accessory. Further measurements are available on our website.

## Colors



Crocodile



Crocodile All-Terrain

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1197.4 Pricing and specifications subject to change without notice.



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

# Crocodile

## Crocodile Frame

*Make a selection of either the standard Crocodile or the Crocodile AT*

### Crocodile Base

Crocodile, White, Size 0 (86800)	\$733
Crocodile, Yellow, Size 1 (86801)	\$808
Crocodile, Orange, Size 2 (86802)	\$808
Crocodile, Metallic Black, Size 2 (86852)	\$808
Crocodile, Anthracite Grey, Size 3 (86803)	\$1333

### Crocodile AT Base

Crocodile AT, Red, Size 1 (86801AT)	\$975
Crocodile AT, Green, Size 2 (86802AT)	\$1025
Crocodile AT, Grey, Size 3 (86803AT)	\$1525

### Anti-Tip Device



#### Swing-Away Anti-Tips for Crocodile

Swing-Away Anti-Tips, for Croc Size 0-2 (86815)	\$175
Swing-Away Anti-Tips, for Croc Size 3 (86815-3S)	\$175



#### Fixed Anti-Tips for Crocodile

Anti-Tip Fixed, Pair, for Croc Size 0-2 (86815-F)	\$84
Anti-Tip Device Fixed, for Croc Size 3 (86815-3F)	\$84

### Swing-Away Anti-Tips for Crocodile AT

Swing-Away Anti-Tips for Croc AT Size 1-2 (868151)	\$225
Swing-Away Anti-Tips for Croc AT Size 3 (868153)	\$245


## Brake



#### Hand Brake

*Eliminates the function of anti-reverse wheel stop*

Hand Brake Kit, for Size 1-2 (86834)	\$262
Hand Brake Kit, for Size 3 (86873)	\$262



#### Deceleration Hand Brake Kit

*Not available with Croc AT. Eliminates the function of anti-reverse wheel stop.*

Deceleration Hand Brake Kit, Sizes 1-2 (86834-DEC)	\$262
Deceleration Hand Brake Kit, Size 3 (86873-DEC)	\$262

### Drum Brake Kit

*Not available with Croc AT*

Drum Brake Kit, Size 1 (86834-1DB)	\$995
Drum Brake Kit, Size 2 (86834-2DB)	\$995
Drum Brake Kit, Size 3 (86873-3DB)	\$995

### Decelerating Hand Brake for Croc AT

Deceleration Hand Brake for Croc AT, Size 1 (868341)	\$295
Deceleration Hand Brake for Croc AT, Size 2 (868342)	\$295
Deceleration Hand Brake for Croc AT, Size 3 (868343)	\$295

### Holder for Hand Brake

*Required with use of Hand Brake and Forearm Support. Allows for mounting of Hand Brake to Forearm Support.*

Small, for Size 1 (86874-1)	\$367
Medium, for Size 2 (86874-2)	\$414
Large, for Size 3 (86874-3)	\$445



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

# Crocodile

## Handles and Grip Bars

### Grip Bar

*Grip Bar is available for Crocodile sizes 0, 1, and 2*



#### Standard Grip Bar

Grip Bar, for Size 0 (86879-0)	\$141
Grip Bar, for Size 1 (86879-1)	\$163
Grip Bar, for Size 2 (86879-2)	\$196



#### Angle Adjustable Grip Bar

*Angle Adjustable Grip Bar requires Surrounding Support or Back and Lateral Supports*

Angle Adjustable Grip Bar, Size 0 (86878-0)	\$235
Angle Adjustable Grip Bar, Size 1 (86878-1)	\$275
Angle Adjustable Grip Bar, Size 2 (86878-2)	\$305

### Attendant Handle



#### Attendant Guide Bar

*Available only for Crocodile sizes 0, 1 and 2*

Attendant Handle (86866-1)	\$249
----------------------------	-------

## Arm and Forearm Options

### Arm Supports



#### Forearm Support with Handgrip

Forearm Support Pair - for Sizes 1 (86858-1)	\$472
Forearm Support Pair - for Size 2 (86859-1)	\$472
Forearm Support Pair - for Size 3 (86860-1)	\$472



#### Handle Knobs

Handle Knobs, Pair (86861)	\$56
----------------------------	------

### Padded Forearm Positioner (Each)

*For use with Forearm Support with Handgrip*

Padded Forearm Positioner Strap Qty 1 (869370)	\$75
Padded Forearm Positioner Strap Qty 2 (869370)	\$150



### Starter Kit for Crocodile

*Available for Crocodile size 1*



Crocodile Starter Kit, for Croc size 1 (86801-0)	\$513
--	-------



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

# Crocodile

## Back and Pelvic and Hip Supports



### Surrounding Supports

Surrounding Support for Size 0 (868610-0S)	\$595
Surrounding Support for Size 1 (868610-1S)	\$595
Surrounding Support for Size 2 (868610-2M)	\$625



### Back Support

Back Support, Adjustable Padded for Size 3 (86886)	\$552
--	-------



### Pelvic Support With Laterals

Pelvic Support With Laterals, for Sizes 0-2 (86808)	\$287
Pelvic Support With Laterals, For Sizes 1-2 (86809)	\$287

### Abdominal Support Belt

*For size 3 Back Support only*

Abdominal Support Belt, for Size 3 (85490-2)	\$134
--	-------



### Lateral Strap

*Used with Pelvic Support with Laterals*

Lateral Strap (SN4016)	\$56
------------------------	------



### Mounting Hardware

*Required for use of Pelvic and Back Support on size 3*

Cross Bar for Mounting Back/Pelvic Supports, for Size 3 (86880)	\$107
---	-------



### Hip Support Pair

Hip Supports, Pair (86819)	\$257
----------------------------	-------



### Pad for Back Support

Pad for Back Support (86816)	\$195
------------------------------	-------

## Seating Options



### Solid Flip Down Seat

Size 0 (86820-01)	\$186
Size 1 (86820)	\$186
<i>Not compatible with Crocodile starter kit</i>	
Size 2 (86822)	\$186
Metallic Black, Size 2 (86823)	\$186
Size 3 (86863)	\$186



### Seat Sling

*Not Available with Size 3*

Seat Sling, for Size 1 (86824)	\$147
Seat Sling, for Size 2 (86825)	\$147



### Lightweight Sling Seat

*Requires Pelvic Support with Laterals and Lateral Strap or Surrounding Supports*

Lightweight Sling Seat, for Size 0-1 (86824-1)	\$147
Lightweight Sling Seat, for Size 2 (86825-1)	\$147

## Accessories



Storage Bag for Sizes 1-2 (86839)	\$78
Storage Bag, for Size 3 (86843)	\$230