

# Convoid Carrot 3 Child Restraint System

## Canadian Retail Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ ATP/RTS: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

<b>Carrot 3 Standard Features:</b>								
<ul style="list-style-type: none"> <li>Recline Adjustment (90° - 180°)</li> <li>Head Support: 8.5" x 10.5" (H x L)</li> <li>Lower Anchors Connector System</li> <li>5-Point Positioning Harness</li> </ul>	<ul style="list-style-type: none"> <li>Trunk Positioning Inserts (Body Balance Pads)</li> <li>Seat Positioning Inserts</li> <li>Hip Positioning Inserts</li> <li>Chest Support: 9" Width</li> </ul>	<ul style="list-style-type: none"> <li>Tether Strap (Black)</li> </ul>						
<p><b>Standard Models:</b></p> <table style="width: 100%; border: none;"> <tr> <td style="width: 60%;"><input type="checkbox"/> Carrot 3 CT3BLACK - CCT3000-100</td> <td style="text-align: right;"><b>Price</b></td> </tr> <tr> <td><input type="checkbox"/> Carrot 3 CT3BLUE - CCT3000-200</td> <td style="text-align: right;"><b>\$3,124</b></td> </tr> <tr> <td><input type="checkbox"/> Carrot 3 CT3RED - CCT3000-300</td> <td style="text-align: right;"><b>\$3,124</b></td> </tr> </table>			<input type="checkbox"/> Carrot 3 CT3BLACK - CCT3000-100	<b>Price</b>	<input type="checkbox"/> Carrot 3 CT3BLUE - CCT3000-200	<b>\$3,124</b>	<input type="checkbox"/> Carrot 3 CT3RED - CCT3000-300	<b>\$3,124</b>
<input type="checkbox"/> Carrot 3 CT3BLACK - CCT3000-100	<b>Price</b>							
<input type="checkbox"/> Carrot 3 CT3BLUE - CCT3000-200	<b>\$3,124</b>							
<input type="checkbox"/> Carrot 3 CT3RED - CCT3000-300	<b>\$3,124</b>							
<b>All prices in CAD</b>								

- Carrot 3 Recommended Item**
- |   |              |
|---|--------------|
| <input type="checkbox"/> Top Securement Strap*        | <b>Price</b> |
| <input type="checkbox"/> Lower Securement Strap*      | <b>\$73</b>  |
| <b>*For vehicles not equipped with anchor points.</b> | <b>\$92</b>  |

Black	Blue	Red
CCT3HSF-10	CCT3HSF-10	CCT3HSF-10
CCT3LFS-10	CCT3LFS-10	CCT3LFS-10

- Trunk Positioning**
- |  |              |
|--|--------------|
| <input type="checkbox"/> Backrest Extension - 2" (Extends back height to 13.8")* | <b>Price</b> |
| <input type="checkbox"/> Backrest Extension - 4" (Extends back height to 15.8")* | <b>\$96</b>  |
| <input type="checkbox"/> Shoulder Protector Wings (Extends back height 6")       | <b>\$127</b> |
| <input type="checkbox"/> Firm Lateral Trunk Supports - Black***                  | <b>\$737</b> |
|  | <b>\$262</b> |
- \*The 2" & 4" Backrest Extensions cannot be used simultaneously.**  
**^Width adjustment from 10" - 12", and 2" of height adjustment.**  
**\*\*Manufacturer installation required.**

Black	Blue	Red
CCT3BE2-10	CCT3BE2-20	CCT3BE2-30
CCT3BE4-10	CCT3BE4-20	CCT3BE4-30
CCT3SP0-10	CCT3SP0-20	CCT3SP0-30
CCT3LS0-10	CCT3LS0-10	CCT3LS0-10

- Pelvic Positioning**
- |   |              |
|---|--------------|
| <input type="checkbox"/> Rodded Seat Extension [Large/Curved (extends seat depth 7")] | <b>Price</b> |
| <input type="checkbox"/> Rodded Seat Extension [Small/90° (extends seat depth 2")]    | <b>\$291</b> |
| <input type="checkbox"/> Seat Extension [2" (Extends seat depth to 12")]              | <b>\$237</b> |
| <input type="checkbox"/> Seat Extension [4" (Extends seat depth to 14")]              | <b>\$131</b> |
| <input type="checkbox"/> Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°)     | <b>\$158</b> |
| <input type="checkbox"/> Pommel - Black* (Medial Thigh Separator)                     | <b>\$158</b> |
| <input type="checkbox"/> Anti-Adduction/Knee Separator Pad                            | <b>\$166</b> |
| <input type="checkbox"/> Pelvis Lock Harness (Small) (20.5" - 33.75")                 | <b>\$211</b> |
| <input type="checkbox"/> Pelvis Lock Harness (Medium) (23" - 39.5")                   | <b>\$350</b> |
| <input type="checkbox"/> Pelvis Lock Harness (Large) (25" - 46")                      | <b>\$350</b> |
- ^Growth specifications with listed extensions:**  
**Seat depth: 9.8" - 22.3" / Backrest height: 11.8" - 22.8"**  
**\*Rodded Seat Extension or Rodded Footrest required. 2" & 4" Seat Extensions cannot be used simultaneously when using Rodded Seat Extension.**  
**v7" Rodded Seat Extension required.**

Black	Blue	Red
CCT3RSL-10	CCT3RSL-20	CCT3RSL-30
CCT3RSS-10	CCT3RSS-20	CCT3RSS-30
CCT3SE2-10	CCT3SE2-20	CCT3SE2-30
CCT3SE4-10	CCT3SE4-20	CCT3SE4-30
CCT3SA0-10	CCT3SA0-20	CCT3SA0-30
CCT3PM0-10	CCT3PM0-10	CCT3PM0-10
CCT3AAK-00	CCT3AAK-00	CCT3AAK-00
CCT3PLS-10	CCT3PLS-10	CCT3PLS-10
CCT3PLM-10	CCT3PLM-10	CCT3PLM-10
CCT3PLL-10	CCT3PLL-10	CCT3PLL-10

- Foot & Leg Positioning**
- |   |              |
|---|--------------|
| <input type="checkbox"/> Rodded Footrest - Small [Gray (Adjusts from 10"-12")]  | <b>Price</b> |
| <input type="checkbox"/> Rodded Footrest - Medium [Gray (Adjusts from 12"-16")] | <b>\$449</b> |
| <input type="checkbox"/> Rodded Footrest - Large [Gray (Adjusts from 14"-18")]  | <b>\$449</b> |
- \*Not available with Rodded Seat Extension**

Black	Blue	Red
CCT3FRS-00	CCT3FRS-00	CCT3FRS-00
CCT3FRM-00	CCT3FRM-00	CCT3FRM-00
CCT3FRL-00	CCT3FRL-00	CCT3FRL-00

- Small Support Kit**
- |  |              |
|--|--------------|
| <input type="checkbox"/> Small Support Kit | <b>Price</b> |
| • Small Head Support: 7.25" x 8.5" (H x L) | <b>\$644</b> |
| • Small Chest Support: Width 6.9"          | <b>\$0</b>   |
|  | <b>\$0</b>   |

Black	Blue	Red
CCT3SSK-10	CCT3SSK-20	CCT3SSK-30

- Accessories**
- |   |                |
|---|----------------|
| <input type="checkbox"/> Upper Extremity Support (Tray) (Crash Tested)            | <b>Price</b>   |
| <input type="checkbox"/> 3D Body Balance Pads L/R Set                             | <b>\$242</b>   |
| <input type="checkbox"/> Air-Mesh Cushions/Pads Set                               | <b>\$211</b>   |
| <input type="checkbox"/> Anti-Escape Vest (Small) (Circum: 24.75" Height: 10.75") | <b>\$385</b>   |
| <input type="checkbox"/> Anti-Escape Vest (Medium) (Circum: 31.5" Height: 12")    | <b>\$224</b>   |
| <input type="checkbox"/> Anti-Escape Vest (Large) (Circum: 37.5" Height: 13")     | <b>\$224</b>   |
| <input type="checkbox"/> Removeable Incontinence Cover                            | <b>\$280</b>   |
| <input type="checkbox"/> Hinge Lock Belt  | <b>\$59</b>    |
| <input type="checkbox"/> FAA Fixing Strap to Aircraft Seat                        | <b>\$280</b>   |
| <input type="checkbox"/> Booster Kit  | <b>\$1,611</b> |

Black	Blue	Red
CCT3ST0-10	CCT3ST0-20	CCT3ST0-30
CCT3BBP-00	CCT3BBP-00	CCT3BBP-00
CCT3AMC-10	CCT3AMC-10	CCT3AMC-10
CCT3AVS-00	CCT3AVS-00	CCT3AVS-00
CCT3AVM-00	CCT3AVM-00	CCT3AVM-00
CCT3AVL-00	CCT3AVL-00	CCT3AVL-00
CCT3RIR-00	CCT3RIR-00	CCT3RIR-00
CCT3HLB-00	CCT3HLB-00	CCT3HLB-00
CCT3FSA-00	CCT3FSA-00	CCT3FSA-00
CCRT3BK-199	CCRT3BK-299	CCRT3BK-399

**Standard Shipping Rates: Ground \$71.95 • For expedited shipping rates, please call 1-844-876-6245 or email [convoidsales.us@etac.com](mailto:convoidsales.us@etac.com)**

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 9/12/2024

OF0202-ENG-REV.10

Convoid  
 2830 California Street  
 Torrance, CA 90503  
 1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)



**Carrot 3 Base Model Specifications:**

Seat Width:	11.5"
Seat Depth Adjustment:	9.8" - 11.3"
Backrest Height Adjustment:	11.8" - 12.8"
Head Support Height:	8.5"
Recline Adjustment:	90 - 180°
Weight:	19.1 lbs.

**User Measurement Guide:**

Minimum User Weight:	30 lbs.
Maximum User Weight:	108 lbs.
Minimum User Height:	37"
Maximum User Height:	60"
Recommended User Age:	Ages 3-15

**Carrot 3 Specifications with Extensions:**

Seat Depth Adjustment:	9.8" - 22.3"
Backrest Height Adjustment:	11.8" - 22.8"
Shoulder Wing Height:	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"

**Base Model Measurements:\*\***

Seat Depth:	9.8" - 11.3"
Backrest Height:	11.8" - 12.8"
Head Support Height:	8.5"

*\*Measurements do not include head support & growth extensions*

*\*\*Base model includes 1.5" of seat depth growth (A) & 1" of backrest growth (B)*

**Growth Extensions:**


- 1a. Rodded seat extension - large
- 1b. Rodded seat extension - small
- 2. Seat extension - 4"
- 3. Seat extension - 2"

**Inches**


- 7"
- 2"
- 4"
- 2"
- 2"
- 4"
- 6"

- 4. Backrest extension - 2"
- 5. Backrest extension - 4"
- 6. Shoulder protector wings

**Note: The same type of 2" & 4" extensions cannot be used simultaneously.**



**Convoid Carrot 3 Seat Depth Growth:**



**Seat Depth Adjustments**

**2" Rodded Seat Depth Extensions**

**2" Rodded + 2" Seat Depth Extensions**

**7" Rodded Seat Depth Extensions**

**7" Rodded + 2" Seat Depth Extensions**

**7" Rodded + 4" Seat Depth Extensions**

Position 1 (inch)	9.8"	11.8"	13.8"	16.8"	18.8	20.8"
Position 2 (inch)	11.3"	13.3"	15.3"	18.3"	20.3	22.3"

**NOTE: 2" & 4" Seat Depth Extensions cannot be used simultaneously when using Rodded Seat Depth Extension**

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166