

# Convaid Carrot 3 Child Restraint System

## Canadian Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_  
Contact: \_\_\_\_\_ ATP/RTS#: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:						
Address:	Address:						
City/State/Zip:	City/State/Zip:						
Phone:	Phone:						
Email:							
<div><b>Carrot 3 Standard Features:</b><ul style="list-style-type: none"><li>• Recline Adjustment (90° - 180°)</li><li>• Head Support: 8.5" x 10.5" (H x L)</li><li>• Lower Anchors Connector System</li><li>• 5-Point Positioning Harness</li><li>• Tether Strap (Black)</li><li>• Trunk Positioning Inserts (Body Balance Pads)</li><li>• Seat Positioning Inserts</li><li>• Hip Positioning Inserts</li><li>• Chest Support: 9" Width</li></ul></div> <div><b>Standard Models:</b><table><tr><td><input type="checkbox"/> Carrot 3 CT3BLACK - CCT3000-100</td><td><b>Price</b></td></tr><tr><td><input type="checkbox"/> Carrot 3 CT3BLUE - CCT3000-200</td><td><b>\$3,368</b></td></tr><tr><td><input type="checkbox"/> Carrot 3 CT3RED - CCT3000-300</td><td><b>\$3,368</b></td></tr></table></div> <div><b>All prices in CAD</b></div>		<input type="checkbox"/> Carrot 3 CT3BLACK - CCT3000-100	<b>Price</b>	<input type="checkbox"/> Carrot 3 CT3BLUE - CCT3000-200	<b>\$3,368</b>	<input type="checkbox"/> Carrot 3 CT3RED - CCT3000-300	<b>\$3,368</b>
<input type="checkbox"/> Carrot 3 CT3BLACK - CCT3000-100	<b>Price</b>						
<input type="checkbox"/> Carrot 3 CT3BLUE - CCT3000-200	<b>\$3,368</b>						
<input type="checkbox"/> Carrot 3 CT3RED - CCT3000-300	<b>\$3,368</b>						

- Carrot 3 Recommended Item**
- ☐ Top Securement Strap\* **Price \$79**
- ☐ Lower Securement Strap\* **\$99**

**\*For vehicles not equipped with anchor points.**

- Trunk Positioning**
- ☐ Backrest Extension - 2" (Extends back height to 13.8")\* **Price \$103**
- ☐ Backrest Extension - 4" (Extends back height to 15.8")\* **\$137**
- ☐ Shoulder Protector Wings (Extends back height 6") **\$795**
- ☐ Firm Lateral Trunk Supports - Black\*\*\* **\$282**

**\*The 2" & 4" Backrest Extensions cannot be used simultaneously.**

**^Width adjustment from 10" - 12", and 2" of height adjustment.**

**\*\*Manufacturer installation required.**

- Pelvic Positioning**
- ☐ Rodded Seat Extension [Large/Curved (extends seat depth 7")]<sup>^</sup> **Price \$313**
- ☐ Rodded Seat Extension [Small/90° (extends seat depth 2")]<sup>^</sup> **\$255**
- ☐ Seat Extension [2" (Extends seat depth to 12")]<sup>^</sup> **\$141**
- ☐ Seat Extension [4" (Extends seat depth to 14")]<sup>^v</sup> **\$171**
- ☐ Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°) **\$171**
- ☐ Pommel - Black\* (Medial Thigh Separator) **\$179**
- ☐ Anti-Adduction/Knee Separator Pad **\$227**
- ☐ Pelvis Lock Harness (Small) (20.5" - 33.75") **\$377**
- ☐ Pelvis Lock Harness (Medium) (23" - 39.5") **\$377**
- ☐ Pelvis Lock Harness (Large) (25" - 46") **\$377**

**^Growth specifications with listed extensions:**

**Seat depth: 9.8" - 22.3" / Backrest height: 11.8" - 22.8".**

**\*Rodded Seat Extension or Rodded Footrest required. 2" & 4" Seat Extensions cannot be used simultaneously when using Rodded Seat Extension.**

**v7" Rodded Seat Extension required.**

- Foot & Leg Positioning**
- ☐ Rodded Footrest - Small [Gray (Adjusts from 10"-12")]<sup>^</sup> **Price \$484**
- ☐ Rodded Footrest - Medium [Gray (Adjusts from 12"-16")]<sup>^</sup> **\$484**
- ☐ Rodded Footrest - Large [Gray (Adjusts from 14"-18")]<sup>^</sup> **\$484**

**\*Not available with Rodded Seat Extension.**

- Small Support Kit**
- ☐ Small Support Kit **Price \$694**
- Small Head Support: 7.25" x 8.5" (H x L)
  - Small Chest Support: Width 6.9"

- Accessories**
- ☐ Upper Extremity Support (Tray) (Crash Tested) **Price \$261**
- ☐ 3D Body Balance Pads L/R Set **\$227**
- ☐ Air-Mesh Cushions/Pads Set **\$415**
- ☐ Anti-Escape Vest (Small) (Circum: 24.75" Height: 10.75") **\$241**
- ☐ Anti-Escape Vest (Medium) (Circum: 31.5" Height: 12") **\$241**
- ☐ Anti-Escape Vest (Large) (Circum: 37.5" Height: 13") **\$241**
- ☐ Removeable Incontinence Cover **\$302**
- ☐ Hinge Lock Belt **\$64**
- ☐ FAA Fixing Strap to Aircraft Seat **\$302**
- ☐ Booster Kit **\$1,736**

Black	Blue	Red
CCT3HSF-10	CCT3HSF-10	CCT3HSF-10
CCT3LFS-10	CCT3LFS-10	CCT3LFS-10

Black	Blue	Red
CCT3BE2-10	CCT3BE2-20	CCT3BE2-30
CCT3BE4-10	CCT3BE4-20	CCT3BE4-30
CCT3SP0-10	CCT3SP0-20	CCT3SP0-30
CCT3LS0-10	CCT3LS0-10	CCT3LS0-10

Black	Blue	Red
CCT3RSL-10	CCT3RSL-20	CCT3RSL-30
CCT3RSS-10	CCT3RSS-20	CCT3RSS-30
CCT3SE2-10	CCT3SE2-20	CCT3SE2-30
CCT3SE4-10	CCT3SE4-20	CCT3SE4-30
CCT3SA0-10	CCT3SA0-20	CCT3SA0-30
CCT3PM0-10	CCT3PM0-10	CCT3PM0-10
CCT3AAK-00	CCT3AAK-00	CCT3AAK-00
CCT3PLS-10	CCT3PLS-10	CCT3PLS-10
CCT3PLM-10	CCT3PLM-10	CCT3PLM-10
CCT3PLL-10	CCT3PLL-10	CCT3PLL-10

Black	Blue	Red
CCT3FRS-00	CCT3FRS-00	CCT3FRS-00
CCT3FRM-00	CCT3FRM-00	CCT3FRM-00
CCT3FRL-00	CCT3FRL-00	CCT3FRL-00

Black	Blue	Red
CCT3SSK-10	CCT3SSK-20	CCT3SSK-30

Black	Blue	Red
CCT3ST0-10	CCT3ST0-20	CCT3ST0-30
CCT3BBP-00	CCT3BBP-00	CCT3BBP-00
CCT3AMC-10	CCT3AMC-10	CCT3AMC-10
CCT3AVS-00	CCT3AVS-00	CCT3AVS-00
CCT3AVM-00	CCT3AVM-00	CCT3AVM-00
CCT3AVL-00	CCT3AVL-00	CCT3AVL-00
CCT3RIR-00	CCT3RIR-00	CCT3RIR-00
CCT3HLB-00	CCT3HLB-00	CCT3HLB-00
CCT3FSA-00	CCT3FSA-00	CCT3FSA-00
CCRT3BK-199	CCRT3BK-299	CCRT3BK-399

**Standard Shipping Rate: Ground \$71.95 • For expedited shipping rates, please call 1-844-876-6245 or email [convaidsales.us@etac.com](mailto:convaidsales.us@etac.com)**

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 2/1/2025

Convaid  
2830 California Street  
Torrance, CA 90503  
1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF0101-ENG-REV.01

**etac**<sup>®</sup>  
Creating Possibilities

Carrot 3 Base Model Specifications:

Seat Width:	11.5"
Seat Depth Adjustment:	9.8" - 11.3"
Backrest Height Adjustment:	11.8" - 12.8"
Head Support Height:	8.5"
Recline Adjustment:	90 - 180°
Weight:	19.1 lbs.

User Measurement Guide:

Minimum User Weight:	30 lbs.
Maximum User Weight:	108 lbs.
Minimum User Height:	37"
Maximum User Height:	60"
Recommended User Age:	Ages 3-15

Carrot 3 Specifications with Extensions:

Seat Depth Adjustment:	9.8" - 22.3"
Backrest Height Adjustment:	11.8" - 22.8"
Shoulder Wing Height:	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"

Base Model Measurements:\*\*

Seat Depth:	9.8" - 11.3"
Backrest Height:	11.8" - 12.8"
Head Support Height:	8.5"

\*Measurements do not include head support & growth extensions  
^Base model includes 1.5" of seat depth growth (A) & 1" of backrest growth (B)

Growth Extensions:

- 1a. Rodded seat extension - large
- 1b. Rodded seat extension - small


- 2.a Seat extension - 4"
- 3.b Seat extension - 2"

- 4.a Backrest extension - 2"
- 5.a Backrest extension - 4"
- 6.a Shoulder protector wings


Note: The same type of 2" & 4" extensions cannot be used simultaneously.

Inches

- 7"
- 2"
- 4"
- 2"
- 2"
- 4"
- 6"



Convoid Carrot 3 Seat Depth Growth:



Position 1 (inch)	9.8"	11.8"	13.8"	16.8"	18.8"	20.8"
Position 2 (inch)	11.3"	13.3"	15.3"	18.3"	20.3"	22.3"

NOTE: 2" & 4" Seat Depth Extensions cannot be used simultaneously when using Rodded Seat Depth Extension

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166