

Convoid Carrot 3 Child Restraint System

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Carrot 3 Standard Features:

- Recline Adjustment (90° - 180°)
- Head Support: 8.5" x 10.5" (H x L)
- Lower Anchors Connector System
- 5-Point Positioning Harness
- Trunk Positioning Inserts (Body Balance Pads)
- Seat Positioning Inserts
- Hip Positioning Inserts
- Chest Support: 9" Width
- Tether Strap (Black)

Standard Models:

	Price
<input type="checkbox"/> Carrot 3 CT3BLACK - CCT3000-100	\$2,386
<input type="checkbox"/> Carrot 3 CT3BLUE - CCT3000-200	\$2,386
<input type="checkbox"/> Carrot 3 CT3RED - CCT3000-300	\$2,386

Carrot 3 Recommended Item

<input type="checkbox"/> Top Securement Strap*	\$56
<input type="checkbox"/> Lower Securement Strap*	\$70

**For vehicles not equipped with anchor points.*

Trunk Positioning

<input type="checkbox"/> Backrest Extension - 2" (Extends back height to 13.8")*	\$73
<input type="checkbox"/> Backrest Extension - 4" (Extends back height to 15.8")*	\$97
<input type="checkbox"/> Shoulder Protector Wings (Extends back height 6")	\$563
<input type="checkbox"/> Firm Lateral Trunk Supports - Black***	\$200

**The 2" & 4" Backrest Extensions cannot be used simultaneously.*

^Width adjustment from 10" - 12", and 2" of height adjustment.

***Manufacturer installation required.*

Pelvic Positioning

<input type="checkbox"/> Rodded Seat Extension [Large/Curved (extends seat depth 7")] [^]	\$222
<input type="checkbox"/> Rodded Seat Extension [Small/90° (extends seat depth 2")] [^]	\$181
<input type="checkbox"/> Seat Extension [2" (Extends seat depth to 12")] ^{^*}	\$100
<input type="checkbox"/> Seat Extension [4" (Extends seat depth to 14")] ^{^*v}	\$121
<input type="checkbox"/> Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°)	\$121
<input type="checkbox"/> Pommel - Black* (Medial Thigh Separator)	\$127
<input type="checkbox"/> Anti-Adduction/Knee Separator Pad	\$161
<input type="checkbox"/> Pelvis Lock Harness (Small) (20.5" - 33.75")	\$267
<input type="checkbox"/> Pelvis Lock Harness (Medium) (23" - 39.5")	\$267
<input type="checkbox"/> Pelvis Lock Harness (Large) (25" - 46")	\$267

^Growth specifications with listed extensions:

Seat depth: 9.8" - 22.3" / Backrest height: 11.8" - 22.8".

**Rodded Seat Extension or Rodded Footrest required. 2" & 4" Seat Extensions cannot be used simultaneously when using Rodded Seat Extension.*

v7" Rodded Seat Extension required.

Foot & Leg Positioning

<input type="checkbox"/> Rodded Footrest - Small [Gray (Adjusts from 10"-12")] [*]	\$343
<input type="checkbox"/> Rodded Footrest - Medium [Gray (Adjusts from 12"-16")] [*]	\$343
<input type="checkbox"/> Rodded Footrest - Large [Gray (Adjusts from 14"-18")] [*]	\$343

**Not available with Rodded Seat Extension.*

Small Support Kit

<input type="checkbox"/> Small Support Kit	\$492
• Small Head Support: 7.25" x 8.5" (H x L)	\$0
• Small Chest Support: Width 6.9"	\$0

Accessories

<input type="checkbox"/> Upper Extremity Support (Tray) (Crash Tested)	\$185
<input type="checkbox"/> 3D Body Balance Pads L/R Set	\$161
<input type="checkbox"/> Air-Mesh Cushions/Pads Set	\$294
<input type="checkbox"/> Anti-Escape Vest (Small) (Circum: 24.75" Height: 10.75")	\$171
<input type="checkbox"/> Anti-Escape Vest (Medium) (Circum: 31.5" Height: 12")	\$171
<input type="checkbox"/> Anti-Escape Vest (Large) (Circum: 37.5" Height: 13")	\$171
<input type="checkbox"/> Removeable Incontinence Cover	\$214
<input type="checkbox"/> Hinge Lock Belt	\$45
<input type="checkbox"/> FAA Fixing Strap to Aircraft Seat	\$214
<input type="checkbox"/> Booster Kit	\$1,230

Black	Blue	Red
CCT3HSF-10	CCT3HSF-10	CCT3HSF-10
CCT3LFS-10	CCT3LFS-10	CCT3LFS-10

Black	Blue	Red
CCT3BE2-10	CCT3BE2-20	CCT3BE2-30
CCT3BE4-10	CCT3BE4-20	CCT3BE4-30
CCT3SP0-10	CCT3SP0-20	CCT3SP0-30
CCT3LS0-10	CCT3LS0-10	CCT3LS0-10

Black	Blue	Red
CCT3RSL-10	CCT3RSL-20	CCT3RSL-30
CCT3RSS-10	CCT3RSS-20	CCT3RSS-30
CCT3SE2-10	CCT3SE2-20	CCT3SE2-30
CCT3SE4-10	CCT3SE4-20	CCT3SE4-30
CCT3SA0-10	CCT3SA0-20	CCT3SA0-30
CCT3PM0-10	CCT3PM0-10	CCT3PM0-10
CCT3AAK-00	CCT3AAK-00	CCT3AAK-00
CCT3PLS-10	CCT3PLS-10	CCT3PLS-10
CCT3PLM-10	CCT3PLM-10	CCT3PLM-10
CCT3PLL-10	CCT3PLL-10	CCT3PLL-10

Black	Blue	Red
CCT3FRS-00	CCT3FRS-00	CCT3FRS-00
CCT3FRM-00	CCT3FRM-00	CCT3FRM-00
CCT3FRL-00	CCT3FRL-00	CCT3FRL-00

Black	Blue	Red
CCT3SSK-10	CCT3SSK-20	CCT3SSK-30

Black	Blue	Red
CCT3ST0-10	CCT3ST0-20	CCT3ST0-30
CCT3BBP-00	CCT3BBP-00	CCT3BBP-00
CCT3AMC-10	CCT3AMC-10	CCT3AMC-10
CCT3AVS-00	CCT3AVS-00	CCT3AVS-00
CCT3AVM-00	CCT3AVM-00	CCT3AVM-00
CCT3AVL-00	CCT3AVL-00	CCT3AVL-00
CCT3RIR-00	CCT3RIR-00	CCT3RIR-00
CCT3HLB-00	CCT3HLB-00	CCT3HLB-00
CCT3FSA-00	CCT3FSA-00	CCT3FSA-00
CCRT3BK-199	CCRT3BK-299	CCRT3BK-399

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.

Convoid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

Effective: 1/10/2024

OF0200-ENG-REV.08



Carrot 3 Base Model Specifications:

Seat Width:	11.5"
Seat Depth Adjustment:	9.8" - 11.3"
Backrest Height Adjustment:	11.8" - 12.8"
Head Support Height:	8.5"
Recline Adjustment:	90 - 180°
Weight:	19.1 lbs.

User Measurement Guide:

Minimum User Weight:	30 lbs.
Maximum User Weight:	108 lbs.
Minimum User Height:	37"
Maximum User Height:	60"
Recommended User Age:	Ages 3-15

Carrot 3 Specifications with Extensions:

Seat Depth Adjustment:	9.8" - 22.3"
Backrest Height Adjustment:	11.8" - 22.8"
Shoulder Wing Height:	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"

Base Model Measurements:*^

Seat Depth:	9.8" - 11.3"
Backrest Height:	11.8" - 12.8"
Head Support Height:	8.5"

*Measurements do not include head support & growth extensions

^Base model includes 1.5" of seat depth growth (A) & 1" of backrest growth (B)

Growth Extensions:

- 1a. Rodded seat extension - large
- 1b. Rodded seat extension - small


- 2. Seat extension - 4"
- 3. Seat extension - 2"

- 4. Backrest extension - 2"
- 5. Backrest extension - 4"
- 6. Shoulder protector wings


Note: The same type of 2" & 4" extensions cannot be used simultaneously.

Inches

- 7"
- 2"
- 4"
- 2"
- 2"
- 4"
- 6"



Convoid Carrot 3 Seat Depth Growth:



Seat Depth Adjustments

2" Rodded Seat Depth Extensions

2" Rodded + 2" Seat Depth Extensions

7" Rodded Seat Depth Extensions

7" Rodded + 2" Seat Depth Extensions

7" Rodded + 4" Seat Depth Extensions

Position 1 (inch)	9.8"	11.8"	13.8"	16.8"	18.8"	20.8"
Position 2 (inch)	11.3"	13.3"	15.3"	18.3"	20.3"	22.3"

NOTE: 2" & 4" Seat Depth Extensions cannot be used simultaneously when using Rodded Seat Depth Extension

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166