

Convoid Coaster - HCPCS E1236

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Coaster Standard Features:

- Three fixed tilt setting options 10°, 20°, or 30° (Adjustable with Tools)
- Four Fixed Recline Setting Options (Adjustable with Tools)
- Tension Adjustable Back Straps
- 2-Point Positioning Belt
- Seat Depth Growth of 2"
- Standard Convoid Upholstery
- Zippered Pockets in Seat & Back Upholstery
- Removable, Height Adjustable Swing-Away Footplates
- One Piece Height Adjustable Push Handle (Inside Mount)

- 80° Front Rigging
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire
- Quick Release Rolko Rear Wheels
- Single Foot Operated Wheel Lock
- Clamshell Compact Fold
- Closure Strap
- Saddle Bags
- Powder Coated Frame
- Made in the USA

Standard Transit Additional Features:

- Transit Anchor Brackets
- Headrest Support - Headrest Extension
- Trunk Positioning H-Harness
- Pelvic Positioning 3-Point Belt
- Foot Positioners - Size Choice Required

Standard Models:

- Coaster CTR12 - C800121
- Coaster CTR14 - C800141
- Coaster CTR16 - C800161
- Coaster CTR18 - C800181

Price

- \$1,966**
- \$1,966**
- \$1,966**
- \$1,966**

Standard Transit Included Models:

- Coaster Transit CTR12T - C800122
- Coaster Transit CTR14T - C800142
- Coaster Transit CTR16T - C800162
- Coaster Transit CTR18T - C800182

Price

- \$2,642**
- \$2,642**
- \$2,642**
- \$2,642**

Select Tilt and Recline Setting

Factory Default Setting if No Selection is Made: 10° Fixed Tilt and 105° Recline

- Fixed Tilt set at 10°
Fixed Recline set at: 100 | 105 | 110 | 115
- Fixed Tilt set at 20°
Fixed Recline set at: 90 | 95 | 100 | 105
- Fixed Tilt set at 30°
Fixed Recline set at: 80 | 85 | 90 | 95

Tilt and Recline Options



Upholstery Options

Colors:

- Sassy Purple
- Cherry Red
- Princess Pink
- Royal Blue
- Forest Green
- Panther Black

- C903851
- C903855
- C903759
- C903850
- C903857
- C903854



Upholstery Color Options

Foam Kit:

- Positioning Wedges (Foam Kit for Zippered Upholstery Pockets)

Lumbar Support:

- Height and Tension Adjustable Back Positioner (Hardware Required)
- Height and Tension Adjustable Back Positioner Hardware

Price
\$100

HCPCS
K0108

\$126
\$254

E0956
E1028

	CTR12	CTR14	CTR16	CTR18
	C905567	C905568	C905569	C905570
	C905536	C905536	C905536	C905536
	C905538	C9055618	C905539	C905619

Seat Depth Extension (Standard 2")

- 3" Seat Depth Extension Tubes

Price
\$100
HCPCS
E1297

C905540	C905540	C905540	C905540
---------	---------	---------	---------

Seat Depth Ranges

Seat Width	Standard 2"	Option +3"
12"	10-12	12-15
14"	11-13	13-16
16"	12-14	14-17
18"	14.5-16.5	16.5-19.5



Frame Options

- One Piece Height Adjustable Push Handle (Inside Mount)
 - Height Adjustable Flip-Up Armrests (Pair) Required for Tray
 - Caster Lock for 2" Wide Front Wheels
 - Curb Tipper (*Elevates front of chair to go up curbs*)***
 - Rear Anti-Tip Tubes
- *****Cannot be combined with Rear Anti-Tip Tubes but can still be ordered.**

Price
N/C
\$244
\$100
\$105
\$123
HCPCS
E0973
K0073
E0971

CTR12	CTR14	CTR16	CTR18
C905620	C905621	C905622	C905623
C905611	C905611	C905611	C905611
C905526	C905526	C905526	C905526
C510737	C510737	C510737	C510737
C905527	C905527	C905527	C905527

Transit Required Items^**

- Transit Anchor Brackets**
- Head Support:**
- Headrest Extension**
- Trunk Positioning:**
- H-Harness Anterior Trunk Support w/ Padded Covers (3-Pt Belt Required)**
- Pelvic Positioning:**
- 3-Point Pelvic Positioning Belt**
- Foot Positioners - Select One Size Below:**
- Small 7" - 10" (Pair)**
- Medium 8" - 11" (Pair)** (Default if no size is selected)
- Large 9" - 12" (Pair)**

Price
\$214
\$138
\$105
\$52
\$167
\$167
\$167
HCPCS
E0966
E0960
E0978

CTR12	CTR14	CTR16	CTR18
C905531	C905531	C905531	C905531
C902863	C900484	C900138	C905571
C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

^Crash Tested - WC19 and FMVSS213 compliant.

****Transit required items are included in base price when selecting Coaster Transit Model**

Transit Option

- Q'Straint Transit Lap Belt

Price
\$198
HCPCS
E0978

CTR12	CTR14	CTR16	CTR18
C903094	C903094	C903094	C903094

Tire Options

- 8" x 2" Front Ribbed, 12.5" Rear Solid Knobby Tires
- 7.5" x 2" Front, 12.5" Rear Pneumatic Tires

Price
\$319
\$319
HCPCS

CTR12	CTR14	CTR16	CTR18
C905535	C905535	C905535	C905535
C905561	C905561	C905561	C905561

Head Supports

- Padded Headwing Head Support
- Occi Headwing Head Support
- Headrest Extension

Price
\$260
\$244
\$138
HCPCS
E0955
E0955
E0966

CTR12	CTR14	CTR16	CTR18
C905572	C905573	C905574	C905575
C905270	C905271	C905271	C905272
C902863	C900484	C900138	C905571

Trunk Positioning

- H-Harness Anterior Trunk Support w/ Padded Covers (3-Pt Belt Required)
- Full Torso Support Vest

Price
\$105
\$282
HCPCS
E0960
K0108

CTR12	CTR14	CTR16	CTR18
C903988	C903988	C903988	C903988
C905598	C905599	C905600	C905601

Pelvic Positioning

- 3-Point Pelvic Positioning Belt

Price
\$52
HCPCS
E0978

CTR12	CTR14	CTR16	CTR18
C904017	C904017	C904017	C904017

Foot & Leg Positioning

- Medial Thigh Support-Abductor
- Lateral Thigh Support-Adductor
- Calf Panel*
- Front Rigging Options: (80° Fixed Front Riggings Standard)**
- Angle Adjustable Elevating Leg Rest (90° - 180°) (Pair)
- Padded Elevating Legrest Cover (Pair)
- Footplate Options:**
- Individual Angle Adjustable Footplates (Pair) (Elevating Leg Rest Required)
 - Angle Adjustable Elevating Leg Rest (90° - 180°) (Pair)
- Angle Adjustable Footplates
- One Piece Angle Adjustable Flip-Up Footplate (-15° - +15°)
- Fully Padded Footbox (One Piece Footplate Required)
- Footplate Securement Strap
- Foot Positioners - Select One Size Below:**
- Small 7" - 10" (Pair)
- Medium 8" - 11" (Pair) (Default if no size is selected)
- Large 9" - 12" (Pair)

Price
\$173
\$173
\$66
\$267
\$100
\$170
\$267
\$170
\$118
\$347
\$33
\$167
\$167
\$167
HCPCS
E0957
E0956
E0995
K0990
K0108
K0040
K0990
K0040
K0040
K0108
K0108

CTR12	CTR14	CTR16	CTR18
C905576	C905577	C905578	C905579
C905580	C905581	C905582	C905583
C904580	C904580	C904580	C904580
C905541	C905542	C905543	C905544
C904619	C904619	C904619	C904619
C905624	C905625	C905626	C905627
C905541	C905542	C905543	C905544
C903445	C903268	C903268	C903269
C905557	C905558	C905559	
C904954	C904955	C904956	
C903584	C903584	C903584	C903584
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

***Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.**

Accessories

- Headrest Cover (Canopy)
- Extended Headrest Cover (Canopy)
- Rain Cover (Canopy Required)
- Mosquito Net (Canopy Required)
- Height Adjustable Upper Extremity Support Surface (Clear Tray)
 - Adj. Removable Hardware (Required for Clear Tray)
- Under Seat Storage Basket
- Medical Necessity Storage Basket
- Utility Bag
- IV Pole (Height Adjustable and Collapsible)^
- Adj. Swing Away Hardware (Required for IV Pole)

Price
\$244
\$267
\$161
\$161
\$119
\$243
\$122
\$315
\$45
\$209
\$147
HCPCS
E0950
E1028

CTR12	CTR14	CTR16	CTR18
C905017	C905018	C905019	C905020
C905602	C905603	C905604	C905605
C905590	C905591	C905592	C905593
C905594	C905595	C905596	C905597
C905547	C905548	C905549	C905550
C905551	C905551	C905551	C905551
C905584	C905585	C905585	C905586
C905587	C905588	C905588	C905589
C904531	C904531	C904531	C904531
C904566	C904566	C904566	C904566
C905612	C905612	C905612	C905612

^^Cannot be used with canopy.

The Coaster transit option includes four manufacturer installed wheelchair transportation anchors. The Coaster transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are CTR12: 66 lbs./30kg; CTR14: 100 lbs./45.5kg; CTR16 and CTR18: 170 lbs./77kg.

The Convoid Coaster was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convoid is unknown, and no warranty or claim about their performance is being, or can be, made.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.

Effective: 2/9/2024

Convoid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF2600-ENG-REV.01

 **etac**[®]
Creating Possibilities