

Convoid Cruiser

Canadian Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Cruiser Standard Features:

- Removable, Height Adjustable Swing-Away Footplates
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire (CX10/12/14/16)
- 8" Front, 11" Rear Solid Tire (CX18)
- Quick Release Rolko Rear Wheels (CX10/12/14/16)
- Foot Operated Wheel Lock (CX10/12/14/16)
- Zippered Pockets in Seat & Back Upholstery (Standard Upholstery: CX10/12/14/16)
- One Piece Height Adjustable Push Handle (CX10/12/14/16)
- Two Piece Push Handles (CX18)
- Tool-less Adjustable Seat Depth Growth of 5"

Standard Models:

- Cruiser CX10 - C903314
- Cruiser CX12 - C902845
- Cruiser CX14 - C900490
- Cruiser CX16 - C900145
- Cruiser CX18 - C902594

Price

- \$2,614**
- \$2,614**
- \$2,614**
- \$2,614**
- \$2,614**

Scout Standard Features:*

- Flip-Up Armrests
- 8" x 2" Front Ribbed, 12 1/2" Rear Pneumatic Knobby Tire
- Attendant Hand Brakes
- Anatomic Back
- Heavy Duty Reinforced Upholstery (Cordura and Textilene)

***Total for Scout Standard Features is \$1,038**

All prices in CAD

Cruiser Scout Models:

- Cruiser Scout CX12 - C904532 **\$3,652**
- Cruiser Scout CX14 - C904533 **\$3,652**
- Cruiser Scout CX16 - C904534 **\$3,652**
- Cruiser Scout CX18 - C904535 **\$3,652**

- | <u>Upholstery Colors</u> | <u>Price</u> |
|---|--------------|
| <input type="checkbox"/> Standard Upholstery | |
| <input type="checkbox"/> Panther Black | \$0 |
| <input type="checkbox"/> Royal Blue | \$0 |
| <input type="checkbox"/> Sassy Purple | \$0 |
| <input type="checkbox"/> Princess Pink | \$0 |
| <input type="checkbox"/> Forest Green | \$0 |
| <input type="checkbox"/> Cherry Red | \$0 |
| <input type="checkbox"/> Textilene Upholstery | |
| <input type="checkbox"/> Navy Blue | \$0 |
| <input type="checkbox"/> Violet Purple | \$0 |
| <input type="checkbox"/> Heavy Duty Reinforced Upholstery - Standard Upholstery ^{^*} | \$372 |
| <input type="checkbox"/> Panther Black | N/C |
| <input type="checkbox"/> Royal Blue | N/C |
| <input type="checkbox"/> Sassy Purple | N/C |
| <input type="checkbox"/> Heavy Duty Reinforced Upholstery - Textilene ^{^*} | \$372 |
| <input type="checkbox"/> Navy Blue | N/C |
| <input type="checkbox"/> Violet Purple | N/C |
| <input type="checkbox"/> Curved Standard Upholstery Style on Size 18 (Default is Square)* | \$585 |
- [^]**Not available with Planar Option.**
^{*}**Not available with Cushion or Seat Depth Extension 3".**

	CX10	CX12	CX14	CX16	CX18
C903854	C903854	C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851	C903851	C903851
C903759	C903759	C903759	C903759	C903759	C903759
C903857	C903857	C903857	C903857	C903857	C903857
C903855	C903855	C903855	C903855	C903855	C903855
C903464	C903464	C903464	C903464	C903464	C903464
C903465	C903465	C903465	C903465	C903465	C903465
C903322	C902853	C900443	C900139	C902601	C902601
C903854	C903854	C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851	C903851	C903851
C904485	C904486	C904487	C904488	C904489	C904489
C903464	C903464	C903464	C903464	C903464	C903464
C903465	C903465	C903465	C903465	C903465	C903465
C903465	C903465	C903465	C903465	C903465	C510445

- | <u>Frame Option</u> | <u>Price</u> |
|--|----------------|
| <input type="checkbox"/> Two Piece Push Handle ¹ | \$0 |
| <input type="checkbox"/> One Piece Push Handle on Size 18 (Default is 2 piece) ¹ | \$1,196 |
| <input type="checkbox"/> Attendant Hand Brake ¹ (Not Available with Two Piece Handle) | \$327 |
| <input type="checkbox"/> Curb Tipper (Elevates front of chair to go up curbs) ^{***1} | \$136 |
| <input type="checkbox"/> Caster Lock for 2" Wide Front Wheels ¹ | \$130 |
| <input type="checkbox"/> Rear Anti-Tip Tubes | \$160 |
- ¹**Manufacturer Installation Required.**
^{***}**Cannot be combined with Rear Anti-Tip Tubes but can still be ordered.**

	CX10	CX12	CX14	CX16	CX18
C904124	C904124	C904125	C904126	C904127	C904461
C904134	C904134	C904134	C904134	C904134	C904923
C904128	C904128	C904128	C904128	C904128	C903916
C902920	C902920	C902920	C902920	C902920	C902920
C904520	C904520	C904540	C904540	C904540	C510840

- | <u>Transit Required Items</u> | <u>Price</u> |
|--|--------------|
| <input type="checkbox"/> Transit Option ¹ | \$278 |
| Head Support - Select One Style Below: | |
| <input type="checkbox"/> Headrest Extension - Curved - Standard Upholstery (Not Available in Planar) | \$179 |
| <input type="checkbox"/> Headrest Extension - Rectangular - Textilene ^{**} | \$179 |
| <input type="checkbox"/> Contoured Firm Headrest (Planar Option Only/Hrdw Req) | \$338 |
| <input type="checkbox"/> Adj. Removable Hardware (Req for Contoured Firm Headrest) | \$342 |
| Trunk Positioning - Select Style Below: | |
| <input type="checkbox"/> 5-Point Positioning/Restraint Harness (CX10/12 only) ^{*1} | \$229 |
| <input type="checkbox"/> H-Harness with Padded Covers (3-Pt Positioning Belt Req) | \$136 |
| <input type="checkbox"/> 3-Pt Positioning Belt | \$67 |
| <input type="checkbox"/> Full Torso Support Vest for Standard Upholstery (Not Available in Planar) | \$366 |
| <input type="checkbox"/> Full Torso Support Vest for Textilene | \$366 |
| Foot Positioners - Select One Size Below: | |
| <input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair) | \$217 |
| <input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair) | \$217 |
| <input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair) | \$217 |
- ¹**Manufacturer Installation Required.**
^{*}**Crash Tested - WC19 and FMVSS213 compliant. Not available with Reinforced Upholstery, Anatomic Back, Cushion, Planar Option, or H-Harness/3-Pt Belt.**
^{**}**Required when ordering the Contoured firm headrest and flat firm backrest extension.**

	CX10	CX12	CX14	CX16	CX18
C903316	C902847	C900491	C900173	C902596	C902596
C903332	C902863	C900484	C900138	C902619	C902619
C904514	C904515	C904516	C904517	C904518	C904518
C904032	C904032	C904033	C904034	C904034	C904034
C904035	C904035	C904035	C904035	C904035	C904035
C903682	C903683	C903988	C903988	C903988	C903988
C903988	C903988	C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017	C904017	C904017
C903330	C902861	C900680	C900339	C902616	C902616
C905436	C905437	C905438	C905439	C902616	C902616
C903415	C903415	C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417	C903417	C903417

Optional Accessories for Transit Vehicles

Q'Straint Transit Lap Belt **Price \$257**

Planar Option*

Planar Seat Modification **Price \$477**
 Planar Back with Attachment (Incl. 2 Backrest Ext. w/ Hrdw.) **Price \$443**
 Cushion - Select One Cushion Below:
 Support - General Use Cushion (Solid Seat Insert Included)¹ **\$334**
 Position - Medial Thigh Support/Anti Thrust (Solid Seat Insert Included)¹ **\$646**
 Align - Lateral Pelvis/Lateral Thigh Support (Solid Seat Insert Included)¹ **\$646**
 Planar Head Support Options:
 Contoured Firm Headrest (Hardware Required) **\$338**
 Adj. Removable Hardware (Req for Contoured Firm Headrest) **\$342**
 Flat Firm Backrest Extension (Hardware Req) **\$316**
 Adj. Removable Hardware (Req for Flat Firm Backrest Ext) **\$310**
 Planar Hard Lateral Supports
 Planar Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)
 Small = 3"x 5" **\$189**
 Medium = 3.75"x 5.5" **\$189**
 Large = 4.75"x 6" **\$189**
 Planar Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)
 Small = 3"x 5" **\$189**
 Medium = 3.75"x 5.5" **\$189**
 Large = 4.75"x 6" **\$189**
 Planar Right Adj. Removable Hardware (Req for Lateral Support) **\$246**
 Planar Left Adj. Removable Hardware (Req for Lateral Support) **\$246**

¹Manufacturer Installation Required.

*Planar Option must be ordered with Planar Seat, Planar Back, Cushion, Planar Head Support, and Textilene Upholstery.

Scout Model Accessories

Scout Curb Tipper¹ **Price \$130**
 Upholstery Style - Choose One
 Heavy-Duty Reinforced Upholstery - Cordura **N/C**
 Heavy-Duty Reinforced Upholstery - Textilene **N/C**
¹Manufacturer Installation Required.

Head Supports

Padded Headwings[^] **Price \$338**
 Occi Headwings[^] **\$316**
 Headrest Extension - Curved - Standard Upholstery[^] **\$179**
 Headrest Extension - Rectangular - Textilene[^] **\$179**
[^]Not available with Planar Option.

Trunk Positioning

H-Harness with Padded Covers (3-Pt Positioning Belt Req) **Price \$136**
 Full Torso Support Vest for Standard Upholstery* **\$366**
 Full Torso Support Vest for Textilene* **\$366**
 Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")[^] **\$530**
[^]Not available with Planar Option.

Pelvic Positioning

3-Pt Positioning Belt **Price \$67**
 Support - General Use Cushion (Solid Seat Insert Included)^{*1} **\$334**
 Position - Medial Thigh Support/Anti Thrust (Solid Seat Insert Included)^{*1} **\$646**
 Align - Lateral Pelvis/Lateral Thigh Support (Solid Seat Insert Included)^{*1} **\$646**
 Incontinence Cushion Liner (Cushion Req)¹ **\$100**
 Seat Depth Extension 3"^{^1} **\$130**
 Seat Depth Extension 3" with Cushion **\$415**
 Standard Upholstery with Cushions (Default is Textilene) **\$585**
¹Manufacturer Installation Required.
[^]Not available with Planar Option.
^{*}Abd/Add Different Design.

Foot & Leg Positioning

Medial Thigh Support-Abductor¹ **Price \$225**
 Lateral Thigh Support-Adductor¹ **\$225**
 Calf Panel* **\$85**
 Angle Adjustable Footplates (pair) **\$221**
 Footplate Securement Strap **\$43**
 Foot Positioners - Small 7" - 10" (pair) **\$217**
 Foot Positioners - Medium 8" - 11" (pair) **\$217**
 Foot Positioners - Large 9" - 12" (pair) **\$217**
¹Manufacturer Installation Required.
^{*}Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.

Tire Options

12.5" Rear Solid Knobby Tire* **Price \$277**
 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire* **\$413**
 8" x 2" Front Ribbed, 12.5" Rear Solid Knobby Tire **\$413**
 8" x 2" Front Ribbed, 12.5" Rear Pneumatic Knobby Tire **\$413**

CX10	CX12	CX14	CX16	CX18
C903094	C903094	C903094	C903094	C903094

CX10	CX12	CX14	CX16	CX18
C904019	C904020	C904021	C904022	C904023
C904024	C904025	C904026	C904027	C904028
C903639	C903642	C903645	C903648	C903651
C903641	C903644	C903647	C903650	C903653
C903640	C903643	C903646	C903649	C903652
C904032	C904032	C904033	C904034	C904034
C904035	C904035	C904035	C904035	C904035
C904089	C904090	C904091	C904092	C904093
C904057	C904057	C904057	C904057	C904057
C902965	C902965	C902965	C902965	C902965
C902969	C902969	C902969	C902969	C902969
C902973	C902973	C902973	C902973	C902973
C904462	C904462	C904462	C904462	C904462
C904463	C904463	C904463	C904463	C904463
C904464	C904464	C904464	C904464	C904464
C904030	C904030	C904030	C904031	C904031
C904030	C904030	C904030	C904031	C904031

CX10	CX12	CX14	CX16	CX18
C903916	C903916	C903916	C903916	C903916
C902853B	C900443B	C900139B	C902601B	C902601B
C904486B	C904487B	C904488B	C904489B	C904489B

CX10	CX12	CX14	CX16	CX18
C903331	C902862	C900681	C900142	C902618
C905270	C905270	C905271	C905271	C905272
C903332	C902863	C900484	C900138	C902619
C904514	C904515	C904516	C904517	C904518

CX10	CX12	CX14	CX16	CX18
C903988	C903988	C903988	C903988	C903988
C903330	C902861	C900680	C900339	C902616
C905436	C905437	C905438	C905439	C902616
C904211	C904211	C904110	C904110	C904110

CX10	CX12	CX14	CX16	CX18
C904017	C904017	C904017	C904017	C904017
C903639	C903642	C903645	C903648	C903651
C903641	C903644	C903647	C903650	C903653
C903640	C903643	C903646	C903649	C903652
C902324	C902324	C902324	C902324	C902324
C903995	C903997	C903999	C904001	
C905606	C905606	C905606	C905606	
C904900	C904900	C904900	C904900	

CX10	CX12	CX14	CX16	CX18
C903326	C902857	C900480	C900137	C902607
C903327	C902858	C900482	C900188	C902608
C904580	C904580	C904580	C904580	C904580
C903445	C903268	C903268	C903269	C903269
C903584	C903584	C903584	C903584	C903584
C903415	C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417	C903417

CX10	CX12	CX14	CX16	CX18
C904117	C904117	C904117	C904117	
C904118	C904118	C904118	C904118	
				C903084
				C903083

<u>Accessories</u>	<u>Price</u>
<input type="checkbox"/> Headrest Cover (Canopy)	\$316
<input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy)	\$346
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray)^	\$154
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$315
<input type="checkbox"/> Under Seat Storage Basket - Fixed Tilt	\$158
<input type="checkbox"/> Utility Bag	\$58
<input type="checkbox"/> Saddle Bags (pair)	\$185
<input type="checkbox"/> Wheelchair Travel Bag	\$387
<input type="checkbox"/> Rain Cover Size 1 (Canopy Req)	\$209
<input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (CX14 Extended Canopy Req)	\$209
<input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req)	\$209
<input type="checkbox"/> Mosquito Net Size 1 (Canopy req)	\$209
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (CX14 Extended Canopy Req)	\$209
<input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req)	\$209

CX10	CX12	CX14	CX16	CX18
C903336	C902867	C900162	C900184	C902623
C904720	C904721	C904722	C904723	C904724
C903078	C903385	C903386	C903387	C903388
C903352	C903353	C903354	C903355	C903356
C904046	C904046	C904047	C904047	C904048
C904531	C904531	C904531	C904531	C904531
C904094	C904094	C904094	C904094	C904094
C902022	C902022	C902022	C902022	C902022
C904939	C904939	C904939		
		C904940	C904940	C904940
			C905283	C905283
C904655	C904655	C904655		
		C904656	C904656	C904656
			C905319	C905319

The Cruiser transit option includes four manufacturer installed wheelchair transportation anchors. The Cruiser transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are CX10 and CX12: 66 lbs./30kg; CX14: 100 lbs./45.5kg; CX16 and CX18: 170 lbs./77kg. The Convaid Cruiser was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convaid is unknown, and no warranty or claim about their performance is being, or can be, made.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Shipping Rates: Ground \$71.95 • For expedited shipping rates, please call 1-844-876-6245 or email convaidsales.us@etac.com
 Prices in CAD & subject to change, ex works Torrance, California.

Convaid
 2830 California Street
 Torrance, CA 90503
 1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

Effective: 2/16/2024

OF1902-ENG-REV.10

