

Convoid Cruiser Scout - HCPCS E1236

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Cruiser Scout Standard Features:

- Removable, Height Adjustable Swing-Away Footplates
- 8" x 2" Front Ribbed, 12 1/2" Rear Pneumatic Knobby Tire
- Foot Operated Wheel Lock (CX10/12/14/16)
- Zippered Pockets in Seat & Back Upholstery (Standard Upholstery Only: CX10/12/14/16)
- One Piece Height Adjustable Push Handle (CX10/12/14/16)
- Tool-less Adjustable Seat Depth Growth of 5"
- Tool-less Upholstery Removal
- 2-Point Positioning Belt
- 30° Fixed Tilt

- Silver Powder Coated Frame
- Self-Tensio®
- Toggle Wheel Locks (CX18)
- Two Piece Push Handles (CX18)
- Tension Adjustable Back
- Flip-Up Armrests
- Attendant Hand Brakes
- Heavy Duty Reinforced Upholstery
- Made in the USA

Cruiser Scout Transit Features:

- Transit Anchor Brackets
- Headrest Extension
- H-Harness w/Padded Covers
- 3-Pt Positioning Belt
- Foot Positioners

Standard Scout Models:

- Cruiser Scout CX12 - C904532
- Cruiser Scout CX14 - C904533
- Cruiser Scout CX16 - C904534
- Cruiser Scout CX18 - C904535

Price

- \$2,815**
- \$2,815**
- \$2,815**
- \$2,815**

Transit Included Scout Models:

- Cruiser Scout CX12T - C905502
- Cruiser Scout CX14T - C905503
- Cruiser Scout CX16T - C905504
- Cruiser Scout CX18T - C905505

Price

- \$3,625**
- \$3,625**
- \$3,625**
- \$3,625**

Upholstery Colors

- Standard Heavy Duty Upholstery
 - Panther Black
 - Royal Blue
 - Sassy Purple
- Textilene Heavy Duty Upholstery
 - Navy Blue
 - Violet Purple
- Curved Standard Upholstery Style on Size 18 (Default is Square)

Price

- \$0**
- \$0**
- \$0**
- \$0**
- \$0**
- \$0**
- \$451**

	CX12	CX14	CX16	CX18
	C903854	C903854	C903854	C903854
	C903850	C903850	C903850	C903850
	C903851	C903851	C903851	C903851
	C903464	C903464	C903464	C903464
	C903465	C903465	C903465	C903465
				C510445

Frame Option

- Two Piece Push Handle¹
- One Piece Push Handle on Size 18 (Default is 2 piece)¹
- Curb Tipper (*Elevates front of chair to go up curbs*)¹
- Caster Lock for 2" Wide Front Wheels¹

Price

- \$0**
- \$923**
- \$105**
- \$100**

HCPCS

K0073

	CX12	CX14	CX16	CX18
	C904125	C904126	C904127	C904461
				C904923
	C903916	C903916	C903916	C903916
	C902920	C902920	C902920	C902920

¹Manufacturer Installation Required.

Transit Required Items

- Transit Option¹
 - Head Support - Select Style Below:
 - Headrest Extension - Curved - Standard Upholstery
 - Headrest Extension - Rectangular - Textilene
 - Trunk Positioning - Select Style Below:
 - H-Harness with Padded Covers (3-Pt Positioning Belt Req)
 - 3-Pt Positioning Belt
 - Full Torso Support Vest for Standard Upholstery
 - Full Torso Support Vest for Textilene
 - Foot Positioners - Select Size Below:
 - Foot Positioners - Small 7" - 10" (pair)
 - Foot Positioners - Medium 8" - 11" (pair)
 - Foot Positioners - Large 9" - 12" (pair)

Price (CX# / CX#T)

- \$214 / \$0**
- \$138 / \$0**
- \$138 / \$0**
- \$105 / \$0**
- \$52 / \$0**
- \$282**
- \$282**
- \$167 / \$0**
- \$167 / \$0**
- \$167 / \$0**

HCPCS

E0966
E0966
E0960
E0978

	CX12	CX14	CX16	CX18
	C902847	C900491	C900173	C902596
	C902863	C900484	C900138	C902619
	C904515	C904516	C904517	C904518
	C903988	C903988	C903988	C903988
	C904017	C904017	C904017	C904017
	C902861	C900680	C900339	C902616
	C905437	C905438	C905439	C902616
	C903415	C903415	C903415	C903415
	C903416	C903416	C903416	C903416
	C903417	C903417	C903417	C903417

¹Manufacturer Installation Required.

Optional Accessories for Transit Vehicles

- Q'Straint Transit Lap Belt

Price

- \$198**

HCPCS

E0978

	CX12	CX14	CX16	CX18
	C903094	C903094	C903094	C903094

Head Supports

- Padded Headwings
- Occi Headwings
- Headrest Extension - Curved - Standard Upholstery
- Headrest Extension - Rectangular - Textilene

Price (CX# / CX#T)

- \$260**
- \$244**
- \$138 / \$0**
- \$138 / \$0**

HCPCS

E0955
E0955
E0966
E0966

	CX12	CX14	CX16	CX18
	C902862	C900681	C900142	C902618
	C905270	C905271	C905271	C905272
	C902863	C900484	C900138	C902619
	C904515	C904516	C904517	C904518

Trunk Positioning

- H-Harness with Padded Covers (3-Pt Positioning Belt Req)
- Full Torso Support Vest for Standard Upholstery
- Full Torso Support Vest for Textilene*
- Reducer Seat Insert (Reduces Seat Width 2" & Depth 1")

Price (CX# / CX#T)

- \$105 / \$0**
- \$282**
- \$282**
- \$409**

HCPCS

E0960

	CX12	CX14	CX16	CX18
	C903988	C903988	C903988	C903988
	C902861	C900680	C900339	C902616
	C905437	C905438	C905439	C902616
	C904211	C904110	C904110	C904110

<u>Pelvic Positioning</u>	<u>Price</u> <u>(CX# / CX#T)</u>	<u>HCPCS</u>	<u>CX12</u>	<u>CX14</u>	<u>CX16</u>	<u>CX18</u>
<input type="checkbox"/> 3-Pt Positioning Belt	\$52 / \$0	E0978	C904017	C904017	C904017	C904017
<input type="checkbox"/> Support - General Use Cushion (Solid Seat Insert Included)*1	\$257	E2601	C903642	C903645	C903648	C903651
<input type="checkbox"/> Position - Medial Thigh Support/Anti Thrust (Solid Seat Insert Included)*1	\$498	E2607	C903644	C903647	C903650	C903653
<input type="checkbox"/> Align - Lateral Pelvis/Lateral Thigh Support (Solid Seat Insert Included)*1	\$498	E2607	C903643	C903646	C903649	C903652
<input type="checkbox"/> Incontinence Cushion Liner (Cushion Req) ¹	\$77	E2619	C902324	C902324	C902324	C902324
<input type="checkbox"/> Custom Seat Depth Extension 3" ^{^1}	\$100	E1297	C903997	C903999	C904001	

'Manufacturer Installation Required.
****Abd/Add Different Design.***

<u>Foot & Leg Positioning</u>	<u>Price</u> <u>(CX# / CX#T)</u>	<u>HCPCS</u>	<u>CX12</u>	<u>CX14</u>	<u>CX16</u>	<u>CX18</u>
<input type="checkbox"/> Medial Thigh Support-Abductor ¹	\$173	E0957	C902857	C900480	C900137	C902607
<input type="checkbox"/> Lateral Thigh Support-Adductor ¹	\$173	E0956	C902858	C900482	C900188	C902608
<input type="checkbox"/> Calf Panel*	\$66	E0995	C904580	C904580	C904580	C904580
<input type="checkbox"/> Angle Adjustable Footplates (pair)	\$170	K0040	C903268	C903268	C903269	C903269
<input type="checkbox"/> Footplate Securement Strap	\$33		C903584	C903584	C903584	C903584
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	\$167 / \$0		C903415	C903415	C903415	C903415
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	\$167 / \$0		C903416	C903416	C903416	C903416
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	\$167 / \$0		C903417	C903417	C903417	C903417

'Manufacturer Installation Required.
****Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.***

<u>Tire Options</u>	<u>Price</u>	<u>HCPCS</u>	<u>CX12</u>	<u>CX14</u>	<u>CX16</u>	<u>CX18</u>
<input type="checkbox"/> 8" x 2" Front Ribbed, 12.5" Rear Solid Knobby Tire	\$124		C903224	C903224	C903224	C903224

<u>Accessories</u>	<u>Price</u>	<u>HCPCS</u>	<u>CX12</u>	<u>CX14</u>	<u>CX16</u>	<u>CX18</u>
<input type="checkbox"/> Headrest Cover (Canopy)	\$244		C902867	C900162	C900184	C902623
<input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy)	\$267		C904721	C904722	C904723	C904724
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray) [^]	\$119	E0950	C903385	C903386	C903387	C903388
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$243	E1028	C903353	C903354	C903355	C903356
<input type="checkbox"/> Under Seat Storage Basket - Fixed Tilt	\$122		C904046	C904047	C904047	C904048
<input type="checkbox"/> Utility Bag	\$45		C904531	C904531	C904531	C904531
<input type="checkbox"/> Saddle Bags (pair)	\$143		C904094	C904094	C904094	C904094
<input type="checkbox"/> Wheelchair Travel Bag	\$298		C902022	C902022	C902022	C902022
<input type="checkbox"/> Rain Cover Size 1 (Canopy Req)	\$161		C904939	C904939		
<input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (CX14 Extended Canopy Req)	\$161			C904940	C904940	C904940
<input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req)	\$161				C905283	C905283
<input type="checkbox"/> Mosquito Net Size 1 (Canopy req)	\$161		C904655	C904655		
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (CX14 Extended Canopy Req)	\$161			C904656	C904656	C904656
<input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req)	\$161				C905319	C905319

The Cruiser transit option includes four manufacturer installed wheelchair transportation anchors. The Cruiser transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are CX12: 66 lbs./30kg; CX14: 100 lbs./45.5kg; CX16 and CX18: 170 lbs./77kg.

The Convaid Cruiser was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convaid is unknown, and no warranty or claim about their performance is being, or can be, made.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convaidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.

Convaid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

Effective: 1/12/2024

OF1910-ENG-REV.01

