

Convaid EZRider - HCPCS E1236

Retail Price List / Order Form

Account #: _____
 Contact: _____ ATP/RTS: _____

Purchase Order #: _____ Date: _____
 Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

EZRider Standard Features:

- Silver Powder Coated Frame
- Removable, Height Adjustable Swing-Away Footplates
- Tool-less Adjustable Seat Depth Growth of 5"
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire (EZ12/14/16)
- 8" Front, 11" Rear Solid Tire (EZ18)
- Quick Release Rolko Rear Wheels (EZ12/14/16)
- One Piece Height Adjustable Push Handle (EZ12/14/16)
- Foot Operated Wheel Lock (EZ12/14/16)
- Zippered Pockets in Seat & Back Upholstery (Standard Upholstery Only: EZ12/14/16)

- Toggle Wheel Locks (EZ18)
- Two Piece Push Handles (EZ18)
- 2-Point Pelvic Belt
- 10° Fixed Tilt
- Self-Tensio®
- Tool-less Upholstery Removal
- Tension Adjustable Back
- Made in the USA

EZRider Transit Features:

- Transit Anchor Brackets
- Headrest Extension
- H-Harness with Padded Covers
- 3-Point Pelvic Belt
- Foot Positioners

Standard Models:

- EZRider EZ12 - C900860
- EZRider EZ14 - C900301
- EZRider EZ16 - C900996
- EZRider EZ18 - C900351

Price
\$1,990
\$1,990
\$1,990
\$1,990

Transit Included Models:

- EZRider Transit EZ12T - C904856
- EZRider Transit EZ14T - C904857
- EZRider Transit EZ16T - C904858
- EZRider Transit EZ18T - C904859

Price
\$2,731
\$2,731
\$2,731
\$2,731

Upholstery Options

- Convaid Upholstery**
 - Panther Black
 - Royal Blue
 - Sassy Purple
 - Princess Pink
 - Forest Green
 - Cherry Red
- Textilene Upholstery - (Required for Planar Option)**
 - Navy Blue
 - Violet Purple
- Heavy Duty Reinforced Upholstery - Convaid Upholstery[^]**
 - Panther Black
 - Royal Blue
 - Sassy Purple
- Heavy Duty Reinforced Upholstery - Textilene Upholstery[^]**
 - Navy Blue
 - Violet Purple

Price

\$0
 \$0
 \$0
 \$0
 \$0
 \$0
 \$0
 \$0
 \$0
 \$0
 \$287
 N/C
 N/C
 N/C
 \$287
 N/C
 N/C

[^]Not available with Planar Option.

Frame Options

- Two Piece Push Handle¹
- Custom One Piece Push Handle on Size 18 (Default is 2 piece)¹
- Height Adjustable Flip-Up Armrests
- Caster Lock for 2" Wide Front Wheels¹
- Rear Anti-Tip Tubes

Price
 \$0
 \$923
 \$244
 \$100
 \$123

¹Manufacturer Installation Required.

Transit Required Items[^]

- Transit Anchor Brackets¹
- Head Support - Select Style Below:**
 - Headrest Extension for Convaid Upholstery (Not available in Planar)
 - Headrest Extension for Textilene Upholstery^{**}
 - Contoured Firm Headrest (Planar Option Only/Hrdw Req)
 - Adj. Removable Hardware (Req for Contoured Firm Headrest)
- Trunk Positioning - Select Style Below:**
 - 5-Point Harness with Padded Covers (EZ12 only)^{*1}
 - H-Harness with Padded Covers (3-Point Pelvic Belt Required)
 - 3-Point Pelvic Belt
 - Full Torso Support Vest for Convaid Upholstery (Not available in Planar)
 - Full Torso Support Vest for Textilene Upholstery
- Foot Positioners - Select Size Below:**
 - Small 7" - 10"[^]
 - Medium 8" - 11" (Default if no size is selected)[^]
 - Large 9" - 12"[^]

Price
 \$214
 \$138
 \$138
 \$260
 \$263
 \$176
 \$105
 \$52
 \$282
 \$282
 \$167
 \$167
 \$167

¹Manufacturer Installation Required.

[^]Transit required items are included in base price when selecting EZ Rider Transit Model

^{*}Crash Tested - WC19 Compliant. Not available with Reinforced Upholstery, Cushion, Planar Option or H-Harness and 3-Point Pelvic Belt.

^{**}Required when ordering the Contoured firm headrest and flat firm backrest extension.

	EZ12	EZ14	EZ16	EZ18
C903854	C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851	C903851
C903759	C903759	C903759	C903759	C903759
C903857	C903857	C903857	C903857	C903857
C903855	C903855	C903855	C903855	C903855
C903464	C903464	C903464	C903464	C903464
C903465	C903465	C903465	C903465	C903465
C900835	C900321	C901031	C900392	C900392
C903854	C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851	C903851
C904494	C904495	C904496	C904497	C904497
C903464	C903464	C903464	C903464	C903464
C903465	C903465	C903465	C903465	C903465

HCPCS

	EZ12	EZ14	EZ16	EZ18
C904510	C904511	C904512	C904513	C904513
C904523	C904523	C904524	C904524	C904924
C902920	C902920	C902920	C902920	C902920
C900377	C900377	C900378	C900378	C900378

HCPCS

	EZ12	EZ14	EZ16	EZ18
C900861	C900535	C900997	C900536	C900536
C902863	C900484	C900138	C902619	C902619
C900850	C900685	C901029	C900390	C900390
C904032	C904033	C904034	C904034	C904034
C904035	C904035	C904035	C904035	C904035
C903702	C903988	C903988	C903988	C903988
C903988	C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017	C904017
C900848	C900338	C901020	C905501	C905501
C905437	C905438	C905439	C900387	C900387
C903415	C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417	C903417

		<u>Transit Option</u>	Price	HCPCS	EZ12	EZ14	EZ16	EZ18
<input type="checkbox"/>	Q'Straint Transit Lap Belt		\$198	E0978	C903094	C903094	C903094	C903094
		<u>Planar Option*</u>	Price	HCPCS	EZ12	EZ14	EZ16	EZ18
<input type="checkbox"/>	Planar Seat Modification ¹		\$368		C904049	C904050	C904051	C904052
<input type="checkbox"/>	Planar Back with Attachment (Incl. 1 Backrest Ext. w/ Hrdw.) ¹		\$341	E2291	C904053	C904054	C904055	C904056
<input type="checkbox"/>	Custom Standard Upholstery with Planar Cushions (Default is Textilene)		\$451		C904900	C904900	C904900	C904900
Cushion - select one cushion below:								
<input type="checkbox"/>	Support - General Use Cushion (Solid Seat Insert Included) ¹		\$257	E2601	C903684	C903687	C903690	C903693
<input type="checkbox"/>	Position - Medial Thigh Support/Anti Thrust (Solid Seat Insert Included) ¹		\$498	E2607	C903686	C903688	C903692	C903695
<input type="checkbox"/>	Align - Lateral Pelvis/Lateral Thigh Support (Solid Seat Insert Included) ¹		\$498	E2607	C903685	C903689	C903691	C903694
Planar Head Support Options - Must Select One								
<input type="checkbox"/>	Contoured Firm Headrest (Hardware Required) ¹		\$260	E0955	C904032	C904033	C904034	C904034
<input type="checkbox"/>	Adj. Removable Hardware (Required for Contoured Headrest) ¹		\$263	E1028	C904035	C904035	C904035	C904035
<input type="checkbox"/>	Flat Firm Backrest Extension (Hardware Required) ¹		\$244	E0955	C903928	C903582	C902151	C903674
<input type="checkbox"/>	Adj. Removable Hardware (Required for Flat Firm Backrest Ext.) ¹		\$239	E1028	C904057	C904057	C904057	C904057
Planar Hard Lateral Supports								
Planar Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)¹								
<input type="checkbox"/>	Small = 3"x 5"		\$146	E0956	C902965	C902965	C902965	C902965
<input type="checkbox"/>	Medium = 3.75"x 5.5"		\$146	E0956	C902969	C902969	C902969	C902969
<input type="checkbox"/>	Large = 4.75"x 6"		\$146	E0956	C902973	C902973	C902973	C902973
Planar Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Req)¹								
<input type="checkbox"/>	Small = 3"x 5"		\$146	E0956	C904462	C904462	C904462	C904462
<input type="checkbox"/>	Medium = 3.75"x 5.5"		\$146	E0956	C904463	C904463	C904463	C904463
<input type="checkbox"/>	Large = 4.75"x 6"		\$146	E0956	C904464	C904464	C904464	C904464
<input type="checkbox"/>	Planar Right Adj. Removable Hardware (Required for Lateral Support) ¹		\$190	E1028	C904030	C904030	C904031	C904031
<input type="checkbox"/>	Planar Left Adj. Removable Hardware (Required for Lateral Support) ¹		\$190	E1028	C904030	C904030	C904031	C904031

¹Manufacturer Installation Required.

*Planar Option must be ordered with Planar Seat, Planar Back, Cushion, Planar Head Support, and Textilene Upholstery.

		<u>Head Supports</u>	Price	HCPCS	EZ12	EZ14	EZ16	EZ18
<input type="checkbox"/>	Padded Headwings ^A		\$260	E0955	C900849	C900682	C900683	C900684
<input type="checkbox"/>	Occi Headwings ^A		\$244	E0955	C905271	C905271	C905271	C905272
<input type="checkbox"/>	Headrest Extension for Convaid Upholstery ^A		\$138	E0966	C902863	C900484	C900138	C902619
<input type="checkbox"/>	Headrest Extension for Textilene Upholstery		\$138	E0966	C900850	C900685	C901029	C900390

^ANot available with Planar Option.

		<u>Trunk Positioning</u>	Price	HCPCS	EZ12	EZ14	EZ16	EZ18
<input type="checkbox"/>	H-Harness with Padded Covers (3-Point Pelvic Belt Required)		\$105	E0960	C903988	C903988	C903988	C903988
<input type="checkbox"/>	Full Torso Support Vest for Convaid Upholstery		\$282	K0108	C900848	C900338	C901020	C905501
<input type="checkbox"/>	Full Torso Support Vest for Textilene Upholstery		\$282	K0108	C905437	C905438	C905439	C900387
<input type="checkbox"/>	Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1") ^A		\$409		C904211	C904110	C904110	C904110

^ANot available with Planar Option.

		<u>Pelvic Positioning</u>	Price	HCPCS	EZ12	EZ14	EZ16	EZ18
<input type="checkbox"/>	3-Point Pelvic Belt		\$52	E0978	C904017	C904017	C904017	C904017
<input type="checkbox"/>	Support Cushion, General Use ^{1*2}		\$257	E2601	C903684	C903687	C903690	C903693
<input type="checkbox"/>	Position Cushion, Medial Thigh Support/Anti Thrust ^{1*2}		\$498	E2607	C903686	C903688	C903692	C903695
<input type="checkbox"/>	Align Cushion, Lateral Pelvis/Lateral Thigh Support ^{1*2}		\$498	E2607	C903685	C903689	C903691	C903694
<input type="checkbox"/>	Incontinence Cushion Liner (Cushion Required) ¹		\$77	E2619	C902324	C902324	C902324	C902324
<input type="checkbox"/>	3" Seat Depth Extension ^{1^A}		\$100	E1297	C904003	C904005	C904007	C904009
<input type="checkbox"/>	Custom 3" Seat Depth Extension with Cushion ^{1^A}		\$320	E1297	C905606	C905606	C905606	C905606

¹Manufacturer Installation Required.

²Solid Seat Insert Included with all Cushions.

^ANot available with Planar Option.

		<u>Foot & Leg Positioning</u>	Price	HCPCS	EZ12	EZ14	EZ16	EZ18
<input type="checkbox"/>	Medial Thigh Support, Abductor		\$173	E0957	C900843	C900329	C901024	C900388
<input type="checkbox"/>	Lateral Thigh Support, Adductor		\$173	E0956	C900844	C900330	C901026	C900389
<input type="checkbox"/>	Calf Panel [*]		\$66	E0995	C904580	C904580	C904580	C904580
<input type="checkbox"/>	Angle Adjustable Footplates		\$170	K0040	C903268	C903268	C903269	C903269
<input type="checkbox"/>	Footplate Securement Strap		\$33		C903584	C903584	C903584	C903584
Foot Positioners - Select Size Below:								
<input type="checkbox"/>	Small 7" - 10"		\$167	K0108	C903415	C903415	C903415	C903415
<input type="checkbox"/>	Medium 8" - 11"		\$167	K0108	C903416	C903416	C903416	C903416
<input type="checkbox"/>	Large 9" - 12"		\$167	K0108	C903417	C903417	C903417	C903417

*Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.

		<u>Tire Options</u>	Price	HCPCS	EZ12	EZ14	EZ16	EZ18
<input type="checkbox"/>	12.5" x 2" Rear Solid Knobby Tires*		\$213		C904117	C904117	C904117	
<input type="checkbox"/>	7.5" x 2" Front, 12.5" x 2" Rear Pneumatic Knobby Tires		\$319		C904118	C904118	C904118	
<input type="checkbox"/>	8" x 2" Front Ribbed, 12.5" x 2" Rear Solid Knobby Tires		\$319					C903084
<input type="checkbox"/>	8" x 2" Front Ribbed, 12.5" x 2" Rear Pneumatic Knobby Tires		\$319					C903083

	<u>Accessories</u>	<u>Price</u>	<u>HCPCS</u>	<u>EZ12</u>	<u>EZ14</u>	<u>EZ16</u>	<u>EZ18</u>
<input type="checkbox"/>	Headrest Cover (Canopy)	\$244	K0108	C900853	C901210	C901205	C901211
<input type="checkbox"/>	Extended Headrest (Canopy)	\$267	K0108	C904721	C904725	C904726	C904727
<input type="checkbox"/>	Upper Extremity Support Surface (Clear Tray)^	\$119	E0950	C904525	C904526	C904527	C904528
<input type="checkbox"/>	<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$243	E1028	C903369	C903369	C903369	C903369
<input type="checkbox"/>	Upper Extremity Support Surface (Foam Tray)	\$162	E0950	C903669	C903670	C903670	C903671
<input type="checkbox"/>	Wheelchair Travel Bag	\$298		C902022	C902022	C902022	C902022
<input type="checkbox"/>	Utility Bag	\$45		C904531	C904531	C904531	C904531
<input type="checkbox"/>	Rain Cover Size 1 (Canopy Required)	\$161		C904939	C904939		
<input type="checkbox"/>	Rain Cover Size 2 (Canopy Required) (EZ14 Extended Canopy Required)	\$161			C904940	C904940	C904940
<input type="checkbox"/>	Rain Cover Size 3 (Extended Canopy Required)	\$161					C905283
<input type="checkbox"/>	Mosquito Net Size 1 (Canopy Required)	\$161		C904655	C904655	C904655	
<input type="checkbox"/>	Mosquito Net Size 2 (Canopy Required) (EZ14 Extended Canopy Required)	\$161			C904656	C904656	C904656
<input type="checkbox"/>	Mosquito Net Size 3 (Extended Canopy Required)	\$161					C905319

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

The EZ Rider transit option includes four manufacturer installed wheelchair transportation anchors. The EZ Rider transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are EZ12: 66 lbs./30kg.; EZ14: 100 lbs./45.5kg; EZ16 and EZ18: 170 lbs./77kg. The Convoid EZ Rider was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convoid is unknown, and no warranty or claim about their performance is being, or can be, made.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.

Effective: 3/22/2024

Convoid
 2830 California Street
 Torrance, CA 90503
 1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF1200-ENG-REV.09

