

Convoid Rodeo Tilt - HCPCS E1232/E1161

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Rodeo Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- One Piece Height Adjustable Push Handle
- One Piece Angle Adjustable Footplate
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire
- Quick Release Rolko Rear Tires
- Single Foot Activated Wheel Lock Bar

- Zippered Pockets in Seat Bottom Upholstery
- 2-Point Pelvic Belt
- Seat Depth Growth of 5"
- Elevating Legrest
- HCPCS code E1161 is only valid on size 16 (RD16)
- Tool-less Upholstery Removal
- Made in the USA

Transit Standard Features:

- Transit Anchor Brackets
- Headrest Extension
- H-Harness with Padded Covers
- 3-Point Pelvic Belt
- Foot Positioners

Standard Models:

- Rodeo RD10 - C904340
- Rodeo RD12 - C902689
- Rodeo RD14 - C902691
- Rodeo RD16 - C903134

Price

- \$2,638
- \$2,638
- \$2,638
- \$2,638

Transit Included Models:

- Rodeo Transit RD10T - C904860
- Rodeo Transit RD12T - C904861
- Rodeo Transit RD14T - C904862
- Rodeo Transit RD16T - C904863

Price

- \$3,449
- \$3,449
- \$3,449
- \$3,449

Upholstery Colors

- Panther Black
- Royal Blue
- Sassy Purple
- Princess Pink
- Forest Green
- Cherry Red

Price

- \$0
- \$0
- \$0
- \$0
- \$0
- \$0

Frame Options

- Recline Option¹
- Height Adjustable Flip-Up Armrests
- Attendant Hand Brakes¹
- Caster Lock for 2" Wide Front Tires¹
- Curb Tipper (Cannot combine with Rear Anti-Tip Tubes)
- Rear Anti-Tip Tubes (Cannot combine with Curb Tipper)

Price

- \$493
- \$244
- \$252
- \$100
- \$105
- \$123

HCPCS

- E1225
- E0973
- K0073
- E0971

¹Manufacturer Installation Required.

Transit Required Items^{A**}

- Transit Anchor Brackets^{1A}
- Headrest Extension^A
- H-Harness with Padded Covers^A
- 3-Point Pelvic Belt^A

Price

- \$214
- \$138
- \$105
- \$52

HCPCS

- K0108
- E0966
- E0960
- E0978

Foot Positioners^{**} - Select Size Below:

- Small 7"-10"
- Medium 8"-11"
- Large 9"-12"

- \$167
- \$167
- \$167

- K0108
- K0108
- K0108

¹Manufacturer Installation Required.

^ATransit required items are included in base price when selecting Rodeo Transit Model.

^{*}Crash Tested - WC19 Compliant.

Transit Options

- Q'Straint Transit Lap Belt

Price

- \$198

HCPCS

- E0978

Head Supports

- Padded Headwings
- Occi Headwings
- Headrest Extension

Price

- \$260
- \$244
- \$138

HCPCS

- E0955
- E0955
- E0966

Trunk Positioning

- H-Harness with Padded Covers (3-Point Pelvic Belt Required)

Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Required)

- X Small= 2.75"x 3.75"
- Small = 3"x 5"
- Medium = 3.75"x 5.5"
- Large = 4.75"x 6"

- \$126
- \$126
- \$126
- \$126

- E0956
- E0956
- E0956
- E0956

Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Required)

- X Small= 2.75"x 3.75"
- Small = 3"x 5"
- Medium = 3.75"x 5.5"
- Large = 4.75"x 6"

- \$126
- \$126
- \$126
- \$126

- E0956
- E0956
- E0956
- E0956

- Right S/A Adj. Hardware for Lateral Trunk Support

- \$254

- E1028

- Left S/A Adj. Hardware for Lateral Trunk Support

- \$254

- E1028

- Full Torso Support Vest

- \$282

- K0108

- Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")

- \$409

^{*}RD10 requires X Small laterals.

RD10	RD12	RD14	RD16
C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851
C903759	C903759	C903759	C903759
C903857	C903857	C903857	C903857
C903855	C903855	C903855	C903855

RD10	RD12	RD14	RD16
C904342	C904343	C904344	C904345
C904346	C904346	C904346	C904346
C904650	C904650	C904650	C904650
C905310	C905310	C905310	C905310
C904649	C904649	C904649	C904649

RD10	RD12	RD14	RD16
C904341	C903167	C903168	C903169
C903332	C902863	C900484	C900138
C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

RD10	RD12	RD14	RD16
C903094	C903094	C903094	C903094

RD10	RD12	RD14	RD16
C904365	C902673	C902653	C903147
C905270	C905270	C905271	C905271
C903332	C902863	C900484	C900138

RD10	RD12	RD14	RD16
C903988	C903988	C903988	C903988
C904385	C904385	C904385	C904385
C904071	C904071	C904071	C904071
C904072	C904072	C904072	C904072
C904073	C904073	C904073	C904073
C904465	C904465	C904465	C904465
C904466	C904466	C904466	C904466
C904467	C904467	C904467	C904467
C904468	C904468	C904468	C904468
C904362	C904362	C904362	C904362
C904363	C904363	C904363	C904363
C904364	C903261	C903260	C903259
C904211	C904211	C904110	C904110

<u>Pelvic Positioning</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 3-Point Pelvic Belt	\$52	E0978
<input type="checkbox"/> Support Cushion, General Use*1	\$257	E2601
<input type="checkbox"/> Position Cushion, Medial Thigh Support/Anti Thrust*1	\$498	E2607
<input type="checkbox"/> Align Cushion, Lateral Pelvis/Lateral Thigh Support*1	\$498	E2607
<input type="checkbox"/> Incontinence Cushion Liner (Cushion Required)*	\$77	E2619

***Manufacturer Installation Required.**
***Solid Seat Insert included with all Cushions.**

<u>Foot & Leg Positioning</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Medial Thigh Support, Abductor	\$173	E0957
<input type="checkbox"/> Lateral Thigh Support, Adductor	\$173	E0956
<input type="checkbox"/> Padded Elevating Legrest Cover	\$100	K0108
<input type="checkbox"/> 2" Legrest Extension*1	\$112	
<input type="checkbox"/> Calf Panel	\$66	E0995
<input type="checkbox"/> Individual Angle Adjustable Footplates*1	\$368	K0040
<input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required)	\$347	K0108
Foot Positioners - Select Size Below:		
<input type="checkbox"/> Small 7"-10"	\$167	K0108
<input type="checkbox"/> Medium 8"-11"	\$167	K0108
<input type="checkbox"/> Large 9"-12"	\$167	K0108

***Manufacturer installation required if model was made prior to June 2019.**
***2" Legrest Extension interferes with caster swivel when seat is in upright position.**
Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.

<u>Tire Options</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 12.5" x 2" Rear Solid Knobby Tires	\$213	
<input type="checkbox"/> 7.5" x 2" Front, 12.5" x 2" Rear Pneumatic Knobby Tires	\$319	

<u>Medical Support Accessories</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Medical Necessity Storage Basket	\$315	
<input type="checkbox"/> IV Pole (height adjustable and collapsible)*	\$204	K0105
<input type="checkbox"/> Adj. Swing Away Hardware (Required for IV Pole)	\$147	E1028
<input type="checkbox"/> Oxygen Tank Bag	\$187	E2208

***Anti-tip Tubes Required**

<u>Accessories</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Headrest Cover (Canopy)	\$244	K0108
<input type="checkbox"/> Extended Headrest(Canopy)	\$267	K0108
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray)*^	\$119	E0950
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$243	E1028
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray)^	\$162	E0950
<input type="checkbox"/> Under Seat Storage Basket	\$122	
<input type="checkbox"/> Utility Bag	\$45	
<input type="checkbox"/> Solid Side Panels	\$226	K0108
<input type="checkbox"/> Mesh Side Panels	\$159	K0108
<input type="checkbox"/> Rain Cover Size 1 (Canopy Required)	\$161	
<input type="checkbox"/> Rain Cover Size 2 (Canopy Required) (RD14 Extended Canopy Required)	\$161	
<input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Required)	\$161	
<input type="checkbox"/> Mosquito Net Size 1 (Canopy Required)	\$161	
<input type="checkbox"/> Mosquito Net Size 2 (Canopy Required) (RD14 Extended Canopy Required)	\$161	
<input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy Required)	\$161	

***Flip-Up Armrests and Hardware required.**
^Cannot be used during Transit.

RD10	RD12	RD14	RD16
C904017	C904017	C904017	C904017
C904347	C904350	C904353	C904356
C904348	C904351	C904354	C904357
C904349	C904352	C904355	C904358
C902324	C902324	C902324	C902324

RD10	RD12	RD14	RD16
C904359	C902667	C902647	C903140
C904360	C902668	C902648	C903141
C904619	C904619	C904619	C904619
C904095	C904095	C904096	C904097
C904361	C904361	C904361	C904361
C904957	C904957	C904958	C904959
C904953	C904954	C904955	C904956
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

RD10	RD12	RD14	RD16
C904117	C904117	C904117	C904117
C904118	C904118	C904118	C904118

RD10	RD12	RD14	RD16
C904476	C904476	C904477	C904477
C904566	C904566	C904566	C904566
C905612	C905612	C905612	C905612
C902870	C902870	C902870	C902870

RD10	RD12	RD14	RD16
C904366	C902679	C902659	C903152
C904728	C904729	C904730	C904731
C904367	C903405	C903401	C903402
C903369	C903369	C903369	C903369
C903669	C903669		
C904529	C904529	C904530	C904530
C904531	C904531	C904531	C904531
C904388	C902672	C902652	C903145
C904389	C902680	C902681	C903146
C904939	C904939	C904939	
		C904940	C904940
			C905283
C904655	C904655	C904655	
		C904656	C904656
			C905319

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo mobility base was tested with Rodeo seating module equipped with head support, H-Harness with padded Covers, 3-point restraint system, foot positioners, and 4 transit anchor brackets.

It is Convaid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specially designed for the wheelchair to be secured in a forward facing configuration and not in any other direction.
 For more information on RESNA WC-19 and WC-20 standards please visit <http://www.rercwts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convaidsales.us@etac.com

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166

Prices in USD & subject to change, ex works Torrance, California.

Effective: 3/22/2024

Convaid
 2830 California Street
 Torrance, CA 90503
 1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF0100-ENG-REV.10

