

Convoid Rodeo Tilt w/out Seating - HCPCS E1234

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Rodeo Tilt w/out Seating Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- Adjustable One-Piece Push Handle
- Elevating Legrest / One-Piece, Angle Adjustable Footplate
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tires

- Quick Release Rolko Rear Wheels
- Single Foot Activated Wheel Lock Bar
- Seat Depth Growth of 5"
- Made in the USA

Transit Standard Features:

- Transit Anchor Brackets
- 2-Point Positioning Belt

Standard Models:

- ☐ Rodeo w/out Seating RD10WOS - C904971
- ☐ Rodeo w/out Seating RD12WOS - C904972
- ☐ Rodeo w/out Seating RD14WOS - C904973
- ☐ Rodeo w/out Seating RD16WOS - C904974

Price

\$2,418
\$2,418
\$2,418
\$2,418

Standard Models:

- ☐ Rodeo Transit w/out Seating RD10TWOS - C904967
- ☐ Rodeo Transit w/out Seating RD12TWOS - C904968
- ☐ Rodeo Transit w/out Seating RD14TWOS - C904969
- ☐ Rodeo Transit w/out Seating RD16TWOS - C904970

Price

\$2,783
\$2,783
\$2,783
\$2,783

Seating Module

- ☐ Depth Adjustable Solid Seat Pan (Manufacturer Installation Required)

Price

\$161

HCPCS

E2293

Frame Option

- ☐ Recline Option (Manufacturer Installation Required)
☐ Height Adjustable Flip-Up Armrests (pair)
☐ Attendant Hand Brakes (pair)¹
☐ Caster Lock for 2" Wide Front Tires¹
☐ Curb Tipper
☐ Rear Anti-Tip Tubes (Pair)

Price

\$493
\$244
\$252
\$100
\$105
\$123

HCPCS

E1225
E0973

K0073

E0971

¹**Manufacturer Installation Required.**

Transit Required Items

- ☐ Transit Option (Manufacturer Installation Required)
☐ 2-Point Positioning Belt

Price (RD# / RD#T)

\$214 / \$0
\$48 / \$0

HCPCS

E0978

Transit Options

- ☐ Q'Straint Transit Lap Belt

Price

\$198

HCPCS

E0978

Foot & Leg Positioning

- ☐ Individual Angle Adjustable Footplates
☐ Full Padded Footbox (One-Piece Footplate Required)
☐ Legrest Extension - 2"*
☐ Foot Positioners - Small 7" - 10" (pair)
☐ Foot Positioners - Medium 8" - 11" (pair)
☐ Foot Positioners - Large 9" - 12" (pair)

Price (RD# / RD#T)

\$368
\$347
\$112
\$167 / \$0
\$167 / \$0
\$167 / \$0

HCPCS

K0040
K0108

K0108
K0108
K0108

***2" Legrest Extension interferes with caster swivel when seat is in upright position.**

Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.

Tire Options

- ☐ 12.5" Rear Solid Knobby Tire
☐ 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire

Price

\$213
\$319

HCPCS

Canopies

- ☐ Headrest Cover (Canopy)
☐ Extended Headrest Cover - No Windows (Canopy)
Canopy Upholstery (Please Choose One Below)
☐ Panther Black
☐ Royal Blue
☐ Sassy Purple
☐ Princess Pink
☐ Forest Green
☐ Cherry Red

Price

\$244
\$267

\$0
\$0
\$0
\$0
\$0
\$0

HCPCS

RD10	RD12	RD14	RD16
C509549	C509550	C509551	C509552

RD10	RD12	RD14	RD16
C904342	C904343	C904344	C904345
C904346	C904346	C904346	C904346
		C904370	C904371
C904650	C904650	C904650	C904650
C905310	C905310	C905310	C905310
C904649	C904649	C904649	C904649

RD10	RD12	RD14	RD16
C904341	C903167	C903168	C903169
C903931	C903931	C904559	C904559

RD10	RD12	RD14	RD16
C903094	C903094	C903094	C903094

RD10	RD12	RD14	RD16
	C904957	C904958	C904959
C904953	C904954	C904955	C904956
C904095	C904095	C904096	C904097
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

RD10	RD12	RD14	RD16
C904117	C904117	C904117	C904117
C904118	C904118	C904118	C904118

RD10	RD12	RD14	RD16
C904366	C902679	C902659	C903152
C904728	C904729	C904730	C904731
C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851
C903759	C903759	C903759	C903759
C903857	C903857	C903857	C903857
C903855	C903855	C903855	C903855

Accessories	Price	HCPCS
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray)^	\$119	E0950
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$243	E1028
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray)^	\$162	E0950
<input type="checkbox"/> Solid Side Panels	\$226	K0108
<input type="checkbox"/> Mesh Side Panels	\$159	K0108
<input type="checkbox"/> Medical Necessity Storage Basket	\$315	
<input type="checkbox"/> Under Seat Storage Basket - Tilt-in-Space	\$128	
<input type="checkbox"/> IV Pole (height adjustable and collapsible)^2	\$209	K0105
<input type="checkbox"/> Adj. Swing Away Hardware (Required for IV Pole)	\$147	K1028
<input type="checkbox"/> Oxygen Tank Bag	\$191	E2208
<input type="checkbox"/> Utility Bag	\$45	
<input type="checkbox"/> Rain Cover Size 1 (Canopy Req)	\$161	
<input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (RD14 Extended Canopy Req)	\$161	
<input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req)	\$161	
<input type="checkbox"/> Mosquito Net Size 1 (Canopy req)	\$161	
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (RD14 Extended Canopy Req)	\$161	
<input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req)	\$161	
^Flip-Up Armrests and Hardware required.		
^2Anti-tip Tubes Required.		

RD10	RD12	RD14	RD16
C904367	C903405	C903401	C903402
C903369	C903369	C903369	C903369
C903669	C903669		
C904388	C902672	C902652	C903145
C904389	C902680	C902681	C903146
C904476	C904476	C904477	C904477
C904529	C904529	C904530	C904530
C904566	C904566	C904566	C904566
C905612	C905612	C905612	C905612
C902870	C902870	C902870	C902870
C904531	C904531	C904531	C904531
C904939	C904939	C904939	
		C904940	C904940
			C905283
C904655	C904655	C904655	
		C904656	C904656
			C905319

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo Mobility Base was tested with the Rodeo Seating Module, equipped with a head support and wheelchair-anchored 2-point pelvic lap belt. Do not use the foot positioners on the Convaid Rodeo during transit in a motor vehicle. The performance and safety of accessories and/or other seating options not tested by Convaid is unknown, and no warranty or claim about their performance is or can being made.

It is Convaid’s position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specifically designed for the wheelchair to be secured in a forward facing configurartion and not in any other direction.
For more information on RESNA WC-19 and WC-20 standards please visit <http://www.ercwts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convaidsales.us@etac.com
Prices in USD & subject to change, ex works Torrance, California.

Convaid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

Effective: 1/10/2024

