

Convoid Rodeo Tilt w/out Seating - HCPCS E1234

Retail Price List / Order Form

Account #: _____

Purchase Order #: _____

Date: _____

Contact: _____ ATP/RTS: _____

Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

Rodeo Tilt w/out Seating Standard Features:

- | | |
|---|---|
| <ul style="list-style-type: none"> • 5° to 45° of Adjustable Tilt • Silver Powder Coated Frame • Adjustable One-Piece Push Handle • Elevating Legrest / One-Piece, Angle Adjustable Footplate • 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tires | <ul style="list-style-type: none"> • Quick Release Rolko Rear Wheels • Single Foot Activated Wheel Lock Bar • Seat Depth Growth of 5" • Made in the USA |
|---|---|

Standard Models:

	<u>Price</u>
<input type="checkbox"/> Rodeo w/out Seating RD10WOS - C904971	\$2,418
<input type="checkbox"/> Rodeo w/out Seating RD12WOS - C904972	\$2,418
<input type="checkbox"/> Rodeo w/out Seating RD14WOS - C904973	\$2,418
<input type="checkbox"/> Rodeo w/out Seating RD16WOS - C904974	\$2,418

<input type="checkbox"/> Seating Module Depth Adjustable Solid Seat Pan (Manufacturer Installation Required)	Price \$157	HCPCS E2293	RD10 C509549	RD12 C509550	RD14 C509551	RD16 C509552
--	------------------------------	------------------------------	------------------------	------------------------	------------------------	------------------------

<input type="checkbox"/> Frame Option Recline Option (Manufacturer Installation Required)	Price \$481	HCPCS E1225	RD10 C904342	RD12 C904343	RD14 C904344	RD16 C904345
<input type="checkbox"/> Height Adjustable Flip-Up Armrests (pair)	\$238	E0973	C904346	C904346	C904346	C904346
<input type="checkbox"/> Curb Tipper	\$102		C905310	C905310	C905310	C905310
<input type="checkbox"/> Rear Anti-Tip Tubes (Pair)	\$120	E0971	C904649	C904649	C904649	C904649

<input type="checkbox"/> Transit Required Items Transit Option (Manufacturer Installation Required)	Price \$209	HCPCS E0978	RD10 C904341	RD12 C903167	RD14 C903168	RD16 C903169
<input type="checkbox"/> 2-Point Positioning Belt	\$47		C903931	C903931	C904559	C904559

<input type="checkbox"/> Tire Options 12.5" Rear Solid Knobby Tire	Price \$208	HCPCS	RD10 C904117	RD12 C904117	RD14 C904117	RD16 C904117
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire	\$311		C904118	C904118	C904118	C904118

<input type="checkbox"/> Foot & Leg Positioning Individual Angle Adjustable Footplates	Price \$186	HCPCS K0040	RD10	RD12 C904957	RD14 C904958	RD16 C904959
<input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required)	\$339	K0108	C904953	C904954	C904955	C904956
<input type="checkbox"/> Legrest Extension - 2"	\$57		C904095	C904095	C904096	C904097
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	\$163	K0108	C903415	C903415	C903415	C903415
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	\$163	K0108	C903416	C903416	C903416	C903416
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	\$163	K0108	C903417	C903417	C903417	C903417

***2" Legrest Extension interferes with caster swivel when seat is in upright position.
Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.**

<input type="checkbox"/> Canopies Headrest Cover (Canopy)	Price \$238	HCPCS	RD10 C904366	RD12 C902679	RD14 C902659	RD16 C903152
<input type="checkbox"/> Extended Headrest Cover - No Windows (Canopy)	\$260		C904728	C904729	C904730	C904731
<input type="checkbox"/> Extended Headrest Cover - with Windows (Canopy)	\$307		C904744	C904745	C904746	C904747
Canopy Upholstery (Please Choose One Below)						
<input type="checkbox"/> Panther Black	\$0		C903854	C903854	C903854	C903854
<input type="checkbox"/> Royal Blue	\$0		C903850	C903850	C903850	C903850
<input type="checkbox"/> Sassy Purple	\$0		C903851	C903851	C903851	C903851
<input type="checkbox"/> Princess Pink	\$0		C903759	C903759	C903759	C903759
<input type="checkbox"/> Forest Green	\$0		C903857	C903857	C903857	C903857
<input type="checkbox"/> Cherry Red	\$0		C903855	C903855	C903855	C903855
<input type="checkbox"/> Chocolate Brown	\$0		C904189	C904189	C904189	C904189
<input type="checkbox"/> Aqua Blue	\$0		C904753	C904753	C904753	C904753
<input type="checkbox"/> Apple Green	\$0		C904754	C904754	C904754	C904754

<input type="checkbox"/> Transit Options Q'Straint Transit Lap Belt	Price \$94	HCPCS E0978	RD10 C903094	RD12 C903094	RD14 C903094	RD16 C903094
<input type="checkbox"/> Sure-Lok Transit Lap Belt	\$94	E0978	C903232	C903232	C903232	C903232
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Small^	\$158	E0950	C903669	C903669		
<input type="checkbox"/> Upper Extremity Support Surface (Foam Tray) Medium^	\$158	E0950			C903670	C903670

^Crash tested Foam Tray with vinyl cover. Flip-Up Armrests required.

Accessories	Price	HCPCS
<input type="checkbox"/> Upper Extremity Support Surface (Clear Tray)^	\$116	E0950
<input type="checkbox"/> Adj. Removable Hardware (Required for Clear Tray)	\$237	E1028
<input type="checkbox"/> Solid Side Panels	\$220	K0108
<input type="checkbox"/> Mesh Side Panels	\$155	K0108
<input type="checkbox"/> Medical Necessity Storage Basket	\$307	
<input type="checkbox"/> Under Seat Storage Basket - Tilt-in-Space	\$78	
<input type="checkbox"/> IV Pole (height adjustable and collapsible) ²	\$157	K0105
<input type="checkbox"/> Oxygen Tank Bag	\$186	E2208
<input type="checkbox"/> Utility Bag	\$40	
<input type="checkbox"/> Attendant Hand Brakes (pair) ¹	\$246	
<input type="checkbox"/> Caster Lock for 2" Wide Front Tires ¹	\$98	K0073
<input type="checkbox"/> Rain Cover Size 1 (Canopy Req)	\$157	
<input type="checkbox"/> Rain Cover Size 2 (Canopy Req) (RD14 Extended Canopy Req)	\$157	
<input type="checkbox"/> Rain Cover Size 3 (Extended Canopy Req)	\$157	
<input type="checkbox"/> Mosquito Net Size 1 (Canopy req)	\$157	
<input type="checkbox"/> Mosquito Net Size 2 (Canopy req) (RD14 Extended Canopy Req)	\$157	
<input type="checkbox"/> Mosquito Net Size 3 (Extended Canopy req)	\$157	

RD10	RD12	RD14	RD16
C904367	C903405	C903401	C903402
C903369	C903369	C903369	C903369
C904388	C902672	C902652	C903145
C904389	C902680	C902681	C903146
C904476	C904476	C904477	C904477
C904529	C904529	C904530	C904530
C904566	C904566	C904566	C904566
C902870	C902870	C902870	C902870
C904531	C904531	C904531	C904531
		C904370	C904371
C904650	C904650	C904650	C904650
C904939	C904939	C904939	
		C904940	C904940
			C905283
C904655	C904655	C904655	
		C904656	C904656
			C905319

¹Manufacturer Installation Required.
[^]Flip-Up Armrests and Hardware required.
²Anti-tip Tubes Required.

The Rodeo Transit option includes four manufacturer installed red wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are RD10T, RD12T 66lbs/30kg, RD14T 100lbs/45kg, and RD16T 140lbs/64kg.

The Convaid Rodeo Mobility Base was tested with the Rodeo Seating Module, equipped with a head support and wheelchair-anchored 2-point pelvic lap belt. Do not use the foot positioners on the Convaid Rodeo during transit in a motor vehicle. The performance and safety of accessories and/or other seating options not tested by Convaid is unknown, and no warranty or claim about their performance is or can being made.

It is Convaid's position that the wheelchair user should be transferred from the wheelchair into appropriate vehicle seating whenever possible. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

The WC19 standard is specifically designed for the wheelchair to be secured in a forward facing configuration and not in any other direction. For more information on RESNA WC-19 and WC-20 standards please visit <http://www.ercwts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convaidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.

Effective: 6/1/2022

Convaid
 2830 California Street
 Torrance, CA 90503
 1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF0109-ENG-REV.06

