

Convoid Trekker

Canadian Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Trekker Standard Features:

- Trekker Folding Mobility Frame Accepts TR12 & TR14 Seating Systems
- Tilt-in-Space (-10° to +45°)
- 180° Reversible Seating Configuration
- 90° Fixed Front Rigging
- 7.5" x 2" Front and 11.5" x 2.5" Rear Solid Tires
- One Piece Height Adjustable Push Handle (Inside Mount)
- Silver Powder Coated Frame
- 90° Fixed Footplate
- Quick Release Rear Wheels
- Standard Recline 80° to 110°
- Depth Adjustable Seat Pan
- Made in the USA

Standard Models:

- TR Seating System - C904768
- TR Mobility Frame - C904661

Price

\$0
\$3,352

All prices in CAD

	<u>Price</u>		<u>TR12</u>		<u>TR14</u>
<input type="checkbox"/> Seat Width 12"	N/C		C905515		
<input type="checkbox"/> Back Canes Only	N/C		C508446		
<input type="checkbox"/> TR12 Height Adjustable Solid Back 9" - 17"	\$735		C904662		
<input type="checkbox"/> Seat Width 14"	N/C				C905516
<input type="checkbox"/> Back Canes Only	N/C				C508447
<input type="checkbox"/> TR14 Height Adjustable Solid Back 17" - 26"	\$735				C904663
Back Additional Recline Option¹					
<input type="checkbox"/> Recline System (80 to 130 degrees) ¹	\$673		C905053		C905053
Snuggle Back Cushions:**					
<input type="checkbox"/> TR12 Short 12" x 18"	N/C		C904696		
<input type="checkbox"/> TR14 Tall 14" x 26"	N/C				C904699
Comfort Side Support Wings					
<input type="checkbox"/> Comfort Side Support Wing - Left	N/C		C904617		C904617
<input type="checkbox"/> Comfort Side Support Wing - Right	N/C		C904618		C904618
Extra Back Cushion Options					
<input type="checkbox"/> Lumbar Support	\$50		C904664		C904664
Seating Depth Position (Select One - See Chart on Bottom of Order Form):					
<input type="checkbox"/> Seat Depth Position 1 (TR12: 5.5"-8"/TR14: 8"-10")	N/C		C905354		C905354
<input type="checkbox"/> Seat Depth Position 2 (TR12: 7.5"-10"/TR14: 10"-12")	N/C		C905355		C905355
<input type="checkbox"/> Seat Depth Position 3 (TR12: 9.5"-12"/TR14: 12"-14")	N/C		C905356		C905356
<input type="checkbox"/> Seat Depth Position 4 (TR12: NA/TR14: 13"-15")	N/C				C905360
Planar Seat Cushion*²					
<input type="checkbox"/> Snuggle Planar Seat Short Cushion - Standard Fabric (TR12: 12x12) / (TR14: 14x14)	\$412		C905412		C904710
<input type="checkbox"/> Snuggle Planar Seat Short Cushion - Spacer Fabric (TR12: 12x12) / (TR14: 14x14)	\$412		C904708		C905425
Contoured Seat Cushion*²					
<input type="checkbox"/> Snuggle Contoured Seat Long Cushion - Standard Fabric (TR12: 12x14) / (TR14: 14x16)	\$454		C905415		C904783
<input type="checkbox"/> Snuggle Contoured Seat Long Cushion - Spacer Fabric (TR12: 12x14) / (TR14: 14x16)	\$454		C905429		C905430
Upholstery Color:					
<input type="checkbox"/> Panther Black	N/C		C903854		C903854
<input type="checkbox"/> Royal Blue	N/C		C903850		C903850
<input type="checkbox"/> Forest Green	N/C		C903857		C903857
<input type="checkbox"/> Sassy Purple	N/C		C903851		C903851
<input type="checkbox"/> Princess Pink	N/C		C903759		C903759
<input type="checkbox"/> Cherry Red	N/C		C903855		C903855

¹Manufacturer Installation Required.

²Fast tab is only added to Cushion if 3-Pt Belt is added. New Cushion is required if 3-Pt Belt is added after-market.

*All seat cushions must have fast tab for aftermarket 3-Pt Belt and H-Harness.

**Back cushion dimensions are seat pan to top of back pan.

	<u>Price</u>		<u>TR12</u>		<u>TR14</u>
<input type="checkbox"/> Height Adjustable Push Handle ¹ (Outside Mount)	N/C		C904688		C904688
<input type="checkbox"/> Caster Locks ¹	\$131		C904689		C904689
<input type="checkbox"/> Rear Anti-Tip Tubes	\$161		C904649		C904649
<input type="checkbox"/> Curb Tipper	\$137		C905310		C905310

¹Manufacturer Installation Required.

	<u>Price</u>		<u>TR12</u>		<u>TR14</u>
<input type="checkbox"/> Transit Anchor Brackets ^{1^*}	\$280		C904576		C904576
<input type="checkbox"/> 2-Point Pelvic Belt	\$63		C905386		C905386
<input type="checkbox"/> Using other manufacturer's Head Support	\$0		C995681		C995681

¹Manufacturer Installation Required.

[^]Transit required items are included in base price when selecting Trekker Transit Model.

*Requires Head Support and Pelvic Belt.

	<u>Price</u>		<u>TR12</u>		<u>TR14</u>
<input type="checkbox"/> Q'Straint Transit Lap Belt	\$259		C903094		C903094

<u>Head Supports (Required for Transit)*</u>	<u>Price</u>
<input type="checkbox"/> Integrated Adj Head Support Assembly (Option only o TR12 Seating System)	\$247
<input type="checkbox"/> R82 Adjustable Width Headrest, Black Cover - Short (Hardware Required)	\$427
<input type="checkbox"/> Turtle Neck Bar System, Short (Required)	\$322
<input type="checkbox"/> Headrest Universal Mount (Required)	\$76
R82 Anatomic Headrest with Black Cover (Hardware Required)	
<input type="checkbox"/> Small	\$357
<input type="checkbox"/> Medium	\$357
<input type="checkbox"/> Turtle Neck Bar System, Short (Required)	\$322
<input type="checkbox"/> Headrest Universal Mount (Required)	\$76

**The Transit model requires head support selection if other manufacturer's head support is not selected under "Transit Required Items."*

TR12	TR14
C904615	
C905287	C905287
C904948	C904948
C904946	C904946
C904950	C904950
C904951	C904951
C904948	C904948
C904946	C904946

<u>Trunk Positioning (Required for Transit)</u>	<u>Price</u>
<input type="checkbox"/> H-Harness w/ Padded Covers^^	\$137
Torso Vest by R82	
<input type="checkbox"/> LT Shoulder to RT Hip: 7.75" - 10.75" (Size XS)	\$280
<input type="checkbox"/> LT Shoulder to RT Hip: 10.75" - 12.5" (Size S)	\$280
<input type="checkbox"/> LT Shoulder to RT Hip: 12.5" - 16" (Size M)	\$280
<input type="checkbox"/> LT Shoulder to RT Hip: 16" - 18.5" (Size L)	\$280
<input type="checkbox"/> LT Shoulder to RT Hip: 18.5" - 20.5" (Size XL)	\$280
<input type="checkbox"/> LT Shoulder to RT Hip: 20.5" - 23.5" (Size XXL)	\$280
<input type="checkbox"/> Camlock Kit for Vest (Required)	N/C

Cross Vest by R82	
<input type="checkbox"/> LT Shoulder to RT Hip: 7.75" - 12.5" (Size XS)	\$266
<input type="checkbox"/> LT Shoulder to RT Hip: 10.25" - 13.5" (Size S)	\$266
<input type="checkbox"/> LT Shoulder to RT Hip: 13.5" - 17.75" (Size M)	\$266
<input type="checkbox"/> LT Shoulder to RT Hip: 17.75" - 20.75" (Size L)	\$266
<input type="checkbox"/> LT Shoulder to RT Hip: 20.75" - 24.5" (Size XL)	\$266
<input type="checkbox"/> Camlock Kit for Vest (Required)	N/C

^^Requires the Height Adjustable Solid Back System and 3-Pt Belt.

TR12	TR14
C904695	C904695
C905476	C905476
C905152	C905152
C905153	C905153
C905154	C905154
C905480	C905480
C905481	C905481
C905394	C905394
C905176	C905176
C905177	C905177
C905178	C905178
C905179	C905179
C905180	C905180
C905394	C905394

<u>Laterals</u>	<u>Price</u>
Fixed Laterals - Width, Height, Adj. Lateral Thoracic Support	
Left Lateral Padding Selection:	
<input type="checkbox"/> Size 1: 4" (D) x 4.5" (H)	\$204
<input type="checkbox"/> Size 2: 5" (D) x 4.5" (H)	\$204
Right Lateral Padding Selection:	
<input type="checkbox"/> Size 1: 4" (D) x 4.5" (H)	\$204
<input type="checkbox"/> Size 2: 5" (D) x 4.5" (H)	\$204
Swing-Away Laterals - Height Adj. Lateral Thoracic Support	
Left Lateral Padding Selection:	
<input type="checkbox"/> Size 1: 4" (D) x 2.75" (H)	\$165
<input type="checkbox"/> Size 2: 5.75" (D) x 4.25" (H)	\$165
<input type="checkbox"/> Swing Away Adj. Hardware (Required for S/A Laterals)	\$333
Right Lateral Padding Selection:	
<input type="checkbox"/> Size 1: 4" (D) x 2.75" (H)	\$165
<input type="checkbox"/> Size 2: 5.75" (D) x 4.25" (H)	\$165
<input type="checkbox"/> Swing Away Adj. Hardware (Required for S/A Laterals)	\$333
<input type="checkbox"/> 1" Thick Padded Lateral Support Cover	\$69

C904760	C904760
C904876	C904876
C904761	C904761
C904886	C904886
C904864	C904864
C904866	C904866
C904884	C904884
C904865	C904865
C904867	C904867
C904885	C904885
C905411	C905411

<u>Pelvic Positioning (Required for Transit)*</u>	<u>Price</u>
<input type="checkbox"/> 2-Point Pelvic Belt	\$63
<input type="checkbox"/> 3-Point Pelvic Belt	\$68
4-Point Pelvic Positioning Belt by R82 (Cannot be paired with 2-Point Pelvic Belt)*	
<input type="checkbox"/> Hip Width: 6.75" - 12.5" (Size XS)	\$189
<input type="checkbox"/> Hip Width: 8.25" - 13.75" (Size S)	\$189
<input type="checkbox"/> Hip Width: 9.75" - 16.25" (Size M)	\$189
<input type="checkbox"/> Hip Width: 12.5" - 17.25" (Size L)	\$189
<input type="checkbox"/> Hip Width: 16.25" - 20.5" (Size XL)	\$189

**The Transit model requires belt selection if 2-Point Belt is not selected under "Transit Required Items."*

TR12	TR14
C905386	C905386
C904017	C904017
C905148	C905148
C905149	C905149
C905150	C905150
C905151	C905151
C905475	C905475

<u>Hip Guides</u>	<u>Price</u>
Hip Guides without Armrest - Height Adj. Lateral Pelvic Support	
Left Size Selection:	
<input type="checkbox"/> Small - 4.5" - 5.5" (H)*	\$168
<input type="checkbox"/> Large - 5.25" - 6.25" (H)*	\$168
Right Size Selection:	
<input type="checkbox"/> Small - 4.5" - 5.5" (H)*	\$168
<input type="checkbox"/> Large - 5.25" - 6.25" (H)*	\$168
Hip Guides with Armrest - Height Adj. Lateral Pelvic Support (Cannot be used with Reversible Tray)	
Left Size Selection:	
<input type="checkbox"/> Small - 5"-6"(H)*	\$204
<input type="checkbox"/> Large - 7"-8" (H)*	\$204
Right Size Selection:	
<input type="checkbox"/> Small - 5"-6"(H)	\$204
<input type="checkbox"/> Large - 7"-8" (H)*	\$204
Medial Knee Abductor (Pommel)	
<input type="checkbox"/> Small (1.5" front width, 3.5" depth, 2" back width, 2.25" height)	\$227
<input type="checkbox"/> Medium (2.6" front width, 4.1" depth, 2" back width, 2.7" height)	\$227

**Height is reduced by 1" when used with Planar Seat Cushion or 2" if used with Contoured Seat Cushion.*

TR12	TR14
C904935	C905042
C905040	C904937
C904936	C905043
C905041	C904938
C904611	C904613
	C904927
C904612	C904614
	C904928
C904557	C904557
C904591	C904591

<u>Foot & Leg Positioning</u>	<u>Price</u>
<input type="checkbox"/> Elevating Leg Rest (90° - 180°) ¹	\$350
<input type="checkbox"/> Padded Elevating Legrest Cover	\$131
<input type="checkbox"/> Calf Panel (E.I.) H: 3.5"	\$86
<input type="checkbox"/> Calf Panel (Standard) H: 6"	\$86
<input type="checkbox"/> Individual Angle Adjustable Footplates (Elevating Leg Rest Required) ¹	\$482
<input type="checkbox"/> One Piece Angle Adjustable Flip Up Footplate (-15° - +15°)	\$155
<input type="checkbox"/> Full Padded Footbox (One-Piece Footplate Required)	\$454
Foot Positioners	
<input type="checkbox"/> Small 7" - 10"	\$219
<input type="checkbox"/> Medium 8" - 11"	\$219
Therafin Shoe Holders w/Padded Straps	
<input type="checkbox"/> Small	\$263
<input type="checkbox"/> Medium	\$263

¹Manufacturer Installation Required.

<u>Tire Options</u>	<u>Price</u>
<input type="checkbox"/> 12.5" x 2" Rear Solid Knobby Tires	\$279
<input type="checkbox"/> 7.5" x 2" Front, 12.5" X 2" Rear Pneumatic Knobby Tires	\$418

<u>Medical Support Accessories</u>	<u>Price</u>
<input type="checkbox"/> Medical Necessity Storage Basket	\$412
<input type="checkbox"/> IV Pole (height adjustable and collapsible) ^o	\$267
<input type="checkbox"/> Adjustable Swing-Away Hardware (Required for IV Pole)	\$192
<input type="checkbox"/> Vertical O2 Tank Holder* ¹	\$245
<input type="checkbox"/> LTV Ventilator Holder ¹	\$630
<input type="checkbox"/> Ventilator/Suction/Hard Tray	\$630
<input type="checkbox"/> Extended Ventilator/Suction/Hard Tray	\$630

¹Manufacturer Installation Required.

^oRequires Rear Anti-Tip Tubes.

*Cannot be used during transit.

<u>Accessories</u>	<u>Price</u>
<input type="checkbox"/> Headrest Cover (Canopy)	\$319
<input type="checkbox"/> Extended Headrest Cover - (Canopy)	\$350
<input type="checkbox"/> Upper Extremity Support Tray (Clear Tray)* [^]	\$156
<input type="checkbox"/> Adj Removable Hardware (Required for Clear Tray)	\$318
<input type="checkbox"/> Padded Tray Cover* (Clear Tray Required)	\$92
<input type="checkbox"/> Under Seat Storage Basket	\$160
<input type="checkbox"/> Utility Bag	\$59
<input type="checkbox"/> Rain Cover (Canopy Required)	\$211
<input type="checkbox"/> Mosquito Net (Canopy Required)	\$211

*Cannot be used during transit.

[^]Not available with arm rests.

<u>R82 High-Low:xo Frame Options</u>	<u>Price</u>
Manual Lifts	
<input type="checkbox"/> High-Low:xo Frame, Manual Lift, Indoor, Size 1	\$2,730
Push Brace for High-Low:xo - Indoor (will interfere with 170° recline)	
<input type="checkbox"/> Push Brace - Fixed Angle	\$224
<input type="checkbox"/> Push Brace, Angle Adjustable	\$539
<input type="checkbox"/> Stingray Seat Module - Interface Adaptor Kit for Trekker Frame*	N/C

*Required with Trekker Frame.

***The Trekker Transit option includes four manufacturer installed red wheelchair transportation anchors. The Trekker Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are TR12T 75lbs/34kg and TR14T 110lbs/50kg.

The Convoid Trekker Mobility Base was tested with the Trekker Seating Module equipped with a head support and wheelchair anchored pelvic lap belt. Do not use the foot positioners on the Convoid Trekker during transit in a motor vehicle. The performance and safety of accessories and/or other seating options not tested by Convoid is unknown, and no warranty or claim about their performance is or can being made.

Furthermore, when being transported in a vehicle it is Convoid's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convoid cannot recommend a transit system.

For more information on RESNA WC-19 and WC-20 standards please visit <http://www.rercwts.pitt.edu/WC19.html>.

Standard Shipping Rates: Ground \$71.95 • For expedited shipping rates, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 9/12/2024

Convoid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com



TR12	TR14
C905046	C905047
C904619	C904619
C905370	C905370
C904580	C904580
C905049	C905050
C904568	C904596
C904955	C904956
C903415	C903415
C903416	C903416
C904679	C904679
C904680	C904680

TR12	TR14
C904117	C904117
C904118	C904118

TR12	TR14
C904606	C904606
C904566	C904566
C905612	C905612
C904565	C904565
C904671	C904671
C904595	C904595
C905416	C905416

TR12	TR14
C904572	C904599
C904690	C904691
C904556	C904590
C905361	C905361
C904586	C904610
C904605	C904605
C904531	C904531
C904790	C904791
C904788	C904789

TR12	TR14
C905366	C905366
C905368	C905368
C905369	C905369
C905363	C905363

Seat Depth Adjustable Positions:

Seat Depth Growth				
Position:	1	2	3	4
TR12:	5.5" - 8"	7.5" - 10"	9.5" - 12"	N/A
TR14:	8" - 10"	10" - 12"	12" - 14"	13" - 15"



1
Most Forward Position



2
One Hole Position Rearward



3
Two Holes Position Rearward



4
Three Holes Position Rearward

Trekker Top Vests Belt Slot (Shoulder) Height from Seat Pan

	Minimum	Maximum
TR12 Integrated Headrest	9"	13"
TR12 Separate Headrest	13"	17"
TR14 Separate Headrest	17"	26"

Seating Module Specifications:**Convaid Trekker Seating Module Specifications**

Seat Width with Accessories	TR 12 (Size 1)	TR 14 (Size 2)
Trunk Width between Fixed Laterals (With Comfort Wings/Without Comfort Wings)	7.5" - 8" / 7.5" - 12"	9.5" - 10" / 9.5" - 14"
Trunk Width between Swing-Away Laterals (With Comfort Wings/Without Comfort Wings)	6" - 8" / 6" - 10"	7" - 9.5" / 7" - 14"
Hip Width between Hip Guides Without Armrests	5.5" - 11.5"	6.5" - 14"
Hip Width between Hip Guides With Armrests	5.5" - 11.5"	6.5" - 14"

Convaid Trekker Seat to Footplate Dimensions

Seat to Footplate Dimensions	TR 12 (Size 1)	TR 14 (Size 2)
Fixed 90° with one piece footplate	5.5" - 11.5"	5.5" - 11.5"
Fixed 90° with one piece footplate (flipped)	3.25" - 9.25"	3.25" - 9.25"
ELR with one piece footplate	4.75" - 10.75"	6.25" - 12.25"
ELR with one piece footplate (flipped)	2" - 8"	4.25" - 10.25"
ELR with individual footplates	5" - 11"	7" - 13"
ELR with individual footplates (flipped)	2" - 5"	4" - 7"