

Convaid Trekker - HCPCS E1234
Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To: _____ Ship To: _____
Address: _____ Address: _____
City/State/Zip: _____ City/State/Zip: _____
Phone: _____ Phone: _____
Email: _____
Trekker Standard Features:
• Trekker Folding Mobility Frame Accepts TR12 & TR14 Seating Systems
• Tilt-in-Space (-10° to +45°)
• 180° Reversible Seating Configuration
• 90° Fixed Front Rigging
• 7.5" x 2" Front and 11.5" x 2.5" Rear Solid Tires
• One Piece Height Adjustable Push Handle (Inside Mount)
Transit Standard Features:
• Transit Anchor Brackets
• 2-Point Pelvic Belt
• Head Support
Standard Models:
[] TR Seating System - C904768
[] TR Mobility Frame - C904661
Price
\$0
\$2,560
Transit Included Models:
[] TR Seating System - C904768
[] TR Mobility Frame Transit - C905362
Price
\$0
\$2,877

Seating System Price HCPCS
[] Seat Width 12" N/C
[] Back Canes Only N/C
[] TR12 Height Adjustable Solid Back 9" - 17" \$561 E2293
[] Seat Width 14" N/C
[] Back Canes Only N/C
[] TR14 Height Adjustable Solid Back 17" -26" \$561 E2293
Back Additional Recline Option
[] Recline Option (80 to 170 degrees)^1 \$514 E1226
Snuggle Back Cushions:**
[] TR12 Short 12" x 18" N/C
[] TR14 Tall 14" x 26" N/C
Comfort Side Support Wings
[] Comfort Side Support Wing - Left N/C
[] Comfort Side Support Wing - Right N/C
Extra Back Cushion Options
[] Lumbar Support \$38 K0108
Seating Depth Position (Select One - See Chart on Bottom of Order Form):
[] Seat Depth Position 1 (TR12: 5.5"-8"/TR14: 8"-10") N/C
[] Seat Depth Position 2 (TR12: 7.5"-10"/TR14: 10"-12") N/C
[] Seat Depth Position 3 (TR12: 9.5"-12"/TR14: 12"-14") N/C
[] Seat Depth Position 4 (TR12: NA/TR14: 13"-15") N/C
Planar Seat Cushion^2
[] Snuggle Planar Seat Short Cushion - Standard Fabric (TR12: 12x12) / (TR14: 14x14) \$315 E2292
[] Snuggle Planar Seat Short Cushion - Spacer Fabric (TR12: 12x12) / (TR14: 14x14) \$315 E2292
Contoured Seat Cushion^2
[] Snuggle Contoured Seat Long Cushion - Standard Fabric (TR12: 12x14) / (TR14: 14x16) \$347 E2294
[] Snuggle Contoured Seat Long Cushion - Spacer Fabric (TR12: 12x14) / (TR14: 14x16) \$347 E2294
Upholstery Color:
[] Panther Black N/C
[] Royal Blue N/C
[] Forest Green N/C
[] Sassy Purple N/C
[] Princess Pink N/C
[] Cherry Red N/C
^1Manufacturer Installation Required.
^2Fast tab is only added to Cushion if 3-Pt Belt is added. New Cushion is required if 3-Pt Belt is added after-market.
*All seat cushions must have fast tab for aftermarket 3-Pt Belt and H-Harness.
**Back cushion dimensions are seat pan to top of back pan.

Table with 2 columns: TR12, TR14. Rows include various model numbers like C905515, C508446, C904662, C905053, C904696, C904617, C904618, C904664, C905354, C905355, C905356, C905360, C905412, C904710, C904708, C905425, C905415, C904783, C905429, C905430, C903854, C903850, C903857, C903851, C903759, C903855.

Frame Options Price HCPCS
[] Height Adjustable Push Handle^1 (Outside Mount) N/C
[] Caster Locks^1 \$100 K0073
[] Rear Anti-Tip Tubes \$123 E0971
[] Curb Tipper \$105
^1Manufacturer Installation Required.

Table with 2 columns: TR12, TR14. Rows include model numbers C904688, C904689, C904649, C905310.

Transit Required Items Price HCPCS
[] Transit Anchor Brackets^1^* \$214 K0108
[] 2-Point Pelvic Belt \$48 E0978
[] Using other manufacturer's Head Support \$0
^1Manufacturer Installation Required.
^Transit required items are included in base price when selecting Trekker Transit Model
*Requires Head Support and Pelvic Belt.

Table with 2 columns: TR12, TR14. Rows include model numbers C904576, C905386, C995681.

Transit Option Price HCPCS
[] Q'Straint Transit Lap Belt \$198 E0978

Table with 2 columns: TR12, TR14. Rows include model number C903094.

	Head Supports (Required for Transit)*	Price	HCPCS
<input type="checkbox"/>	Integrated Adj Head Support Assembly (Option only on TR12 Seating System)	\$189	E0955
<input type="checkbox"/>	R82 Adjustable Width Headrest, Black Cover - Short (Hardware Required)	\$326	E0955
<input type="checkbox"/>	Turtle Neck Bar System, Short (Required)	\$246	E1028
<input type="checkbox"/>	Headrest Universal Mount (Required)	\$58	K0108
	R82 Anatomic Headrest with Black Cover (Hardware Required)		
<input type="checkbox"/>	Small	\$273	E0955
<input type="checkbox"/>	Medium	\$273	E0955
<input type="checkbox"/>	Turtle Neck Bar System, Short (Required)	\$246	E1028
<input type="checkbox"/>	Headrest Universal Mount (Required)	\$58	K0108

*The Transit model requires head support selection if other manufacturer's head support is not selected under "Transit Required Items."

TR12	TR14
C904615	
C905287	C905287
C904948	C904948
C904946	C904946
C904950	C904950
C904951	C904951
C904948	C904948
C904946	C904946

	Trunk Positioning	Price	HCPCS
<input type="checkbox"/>	H-Harness w/ Padded Covers^^	\$105	E0960
	Torso Vest by R82		
<input type="checkbox"/>	LT Shoulder to RT Hip: 7.75" - 10.75" (Size XS)	\$214	E0960
<input type="checkbox"/>	LT Shoulder to RT Hip: 10.75" - 12.5" (Size S)	\$214	E0960
<input type="checkbox"/>	LT Shoulder to RT Hip: 12.5" - 16" (Size M)	\$214	E0960
<input type="checkbox"/>	LT Shoulder to RT Hip: 16" - 18.5" (Size L)	\$214	E0960
<input type="checkbox"/>	LT Shoulder to RT Hip: 18.5" - 20.5" (Size XL)	\$214	E0960
<input type="checkbox"/>	LT Shoulder to RT Hip: 20.5" - 23.5" (Size XXL)	\$214	E0960
<input type="checkbox"/>	Camlock Kit for Vest (Required)	N/C	
	Cross Vest by R82		
<input type="checkbox"/>	LT Shoulder to RT Hip: 7.75" - 12.5" (Size XS)	\$203	E0960
<input type="checkbox"/>	LT Shoulder to RT Hip: 10.25" - 13.5" (Size S)	\$203	E0960
<input type="checkbox"/>	LT Shoulder to RT Hip: 13.5" - 17.75" (Size M)	\$203	E0960
<input type="checkbox"/>	LT Shoulder to RT Hip: 17.75" - 20.75" (Size L)	\$203	E0960
<input type="checkbox"/>	LT Shoulder to RT Hip: 20.75" - 24.5" (Size XL)	\$203	E0960
<input type="checkbox"/>	Camlock Kit for Vest (Required)	N/C	

^^^Requires the Height Adjustable Solid Back System and 3-Point Pelvic Belt.

TR12	TR14
C904695	C904695
C905476	C905476
C905152	C905152
C905153	C905153
C905154	C905154
C905480	C905480
C905481	C905481
C905394	C905394
C905176	C905176
C905177	C905177
C905178	C905178
C905179	C905179
C905180	C905180
C905394	C905394

	Laterals	Price	HCPCS
	Fixed Laterals - Width, Height, Adj. Lateral Thoracic Support		
	Left Lateral Padding Selection:		
<input type="checkbox"/>	Size 1: 4" (D) x 4.5" (H)	\$156	E0956
<input type="checkbox"/>	Size 2: 5" (D) x 4.5" (H)	\$156	E0956
	Right Lateral Padding Selection:		
<input type="checkbox"/>	Size 1: 4" (D) x 4.5" (H)	\$156	E0956
<input type="checkbox"/>	Size 2: 5" (D) x 4.5" (H)	\$156	E0956
	Swing-Away Laterals - Height Adj. Lateral Thoracic Support		
	Left Lateral Padding Selection:		
<input type="checkbox"/>	Size 1: 4" (D) x 2.75" (H)	\$126	E0956
<input type="checkbox"/>	Size 2: 5.75" (D) x 4.25" (H)	\$126	E0956
<input type="checkbox"/>	Swing Away Adj. Hardware (Required for S/A Laterals)	\$254	E1028
	Right Lateral Padding Selection:		
<input type="checkbox"/>	Size 1: 4" (D) x 2.75" (H)	\$126	E0956
<input type="checkbox"/>	Size 2: 5.75" (D) x 4.25" (H)	\$126	E0956
<input type="checkbox"/>	Swing Away Adj. Hardware (Required for S/A Laterals)	\$254	E1028
<input type="checkbox"/>	1" Thick Padded Lateral Support Cover	\$53	E0960

C904760	C904760
C904876	C904876
C904761	C904761
C904886	C904886
C904864	C904864
C904866	C904866
C904884	C904884
C904865	C904865
C904867	C904867
C904885	C904885
C905411	C905411

	Pelvic Positioning (Required for Transit)*	Price	HCPCS
<input type="checkbox"/>	2-Point Pelvic Belt	\$48	E0978
<input type="checkbox"/>	3-Point Pelvic Belt	\$52	E0978
	4-Point Pelvic Belt by R82 (Cannot be paired with 2-Point Pelvic Belt)*		
<input type="checkbox"/>	Hip Width: 6.75" - 12.5" (Size XS)	\$144	K0108
<input type="checkbox"/>	Hip Width: 8.25" - 13.75" (Size S)	\$144	K0108
<input type="checkbox"/>	Hip Width: 9.75" - 16.25" (Size M)	\$144	K0108
<input type="checkbox"/>	Hip Width: 12.5" - 17.25" (Size L)	\$144	K0108
<input type="checkbox"/>	Hip Width: 16.25" - 20.5" (Size XL)	\$144	K0108

*The Transit model requires belt selection if 2-Point Belt is not selected under "Transit Required Items."

TR12	TR14
C905386	C905386
C904017	C904017
C905148	C905148
C905149	C905149
C905150	C905150
C905151	C905151
C905475	C905475

	Hip Guides	Price	HCPCS
	Hip Guides without Armrest - Height Adj. Lateral Pelvic Support		
	Left Size Selection:		
<input type="checkbox"/>	Small - 4.5" - 5.5" (H)*	\$128	E0956
<input type="checkbox"/>	Large - 5.25" - 6.25" (H)*	\$128	E0956
	Right Size Selection:		
<input type="checkbox"/>	Small - 4.5" - 5.5" (H)*	\$128	E0956
<input type="checkbox"/>	Large - 5.25" - 6.25" (H)*	\$128	E0956
	Hip Guides with Armrest - Height Adj. Lateral Pelvic Support (Cannot be used with Reversible Tray)		
	Left Size Selection:		
<input type="checkbox"/>	Small - 5"-6"(H)*	\$156	E0956
<input type="checkbox"/>	Large - 7"-8" (H)*	\$156	E0956
	Right Size Selection:		
<input type="checkbox"/>	Small - 5"-6"(H)	\$156	E0956
<input type="checkbox"/>	Large - 7"-8" (H)*	\$156	E0956
	Medial Knee Abductor (Pommel)		
<input type="checkbox"/>	Small (1.5" front width, 3.5" depth, 2" back width, 2.25" height)	\$173	E0957
<input type="checkbox"/>	Medium (2.6" front width, 4.1" depth, 2" back width, 2.7" height)	\$173	E0957

*Height is reduced by 1" when used with Planar Seat Cushion or 2" if used with Contoured Seat Cushion.

TR12	TR14
C904935	C905042
C905040	C904937
C904936	C905043
C905041	C904938
C904611	C904613
	C904927
C904612	C904614
	C904928
C904557	C904557
C904591	C904591

	Foot & Leg Positioning	Price	HCPCS
<input type="checkbox"/>	Elevating Leg Rest (90° - 180°) ¹	\$267	K0990
<input type="checkbox"/>	Padded Elevating Legrest Cover	\$100	K0108
<input type="checkbox"/>	Calf Panel (E.I.) H: 3.5"	\$66	E0995
<input type="checkbox"/>	Calf Panel (Standard) H: 6"	\$66	E0995
<input type="checkbox"/>	Individual Angle Adjustable Footplates (Elevating Leg Rest Required) ¹	\$368	K0040
<input type="checkbox"/>	One Piece Angle Adjustable Flip Up Footplate (-15° - +15°)	\$118	K0040
<input type="checkbox"/>	Full Padded Footbox (One-Piece Footplate Required)	\$347	K0108
Foot Positioners			
<input type="checkbox"/>	Small 7" - 10"	\$167	K0108
<input type="checkbox"/>	Medium 8" - 11"	\$167	K0108
Therafin Shoe Holders w/Padded Straps			
<input type="checkbox"/>	Small	\$201	K0108
<input type="checkbox"/>	Medium	\$201	K0108

¹*Manufacturer Installation Required.*

	Tire Options	Price	HCPCS
<input type="checkbox"/>	12.5" x 2" Rear Solid Knobby Tires	\$213	
<input type="checkbox"/>	7.5" x 2" Front, 12.5" X 2" Rear Pneumatic Knobby Tires	\$319	

	Medical Support Accessories	Price	HCPCS
<input type="checkbox"/>	Medical Necessity Storage Basket	\$315	
<input type="checkbox"/>	IV Pole (Height Adjustable and Collapsible) [°]	\$204	K0105
<input type="checkbox"/>	Adj. Swing-Away Hardware (Required for IV Pole)	\$147	K1028
<input type="checkbox"/>	Vertical O2 Tank Holder ^{*1}	\$187	E2208
<input type="checkbox"/>	LTV Ventilator Holder ¹	\$481	E1029
<input type="checkbox"/>	Ventilator/Suction/Hard Tray	\$481	E1029
<input type="checkbox"/>	Extended Ventilator/Suction/Hard Tray	\$481	E1029

¹*Manufacturer Installation Required.*

[°]*Requires Rear Anti-Tip Tubes.*

**Cannot be used during transit.*

	Accessories	Price	HCPCS
<input type="checkbox"/>	Headrest Cover (Canopy)	\$244	K0108
<input type="checkbox"/>	Extended Headrest (Canopy)	\$267	K0108
<input type="checkbox"/>	Upper Extremity Support Tray (Clear Tray) ^{*^}	\$119	E0950
<input type="checkbox"/>	Adj Removable Hardware (Required for Clear Tray)	\$243	E1028
<input type="checkbox"/>	Padded Tray Cover* (Clear Tray Required)	\$70	K0108
<input type="checkbox"/>	Under Seat Storage Basket	\$122	
<input type="checkbox"/>	Utility Bag	\$45	
<input type="checkbox"/>	Rain Cover (Canopy Required)	\$161	
<input type="checkbox"/>	Mosquito Net (Canopy Required)	\$161	

**Cannot be used during transit.*

[^]*Not available with Hip Guides with Armrest.*

	R82 High-Low:xo Frame Options	Price	HCPCS
Manual Lifts			
<input type="checkbox"/>	High-Low:xo Frame, Manual Lift, Indoor, Size 1	\$2,085	
Push Brace for High-Low:xo - Indoor (will interfere with 170° recline)			
<input type="checkbox"/>	Push Brace - Fixed Angle	\$171	
<input type="checkbox"/>	Push Brace, Angle Adjustable	\$412	
<input type="checkbox"/>	Stingray Seat Module - Interface Adaptor Kit for Trekker Frame*	N/C	

**Required with Trekker Frame.*

***The Trekker Transit option includes four manufacturer installed red wheelchair transportation anchors. The Trekker Transit model has been crash tested and conforms to RESNA WC19 standards. Anthropomorphic Test Device weights are TR12T 75lbs/34kg and TR14T 110lbs/50kg.

The Convaid Trekker Mobility Base was tested with the Trekker Seating Module equipped with a head support and wheelchair anchored pelvic lap belt. Do not use the foot positioners on the Convaid Trekker during transit in a motor vehicle. The performance and safety of accessories and/or other seating options not tested by Convaid is unknown, and no warranty or claim about their performance is or can being made.

Furthermore, when being transported in a vehicle it is Convaid's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

For more information on RESNA WC-19 and WC-20 standards please visit <http://www.rercwts.pitt.edu/WC19.html>.

TR12	TR14
C905046	C905047
C904619	C904619
C905370	C905370
C904580	C904580
C905049	C905050
C904568	C904596
C904955	C904956
C903415	C903415
C903416	C903416
C904679	C904679
C904680	C904680

TR12	TR14
C904117	C904117
C904118	C904118

TR12	TR14
C904606	C904606
C904566	C904566
C905612	C905612
C904565	C904565
C904671	C904671
C904595	C904595
C905416	C905416

TR12	TR14
C904572	C904599
C904690	C904691
C904556	C904590
C905361	C905361
C904586	C904610
C904605	C904605
C904531	C904531
C904790	C904791
C904788	C904789

TR12	TR14
C905366	C905366
C905368	C905368
C905369	C905369
C905363	C905363

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245
Prices in USD & subject to change, ex works Torrance, California.

Convaid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com


Effective: 9/12/2024

OF1800-ENG-REV.11



Seat Depth Adjustable Positions:

Seat Depth Growth				
Postion:	1	2	3	4
TR12:	5.5" - 8"	7.5" - 10"	9.5" - 12"	N/A
TR14:	8" - 10"	10" - 12"	12" - 14"	13" - 15"



Trekker Top Vests Belt Slot (Shoulder) Height from Seat Pan

	Minimum	Maximum
TR12 Integrated Headrest	9"	13"
TR12 Separate Headrest	13"	17"
TR14 Separate Headrest	17"	26"

Seating Module Specifications:**Convaid Trekker Seating Module Specifications**

Seat Width with Accessories	TR 12 (Size 1)	TR 14 (Size 2)
Trunk Width between Fixed Laterals (With Comfort Wings/Without Comfort Wings)	7.5" - 8" / 7.5" - 12"	9.5" - 10" / 9.5" - 14"
Trunk Width between Swing-Away Laterals (With Comfort Wings/Without Comfort Wings)	6" - 8" / 6" - 10"	7" - 9.5" / 7" - 14"
Hip Width between Hip Guides Without Armrests	5.5" - 11.5"	6.5" - 14"
Hip Width between Hip Guides With Armrests	5.5" - 11.5"	6.5" - 14"

Convaid Trekker Seat to Footplate Dimensions

Seat to Footplate Dimensions	TR 12 (Size 1)	TR 14 (Size 2)
Fixed 90° with one piece footplate	5.5" - 11.5"	5.5" - 11.5"
Fixed 90° with one piece footplate (flipped)	3.25" - 9.25"	3.25" - 9.25"
ELR with one piece footplate	4.75" - 10.75"	6.25" - 12.25"
ELR with one piece footplate (flipped)	2" - 8"	4.25" - 10.25"
ELR with individual footplates	5" - 11"	7" - 13"
ELR with individual footplates (flipped)	2" - 5"	4" - 7"