

Company:

Account Number:

Quote:

ATP:

Marked For:

Associated Quote/Order#:

Ship To:

Address:

City:

Zip:

Email/Fax:


Date:

State:

Phone:

Effective: 8/24/2026

x:panda Shape - Advanced

	Infant	1	2
 Effective seat width	8.5" 220 mm	8.5"/10.25"/11.75" 220/260/300 mm	11"/12.5"/14.25"/15.75" 280/320/360/400 mm
 Effective seat depth	4.75-6.25" 120-160 mm	6.25-12.5" 160-320 mm	9.5-16.5" 240-420 mm
 Backrest height	10.75-14.25" 270-360 mm	11-14.5"/12.5-16.5"/14.25-18.5" 280-370/320-420/360-470mm	13.5-17.25"/15-19.25"/17-22" 340-440/380-490/430-560 mm
 Max user mass	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
 Max user mass in transportation	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
 Backrest angle	-10° - 30°	-10° - 30°	-10° - 30°
 Back width, top of backrest	7.5" 190 mm	7"/8.5"/10.25"/11.75"/13.5" 180/220/260/300/340 mm	7"/8.5"/10.25"/11.75"/13.5"/15" 180/220/260/300/340/380 mm
 Overall width	9.5" 240 mm	11.75"/15.75"/19.75" 300/400/500 mm	13.75"/18"/22"/26" 350/460/560/660 mm
 Mass device	13.2 lb 6 kg	14.3 lb 6.5 kg	24.3 lb 11 kg

Measurements taken with both legs in max abduction. Size 1-3 measurements are based on Advanced Seat.

Shape Base Options

Please see measurement guide for growth information. A selection must be made in each section below.

Seat Base

Includes base and plate.

Seat Base Size

Size 1 Advanced Seat (E2292)

\$824

Seat Base Size

Size 2 Advanced Seat (E2292)

\$824

Back Recline

Required Selection

Reclining Dynamic Back (E2398)

\$723

Reclining Adjustable Dynamic Back (E2398)

\$723



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Cushions and Covers

Seat Cushion

Seat Cushion - Split

Required when Quick Release Pommel is selected

Split Cushion, Size 1- W: 8.75" x L: 19.75" (901113)	N/C
Split Cushion, Size 1- W: 10.25" x L: 19.75" (901114)	N/C
Split Cushion, Size 1- W: 11.75" x L: 19.75" (901115)	N/C
Split Cushion, Size 2- W: 12.5" x L: 25.5" (901117)	N/C
Split Cushion, Size 2- W: 14.25" x L: 25.5" (901118)	N/C
Split Cushion, Size 2- W: 15.75" x L: 25.5" (901119)	N/C

Color

Required selection for cover color

Red (-50)	N/C
Grey (-60)	N/C
Black Wipeable (-70)	\$21

Upper Back Size

Required Selection. Select Upper Back Base then select Upper and Lower Back Cushion and Covers. Includes quick release detachment (#9050374).



Size 1 Upper Back

Size 1 Base Upper - H: 7" (905073) (E2291)	\$580
Total back height range 12.25" - 16.25"	
Size 1 Base Upper - H: 8.75" (905074) (E2291)	\$580
Total back height range 13.75" - 18.50"	

Size 2 Upper Back

Size 2 Base Upper - H: 7" (905073) (E2291)	\$580
Total back height range 13.5" - 17.25"	
Size 2 Base Upper - H: 8.75" (905074) (E2291)	\$580
Total back height range 15" - 19.75"	
Size 2 Base Upper - H: 10.75" (905075) (E2291)	\$580
Total back height range 17" - 21.75"	

Upper Back Cushion and Cover

Upper back cushion height must match Upper Back Base height

Upper Back Cushion

Upper Back Cushion - W: 7" x H: 7" (902103)	N/C
Upper Back Cushion - W: 8.75" x H: 7" (902104)	N/C
Upper Back Cushion - W: 10.25" x H: 7" (902105)	N/C
Upper Back Cushion - W: 8.75" x H: 8.75" (902106)	N/C
Upper Back Cushion - W: 10.25" x H: 8.75" (902107)	N/C
Upper Back Cushion - W: 11.75" x H: 8.75" (902108)	N/C
Upper Back Cushion - W: 13.5" x H: 8.75" (902109)	N/C
Upper Back Cushion - W: 10.25" x H: 10.75" (902110)	N/C
Upper Back Cushion - W: 11.75" x H: 10.75" (902111)	N/C
Upper Back Cushion - W: 13.5" x H: 10.75" (902112)	N/C
Upper Back Cushion - W: 15" x H: 10.75" (902113)	N/C

Upper Back Cover

Selection required

Red (-50)	N/C
Grey (-60)	N/C
Black Wipeable (-70)	\$21

Lower Back Cushion and Cover

The width measurements of the cushion are not actual measurements of the cushion. The number matches the lower back of the seating system. If the seating system is configured without hip and pelvis support, consider going up a size in width when choosing the lower back cushion.

Lower Back Cushion

Shape Size 1 Lower Back Cushion - W: 8.75" x H: 4.75" (902120)	N/C
Shape Size 1 Lower Back Cushion - W: 10.25" x H: 4.75" (902121)	N/C
Shape Size 1 Lower Back Cushion - W: 11.75" x H: 4.75" (902122)	N/C
Shape Size 2 Lower Back Cushion - W: 11" x H: 5.5" (902123)	N/C
Shape Size 2 Lower Back Cushion - W: 12.5" x H: 5.5" (902124)	N/C
Shape Size 2 Lower Back Cushion - W: 14.25" x H: 5.5" (902125)	N/C
Shape Size 2 Lower Back Cushion - W: 15.75" x H: 5.5" (902126)	N/C

Lower Back Cover

Selection Required

Red (-50)	N/C
Grey (-60)	N/C
Black Wipeable (-70)	\$21

Daily Adjustment Kit for Upper Back

Allows for toolless adjustment of back height. Ideal for frequent adjustment.



Daily Adjustment Kit - Upper Back (905081)	\$21
--	------



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Spacer Blocks

Spacer Block Holder without Cushion

Cushion width must match spacer block holder width

Spacer Block Holder, Size 1 - Back W: 7", 8.75", 10.25"	\$183
Spacer Block Holder, Size 2 - Back W: 8.75", 10.25", 11.75", 13.5"	\$183

Spacer Block Cushion Cover

Black (-80)	N/C
Black Wipeable (-70)	N/C

Spacer Block Cushion

Spacer Block Cushion Size 1 - W: 7" x H: 2.25" (903011)	N/C
Spacer Block Cushion Size 2 - W: 8.75" x H: 2.25" (903012)	N/C
Spacer Block Cushion Size 3 - W: 10.25" x H: 2.25" (903013)	N/C
Spacer Block Cushion Size 4 - W: 11.75" x H: 2.25" (903014)	N/C
Spacer Block Cushion Size 5 - W: 13.5" x H: 2.25" (903015)	N/C



Headrest and Hardware Options

Headrest Style

Includes Cushion. Select desired headrest in chart and choose the corresponding required hardware below.

Style	Size	Price	Type
Curved	Small	\$220	Curved, Small - Black Wipeable - W: 7", H: 5" (E0955)
			Curved, Small - Black - W: 7", H: 5" (E0955)
	Medium	\$240	Curved, Medium - Black Wipeable W: 10.2", H: 6.3" (E0955)
			Curved, Medium - Black W: 10.2", H: 6.3" (E0955)
Anatomic	Small	\$271	Anatomic, Small - Black Wipeable A: 5.5" x B: 5.5" x C: 9.5" (E0955)
			Anatomic, Small - Black A: 5.5" x B: 5.5" x C: 9.5" (E0955)
	Medium	\$281	Anatomic, Medium - Black Wipeable A: 6" x B: 5.5" x C: 10" (E0955)
			Anatomic, Medium - Black A: 6" x B: 5.5" x C: 10" (E0955)
	Large	\$291	Anatomic, Large - Black Wipeable A: 6.5" x B: 5.5" x C: 12.5" (E0955)
			Anatomic, Large - Black A: 6.5" x B: 5.5" x C: 12.5" (E0955)
Adjustable		\$332	Adjustable with Large Laterals - Black A: 4.75" x B: 6 C: max 9.5" (E0955)
		\$440	Adjustable with Large Laterals and Speakers - Black A: 4.75" x B: 6 C: max 9.5" (E0955)
		\$411	Adjustable with Four Laterals - Black A: 4.75" x B: max 9.5" C: max 8" (E0955)

Mounting Hardware for Curved and Anatomic Headrest

Required with Curved and Anatomic Head support

Headrest Mounting Hardware

Mounting Brackets, Standard (E1033) 5.5" of forward/backward adjustment	\$245
Mounting Brackets, Extended (E1033) 8.75" of forward/backward adjustment	\$255
Mounting Brackets, 3D Adjustment (E1033) 8.75" of forward/backward adjustment	\$266

Mounting Hardware for Adjustable Headrest

Required with Adjustable Head Support

Short Mounting Hardware (E1033) Suitable for upper back heights of 7", 8.75", 10.75" and 13.5"	\$250
Long Mounting Hardware (E1033) Suitable for upper back heights of 10.75" and 13.5"	\$250

Headrest Options



Daily Adjustment Kits

*Creates toolless adjustability, For angle adjustment of head support.
For use with Anatomic and Curved Headrests only.*

Daily Adjustment Kit Standard - Set of Two (905034)	\$32
Daily Adjustment Kit Extended - Set of Three (905035)	\$47

Additional Headrest Options

Headband System For use with Anatomic Headrest	\$250
Universal Headrest Mount Required if using a third party headrest	\$67



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Arm Supports

Arm Supports

Arm Rests

Cover not included, select separately

Arm Rests, Fixed W: 3" L: 10.25" H: 4" - 6.5" (8910616-1) (E0973)	\$339
Arm Rest, Fixed W: 3" L: 12.25" H: 4" - 8.75" (8910617-1) (E0973)	\$339

Arm Rest Cover

Black	\$79
Black Wipeable	\$79

Daily Adjustment Kit

For toolless height adjustment of arm supports

Height Adjustment Kit of Arm Supports (8910615-11)	\$21
--	------



Fixed Upper Arm Support

Includes Black Wipeable cover. Not compatible with Swing-Away Lateral Supports.

Fixed Upper Arm Support, Size 1 - L: 5.5" H: 7" (905017-70) (K0108)	\$260
<i>For size 1 Base</i>	
Fixed Upper Arm Support Size 2 - L: 7.5" H: 10.75" (905018-70) (K0108)	\$260
<i>For size 2 Base</i>	

Upper Arm Supports - Adjustable

Not compatible with Swing-Away Lateral Supports

Adjustable Upper Arm Support, Size 1 - L: 5.5" H: 7" (905019-70) (K0108)	\$349
<i>For Size 1 Base</i>	
Adjustable Upper Arm Support, Size 2 - L: 7.5" H: 10.75" (905020-70) (K0108)	\$349
<i>For size 2 Base</i>	

Knee and Postural Supports

Abduction and Adduction

Inner Knee Support



Inner Knee Support For Size 1 and 2 Black - L: 7.5", H: 7.5" (905001-80) (E0957)	\$302
Inner Knee Support For Size 1 and 2 Black Wipeable - L: 7.5", H: 7.5" (905001-70) (E0957)	\$302
Inner Knee Support For Size 2 Black - L: 9.5", H: 9.5" (905002-80) (E0957)	\$302
Inner Knee Support For Size 2 Black Wipeable - L: 9.5", H: 9.5" (905002-70) (E0957)	\$302

Quick Release Pommel

Quick Release Pommel For Size 1 Small, Black - W: 2.25", L: 6.25", H: 7.5" (905077-80) (E0957)	\$364
Quick Release Pommel For Size 1 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" (905077-70) (E0957)	\$364
Quick Release Pommel For Size 2 Small, Black - W: 2.25", L: 6.25", H: 7.5" (905078-80) (E0957)	\$364
Quick Release Pommel For Size 2 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" (905078-70) (E0957)	\$364
Quick Release Pommel For Size 2 Large, Black - W: 3.25", L: 7.5", H: 9" (905082-80) (E0957)	\$364
Quick Release Pommel For Size 2 Large, Black Wipeable - W: 3.25", L: 7.5", H: 9" (905082-70) (E0957)	\$364

Outer Knee Support



Outer Knee Support, Standard Set for Size 1, Black - L: 7.5", H: 7.5" (905008-80) (E0956)	\$219
Outer Knee Support, Standard Set for Size 1, Black Wipeable - L: 7", H: 7.5" (905008-70) (E0956)	\$219
Outer Knee Support, Standard Set for Size 2, Black - L: 9.5", H: 9.5" (905009-80) (E0956)	\$219
Outer Knee Support, Standard Set for Size 2, Black Wipeable - L: 9.5", H: 9.5" (905009-70) (E0956)	\$219
Outer Knee Support, Flip Down Set for Size 2, Black - L: 7.75", H: 7.75" (905010-80) (E0956)	\$260
Outer Knee Support, Flip Down Set for Size 2, Black Wipeable - L: 7.75", H: 7.75" (905010-70) (E0956)	\$260



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

x:panda shape advanced seat

Postural Supports

Width Adjustable Pelvis Supports

Mounted from the back



Width Adjustable Pelvis Supports, Small Black - L: 4", H: 4.25" (905041-80) (E0956)	\$245
Width Adjustable Pelvis Supports, Small Black Wipeable - L: 4", H: 4.25" (905041-70) (E0956)	\$245
Width Adjustable Pelvis Supports, Large Black - L: 5.5", H: 5" (905042-80) (E0956)	\$245
Width Adjustable Pelvis Supports, Large Black Wipeable - L: 5.5", H: 5" (905042-70) (E0956)	\$245

Width Adjustable Hip Supports

Mounted from the sides



Width Adjustable Hip Supports, Size 1, Black - L: 4" H: 4.25" (905044-80) (E0956)	\$349
Width Adjustable Hip Supports, Size 1, Black Wipeable - L: 4" H: 4.25" (905044-70) (E0956)	\$349
Width Adjustable Hip Supports, Size 2, Black - L: 5.5" H: 5" (905045-80) (E0956)	\$349
Width Adjustable Hip Supports, Size 2, Black Wipeable - L: 5.5", H: 5" (905045-70) (E0956)	\$349

Adjustable Anti-Thrust Wedges

Adjustable Anti-Thrust Wedges, Short - W: 4.5", L: 2.75-4.25", H: .75" (905052) (K0108)	\$115
<i>For Size 1</i>	

Adjustable Anti-Thrust Wedges, Medium - W: 5.25", L: 4.25-6.75", H: 1.25" (905053) (K0108)	\$115
<i>For Size 1</i>	

Adjustable Anti-Thrust Wedges, Long - W: 5.25", L: 4.75-7.75", H: 1.75" (905054) (K0108)	\$115
<i>For Size 1</i>	

Adjustable Anti-Thrust Wedges, Medium - W: 6", L: 5 -8.75", H: 1.25" (905055) (K0108)	\$115
<i>For Size 2</i>	

Adjustable Anti-Thrust Wedges, Long - W: 6", L: 5.50 -9.5", H: 2" (905056) (K0108)	\$115
<i>For Size 2</i>	

Trays and Accessories



Shape 5 Tray with Edge Size 1 (8910725-1) (E0950)	\$334
<i>For Size 1. Fits Size 2 widths 11", 12.5", 14.25". Double Rod, Grey ABS.</i>	



Shape 6 Tray with Edge Size 2 (8910725-2) (E0950)	\$334
<i>For Size 2. Double rod, Grey ABS.</i>	



Shape 8 Tray with Edge Size 1 (8976350) (E0950)	\$334
---	-------

Padding for Size 1 Tray, Shape 5 Plate, Over and Under (8910725-10) (E1028)	\$443
Padding for Size 1 Tray, Top, Shape 5 Plate, Over (8910725-10T) (E1028)	\$312
Padding for Size 2 Tray, Shape 6 Plate, Over and Under (8910725-20) (E1028)	\$443
Padding for Size 2 Tray, Top, Shape 6 Plate, Over (8910725-20T) (E1028)	\$312
Arm Pads for Tray, for Shape 5 and 6 (8910730)	\$315
Toy Bar (8910740)	\$250

Accessory Holders and Brackets

Remote Control Holder for Size 1 and 2 - Arm Support Mounted (905030-1)	\$68
Remote Control Holder for Size 3 - Arm Support Mounted (905030-2)	\$68
Remote Control Holder - Back Mounted (905030)	\$131

Mounting Hardware

Bracket for Communication Aid, Size 1 and 2 - Arm Support Mounted (905033-1)	\$82
Bracket for Communication Aid, Size 3 - Arm Support Mounted (905033-2)	\$82

Canopy

Canopy does not come with the brackets. Both the Canopy and brackets must be quoted when mounting to shape.

Black Canopy - Size 1 (8740110)	\$339
Black Canopy - Size 2 (8740120)	\$339

Canopy Brackets

Canopy Brackets with Push Handle for Size 1 (905039-1)	\$203
Canopy Brackets with Push Handle for Size 2 (905039-2)	\$203



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Vest and Belts

More options can be found at etac.com

Vests



Torso Vests

Torso Vest Small - S-P: 10.75" - 12.5" (85705-S) (E0960)	\$232
Torso Vest Medium - S-P: 12.5" - 16" (85705-M) (E0960)	\$232
Torso Vest Large - S-P: 16" - 18.5" (85705-L) (E0960)	\$232
Torso Vest X-Large - S-P: 18.5" - 20.5" (85705-XL) (E0960)	\$232



Cross Vests

Cross Vest X-Small - S-P: 7.75" - 12.5" (85720-XS) (E0960)	\$220
Cross Vest Small - S-P: 10.25" - 13.5" (85720-S) (E0960)	\$220
Cross Vest Medium - S-P: 13.5" - 17.75" (85720-M) (E0960)	\$220
Cross Vest Large - S-P: 17.75" - 20.75" (85720-L) (E0960)	\$220
Cross Vest X-Large - S-P: 20.75" - 24.5" (85720-XL) (E0960)	\$220



Fixlocks Bracket

Higher Brackets (8910260)	\$77
<i>Bracket for higher strap guidance. Includes set of fixlocks. Requires the wider brackets as well.</i>	
Wider Brackets (905068)	\$57

Pelvic Belts and Harnesses



4-Point Belts

4-Point Padded Hip Belts 1" X-Small - H-H: 6.75" - 12.5" (85725-XS) (E0960)	\$156
4-Point Padded Hip Belts 1" Small - H-H: 8.25" - 13.75" (85725-S) (E0960)	\$156
4-Point Padded Hip Belts 1" Medium - H-H: 9.75" - 16.25" (85725-M) (E0960)	\$156
4-Point Padded Hip Belts 1" Large - H-H: 12.5" - 17.25" (85725-L) (E0960)	\$156



Pelvic Harness

Pelvic Harness, 0 - A:6.5" x B:7.5" (89108-0)	\$261
Pelvic Harness, 1 - A:8.5" x B:8.75" (89108-1)	\$261
Pelvic Harness, 2 - A:10.5" x B:9.75" (89108-2)	\$261



2-Point Pelvic Belt

2-Point Pelvic Belt Black - L: 13-14.5" (85424-1) (E0978)	\$98
---	------



Pelvis Belt

Pelvis Belt, X-Small - H-H: 6.25" - 7.75" (85735-XS) (E0978)	\$175
Pelvis Belt, Small - H-H: 7.75" - 9.5" (85735-S) (E0978)	\$175
Pelvis Belt, Medium - H-H: 9.5" - 11" (85735-M) (E0978)	\$175

Lateral Support Cushion and Fittings

Lateral Supports

If choosing Laterals, both the Left and Right cushions must be selected.

Left Cushion

Left Rectangular Small - L: 4.25" x H: 3.25" (E0956)	\$115
Left Rectangular Medium - L: 6" x H: 3.5" (E0956)	\$115
Left Rectangular Large - L: 7.5" x H: 3.5" (E0956)	\$115

Right Cushion

Right Rectangular Small - L: 4.25" x H: 3.25" (E0956)	\$115
Right Rectangular Medium - L: 6" x H: 3.5" (E0956)	\$115
Right Rectangular Large - L: 7.5" x H: 3.5" (E0956)	\$115

Type of Fitting

Selection required, fittings are sold as a pair.

Fixed Fitting (E1034)	\$172
Swing-Away and Angle Adjustable Fitting (E1034)	\$260



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

x:panda shape advanced seat

Footrests for Mobility Bases

Foot plate(s) attach to each Base. If multiple bases are requested, add additional footrest for each. Choose appropriate base(s) below.

Footrests

Select footrest type and desired mounting bar

Footrest Type

Angle and Height Adjustable Foot Plate, Small, 12" x 7.5" (95145-1) (K0040)	\$325
Angle and Height Adjustable Foot Plate, Medium, 14" x 9" (95145-2)(K0040)	\$325
Angle and Height Adjustable Foot Plate, Large, 16" x 9.5" (95145-3) (K0040)	\$325

Footrest Mounting Bar

Required for Angle and Height Adjustable Footplates

Short Mounting Bar - 6-10.5" (8910821-1)	N/C
Medium Mounting Bar - 6-14.5" (8910821-2)	N/C
Long Mounting Bar - 6-17.5" (8910821-3)	N/C

Width Extensions

Available only with Multi-Adjustable Foot Supports

Width Extension - Multi Adjustable (8910855)	\$82
--	------

Individual Multi-Adjustable Foot Supports

Mounting hardware included. Do not quote mounting bar. Not available with Foot Box.

Individual Multi-Adjust Foot Supports, Small - 6" - 10.5" (8910852-1)	\$495
Individual Multi-Adjust Foot Supports, Medium - 6" - 14.5" (8910852-2)	\$495
Individual Multi-Adjust Foot Supports, Large - 6" - 17.5" (8910852-3)	\$495

Footrest Accessories

Neoprene Sandals - Small (85696-1)	\$165
Neoprene Sandals - Medium (85696-2)	\$165
Formed Shoeholder - Small (86147)	\$195
Formed Shoeholder - Medium (86148)	\$195
Footstraps w/ fix lock (85096) (E0951)	\$43
Ankle Straps, X-Small - C: 6.25"-9" (85745-XS)	\$197
Ankle Straps, Medium - C: 9"-12.25" (85745-M)	\$197
Ankle Straps, X-Large - C: 11.5"-15" (85745-XL)	\$197

Foot Box Kits

Select desired size for Foot Box and corresponding Foot Plate Padding. Strap for Foot Box is required with all sizes. Not available with Multi Adjustable Footrests.

Foot Box Kit

Foot Box Kit, Size 1 Plate (95145-1BOX)	\$240
Foot Box Kit, Size 2 Plate (95145-2BOX)	\$255
Foot Box Kit, Size 3 Plate (95145-3BOX)	\$275

Strap for Foot Box

Strap for Foot Box (95145-20)	\$98
-------------------------------	------

Foot Plate Padding

Foot Plate Padding, Size 1 (95145-1PAD)	\$120
Foot Plate Padding, Size 2 (95145-2PAD)	\$165
Foot Plate Padding, Size 3 (95145-3PAD)	\$205

Padding for Foot Plate Rod

Padding for Foot Plate Rod (95145-150)	\$64
--	------



Calf Support for Centrally Mounted Foot Plate (85131) \$245
Attaches to Mounting Bar

Chilla Mobility Base

Chilla Frame

Chilla Folding Frame

Shape Adapter required

Chilla, Size 2 Frame Solid Adjustable Suspension (881121110) E1234	\$2630
--	--------

Shape Adapter

x:panda Shape Adapter Plate (881091) K0108	\$875
--	-------

Chilla Frame Accessories

Chilla Transit Fitting for Shape Adapter (881089)	\$224	O2 Holder (881096)	\$175
Under Storage Basket (881092)	N/C	IV Pole (881097)	\$475
Vent Tray (881095)	\$330		



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

High-Low:xo Base Options

High-Low:xo

Manual Lifts

Manual Lift - Indoor Size 1 (K0108)	\$1976
Manual Lift - Indoor, Size 2 (K0108)	\$2121

Outdoor High-low:xo

Outdoor - Manual Lift, Size 2 (89120152)	\$3329
Outdoor- power lift, Size 2	\$4266

Power Lifts

Indoor- Power Lift, Size 1 (K0108)	\$3175
Indoor- Power Lift, Indoor, Size 2 (K0108)	\$3409

High-Low:xo Options

Transit - Outdoor High-low

Transit (890005)	\$260
------------------	-------

Anti-Tips on the High-low

Anti-tips - Indoor (890007-1)	\$223
Anti-tips - Outdoor (890008)	\$240

Push Brace for High-low

Required selection

Fixed Angle Push Brace (890009-04)	\$198
Angle Adjustable Push Brace (8910555)	\$419

Multi Frame Mobility Base and Options

Adaptor plate for x:panda/Shape required

Multi Frame Base Options

Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Standard Wheel Locks	\$2995
Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Drum Brake	\$3595

Anti- Tip Option

Anti-tips	\$165
-----------	-------

Transit Option

Transit for Multi Frame	\$240
-------------------------	-------

Multi Frame Accessories

Adapter for x:panda on Multi Frame - for x:panda sz 1 & 2	\$210	Vent Tray	\$525
O2 Holder	\$395		

Ki Mobility Bases and Options

Ki Mobility Bases

Separate Order form for base must be filled out. Order forms are located at Kimobility.com. Footrests are not provided with base. If footrests are desired, please make selection at Footrests for Mobility Bases section above.

Tyke x:panda

Prices will be depicted on the Kimobility.com order form. Must be configured separately.

Tyke for x:panda	N/C
Tyke size 1 is only Available for Infant Package	
Tyke size 2 is only available for Shape Simple/Advanced size 1	

Little Wave Flip for x:panda

Prices will be depicted on the Kimobility.com order form. Must be configured separately.

Flip for x:panda	N/C
------------------	-----