



by Etac
Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

x:panda shape - Simple Seat

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 8/24/2026

x:panda Shape - Simple Seat

	Infant	1	2
 Effective seat width	8.5" 220 mm	8.5"/10.25"/11.75" 220/260/300 mm	11"/12.5"/14.25"/15.75" 280/320/360/400 mm
 Effective seat depth	4.75-6.25" 120-160 mm	6.25-12.5" 160-320 mm	9.5-16.5" 240-420 mm
 Backrest height	10.75-14.25" 270-360 mm	11-14.5"/12.5-16.5"/14.25-18.5" 280-370/320-420/360-470mm	13.5-17.25"/15-19.25"/17-22" 340-440/380-490/430-560 mm
 Max user mass	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
 Max user mass in transportation	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
 Backrest angle	-10° - 30°	-10° - 30°	-10° - 30°
 Back width, top of backrest	7.5" 190 mm	7"/8.5"/10.25"/11.75"/13.5" 180/220/260/300/340 mm	7"/8.5"/10.25"/11.75"/13.5"/15" 180/220/260/300/340/380 mm
 Overall width	9.5" 240 mm	11.75"/15.75"/19.75" 300/400/500 mm	13.75"/18"/22"/26" 350/460/560/660 mm
 Mass device	13.2 lb 6 kg	14.3 lb 6.5 kg	24.3 lb 11 kg

Measurements taken with both legs in max abduction. Size 1-3 measurements are based on Advanced Seat.

Shape Base Options

Please see measurement guide for growth information. A selection must be made in each section below.

Shape Base Options

Selection required

Seat Base, Size 1

Size 1 Simple Seat (E2292) \$611

Seat Base 2

Size 2 Simple Seat (E2292) \$611

Seat Plate

Selection required

Size 1 Seat Plate W: 10 LB: 8.75 N/C
Size 1, Seat Plate W: 11.75" LB 11.75" N/C

Seat Plate

Required Selection

Size 2, Seat Plate W: 12.5" LB: 8.75 N/C
Size 2, Seat Plate W: 15.75" LB: 11.75 N/C

Recline Options

Back Recline

Selection Required

Adjustable Stepless Recline, Toolless (K1014) \$443
Adjustable Stepless Recline with Dynamic, Lock of Dynamic, Toolless (E2398) \$723

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0627.6 Pricing and specifications subject to change without notice.



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

x:panda shape - Simple Seat

Cushions and Covers

Seat Cushions

Seat Cushion Size - One Piece

One Piece Cushion, Size 1 - W: 9.75" x L: 19.75" (901101)	N/C
One Piece Cushion, Size 1 - W: 11.75" x L: 19.75" (901103)	N/C
One Piece Cushion, Size 2 - W: 12.50" x L: 25.5" (901105)	N/C
One Piece Cushion, Size 2 - W: 15.75" x L: 25.5" (901107)	N/C

Seat Cushion Size - Split

Required when Quick Release Pommel is selected

Split Cushion, Size 1 - W: 10.25" x L: 19.75" (901114)	N/C
Split Cushion, Size 1 - W: 11.75" x L: 19.75" (901115)	N/C
Split Cushion, Size 2 - W: 12.5" x L: 25.5" (901117)	N/C
Split Cushion, Size 2 - W: 15.75" x L: 25.5" (901119)	N/C

Seat Cover

Required selection for cover color

Red (-50)	N/C
Grey (-60)	N/C
Black Wipeable (-70)	\$21

Back and Cushions

Upper Back Size

Select upper back base then select Upper and Lower back cushion and cover. Includes quick release detachment (#9050374)



Size 1 Base

Size 1 Base Upper - H: 7" (905073) (E2291)	\$580
<i>Height range: 12.25" - 16.25"</i>	
Size 1 Base Upper - H: 8.75" (905074) (E2291)	\$580
<i>Height range: 13.75" - 18.50"</i>	

Size 2 Base

Size 2 Base Upper - H: 7" (905073) (E2291)	\$580
<i>Height range: 13.5" - 17.25"</i>	
Size 2 Base Upper - H: 8.75" (905074) (E2291)	\$580
<i>Height range: 15" - 19.25"</i>	
Size 2 Base Upper - H: 10.75" (905075) (E2291)	\$580
<i>Height range: 17" - 21.75"</i>	

Upper Back Cushion and Cover

Upper back cushion height must match Upper Back Base height

Upper Back Cushion Size

Upper Back Cushion-W: 7"xH: 5.5" (902101)	N/C
Upper Back Cushion-W: 7"xH: 7" (902103)	N/C
Upper Back Cushion-W: 8.75"xH: 7" (902104)	N/C
Upper Back Cushion-W: 10.25"xH: 7" (902105)	N/C
Upper Back Cushion-W: 8.75"xH: 8.75" (902106)	N/C
Upper Back Cushion-W: 10.25"xH: 8.75" (902107)	N/C
Upper Back Cushion-W: 11.75"xH: 8.75" (902108)	N/C
Upper Back Cushion-W: 13.25"xH: 8.75" (902109)	N/C
Upper Back Cushion-W: 10.25"xH: 10.75" (902110)	N/C
Upper Back Cushion-W: 11.75"xH: 10.75" (902111)	N/C
Upper Back Cushion-W: 13.75"xH: 10.75" (902112)	N/C
Upper Back Cushion-W: 15"xH: 10.75" (902113)	N/C

Upper Back Cushion Cover

Selected required

Red (-50)	N/C
Grey (-60)	N/C
Black Wipeable (-70)	\$21

Lower Back Cushion and Cover

The width measurements of the cushion are not actual measurements of the cushion. The number matches the lower back of the seating system. If the seating system is configured without hip and pelvis support, consider going up a size in width when choosing the lower back cushion.

Lower Back Cushion Size

Shape Size 1 Lower Back Cushion - W: 8.75" x H: 4.75" (902120)	N/C
Shape Size 1 Lower Back Cushion - W: 11.75" x H: 4.75" (902121)	N/C
Shape Size 1 Lower Back Cushion - W: 11" x H: 5.5" (902122)	N/C
Shape Size 2 Lower Back Cushion - W: 11" x H: 5.5" (902123)	N/C
Shape Size 2 Lower Back Cushion - W: 12.5" x H: 5.5" (902124)	N/C
Shape Size 2 Lower Back Cushion - W: 14.25" x H: 5.5" (902125)	N/C
Shape Size 2 Lower Back Cushion - W: 15.75" x H: 5.5" (902126)	N/C

Lower Back Cushion Cover

Selection required

Red (-50)	N/C
Grey (-60)	N/C
Black Wipeable (-70)	\$21



by Etac
Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

x:panda shape - Simple Seat

Back Adjustments



Daily Adjustment Kit

Allows for toolless adjustable of back height. Ideal for frequent adjustment.

Daily Adjustment Kit - Upper Back (905081) \$21

Spacer Blocks

Standard Spacer Block Holder without Cushion

Cushion widths must match Spacer Block Holder width

Spacer Block Holder, Size 1 - Back widths 7", 8.75", 10.25" (903100-1) \$183
Spacer Block Holder, Size 2 - Back widths 8.75", 10.25", 11.75", 13.5" (903100-2) \$183

Standard Spacer Block Cushion Cover

Black (-80) N/C
Black Wipeable (-70) N/C

Spacer Block Cushion

Spacer Block Cushion Size 1 - W: 7" x H: 2.25" (903011) N/C
Spacer Block Cushion Size 2 - W: 8.75" x H: 2.25" (903012) N/C
Spacer Block Cushion Size 3 - W: 10.25" x H: 2.25" (903013) N/C
Spacer Block Cushion Size 4 - W: 11.75" x H: 2.25" (903014) N/C
Spacer Block Cushion Size 5 - W: 13.25" x H: 2.25" (903015) N/C



Headrest and Hardware Options

Headrest and Hardware Options

Includes Cushion. Select desired headrest in chart and choose the corresponding required hardware below.

Style	Size	Price	Type
Curved	Small	\$220	Curved, Small - Black Wipeable - W: 7", H: 5" (E0955)
			Curved, Small - Black - W: 7", H: 5" (E0955)
	Medium	\$240	Curved, Medium - Black Wipeable W: 10.2", H: 6.3" (E0955)
			Curved, Medium - Black W: 10.2", H: 6.3" (E0955)
Anatomic	Small	\$271	Anatomic, Small - Black Wipeable A: 5.5" x B: 5.5" x C: 9.5" (E0955)
			Anatomic, Small - Black A: 5.5" x B: 5.5" x C: 9.5" (E0955)
	Medium	\$281	Anatomic, Medium - Black Wipeable A: 6" x B: 5.5" x C: 10" (E0955)
			Anatomic, Medium - Black A: 6" x B: 5.5" x C: 10" (E0955)
	Large	\$291	Anatomic, Large - Black Wipeable A: 6.5" x B: 5.5" C: 12.5" (E0955)
			Anatomic, Large - Black A: 6.5" x B: 5.5" C: 12.5" (E0955)
Adjustable		\$332	Adjustable w/ Large Laterals - Black A: 4.75" x B: 6 x C: max 9.5" (E0955)
		\$440	Adjustable with Large Laterals and Speakers - Black A: 4.75" x B: 6: max 9.5" (E0955)
		\$411	Adjustable with Four Laterals - Black A: 4.75" x B: Max 9.5" C: max 8" (E0955)

Mounting Hardware for Curved and Anatomic Headrest

Headrest Mounting Hardware

Mounting Brackets, Standard (E1033) \$245
Mounting Brackets, Extended (E1033) \$255
Mounting Brackets, 3D Adjustment (E1033) \$266

Mounting Hardware for Adjustable Headrest Style

Short Mounting Hardware (E1033) \$250
Long Mounting Hardware (E1033) \$250

Headrest Options



Daily Adjustment Kits

Required with Curved and Anatomic Head support

Daily Adjustment Kit Standard - Set of Two (905034) \$32
Daily Adjustment Kit Extended - Set of Three (905035) \$47

Additional Headrest Options

Headband System \$250
Universal Headrest Mount \$67

Arm Supports

Arm Rests

Arm Rests

Cover not included, select separately

Arm Rest, Fixed W: 3" L: 10.25" H: 4" - 6.5" (8910616-1) (E0973)	\$339
Arm Rests, Fixed W: 3" L: 12.25" H: 4" - 8.75" (8910617-1) (E0973)	\$339

Arm Rest Covers

Black (-80)	\$79
Black Wipeable (-70)	\$79

Daily Adjustment Kit

For toolless height adjustment of arm supports

Height Adjustment Kit of Arm Supports (8910615-11)	\$21
--	------

Upper Arm Supports - Fixed

Includes Black Wipeable cover. Not compatible with Swing-Away Lateral Supports.



Upper Arm Supports, Fixed, Size 1 - L: 5.5" H: 7" (905017-70) (K0108)	\$260
<i>For size 1 Base</i>	
Upper Arm Supports, Fixed Size 2 - L: 7.5" H: 10.75" (905018-70) (K0108)	\$260
<i>For Size 2 Base</i>	

Upper Arm Supports - Adjustable

Not compatible with Swing-Away Lateral Supports

Upper Arm Supports, Adjustable, Size 1 - L: 5.5" H: 7" (905019-70) (K0108)	\$349
<i>For Size 1 Base</i>	
Upper Arm Supports, Adjustable, Size 2 - L: 7.5" H: 10.75" (905020-70) (K0108)	\$349
<i>For size 2 Base</i>	

Knee and Postural Supports

Abduction and Adduction

Outer Knee Support



Outer Knee Support, Standard Set for Size 1, Black - L: 6", H: 3.25" (905005-70) (E0956)	\$219
Outer Knee Support, Standard Set for Size 1, Black Wipeable - L: 6", H: 3.25" (905005-80) (E0956)	\$219
Outer Knee Support, Standard Set for Size 2, Black - L: 7.75", H: 3.5" (905006-80) (E0956)	\$219
Outer Knee Support, Standard Set for Size 2, Black Wipeable - L: 7.75", H: 3.5" (905006-70) (E0956)	\$219

Quick Release Pommel

Quick Release Pommel For Size 1 Small, Black - W: 2.25", L: 6.25", H: 7.5" (905077-80) (E0957)	\$364
Quick Release Pommel For Size 1 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" (905077-70) (E0957)	\$364
Quick Release Pommel For Size 2 Small, Black - W: 2.25", L: 6.25", H: 7.5" (905082-80) (E0957)	\$364
Quick Release Pommel For Size 2 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" (905082-70) (E0957)	\$364
Quick Release Pommel For Size 2 Large, Black - W: 3.25", L: 7.5", H: 9" (905078-80) (E0957)	\$364
Quick Release Pommel For Size 2 Large, Black Wipeable - W: 3.25", L: 7.5", H: 9" (905078-70) (E0957)	\$364

Postural Supports

Width Adjustable Pelvis Supports

Mounted from the back



Width Adjustable Pelvis Supports, Small, Black - L: 4", H: 4.25" (905041-80) (E0956)	\$245
Width Adjustable Pelvis Supports, Small, Black Wipeable - L: 4", H: 4.25" (905041-70) (E0956)	\$245
Width Adjustable Pelvis Supports, Large, Black - L: 5.5", H: 5" (905042-80) (E0956)	\$245
Width Adjustable Pelvis Supports, Large, Black Wipeable - L: 5.5", H: 5" (905042-70) (E0956)	\$245

Width Adjustable Hip and Pelvis Supports

Mounted from the sides



Width Adjustable Hip Supports, Size 1, Black - L: 4" H: 4.25" (905044-80) (E0956)	\$349
Width Adjustable Hip Supports, Size 1, Black Wipeable - L: 4" H: 4.25" (905044-70) (E0956)	\$349
Width Adjustable Hip Supports, Size 2, Black - L: 5.5" H: 5" (905045-80) (E0956)	\$349
Width Adjustable Hip Supports, Size 2, Black Wipeable - L: 5.5", H: 5" (905045-70) (E0956)	\$349

Trays and Accessories

Trays and Activity Centers

Tray



Shape 5 Tray with Edge Size 1 (8910725-1) (E0950)	\$334
<i>For Size 1. Fits Size 2 widths 11", 12.5", 14.25". Double Rod, Grey ABS.</i>	
Shape 6 Tray with Edge Size 2 (8910725-2) (E0950)	\$334
<i>For Size 2. Double rod, Grey ABS.</i>	



Padding for Size 1 Tray, Shape 5 Plate, Over and Under (8910725-10) (E1028)	\$443
Padding for Size 1 Tray, Top, Shape 5 Plate, Over (8910725-10T) (E1028)	\$312
Padding for Size 2 Tray, Shape 6 Plate, Over and Under (8910725-20) (E1028)	\$443
Padding for Size 2 Tray, Top, Shape 6 Plate, Over (8910725-20T) (E1028)	\$312



by Etac

Phone (877) 227-0097 | Fax (704) 882-0751

sales.us@etac.com

x:panda shape - Simple Seat

Canopy and Mounting Hardware

Canopy does not come with brackets. Both the canopy and bracket must be quoted when mounting to the shape.



Canopy Without Brackets

Black Canopy - Size 1 (8740110)	\$339
Black Canopy - Size 2 (8740120)	\$339



For Shape 5 and 6

Additional Tray Accessories

Toy Bar (8910740)	\$250
<i>For shape 5 and 6</i>	

Arm Pads for Tray (8910730)	\$315
-----------------------------	-------



Canopy Brackets

Canopy Brackets with Push Handle for Size 1 (905039-1)	\$203
Canopy Brackets with Push Handle for Size 2 (905039-2)	\$203

Accessory Holders and Brackets

Remote Control Holder Size 1 - Arm Support Mounted (905030-1)	\$68
Remote Control Holder Size 2 - Arm Support Mounted (905030-2)	\$68
Remote Control Holder - Back Mounted (905030)	\$131
Bracket for Communication Aid - Arm Support Mounted (905033-1)	\$82

Vests and Belts

More options can be found at etac.com

Vests

Torso Vests



Torso Vest Small - S-P: 10.75" - 12.5" (85705-S) (E0960)	\$232
Torso Vest Medium - S-P: 12.5" - 16" (85705-M) (E0960)	\$232
Torso Vest Large - S-P: 16" - 18.5" (85705-L) (E0960)	\$232
Torso Vest X-Large - S-P: 18.5" - 20.5" (85705-XL) (E0960)	\$232

Cross Vests



Cross Vest X-Small - S-P: 7.75" - 12.5" (85720-XS) (E0960)	\$220
Cross Vest Small - S-P 10.25" - 13.5" (85720-S) (E0960)	\$220
Cross Vest Medium - S-P: 13.5" - 17.75" (85720-M) (E0960)	\$220
Cross Vest Large - S-P: 17.75" - 20.75" (85720-L) (E0960)	\$220
Cross Vest X-Large - S-P: 20.75" - 24.5" (85720-XL) (E0960)	\$220

Fixlocks Bracket



Higher Brackets (8910260)	\$77
<i>Bracket for higher strap guidance. Includes set of fixlocks. Requires the wider brackets as well.</i>	



Wider Brackets (905068)	\$57
-------------------------	------

Hip Belts

4-Point Belts

4-Point Padded Hip Belts 1" X-Small - H-H: 6.75" - 12.5" (85725-XS) (E0960)	\$156
4-Point Padded Hip Belts 1" Small - H-H: 8.25" - 13.75" (85725-S) (E0960)	\$156
4-Point Padded Hip Belts 1" Medium - H-H: 9.75" - 16.25" (85725-M) (E0960)	\$156
4-Point Padded Hip Belts 1" Large - H-H: 12.5" - 17.25" (85725-L) (E0960)	\$156

Pelvic Harness

Pelvic Harness, Small - A:6.5" x B:7.5" (89108-0)	\$261
Pelvic Harness, Medium - A:8.5" x B:8.75" (89108-1)	\$261
Pelvic Harness, Large - A:10.5" x B:9.75" (89108-2)	\$261
Pelvic Harness, XL, A: 13.5" x B: 12.5" (89108-3)	\$284

Pelvis Belt

Pelvis Belt, X-Small - H-H: 6.25" - 7.75" (85735-XS) (E0978)	\$175
Pelvis Belt, Small - H-H: 7.75" - 9.5" (85735-S) (E0978)	\$175
Pelvis Belt, Medium - H-H: 9.5" - 11" (85735-M) (E0978)	\$175

2-Point Pelvic Belt

2-Point Pelvic Belt Black - L: 13-14.5" (85424-1) (E0978)	\$98
---	------

Lateral Support Cushion and Fittings

if choosing Laterals, both the left and right cushions must be selected

Lateral Supports

Left Cushion

Left Rectangular Small - L: 4.25" x H: 3.25" (E0956)	\$115
Left Rectangular Medium - L: 6" x H: 3.5" (E0956)	\$115
Left Rectangular Large - L: 7.5" x H: 3.5" (E0956)	\$115

Right Lateral Supports

Right Cushion

Right Rectangular Small - L: 4.25" x H: 3.25" (E0956)	\$115
Right Rectangular Medium - L: 6" x H: 3.5" (E0956)	\$115
Right Rectangular Large - L: 7.5" x H: 3.5" (E0956)	\$115

Fitting for Lateral Supports

Selection required, fittings are sold as a pair

Type of Fitting

Fixed Fitting (E1028)	\$172
Swing-Away and Angle Adjustable Fitting (E1028)	\$260

Footrests for Bases

Footrests attach to each Base. If multiple bases are requested, add additional footrest for each. Choose appropriate base(s) below.

Footrest

Select footrest type and desired mounting bar. Footrests attach to each Base. If multiple bases are requested, add additional footrest for each. Choose appropriate base(s) below.

Footrest - Angle and Height Adjustable

Angle and Height Adjustable Foot Plate, Small, 12" x 7.5" (95145-1) (K0040)	\$325
Angle and Height Adjustable Foot Plate, Medium, 14" x 9" (95145-2) (K0040)	\$325
Angle and Height Adjustable Foot Plate, Large, 16" x 9.5" (95145-3) (K0040)	\$325

Footrest Mounting Bar

Required for Angle and Height Adjustable Foot Plates

Short Footrest Mounting Bar - 6-10.5"	N/C
Medium Footrest Mounting Bar - 6-14.5"	N/C
Long Footrest Mounting Bar - 6-17.5"	N/C

Width Extensions

Available only with Multi-Adjustable Foot Supports

Width Extension - Multi Adjustable (8910855)	\$82
--	------

Individual Multi-Adjustable Foot Supports

Mounting hardware comes standard. Not compatible with the Footbox.

Individual Multi-Adjust Foot Supports, Small - 6" - 10.5" (8910852-1)	\$495
Individual Multi-Adjust Foot Supports, Medium - 6" - 14.5" (8910852-2)	\$495
Individual Multi-Adjust Foot Supports, Large - 6" - 17.5" (8910852-3)	\$495

Leg and foot supports

Select desired size for Foot Box and corresponding Foot Plate Padding. Strap for Foot Box is required with all sizes. Not available with Multi-Adjustable Footrest.

Foot Box Kit

Required selection

Foot Box Kit, Size 1 Plate (95145-1BOX)	\$240
Foot Box Kit, Size 2 Plate (95145-2BOX)	\$255
Foot Box Kit, Size 3 Plate (95145-3BOX)	\$275

Strap for Foot Box

Required with all sizes

Strap for Foot Box (95145-21)	\$98
-------------------------------	------

Foot Plate Padding

Foot Plate Padding, Size 1 (95145-1PAD)	\$120
Foot Plate Padding, Size 2 (95145-2PAD)	\$165
Foot Plate Padding, Size 3 (95145-3PAD)	\$205

Padding for Foot Plate Rod

Required selection

Padding for Foot Plate Rod (95145-150)	\$64
--	------



Calf Support for Centrally Mounted Foot Plate (85131)
Attaches to Mounting Bar \$245



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

x:panda shape - Simple Seat

Chilla Base

Chilla Frame	Chilla Frame Accessories
<i>Adaptor for x:panda Shape required</i>	
Chilla Folding Frame	
Chilla, Size 2 Frame Solid Adjustable Suspension (881121110) E1234	\$224
	Under Storage Basket (881092) N/C
	Vent Tray (881095) \$330
	O2 Holder (881096) \$175
	IV Pole (881097) \$475
Shape Adapter	
<i>Required to mount to Chilla Base</i>	
x:panda Shape Adapter Plate (881091) K0108	\$875

High-Low:xo Base and Options

High-Low:xo Base	
Manual Lifts	Power Lifts
<i>Manual High-Low bases require a Push Brace selection</i>	
High-low:xo Frame Manual Lift, Indoor, Size 1	\$1976
High-low:xo Frame Manual Lift, Indoor, Size 2	\$2121
Outdoor High-Low:xo	
High-low:xo Frame Outdoor, Size 2, Manual Lift	\$3329
High-low:xo Frame Outdoor, Size 2, Power Lift	\$4266
High-Low:xo Options	
Transit - Outdoor High-Low:xo	Anti-Tips on the High-Low:xo
Transit (890005) (E1022)	\$260
	Anti-Tips - Indoor Base (890007-1) \$223
	Anti-Tips - Outdoor Base (890008) \$240
Push Brace	
<i>Available only with Manual Indoor Lifts</i>	
Fixed Angle Push Brace (890009-04)	\$198
Angle Adjustable Push Brace (8910555)	\$419

Multi Frame Mobility Base and Options

<i>Adaptor plate for x:panda Shape required</i>	
Multi Frame	Base options
Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Standard Wheel Locks	\$2995
Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Drum Brake	\$3595
Adaptor for x:panda on Multi Frame - for x:panda Shape 1 & 2	\$210
	Anti-tips \$165
	Transit for Multi Frame \$240
	Vent Tray \$525
	O2 Holder \$395



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Ki Mobility Bases

Separate order form for base must be filled out. Order form located at Kimobility.com. Footrests are not provide with base. If footrests are desired, make a selection in the Footrest for Base, section.

Ki Mobility Bases

Price for base calculated with Ki Mobility's order form. please refer to Tyke order form to verify compatibility.

Flip for x:panda	N/C
Shape size 2 only for Flip	
Tyke for x:panda	N/C
Tyke size 1 is only available for Infant Package, Tyke size 2 is only available for Shape Simple/Advanced size 1	